



Brand	Compatible Models	P/N	Description	Male Connector	Female Connector
Masimo		U719-2AL	Spo2 extension cable L=2.4M, LNCS-DB9F require Masimo LNCS sensor		
		U719-2AS	Spo2 extension cable L=10cm, LNCS-DB9F require Masimo LNCS sensor		
		U719-2CL	Spo2 extension cable L=2.4M, LNCS-LNOP require Masimo LNOP sensor		
		U719-2CS	Spo2 extension cable L=10cm, LNCS-6J require Masimo LNOP sensor		
		U719-4AL	Spo2 extension cable L=2.4M, LNOP8pin- DB9F require Masimo LNCS sensor		
		U719-4CL	Spo2 extension cable L=2.4M, LNOP8pin-6J require Masimo LNOP sensor		
		U719-5JL	Spo2 extension cable L=2.4M, DB15pin- DB15F require Masimo sensor		
	All Masimo monitors; MDE: Escort II OPT30 Prism; GE Marquette: Dinamap ProCare with Masimo Pro1000, ProSeries V2, V300; Invivo T8; Schiller: Zoll; Datascope: Passport Trio, Passport 2, Passport XG and Spectrum; Bitmos: sat 800, sat 801, sat 805, sat 816, satPuls, Philips Respironics Alice 5 Diagnostic Sleep System, Medtronic Physio Control pak12, 20.	U719-6AL	Spo2 extension cable L=2.4M, 14pin-DB9F require Masimo LNCS sensor		
		U719-6CL	Spo2 extension cable L=2.4M, 14pin-6J require Masimo LNOP sensor		
MEK	MP100, MP110 MP500, MP500, MP600, MP1000	U713-1AL	Spo2 extension cable L=2.4M, DB6pin-DB9F require MEK sensor		