









Bionet	B-SPCBL-N; B-SPCBL- NP	ADP-10027A	BM3 plus, use with Bionet short spo2 sensor, such as, RSP-10066		
	B-SPCBL-O	ADP-10027B	BM3 + (old style connector) use with Bionet short spo2 sensor, such as, RSP-10066P		lemo 7 p-DB9F
Hellige	30344358	ADP-10028B	Servomed 4. Generation SCP 844, 851, 852; SCS 804; SMK151, 153, 154, 155-1, 156, 157, SMK181, 182, SMV 104T, 107, 111, 168; SMS316; SMC108; SMV104, Servomed N4, Servocard, Serie Cardioserv, Vicom SM, SMU511, 611, 612, 630 Serie SMK-200, SMK210, 211, 230. 231. 240, 250, SMP 300, 310, Pepita 201, 202, 202R, 203, Pepita 203, R204, 204R, 205, Pepita 205R, use with Nellcor spo2 sensor		
Primedica		ADP-10029A	use with Nellcor spo2 sensor, or RSP-10041 etc.		8 pins dual notch
		ADP-10029B	use with Nellcor spo2 sensor, or RSP-10041 etc.		7 pins single notch
Edan	M13-36091	ADP-10030A	EDAN: M9, M9A, M8, M8A, M8B, M3, M3B, WANDONG: BYS-8000, MINDRAY: PM-8000 Use with RSP-10001A spo2 sensor, BCI sensors, etc		6 pins Lemo connector to DB9 female
Edan		ADP-10030B	eDAN oXIMAX		8 pins Lemo connector to DB9 female
Choice-med		ADP-10032	Choice:MMED6000D/D; CSM-2000; Goldway: UT4000A, UT4000C; KP-900; use with BCI spo2 sensor; or RSP-10001A		5 pins Redel connector to DB9 female connector;
Goldway UT4000		ADP-10033	Compatible with: UT4000A, UT4000C PACE TECH; use with nellcor spo2 sensor		5 pins Redel connector to DB9 female