BIOHIT *Helicobacter pylori* UFT300 quick test

Ultra-fast H. pylori detection from a biopsy

The Biohit *H. pylori* UFT300 is a true quick test for the detection of *H. pylori* from a biopsy specimen. The biopsy taken during gastroscopy can be tested immediately to diagnose *H. pylori* infection or to determine the success of eradication therapy. The test results are ready in only 5 minutes enabling diagnosis and reporting at the same time. This saves the patient from an unnecessary visit to the doctor for hearing the test results. The Biohit *H. pylori* UFT300 quick test has excellent sensitivity and specificity which makes it a highly reliable and accurate tool for diagnostics.

One-step test procedure

The Biohit *H. pylori* UFT300 quick test is easy and effortless to use. The biopsy specimen is placed into the test tube and mixed with the test reagent. Clear colour change indicates the presence of in the specimen within 5 minutes. The closed test system of the UFT300 quick test makes the test procedure safe and hygienic for the user, and the interpretation of the results doesn't require specialist training.





Helicobacter pylori UFT300 quick test

- Ready-to-use test kit
- Results ready in 5 min (both positive and negative)
- Storage in room temperature
- Sensitivity 94,5%, Specificity 100%
- Testing and reporting during one appointment

Helicobacter pylori Quick test		
Cat. No.	Product	Qty
602019	UFT300 quick test kit	50 tests
602021	UFT300 quick test kit	100 tests
602017	Positive control	1

In the United States and Japan for research use only.