

Vliwasoft® Nonwoven Wound Dressing, sterile

REF 12076, 12077, 12082, 12083, 12065, 12066, 12071,
 12079, 12085, 12068, 12069, 12073, 12080, 12086
 98075, 98076, 80949, 80960, 80966, 80950, 80958, 80965

January 2020

1. Composition of the product

Vliwasoft® Nonwoven Wound Dressing, sterile consist of:

- viscose
- polyester

This product data sheet covers the following items:

Vliwasoft® Non-woven Wound Dressing, sterile

REF		Size	Packaging / SC
12076	4-ply	5 x 5 cm	2-fold sealed
12077			5-fold sealed
12082		7,5 x 7,5 cm	2-fold sealed
12083			5-fold sealed
12065		10 x 10 cm	2-fold sealed
12066			5-fold sealed
12071		10 x 20 cm	2-fold sealed
12079	6-ply	5 x 5 cm	2-fold sealed
12085		7,5 x 7,5 cm	2-fold sealed
12068		10 x 10 cm	2-fold sealed
12069			5-fold sealed
12073		10 x 20 cm	2-fold sealed
12080*		5 x 5 cm	5-fold sealed
12086*		7,5 x 7,5 cm	5-fold sealed
98075*	4-ply	7,5 x 7,5 cm	á 10 Pcs.
98076*		10 x 10 cm	á 10 Pcs.
80949*		7,5 x 7,5 cm	10 x 2 Pcs.
80960*			25 x 2 Pcs.
80966*			50 x 2 Pcs.
80950*		10 x 10 cm	10 x 2 Pcs.
80958*			25 x 2 Pcs.
80965*			50 x 2 Pcs.

* These devices are regional items. They are – with the exception of the packaging size and printing – identically equal to the international products.

2. Packaging, structure and composition

2.1 Unit container / UC

1 four-edge sealed pouch
made of cellulose, polypropylene, polyethylene

REF	pcs.
12065, 12068, 12076, 12079, 12082, 12085, 80949, 80950, 80958, 80960, 80965, 80966	2

1 deep draw package
consisting of cellulose, polyamide und polypropylene

REF	pcs.
12066, 12069, 12077, 12080, 12083, 12086,	5
98075	10

2.2 Shelf container / SC

1 folding box, made of cellulose

REF	pcs.
80949, 80950	20
80958, 80960	50
12068, 12068, 12073, 12079, 12080, 12085, 12086, 80965, 80966	100
12065, 12071, 12076, 12082	150
12066, 12077, 12083, 98075	250

2.3 Transit container /TC

1 corrugated card box, made of cellulose

REF	Stück
80949, 80950, 80958, 80960, 80965, 80966	800
12073, 12073	1000
12071	1500
12068, 12069, 12079, 12080, 12085, 12086	2000
12066	2500
12065, 12076, 12082	3000
12083	4500
12077	6750

3. **Manufacture**

Vliwasoft® Nonwoven Wound Dressing, sterile is produced according to specification in hygienic conditions and packed as described in their relevant packaging specification.

4. **Description**

Vliwasoft® Nonwoven Wound Dressing, sterile is a white, soft, absorbent universal compress.

5. **Properties**

Vliwasoft® Nonwoven Wound Dressing, sterile is made from fine-pored, non-fraying nonwoven material, for a gentle wound cleansing and has good absorption properties. The non-fraying nonwoven material is highly absorbent and is appreciated for its low-adherent properties.

Indications

- for wounds with low to moderate amounts of exudate
- as a wound dressing
- for wound cleansing
- as a carrier material for ointments
- as a padding material between fingers and toes

6. **Intended purpose**

Vliwasoft® Nonwoven Wound Dressing is intended to be used for cleansing and for wound care.

7. **Medical device classification**

Vliwasoft® Nonwoven Wound Dressing, sterile is a medical device of class Is according to rule 4. (Council Directive 93/42/CEE and Regulation 2017/745)

8. **Biological evaluation and biocompatibility (DIN EN ISO 10993)**

The starting materials used in the manufacture of Vliwasoft® Nonwoven Wound Dressing, sterile are safe if the product is used appropriately.

The purpose of this documentation and the statements made therein is to show that there is no risk involved in the use of the medical device Vliwasoft® Nonwoven Wound Dressing and that they are designed, manufactured and packaged in such a way that they will not compromise the clinical condition or the safety of patients, or the safety and health of users and other persons when used under the conditions and for the purposes intended.

9. Requirements regarding the safety/harmlessness of the constituents

There is currently no information available on critical constituents.

10. Stability

Stored appropriately, dry, clean and protected from sunlight – the Vliwasoft® Nonwoven Wound Dressing, sterile have a shelf life of 5 years from manufacturing date.

11. Disposal

The user is advised to observe current legislation, norms and guidelines, regulating the disposal of medical refuse.

Packaging materials must also be disposed of in compliance with applicable national requirements.

Lohmann & Rauscher International GmbH & Co. KG
D-56579 Rengsdorf
signed by

Gerd Pahlisch
Director Global Portfolio Management