

# PRE-TREATMENT

## WIP'ANIOS EXCEL

- » 100% Biodegradable viscose wipes.
- » Impregnating solution formulated without fragrance nor alcohol.
- » Broad virucidal spectrum.



### WALL HOLDER



MORE DETAILS IN THE DEDICATED BROCHURE

#### INDICATIONS

Surface treatment (e.g. external structures of equipment: respirators, dialysis generators, etc.), non-immersible and non-invasive\* medical devices, and wiping of endoscopy sheaths in preliminary treatment, prior to immersion for cleaning and disinfection. \*Examples: conventional ultrasound probes, pressure sensors, glucose testers, tensiometres reflex hammers, etc.

#### COMPOSITION

Wipes impregnated with a cleaning and disinfecting solution: Didecylmethylammonium chloride, non-ionic surfactants, sequestering agent, excipients.

#### PACKAGING - REF.

Wip'Anios Excel (180 x 200 mm)	
6 packs, 100 wipes each.....	2446.655
12 packs, 50 wipes each.....	2446.424
Wall holder .....	425.185

#### INSTRUCTIONS FOR USE



SINGLE USE  
WIPES

Active from 2 min.

#### THE IMPREGNATION SOLUTION IS:

ACTIVE AGAINST	ACCORDING TO STANDARDS
BACTERIA	EN 13727, EN 1276, EN 13697, EN 16615 Active against MRB and XRB according to EN 13727
YEASTS	EN 13624, EN 16615, EN 13697
MOULDS	EN 13624, EN 13697
VIRUSES	Adenovirus, Norovirus MNV according to EN 14476 HIV, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus, HSV, Rotavirus, VRS, Polyomavirus SV 40, Coronavirus

