

Technical Data

Mineral Modified Glutamate Agar Base (Twin Pack)

M643I

Mineral Modified Glutamate Agar Base is used for the enumeration of *Escherichia coli* from meat and meat products. The composition and performance criteria of this medium are as per the specifications laid down in ISO 1996, Draft ISO 1988, Draft ISO/DIS 6391.

Composition**

| Ingredients | Gms / Litre |
|---|---------------|
| Part A | - |
| Lactose | 10.000 |
| Dipotassium phosphate | 0.900 |
| Sodium formate | 0.250 |
| L-Cystine | 0.020 |
| L-Aspartic acid | 0.024 |
| L-Arginine | 0.020 |
| Thiamine | 0.001 |
| Nicotinic acid | 0.001 |
| Pantothenic acid | 0.001 |
| Ferric ammonium citrate | 0.010 |
| Calcium chloride, 2H2O | 0.010 |
| Magnesium sulphate, 7H2O | 0.100 |
| Agar | 15.000 |
| Part B | - |
| Sodium glutamate | 6.350 |
| Final pH (at 25°C) | 6.7 ± 0.2 |
| **Formula adjusted, standardized to suit performance parameters | |

Directions

Suspend 26.29 grams of dehydrated Part A and 6.35 grams of Part B in 1000 ml distilled water containing 2.5 grams ammonium chloride. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 10 lbs (115°C) for 10 minutes. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Folpmers (1) described a glutamic acid based chemically defined medium for the enumeration of coliform bacteria from water and wastewater. This glutamate-containing medium was later modified by Gray (2), by the addition of lactose, which gave less false positive results when compared to MacConkey Broth.

Mineral Modified Glutamate Agar Base is recommended for the enumeration of *Escherichia coli* from meat and meat products by colony count at 44°C by membrane filtration technique. Anderson and Baird-Parker (3) described a direct plate method for the rapid enumeration of *E. coli* in foods and this method was modified by a resuscitation procedure using Mineral Modified Glutamate Agar with successful recovery of damaged cells from frozen, dried, heat processed or low pH foods.

This media contain a variety of nutrients including salts, amino acids and vitamins. Lactose is the fermentable carbohydrate. Because of the nutrients, this media is superior for enumerating coliforms in water and wastewater as it satisfies most of the nutritional requirements of coliforms.

Quality Control

Appearance

Part A : Off-white to light yellow homogeneous free flowing powder Part B : Colourless to White needles

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Yellow coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of the medium (2.634% w/v Part A + 0.635% w/v Part B + 0.25% ammonium chloride) is at 25°C. pH : 6.7 ± 0.2 pH

6.50-6.90

Cultural Response

Cultural characteristics observed after an incubation at 44°C for 24 hours.

Cultural Response

| Organism | Inoculum (CFU) | Growth | Recovery |
|-----------------------|-------------------|-----------|----------|
| Cultural Response | | | |
| Escherichia coli ATCC | 50-100 | luxuriant | >=50% |
| 25922 | | | |

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1.Folpmers T., 1948, Ant. V. Leeuwenhoek, J. Microbiol. Serol., 14:58.

2.Gray R.D., 1959, J. Hyg. Camb., 57:249.

3. Anderson J. M., and Baird Parker A. C., 1975, J. Appl. Bacteriol., 39:111.

4. International Organization for Standardization (ISO), 1988, Draft, ISO/DIS 6391.

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com