

Product Catalogue

Vision

Serve and care for more healthcare professionals and patients.

Mission

Become the most preferred company by providing high quality medical technologies at affordable prices to meet the needs of patients, healthcare professionals and providers in the global markets whilst serving and meeting our Corporate Social Responsibilities.

SERVING AND CARING FOR MORE PATIENTS

Advanced MedTech Solutions (AMS) is dedicated to advancing patient care by providing high quality medical device technologies to serve the needs of more patients and healthcare professionals.

Global healthcare needs are increasing and the collective capacity to answer them is challenged by meeting the demands of an aging population, the growing burden of chronic disease, and complex healthcare policies, it's incredibly clear the pace of change has never been greater.

The desire to make a difference in patients' lives led to the genesis of AMS.

EMPOWERING HEALTHCARE PROFESSIONALS

Headquartered in Vadodara, India, the company was established by an elite team of respected industry leaders and visionaries who understand the challenges of the current healthcare environment.

Redefining and raising the bar on standards of quality, care and trust. Our unwavering mission is to improve patient care through meaningful solutions.

LIFE ENHANCING INNOVATION

Our state-of-the-art facility is located in the industrial center of Vadodara, Gujarat. This eco-friendly manufacturing plant is singular in design and distinguished in efficiency. Designed by renowned Spanish architect Ricardo Bofill, it is equipped to manufacture a cutting-edge class of medical devices that meet the most stringent global standards.

It serves as a technological and healthcare ecosystem focused on the medical science fields of Advanced Surgical & Interventional Cardiology Solutions.

ADVANCING THE SPECTRUM OF HEALTHCARE

With a dedication to excellence, we are committed to optimizing patient outcomes through quality, affordability, and accessibility – the heart of the AMS mission.

Reaching out to healthcare professionals and providers around the world, Advanced MedTech Solutions offers a compelling choice with a corporate and social responsibility founded on advancing the way we care for generations to come.



State of the Art Excellence

Commitment to People



Commitment to Quality




Commitment to Professional Development



Commitment to Innovation





Touching
2 Lives
Every Second

Geographical Footprint

4

Corporate

1

Manufacturing

5

Therapies

161

Countries



Corporate HQ & Manufacturing

P-21-22,25-26 & 34-35, GIDC Manjusar,
Tal-Savli, Vadodara - 391775 Gujarat, India
T: +91 2667671302



Holding Company

Somerset Road, #07-08 Triple One Somerset,
Singapore - 238164
T:(65) 636910 80/89



Corporate Office

305, Ceejay House, F-Block, Shiv Sagar Estate,
Dr. Annie Besant Road, Worli, Mumbai - 400018,
India | T:+91 22 24946100/01



Regional Office


Sm0132, 1st floor, DLF Star Mall, Sector 30,
NH-48, Opp: 32nd Avenue, Gurugram,
Haryana 122001, India

A New Generation of
NEEDLES


Intelligent Geometry

High Precision

A New Generation of Needle Technology




Needle Anatomy




- a : Needle chord length
- b : Needle point
- c : Swage
- d : Needle radius
- e : Needle length
- f : Needle body

Needle Curvature





Minimal Tissue Trauma Enables Faster Recovery



* Certain features are common in specific needles

Needle Point



Straight



Straight needles can be used for suturing without requiring the use of a needle holder, as long as the suturing area offers good visibility

Blunt Point



Needle point is sharp enough to penetrate fascia and muscle but not skin, dissects the friable tissue rather than cutting it

Conventional



It has a triangular cross section with the apex of the triangle on the inside of the needle curvature, very sharp and cuts rather than dilate

Cutting

Reverse



The body of this needle is triangular in cross section, having the apex cutting edge on the outside of the needle curvature


Cutting

J Shape



These needle designs were developed to facilitate controlled depth of penetration such as vesicourethral repair, in ophthalmological applications (cornea), etc.

Round Body



Round Bodied needles are designed to separate tissue fibres rather than cutting them

Taper



This needle combines the initial penetration of a cutting needle with the minimised trauma of a round bodied needle

Cutting

Needle Point



Trocar Point

The needle is of heavy diameter and design of the cutting head ensures powerful penetration, even when deep in dense tissue



Black

These needles have been designed to give outstanding visibility against tissue and where blood is present in the operating field



Taper

Point


Plus

A modified point profile, the tapered cross section immediately behind the tip has been flattened to an oval shape rather than a conventional round shape



Spatula

Spatula Needle has a square body to produce superb penetration characteristics



Prime

The cross sectional geometry of the needle tip reduces the angles of the cutting edges which gives improved penetration and control



Advanced


Coronary

Calcified

(ACC)

The unique point design of the CC Needle provides significantly improved penetration properties for the Cardiac/ Vascular

ADVATRAY™ - the confidence to suture with ease



ADVATRAY™ an advanced dispensing system designed to deliver ease of suture dispensing with minimal memory in the suture strand

Impeccable suturing experience for surgeons | Eliminates the need of suture lubrication



 Advanced Wound Care
where every suture is a seal of confidence

Absorption Profile

ADVACRYL™ Braided Coated Polyglactin 910 (Undyed & Violet)



ADVACRYL PLUS™ Antibacterial with Triclosan Braided Coated Polyglactin 910 (Undyed & Violet)




ADVACAT™ Purified Collagen (Chromic & Plain)



ADVACRYL RAPID™ Fast Absorbing Polyglactin-910 (Undyed)



ADVAMRYL Monofilament Poliglecaprone 25 (Undyed & Violet)



ADVAPD™ Monofilament Polydioxanone Violet



ADVALENE™ Monofilament Polypropylene Blue



ADVALON™ Monofilament Polyamide Black



ADVASYL™ Braided Silk Black



ADVABOND™ Polyester Green & White (Undyed)



Needle Coding

REVERSE CUTTING	RC	1/2 CIRCLE	1/2 CIR
TAPER CUT	TC	3/8 CIRCLE	3/8 CIR
TAPER POINT	TP	5/8 CIRCLE	5/8 CIR
SPATULATED MICRO POINT	SPM	1/4 CIRCLE	1/4 CIR
DIAMOND POINT	DIAP	STRAIGHT	ST
TROCAR POINT	TRO	COMPOUND CURVE / BI-CURVE	BI-CURVE
CONVENTIONAL CUTTING	CT	J-TYPE	J Type
ADVANCED CALCIFIED CORONARY	ACC	HEAVY	HV
LONG TAPER POINT	LTP	ADVABLACK	BLK
LONG TAPER POINT HEAVY	LTPH	TAPER PLUS NEEDLE	ADVAU
SLIM BLADE	SLB		
CARDIOVASCULAR	CV		

Abbreviations

DOUBLE ARM	DA	PLEDGET 8 X 3 X 1.5	PL8
SINGLE ARM	N	CARDIOVASCULAR	V
MULTI STRAND	M,M2,M3	WHITE	W
ADVABLACK	B	UNDYED	U
SHORT LENGTH	S,S1,S2	1 to 2 mm Needle & 10 cm suture length difference	ADVAT
LONG LENGTH	L,L1,L2	ADVALOOP	ADVALP
ANTIBACTERIAL SUTURE	P	ADVALENE LONG PACK CODES	ADVAL
PLEDGET 3 X 3 X 1.5	PL3	OVERSIZE SUTURE CODES	ADVAO
PLEDGET 6 X 3 X 1.5	PL6		

Where every suture is a seal of confidence



Suture Coding

ADVAPD	PDS	UNDYED	UND
ADVACRYL™	PGLA910	BLUE	BLUE
ADVAMRYL™	PGCL25	BLACK	BLK
ADVACRYL PLUS™	PGLA910 +	SILVER	SIL
ADVACRYL RAPID™	PGLA910 RPD	GREEN	GRN
ADVALENE™	POLY.P	WHITE	WHT
ADVASYL™	SILK	PLAIN	PLAIN
ADVALON™	POLYAMIDE	CHROMIC	CRMC
ADVABOND™	BOND	BLUE	BLUE
ADVACAT™	GUT	GREEN & WHITE	GRN & WHT
ADVASTEEL™	STEEL	ADVAMRYL™	MRL
ADVAMESH™	PP MESH	ADVAPD LOOP	LP
ADVAPLEDGET	PLDGT		
ADVAWAX™	WAX		
ADVAPACER™	PCW		
VIOLET	VIO		

ADVACRYL™ Braided Coated Polyglactin 910 (undyed & Violet)



Consistent performance for confidence in a variety of soft tissue approximation



With great versatility comes the confidence to meet specific needs

Product Description

Composition	90% Glycolide, 10% L-lactide		
Coating	Copolymer of Poly (glycolide-co-lactide) and Calcium Stearate		
Available Sizes	U.S.P. 2 to 6-0		
Colour	Undyed and Violet		
Absorption Profile	Essentially complete upto 80 days		
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Taper Cut, Trocar, Reverse Cutting, Blunt Point, Oval Round Bodied, Diamond Point & Spatula Micro Point needles		

Tensile Strength

Days of Implantation	14 Days	21 Days	28 Days
Tensile Strength Retention	75%	55%	20%

ADVANTAGES:

- Excellent knot security & lesser knot slippage
- Provides excellent support for range of critical wound healing periods
- Tighter braided construction

ADVACRYL™ Braided Coated Polyglactin - 910 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
1	ADVA2347WN180	VIOLET	1	180	-	Non Needed	12
2	ADVA2607M	VIOLET	1	(10 X 45)	-	Non Needed	12
3	ADVA2240	VIOLET	0	(6 X 45)	-	Non Needed	12
4	ADVA0416	VIOLET	0	(3 X 45)	-	Non Needed	12
5	ADVA2346WN180	VIOLET	0	180	-	Non Needed	12
6	ADVA2317WN180	VIOLET	2-0	180	-	Non Needed	12
7	ADVA0415	VIOLET	2-0	(3 X 45)	-	Non Needed	12
8	ADVA0414	VIOLET	3-0	(3 X 45)	-	Non Needed	12
9	ADVA2604MS	VIOLET	4-0	(3 X 45)	-	Non Needed	12
10	ADVA2352L	VIOLET	2	110	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
11	ADVA2353	VIOLET	2	90	50	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
12	ADVA2478	VIOLET	2	90	40	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
13	ADVA2351	VIOLET	2	90	40	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
14	ADVA2452	VIOLET	2	90	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
15	ADVA2360L	VIOLET	1	120	40	1/2 CIRCLE, ROUND BODIED	12
16	ADVA2350	VIOLET	1	120	36	1/2 CIRCLE, ROUND BODIED	12
17	ADVA2448L	VIOLET	1	110	50	1/2 CIRCLE, ROUND BODIED, HEAVY	12
18	ADVA2347VM	VIOLET	1	110	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
19	ADVA2448	VIOLET	1	100	50	1/2 CIRCLE, ROUND BODIED, HEAVY	12
20	ADVA2347PB	VIOLET	1	100	45	1/2 CIRCLE, ROUND BODIED, HEAVY	12
21	ADVA2438	VIOLET	1	90	50	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12

ADVACRYL™ Braided Coated Polyglactin - 910 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
22	ADVA2647	VIOLET	1	90	45	1/2 CIRCLE, ROUND BODIED, HEAVY	12
23	ADVA2421	VIOLET	1	90	40	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
24	ADVA2421G	VIOLET	1	90	40	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
25	ADVA2360	VIOLET	1	90	40	1/2 CIRCLE, ROUND BODIED	12
26	ADVA2347	VIOLET	1	90	40	1/2 CIRCLE, ROUND BODIED / TAPER POINT, HEAVY	12
27	ADVA2344	UNDYED	1	90	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
28	ADVA2519	VIOLET	1	90	40	1/2 CIRCLE, TAPER CUT	12
29	ADVA2364	VIOLET	1	90	40	1/2 CIRCLE, TAPER CUT, HEAVY	12
30	ADVA2422	VIOLET	1	90	40	3/8 CIRCLE, REVERSE CUTTING	12
31	ADVA945	UNDYED	2-0	90	36	1/2 CIRCLE, ROUND BODIED	12
32	ADVA2319	VIOLET	1	90	36	1/2 CIRCLE, TAPER CUT	12
33	ADVA9335	VIOLET	1	75	31	J Type, TAPER CUT, HEAVY,	12
34	ADVA2350S1	VIOLET	1	75	36	1/2 CIRCLE, ROUND BODIED	12
35	ADVA2359	VIOLET	1	70	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
36	ADVA2363	VIOLET	1	70	40	1/2 CIRCLE, TAPER CUT, HEAVY	12
37	ADVA2826E	VIOLET	1	70	23	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
38	ADVA2421S	VIOLET	1	45	40	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
39	ADVA2347S	VIOLET	1	45	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
40	ADVA2511	VIOLET	1	45	30	1/2 CIRCLE, TAPER CUT, HEAVY	12
41	ADVA2827	VIOLET	1	45	26	5/8 CIRCLE, REVERSE CUTTING	12
21	ADVA2438	VIOLET	1	90	50	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12


ADVACRYL™ Braided Coated Polyglactin - 910 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
43	ADVA2826	VIOLET	1	35	23	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
44	ADVA2346MM	VIOLET	0	110	40	1/2 CIRCLE, ROUND BODIED	12
45	ADVA2518L	VIOLET	0	110	40	1/2 CIRCLE, TAPER CUT	12
46	ADVA2346	VIOLET	0	90	40	1/2 CIRCLE, ROUND BODIED	12
47	ADVA2358	VIOLET	0	90	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
48	ADVA2518	VIOLET	0	90	40	1/2 CIRCLE, TAPER CUT	12
49	ADVA2533	VIOLET	0	90	40	3/8 CIRCLE, REVERSE CUTTING	12
50	ADVA2534	VIOLET	0	90	36	1/2 CIRCLE, REVERSE CUTTING	12
51	ADVA946	UNDYED	0	90	36	1/2 CIRCLE, ROUND BODIED	12
52	ADVA2517	VIOLET	0	90	36	1/2 CIRCLE, TAPER CUT	12
53	ADVA2338	VIOLET	0	90	30	1/2 CIRCLE, ROUND BODIED	12
54	ADVA0884	UNDYED	0	75	36	1/2 CIRCLE, ROUND BODIED	12
55	ADVA9901	UNDYED	0	75	40	1/2 CIRCLE, TAPER CUT	12
56	ADVA2656	VIOLET	0	75	45	1/2 CIRCLE, ROUND BODIED HEAVY	12
57	ADVA2342	VIOLET	0	70	30	1/2 CIRCLE, ROUND BODIED	12
58	ADVA2346S	VIOLET	0	45	40	1/2 CIRCLE, ROUND BODIED	12
59	ADVA1045	VIOLET	0	45	40	1/2 CIRCLE, TAPER CUT	12
60	ADVA2511SB	VIOLET	0	45	26	5/8 CIRCLE, REVERSE CUTTING	12
61	ADVA2510SB	VIOLET	0	45	23	1/2 CIRCLE, REVERSE CUTTING	12
62	ADVA2825	VIOLET	0	45	23	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
63	ADVA2534S	VIOLET	0	45	36	1/2 CIRCLE, REVERSE CUTTING	12






ADVACRYL™ Braided Coated Polyglactin - 910 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
64	ADVA2382	VIOLET	2-0	90	40	1/2 CIRCLE, REVERSE CUTTING	12
65	ADVA2345	VIOLET	2-0	90	40	1/2 CIRCLE, ROUND BODIED	12
66	ADVA2357	VIOLET	2-0	90	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
67	ADVA2280	VIOLET	2-0	75	36	1/2 CIRCLE, ROUND BODIED	12
68	ADVA945V	UNDYED	2-0	90	36	1/2 CIRCLE, TAPER CUT	12
69	ADVA2262L	UNDYED	2-0	90	35	1/2 CIRCLE, CUTTING	12
70	ADVA2317	VIOLET	2-0	90	30	1/2 CIRCLE, ROUND BODIED	12
71	ADVA2316	VIOLET	2-0	90	30	3/8 CIRCLE, REVERSE CUTTING	12
72	ADVA2356	VIOLET	2-0	90	26	1/2 CIRCLE, ROUND BODIED, HEAVY	12
73	ADVA2404	UNDYED	2-0	90	26	3/8 CIRCLE, CUTTING	12
74	ADVAT602	VIOLET	2-0	75	27	5/8 CIRCLE, TAPER POINT	12
75	ADVA9350	VIOLET	2-0	75	26	1/2 CIRCLE, TAPER CUT	12
76	ADVA9828	VIOLET	2-0	75	22	1/2 CIRCLE, TAPER CUT	12
77	ADVA2341	VIOLET	2-0	70	30	1/2 CIRCLE, ROUND BODIED	12
78	ADVA2122	VIOLET	2-0	70	26	1/2 CIRCLE, ROUND BODIED	12
79	ADVA602	VIOLET	2-0	70	26	5/8 CIRCLE, ROUND BODIED	12
80	ADVA2382S	VIOLET	2-0	45	40	1/2 CIRCLE, REVERSE CUTTING	12
81	ADVA2317S	VIOLET	2-0	45	30	1/2 CIRCLE, ROUND BODIED	12
82	ADVA2511A	VIOLET	2-0	45	30	1/2 CIRCLE, TAPER CUT	12
83	ADVA2694	VIOLET	2-0	45	26	5/8 CIRCLE, REVERSE CUTTING	12
84	ADVA2356S	VIOLET	2-0	35	26	1/2 CIRCLE, ROUND BODIED	12










ADVACRYL™ Braided Coated Polyglactin - 910 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
85	ADVA2511AA	VIOLET	2-0	20	30	1/2 CIRCLE, TAPER CUT	12
86	ADVA2516	VIOLET	3-0	90	36	1/2 CIRCLE, TAPER CUT	12
87	ADVA2528	VIOLET	3-0	90	25	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
88	ADVA2328	UNDYED	3-0	90	25	3/8 CIRCLE, CUTTING, PRIME	12
89	ADVA2241S	VIOLET	3-0	70	22	1/2 CIRCLE, ROUND BODIED	12
90	ADVA2402	UNDYED	3-0	90	19	3/8 CIRCLE, CUTTING, PRIME	12
91	ADVA2629	VIOLET	3-0	75	24	3/8 CIRCLE, CUTTING, PRIME	12
92	ADVA322	VIOLET	3-0	70	36	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
93	ADVA2123	VIOLET	3-0	70	26	1/2 CIRCLE, OVAL ROUND BODIED	12
94	ADVA2407	UNDYED	3-0	70	26	3/8 CIRCLE, REVERSE CUTTING	12
95	ADVA2407SP	UNDYED	3-0	70	25	3/8 CIRCLE, REVERSE CUTTING	12
96	ADVA2936	UNDYED	3-0	70	24	3/8 CIRCLE, CUTTING	12
97	ADVA936	UNDYED	3-0	70	24	3/8 CIRCLE, REVERSE CUTTING	12
98	ADVA603DA	VIOLET	3-0	70	26	 5/8 CIRCLE, ROUND BODIED	12
99	ADVA2402S	UNDYED	3-0	70	19	3/8 CIRCLE, CUTTING, PRIME	12
100	ADVA2427	VIOLET	3-0	70	12	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
101	ADVA683	UNDYED	3-0	45	24	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
102	ADVA2401	VIOLET	3-0	45	22	3/8 CIRCLE, CUTTING	12
103	ADVA2837	VIOLET	3-0	45	20	3/8 CIRCLE, REVERSE CUTTING	12
104	ADVA2515	VIOLET	3-0	45	17	1/2 CIRCLE, TAPER CUT	12
105	ADVA2471	VIOLET	3-0	45	16	3/8 CIRCLE, CUTTING, PRIME	12

ADVACRYL™ Braided Coated Polyglactin - 910 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
106	ADVA2495	VIOLET	3-0	45	12	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
107	ADVA9106	VIOLET	4-0	75	17	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
108	ADVA2304	VIOLET	4-0	70	20	1/2 CIRCLE, ROUND BODIED	12
109	ADVA2305	VIOLET	4-0	70	16	1/2 CIRCLE, ROUND BODIED	12
110	ADVA2443	UNDYED	4-0	70	16	3/8 CIRCLE, CUTTING, PRIME	12
111	ADVA2304S	VIOLET	4-0	45	20	1/2 CIRCLE, ROUND BODIED	12
112	ADVA2443S	UNDYED	4-0	45	16	3/8 CIRCLE, CUTTING, PRIME	12
113	ADVA2465	VIOLET	4-0	45	16	5/8 CIRCLE, CUTTING	12
114	ADVA2494	VIOLET	4-0	45	12	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
115	ADVA2395	VIOLET	5-0	70	10	3/8 CIRCLE, ROUND BODIED, LONG TAPER	12
116	ADVA2303	VIOLET	5-0	45	16	1/2 CIRCLE, OVAL ROUND BODIED	12
117	ADVA2403	VIOLET	5-0	45	16	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
118	ADVA2442	VIOLET	5-0	45	16	3/8 CIRCLE, CUTTING, PRIME	12
119	ADVA2872	UNDYED	5-0	45	13	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
120	ADVA2493	VIOLET	5-0	45	12	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
121	ADVA511G	UNDYED	5-0	45	10	1/2 CIRCLE, REVERSE CUTTING	12
122	ADVA2347DA200	VIOLET	1	200	40	 1/2 CIRCLE, ROUND BODIED & TAPER CUT HEAVY	12
123	ADVA2347DA180H	VIOLET	1	180	40	 1/2 CIRCLE, ROUND BODIED, HEAVY	36
124	ADVA2347DA180	VIOLET	1	180	40	 1/2 CIRCLE, ROUND BODIED, HEAVY	12
125	ADVA2344DA180	VIOLET	1	180	40	 1/2 CIRCLE, ROUND BODIED, HEAVY	12
126	ADVA2347DA160	VIOLET	1	160	40	 1/2 CIRCLE, ROUND BODIED, HEAVY	12

ADVACRYL™ Braided Coated Polyglactin - 910 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
127	ADVA2347DA140	VIOLET	1	140	40	 1/2 CIRCLE, ROUND BODIED, HEAVY	12
128	ADVA2360L1DA	VIOLET	1	180	40	 1/2 CIRCLE, TAPER POINT	12
129	ADVA2346DA180	VIOLET	0	180	40	 1/2 CIRCLE, ROUND BODIED	12
130	ADVA2345DA180	VIOLET	2-0	180	40	 1/2 CIRCLE, TAPER POINT	12
131	ADVA2318	VIOLET	2-0	90	30	 1/2 CIRCLE, ROUND BODIED	12
132	ADVA2436	VIOLET	3-0	90	26	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
133	ADVA2354	VIOLET	4-0	90	16	 3/8 CIRCLE CUTTING PRIME & 1/2 CIRCLE, TAPER POINT	12
134	ADVA2670	VIOLET	6-0	45	8	 1/4 CIRCLE, SPATULA	12
135	ADVA2949	VIOLET	2	75	65	1/2 CIRCLE TAPER POINT	12
136	ADVAY402	VIOLET	9-0	13	3.8	 3/8 CIRCLE TAPER POINT	12
137	ADVA2606	VIOLET	1	5 X 70	-	Non Needed	12
138	ADVA2196M2	VIOLET	2-0	2 X 70	-	Non Needed	12
139	ADVA2229SM2	VIOLET	3-0	2 X 70	-	Non Needed	12
140	ADVA901	VIOLET	4-0	2 X 70	-	Non Needed	12
141	ADVA901M	VIOLET	4-0	5 X 70	-	Non Needed	12
142	ADVA2123	VIOLET	3-0	70	26	1/2 CIRCLE, TAPER POINT	12
143	ADVA2877	VIOLET	4-0	90	19	3/8 CIRCLE, CUTTING	12
144	ADVA2981L	VIOLET	3-0	75	20	3/8 CIRCLE, CUTTING	12
145	ADVA2340	UNDYED	1	90	48	1/2 CIRCLE, TAPER POINT	12
146	ADVA2283L	VIOLET	0	90	48	1/2 CIRCLE, TAPER POINT	12
147	ADVAL122	VIOLET	3-0	75	26	1/2 CIRCLE, TAPER POINT, LONG TAPER	12

ADVACRYL™ Braided Coated Polyglactin - 910 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
148	ADVA2254	VIOLET	3-0	75	20	1/2 CIRCLE, TAPER POINT, HEAVY	12
149	ADVA2193	VIOLET	0	75	30	3/8 CIRCLE, CUTTING	12
150	ADVA2230	VIOLET	4-0	6 X 45	-	NON NEEDLE	12
151	ADVAT2293	VIOLET	4-0	75	26	1/2 CIRCLE, TAPER POINT	12
152	ADVA9150	VIOLET	2-0	75	40	1/2 CIRCLE, TAPER POINT	12
153	ADVA2117S	VIOLET	3-0	75	26	1/2 CIRCLE, TAPER POINT	12
154	ADVA2610	VIOLET	1	75	40	1/2 CIRCLE, TAPER POINT	12
155	ADVA2300S	VIOLET	1	75	37	1/2 CIRCLE, TAPER POINT, HEAVY	12
156	ADVA2231S	VIOLET	1	75	65	1/2 CIRCLE, TAPER POINT	12
157	ADVA2283	VIOLET	0	75	48	1/2 CIRCLE, TAPER POINT	12
158	ADVA2650S	VIOLET	1	75	48	1/2 CIRCLE, TAPER POINT, HEAVY	12
159	ADVA2234	VIOLET	2-0	75	48	1/2 CIRCLE, TAPER POINT	12
160	ADVA2639	VIOLET	0	75	40	1/2 CIRCLE, TAPER POINT	12
161	ADVA2612S1	VIOLET	1	75	40	1/2 CIRCLE, CUTTING	12
162	ADVA9106N	VIOLET	4-0	75	17	1/2 CIRCLE, TAPER POINT	12
163	ADVA2229SM2	VIOLET	3-0	2 X 70	-	Non Needle	12
164	ADVA2356	VIOLET	2-0	90	26	1/2 CIRCLE, TAPER POINT, HEAVY	12
165	ADVA2607LM2	VIOLET	1	2 X 70	0	Non Needle	12
166	ADVA2347	VIOLET	1	90	40	1/2 CIRCLE, TAPER POINT, HEAVY	12
167	ADVA29252	VIOLET	2	75	48	1/2 CIRCLE, TAPER POINT	12

ADVAMRYL™ Monofilament Poliglecaprone 25 (Undyed & Violet)



Supreme handling performance for confidence in subcuticular skin closure



Predictable strength for superior support

Product Description

Composition	Glycolide and ε-Caprolactone
Coating	None
Available Sizes	U.S.P. 0 to 5-0
Colour	Undyed and Violet
Absorption Profile	Essentially complete upto 90 days
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Taper Cut, Trocar, Reverse Cutting, Taper Cut, Blunt Point, Oval Round Bodied, Diamond Point needles

Tensile Strength

Days of Implantation	7 Days (1 to 5-0)	14 Days (1 to 3-0,5-0)	14 Days (4-0)
Tensile Strength Retention DYED	60%	40%	40%
Tensile Strength Retention UNDYED	55%	20%	15%

ADVANTAGES:

- Causes minimal tissue reaction
- High initial tensile strength
- Exceptional handling characteristics
- Lack of stiffness & memory


ADVAMRYL™ Monofilament Poliglecaprone 25 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
1	ADVA3442	VIOLET	0	70	31	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
2	ADVA1642	UNDYED	0	70	30	1/2 CIRCLE, ROUND BODIED	12
3	ADVA1742	VIOLET	0	70	30	1/2 CIRCLE, ROUND BODIED	12
4	ADVA945P	UNDYED	2-0	90	36	1/2 CIRCLE, ROUND BODIED	12
5	ADVA1666	UNDYED	2-0	70	36	1/2 CIRCLE, ROUND BODIED	12
6	ADVA3441	VIOLET	2-0	70	31	1/2 CIRCLE, ROUND BODIED	12
7	ADVA1641	UNDYED	2-0	70	30	1/2 CIRCLE, ROUND BODIED	12
8	ADVA1741	VIOLET	2-0	70	30	1/2 CIRCLE, ROUND BODIED	12
9	ADVA1602	UNDYED	2-0	70	26	5/8 CIRCLE, ROUND BODIED	12
10	ADVA1665	UNDYED	2-0	70	26	1/2 CIRCLE, OVAL ROUND BODIED	12
11	ADVA3665B	VIOLET	2-0	70	26	1/2 CIRCLE, TAPER POINT, ADVABLACK	12
12	ADVA1666S	UNDYED	2-0	45	36	1/2 CIRCLE, ROUND BODIED	12
13	ADVA3665	VIOLET	3-0	75	31	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
14	ADVA3650	UNDYED	3-0	70	60	STRAIGHT, CUTTING	12
15	ADVA3326M	UNDYED	3-0	70	26	3/8 CIRCLE, REVERSE CUTTING	12
16	ADVA1664	UNDYED	3-0	70	26	1/2 CIRCLE, OVAL ROUND BODIED	12
17	ADVA3664	VIOLET	3-0	70	26	1/2 CIRCLE, TAPER POINT, ADVABLACK	12
18	ADVA1326	UNDYED	3-0	70	25	3/8 CIRCLE, CUTTING, PRIME	12
19	ADVA1326MR	VIOLET	3-0	70	25	3/8 CIRCLE, CUTTING, PRIME	12
20	ADVA1663	UNDYED	3-0	70	22	1/2 CIRCLE, OVAL ROUND BODIED	12
21	ADVA1737	VIOLET	3-0	70	20	1/2 CIRCLE, ROUND BODIED	12


ADVAMRYL™ Monofilament Poliglecaprone 25 (Undyed & Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
22	ADVA1326S	UNDYED	3-0	45	25	3/8 CIRCLE, CUTTING, PRIME	12
23	ADVA3326MS	VIOLET	3-0	45	25	3/8 CIRCLE, REVERSE CUTTING	12
24	ADVA1205	UNDYED	4-0	70	16	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
25	ADVA1648	UNDYED	4-0	70	16	3/8 CIRCLE, ROUND BODIED, LONG TAPER	12
26	ADVA1205MR	VIOLET	4-0	70	16	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
27	ADVA1748	VIOLET	4-0	70	16	3/8 CIRCLE, ROUND BODIED, LONG TAPER	12
28	ADVA3205	UNDYED	4-0	45	16	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
29	ADVA845G	UNDYED	4-0	45	16	3/8 CIRCLE, CUTTING, PRIME	12
30	ADVA3204	UNDYED	5-0	45	16	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
31	ADVA844G	UNDYED	5-0	45	16	3/8 CIRCLE, CUTTING, PRIME	12
32	ADVA3203	UNDYED	5-0	45	13	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
33	ADVA493G	UNDYED	5-0	45	13	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
34	ADVA1004	UNDYED	3-0	70	24	3/8 CIRCLE, CUTTING	12
35	ADVAN1048	UNDYED	4-0	70	19	3/8 CIRCLE, CUTTING	12

ADVACRYL RAPID™ Fast Absorbing Polyglactin-910 (Undyed)



Accelerated performance for confidence in delicate tissue approximation



Swift confidence for superior predictability

Product Description

Composition	90% Glycolide, 10% L-lactide	
Coating	Poly (glycolide-co-lactide) (Glacomer 37)	
Available Sizes	U.S.P. 1 to 5-0	
Colour	Undyed	
Absorption Profile	Essentially complete upto 42 days	
Needle Availability	A wide variety of Round bodied (Taper Point), Oval Round bodied, Ball Point, Cutting, Taper Cut, Reverse Cutting, Blunt Point	




Tensile Strength

Days of Implantation	5 Days	10 - 14 Days
Tensile Strength Retention	42%	0%

ADVANTAGES:

- High initial tensile strength
- Predictable & constant degradation rate
- Excellent knot security
- Quick mass absorption

ADVACRYL RAPID™ Fast Absorbing Polyglactin-910 (Undyed)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
1	ADVA9964	UNDYED	1	90	36	1/2 CIRCLE, TAPER CUT	12
2	ADVA2763	UNDYED	0	110	36	1/2 CIRCLE, TAPER CUT	12
3	ADVA2763SP	UNDYED	0	100	40	1/2 CIRCLE, TAPER CUT	12
4	ADVA2763RB	UNDYED	0	100	36	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
5	ADVA2763FS	UNDYED	0	100	36	1/2 CIRCLE, TAPER CUT	12
6	ADVA9963	UNDYED	0	90	36	1/2 CIRCLE, TAPER CUT	12
7	ADVA2763S	UNDYED	0	70	36	1/2 CIRCLE, TAPER CUT	12
8	ADVA2721	UNDYED	0	60	35	3/8 CIRCLE, REVERSE CUTTING	12
9	ADVA2777	UNDYED	2-0	140	36	 1/2 CIRCLE, ROUND BODIED & 1/2 CIRCLE REVERSE CUTTING	12
10	ADVA2777TR	UNDYED	2-0	140	36	 1/2 CIRCLE, ROUND BODIED & 1/2 CIRCLE REVERSE CUTTING	12
11	ADVA2762ER	UNDYED	2-0	110	35	1/2 CIRCLE, TAPER CUT	12
12	ADVA2761SP	UNDYED	2-0	100	40	1/2 CIRCLE, TAPER CUT	12
13	ADVA2762RB	UNDYED	2-0	100	40	1/2 CIRCLE, ROUND BODIED	12
14	ADVA2761	UNDYED	2-0	100	36	1/2 CIRCLE, TAPER CUT	12
15	ADVA9962	UNDYED	2-0	90	36	1/2 CIRCLE, TAPER CUT	12
16	ADVA2762	UNDYED	2-0	90	35	1/2 CIRCLE, TAPER CUT	12
17	ADVA2727LDA	UNDYED	2-0	90	30	 1/2 CIRCLE, CUTTING	12
18	ADVA9933	UNDYED	2-0	75	26	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
19	ADVA2720	UNDYED	2-0	60	36	3/8 CIRCLE, REVERSE CUTTING	12
20	ADVA2762S	UNDYED	2-0	30	35	1/2 CIRCLE, TAPER CUT	12
21	ADVA9932	UNDYED	3-0	75	26	3/8 CIRCLE, REVERSE CUTTING, PRIME	12

ADVACRYL RAPID™ Fast Absorbing Polyglactin-910 (Undyed)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
22	ADVA9919	UNDYED	3-0	75	16	3/8 CIRCLE, CUTTING, PRIME	12
23	ADVA2732	UNDYED	3-0	70	26	3/8 CIRCLE, REVERSE CUTTING	12
24	ADVA2764	UNDYED	3-0	70	26	1/2 CIRCLE, ROUND BODIED	12
25	ADVA2735	UNDYED	3-0	70	22	1/2 CIRCLE, CUTTING	12
26	ADVA2732S	UNDYED	3-0	45	26	3/8 CIRCLE, REVERSE CUTTING	12
27	ADVA9935	UNDYED	3-0	45	22	1/2 CIRCLE, CUTTING	12
28	ADVA2719	UNDYED	3-0	45	16	3/8 CIRCLE, CUTTING	12
29	ADVA9918	UNDYED	4-0	75	16	3/8 CIRCLE, CUTTING, PRIME	12
30	ADVA2718	UNDYED	4-0	45	16	3/8 CIRCLE, CUTTING, PRIME	12
31	ADVA2718RB	UNDYED	4-0	45	16	3/8 CIRCLE, ROUND BODIED, LONG TAPER	12
32	ADVA9969	UNDYED	5-0	75	17	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
33	ADVA2772	UNDYED	5-0	75	16	3/8 CIRCLE, REVERSE CUTTING	12
34	ADVAT9915N	UNDYED	5-0	45	12	3/8 CIRCLE, CUTTING, PRIME	12
35	ADVA2660R	UNDYED	4-0	75	19	3/8 CIRCLE, CUTTING	12

ADVAPDTM Monofilament Polydioxanone (Violet)



Enduring performance for confidence in extended wound support



The confidence to handle complications with ease

Product Description

Composition	Polydioxanone
Coating	None
Available Sizes	U.S.P. 1 to 7-0
Colour	Violet
Absorption Profile	Essentially complete between 180 to 220 days
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Taper Cut, Trocar, Reverse Cutting, Blunt Point, Oval Round Bodied, Diamond Point & Spatula Micro Point needles

Tensile Strength

Days of Implantation	14 Days	28 Days	42 Days
Tensile Strength Retention	80%	70%	45%

ADVANTAGES:

- High initial tensile strength
- Predictable absorption
- Prolonged wound support
- Minimal tissue reaction
- Good handling characteristics
- Smooth passage through tissues

ADVAPD™ Monofilament Polydioxanone (Violet)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
1	ADVA9013	VIOLET	4-0	90	26	 1/2 CIRCLE, TAPER POINT, DA	12
2	ADVAS9032	VIOLET	1	90	17	1/2 CIRCLE, TAPER POINT	12
3	ADVA9108	VIOLET	2-0	90	17	1/2 CIRCLE, TAPER POINT	12
4	ADVA9053	VIOLET	2-0	75	30	1/2 CIRCLE, TAPER POINT	12
5	ADVA9336N	VIOLET	3-0	75	20	1/2 CIRCLE, TAPER POINT	12
6	ADVA9614	VIOLET	3-0	75	26	3/8 CIRCLE, CUTTING	12
7	ADVA9027L1	VIOLET	3-0	90	26	1/2 CIRCLE, TAPER POINT	12
8	ADVA9004L1	VIOLET	3-0	75	30	1/2 CIRCLE, TAPER POINT	12
9	ADVA9034LDA	VIOLET	4-0	90	20	 1/2 CIRCLE, TAPER POINT, DA	12
10	ADVAT9873L1	VIOLET	4-0	75	15	3/8 CIRCLE, REVERSE CUTTING	12
11	ADVA9049	VIOLET	4-0	90	17	1/2 CIRCLE, TAPER POINT	12
12	ADVAS9033	VIOLET	4-0	75	24	1/2 CIRCLE, REVERSE CUTTING	12
13	ADVA9029L	VIOLET	5-0	75	13	1/2 CIRCLE, TAPER POINT	12
14	ADVA9872L	VIOLET	5-0	75	16	3/8 CIRCLE, CUTTING	12
15	ADVA9020L	VIOLET	1	75	30	1/2 CIRCLE, TAPER POINT	12
16	ADVA9076	VIOLET	4-0	75	24	3/8 CIRCLE, REVERSE CUTTING	12
17	ADVAS9049	VIOLET	5-0	75	17	3/8 CIRCLE, CUTTING	12
18	ADVAS9050	VIOLET	3-0	75	25	3/8 CIRCLE, CUTTING	12

ADVALENE™ Monofilament Polypropylene (Blue)

Reliable performance for confidence in contaminated wound closure



With exceptional strength comes confidence that lasts








Product Description

Composition	Polypropylene
Coating	None
Available Sizes	U.S.P. 1 to 8-0
Colour	Blue
Absorption Profile	Non-Absorbable
Tensile Strength	Permanent
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Taper Cut, Trocar, Reverse Cutting, Blunt Point, Oval Round Bodied, Diamond Point needles, ACC Needles & Black Needles














ADVANTAGES:

- Excellent knot security
- Causes minimal tissue reaction
- Non - thrombogenic
- Least potentiate infection in contaminated wounds
- Advanced dispensing system
- Control linear elongation
- Unsurpassed strength
- Minimal suture breakage & fraying















ADVALENE™ Monofilament Polypropylene (Blue)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
1	ADVA834	BLUE	1	150	48	1/2 CIRCLE, ROUND BODIED, HEAVY, LOOP	12
2	ADVA883	BLUE	1	100	45	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
3	ADVA896	BLUE	1	100	30	1/2 CIRCLE, TAPER CUT, HEAVY	12
4	ADVA843	BLUE	1	70	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
5	ADVA840	BLUE	1	70	30	1/2 CIRCLE, ROUND BODIED, HEAVY	12
6	ADVA833	BLUE	0	150	40	1/2 CIRCLE, ROUND BODIED, LOOP	12
7	ADVA846	BLUE	0	90	30	1/2 CIRCLE, ROUND BODIED	12
8	ADVA894	BLUE	0	90	30	1/2 CIRCLE, TAPER CUT, HEAVY	12
9	ADVA824	BLUE	0	90	26	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
10	ADVA830	BLUE	0	70	30	1/2 CIRCLE, ROUND BODIED, HEAVY	12
11	ADVA842	BLUE	0	70	30	1/2 CIRCLE, ROUND BODIED	12
12	ADVA805	BLUE	0	70	25	1/2 CIRCLE, TAPER CUT	12
13	ADVA850	BLUE	2-0	90	36	1/2 CIRCLE, TAPER CUT	12
14	ADVA8526	BLUE	2-0	90	31	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
15	ADVAS8029	BLUE	2-0	75	20	1/2 CIRCLE, ROUND BODIED	12
16	ADVA806	BLUE	2-0	90	26	1/2 CIRCLE, TAPER CUT	12
17	ADVA8977	BLUE	2-0	90	26	 1/2 CIRCLE, TAPER CUT	12
18	ADVAL8977	BLUE	2-0	90	26	 1/2 CIRCLE, TAPER CUT, CV	12
19	ADVA877	BLUE	2-0	90	25	 1/2 CIRCLE, TAPER CUT	12
20	ADVA8937	BLUE	2-0	90	17	 1/2 CIRCLE, TAPER CUT, CV	12
21	ADVA837	BLUE	2-0	90	17	 1/2 CIRCLE, TAPER CUT, CV	12













ADVALENE™ Monofilament Polypropylene (Blue)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
43	ADVAL8021DA	BLUE	3-0	75	22	 1/2 CIRCLE, TAPER POINT, CV	12
44	ADVA8139L	BLUE	3-0	75	19	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
45	ADVA838	BLUE	3-0	70	25	1/2 CIRCLE, ROUND BODIED	12
46	ADVA802	BLUE	3-0	70	25	3/8 CIRCLE, CUTTING	12
47	ADVA822	BLUE	3-0	70	25	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
48	ADVA018	BLUE	3-0	70	22	3/8 CIRCLE, CUTTING	12
49	ADVA8790L	BLUE	3-0	75	30	3/8 CIRCLE, ROUND BODIED	12
50	ADVA897	BLUE	3-0	70	17	 1/2 CIRCLE, TAPER CUT	12
51	ADVA821S	BLUE	3-0	40	25	3/8 CIRCLE, CUTTING	12
52	ADVA818	BLUE	3-0	60	13	 3/8 CIRCLE, ROUND BODIED, LONG TAPER	12
53	ADVA820	BLUE	3-0	45	25	3/8 CIRCLE, CUTTING, PRIME	12
54	ADVA8140	BLUE	3-0	45	24	3/8 CIRCLE, REVERSE CUTTING	12
55	ADVA8521	BLUE	4-0	90	26	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
56	ADVA8521T	BLUE	4-0	90	26	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	24
57	ADVA8845	BLUE	4-0	90	26	 1/2 CIRCLE, CUTTING	12
58	ADVAL8521	BLUE	4-0	90	26	 1/2 CIRCLE, TAPER POINT, LONG TAPER CV	12
59	ADVAL8845	BLUE	4-0	90	26	 1/2 CIRCLE, CUTTING, CV	12
60	ADVA845	BLUE	4-0	90	25	 1/2 CIRCLE, TAPER CUT, CV	12
61	ADVA867	BLUE	4-0	90	25	 1/2 CIRCLE, TAPER CUT	12
62	ADVAL8015	BLUE	4-0	90	22	 1/2 CIRCLE, TAPER POINT, CV	12
63	ADVA8761	BLUE	4-0	90	20	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12


















ADVALENE™ Monofilament Polypropylene (Blue)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
22	ADVAL8937	BLUE	2-0	90	17	 1/2 CIRCLE, TAPER CUT, CV	12
23	ADVA863	BLUE	2-0	70	45	3/8 CIRCLE, REVERSE CUTTING	12
24	ADVA841	BLUE	2-0	70	30	1/2 CIRCLE, ROUND BODIED	12
25	ADVA807	BLUE	2-0	70	26	1/2 CIRCLE, TAPER CUT	12
26	ADVA808	BLUE	2-0	70	25	 1/2 CIRCLE, TAPER CUT	12
27	ADVA898	BLUE	2-0	70	17	 1/2 CIRCLE, TAPER CUT	12
28	ADVA8525	BLUE	3-0	90	31	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
29	ADVA8522	BLUE	3-0	90	26	 1/2 CIRCLE, ROUND BODIED	12
30	ADVA8976	BLUE	3-0	90	26	 1/2 CIRCLE, TAPER CUT	12
31	ADVAL8522	BLUE	3-0	90	26	 1/2 CIRCLE, TAPER POINT, CV	12
32	ADVA800	BLUE	3-0	90	25	3/8 CIRCLE, CUTTING, PRIME	12
33	ADVA825	BLUE	3-0	90	25	1/2 CIRCLE, ROUND BODIED	12
34	ADVA852	BLUE	3-0	90	25	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
35	ADVA836	BLUE	3-0	90	17	 1/2 CIRCLE, TAPER CUT	12
36	ADVA8936	BLUE	3-0	90	17	 1/2 CIRCLE, TAPER CUT, CV	12
37	ADVA8558	BLUE	3-0	90	17	 1/2 CIRCLE, ROUND BODIED	12
38	ADVAL8936	BLUE	3-0	90	17	 1/2 CIRCLE, TAPER CUT, CV	12
39	ADVAL8558	BLUE	3-0	90	17	 1/2 CIRCLE, TAPER POINT, CV	12
40	ADVA8667	BLUE	3-0	90	16	 1/2 CIRCLE, ROUND BODIED	12
41	ADVA860	BLUE	3-0	75	60	STRAIGHT, CUTTING	12
63	ADVA8761	BLUE	4-0	90	20	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12



















ADVALENE™ Monofilament Polypropylene (Blue)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
64	ADVA857	BLUE	4-0	90	17	1/2 CIRCLE, TAPER CUT	12
65	ADVA853	BLUE	4-0	90	17	 1/2 CIRCLE, TAPER CUT	12
66	ADVA8557	BLUE	4-0	90	17	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
67	ADVA8935	BLUE	4-0	90	17	 1/2 CIRCLE, TAPER CUT	12
68	ADVAL8935	BLUE	4-0	90	17	 1/2 CIRCLE, TAPER CUT, CV	12
69	ADVAL8557	BLUE	4-0	90	17	 1/2 CIRCLE, TAPER POINT, LONG TAPER CV	12
70	ADVA848	BLUE	4-0	90	16	1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
71	ADVA8526DA	BLUE	4-0	75	22	 1/2 CIRCLE, ROUND BODIED	12
72	ADVA8665	BLUE	4-0	75	13	 1/2 CIRCLE, ROUND BODIED	12
73	ADVAL8665	BLUE	4-0	75	13	 1/2 CIRCLE, TAPER POINT, CV	12
74	ADVA835	BLUE	4-0	70	17	 1/2 CIRCLE, TAPER CUT, CV	12
75	ADVA899	BLUE	4-0	70	17	1/2 CIRCLE, TAPER CUT	12
76	ADVA849	BLUE	4-0	70	16	1/2 CIRCLE, ROUND BODIED, CV	12
77	ADVA851	BLUE	4-0	70	16	3/8 CIRCLE, CUTTING, PRIME	12
78	ADVA870	BLUE	4-0	70	16	3/8 CIRCLE, CUTTING, PRIME	12
79	ADVA870DA	BLUE	4-0	70	16	 3/8 CIRCLE, CUTTING, PRIME	12
80	ADVA8127	BLUE	4-0	70	12	3/8 CIRCLE, REVERSE CUTTING	12
81	ADVA816	BLUE	4-0	60	16	3/8 CIRCLE, CUTTING, PRIME	12
82	ADVA8204	BLUE	4-0	60	13	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
83	ADVAL8204	BLUE	4-0	60	13	 1/2 CIRCLE, TAPER POINT, LONG TAPER CV	12
84	ADVA815	BLUE	4-0	45	16	3/8 CIRCLE, CUTTING	12



















ADVALENE™ Monofilament Polypropylene (Blue)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
85	ADVA879	BLUE	4-0	45	12	3/8 CIRCLE, REVERSE CUTTING	12
86	ADVAN8556	BLUE	5-0	90	17	 1/2 Circle, ROUND BODIED	12
87	ADVA856	BLUE	5-0	90	16	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
88	ADVA8556VPL3	BLUE	5-0	90	16	 1/2 CIRCLE, TAPER POINT, ADVAPASS	12
89	ADVA8556V	BLUE	5-0	90	16	 1/2 CIRCLE, TAPER POINT, CV	12
90	ADVA8556BV	BLUE	5-0	90	16	 1/2 CIRCLE, TAPER POINT, ADVABLACK CV	12
91	ADVA8321NLV	BLUE	5-0	90	13	 3/8 CIRCLE, TAPER POINT, CV	12
92	ADVA889	BLUE	5-0	90	12	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
93	ADVA810	BLUE	5-0	75	13	 1/2 CIRCLE, ROUND BODIED, LONG TAPER	12
94	ADVA8721	BLUE	5-0	75	13	 3/8 CIRCLE, CUTTING, 345 Microns	12
95	ADVA8321BV	BLUE	5-0	75	13	 3/8 CIRCLE, TAPER POINT, ADVABLACK & ADVAPASS	12
96	ADVA8710V	BLUE	5-0	75	13	 1/2 CIRCLE, TAPER POINT, ADVAPASS	12
97	ADVA8710VPL3	BLUE	5-0	75	13	 1/2 CIRCLE, TAPER POINT, ADVAPASS	12
98	ADVA8710BV	BLUE	5-0	75	13	 1/2 CIRCLE, TAPER POINT, ADVABLACK CV	12
99	ADVA8721V	BLUE	5-0	75	13	 3/8 CIRCLE, TAPER CUT, CV	12
100	ADVA8610LBV	BLUE	5-0	75	11	 1/2 CIRCLE, TAPER POINT, ADVABLACK CV	12
101	ADVA8810V	BLUE	5-0	75	10	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
102	ADVA892	BLUE	5-0	70	20	 3/8 CIRCLE, ROUND BODIED, LONG TAPER	12
103	ADVA881	BLUE	5-0	70	16	3/8 CIRCLE, ROUND BODIED, CV	12
104	ADVA882	BLUE	5-0	70	16	 3/8 CIRCLE, ROUND BODIED, LONG TAPER	12
105	ADVA826	BLUE	5-0	70	15	3/8 CIRCLE, CUTTING, LONG TAPER	12




ADVALENE™ Monofilament Polypropylene (Blue)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
106	ADVA862	BLUE	5-0	70	13	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
107	ADVA880	BLUE	5-0	70	12	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
108	ADVA827	BLUE	5-0	70	12	 3/8 CIRCLE, ROUND BODIED, LONG TAPER	12
109	ADVA8710SV	BLUE	5-0	60	13	 1/2 CIRCLE, TAPER POINT, CV	12
110	ADVA8610V	BLUE	5-0	60	10	 1/2 CIRCLE, TAPER POINT, CV	12
111	ADVA8307L1BV	BLUE	6-0	110	13	 3/8 CIRCLE, TAPER POINT, ADVABLACK & ADVAPASS	12
112	ADVA8726L1V	BLUE	6-0	90	13	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
113	ADVA8307BV	BLUE	6-0	75	13	 3/8 CIRCLE, TAPER POINT, ADVABLACK & ADVAPASS	12
114	ADVA8726LV	BLUE	6-0	75	13	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
115	ADVA8804BV	BLUE	6-0	75	10	 3/8 CIRCLE, TAPER POINT, ADVABLACK CV	12
116	ADVA8807BLV	BLUE	6-0	75	9.3	 3/8 CIRCLE, TAPER CUT, ADVABLACK CV	12
117	ADVA8403LV	BLUE	6-0	75	9	 3/8 CIRCLE, TAPER POINT, CV	12
118	ADVA8906V	BLUE	6-0	75	9	 1/2 CIRCLE, TAPER POINT, CV	12
119	ADVA8726V	BLUE	6-0	60	13	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
120	ADVA8802BV	BLUE	6-0	60	11	 3/8 CIRCLE, CUTTING, ADVABLACK	12
121	ADVA8597SP	BLUE	6-0	60	10	 3/8 Circle, ROUND BODIED , 282 Microns	12
122	ADVA8597V	BLUE	6-0	60	10	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
123	ADVA8597BV	BLUE	6-0	60	10	 3/8 CIRCLE, TAPER POINT, ADVABLACK	12
124	ADVA8207VPL3	BLUE	6-0	60	10	 1/2 CIRCLE, TAPER POINT, CV	12
125	ADVA8207V	BLUE	6-0	60	10	 1/2 CIRCLE, TAPER POINT, CV	12
126	ADVA878V	BLUE	6-0	45	12	3/8 CIRCLE, REVERSE CUTTING, CV	12

ADVALENE™ Monofilament Polypropylene (Blue)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
127	ADVA878V	BLUE	6-0	45	12	3/8 CIRCLE, REVERSE CUTTING, CV	12
128	ADVA8121NBV	BLUE	7-0	90	9	 3/8 CIRCLE, TAPER CUT, ADVABLACK CV	12
129	ADVA8735LV	BLUE	7-0	90	8	 3/8 CIRCLE, TAPER POINT, CV	12
130	ADVA2078DAV	BLUE	7-0	75	8	 1/2 CIRCLE, TAPER POINT, CV	12
131	ADVA8735L1V	BLUE	7-0	75	8	 3/8 CIRCLE, TAPER POINT, CV	12
132	ADVA0804LBV	BLUE	7-0	75	8	 3/8 CIRCLE, TAPER POINT, ADVABLACK CV	12
133	ADVA8747V	BLUE	7-0	60	9.5	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
134	ADVA8304V	BLUE	7-0	60	9.5	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
135	ADVA8088V	BLUE	7-0	60	9.3	 3/8 CIRCLE, TAPER CUT	12
136	ADVA8702V	BLUE	7-0	60	9	 3/8 CIRCLE, TAPER POINT, CV	12
137	ADVA8735V	BLUE	7-0	60	8	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
138	ADVA0804BV	BLUE	7-0	60	8	 3/8 CIRCLE, TAPER POINT, ADVABLACK CV	12
139	ADVA8186V	BLUE	7-0	60	6.5	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
140	ADVA0813LV	BLUE	8-0	90	9	 3/8 CIRCLE, TAPER POINT, CV	12
141	ADVA8732LV	BLUE	8-0	90	6.5	 3/8 CIRCLE, TAPER POINT, CV	12
142	ADVA0813V	BLUE	8-0	60	9	 3/8 CIRCLE, TAPER POINT, CV	12
143	ADVA8741V	BLUE	8-0	60	8	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
144	ADVA08741V	BLUE	8-0	60	8	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
145	ADVA8732V	BLUE	8-0	60	6.5	 3/8 CIRCLE, TAPER POINT, ADVAPASS	12
146	ADVA8411	BLUE	2-0	75	26	1/2 CIRCLE, TAPER POINT	12
147	ADVA8043	BLUE	4-0	75	26	1/2 CIRCLE, TAPER POINT	12

ADVALENE™ Monofilament Polypropylene (Blue)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
148	ADVA846SP	BLUE	0	75	37	1/2 CIRCLE, TAPER POINT	12
149	ADVA8890V	BLUE	5-0	75	13	3/8 CIRCLE, TAPER POINT	12
150	ADVA8201	BLUE	4-0	75	18	 1/2 CIRCLE, TAPER POINT	12
151	ADVA8125	BLUE	4-0	70	25	 1/2 CIRCLE, TAPER POINT	12
152	ADVA8771L	BLUE	3-0	90	22	 1/2 CIRCLE, TAPER POINT	12

ADVALON™ Monofilament Polyamide (Black)



Ideal selection for confidence in skin closure

Appropriate choice for dermal closure and retention suturing

Product Description

Composition	Polyamide
Coating	None
Available Sizes	U.S.P. 1 to 10-0
Colour	Black
Absorption Profile	Non-Absorbable
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Taper Cut, Trocar, Reverse Cutting, Blunt Point, Oval Round Bodied, Diamond Point & Spatula Micro Point needles

Tensile Strength

Non Absorbable, however progressive hydrolysis of the suture in-vivo may result in gradual loss of its tensile strength over time

ADVANTAGES:

- Flexible, easy to handle and tie
- Smooth passage through tissue
- Excellent histocompatibility






ADVALON™ Monofilament Polyamide (Black)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
1	ADVA905	BLACK	1	(2 X 76)	-	Non Needed	12
2	ADVA3348	BLACK	1	150	48	1/2 CIRCLE, ROUND BODIED, HEAVY, LOOP	12
3	ADVA3362	BLACK	1	100	50	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
4	ADVA3338	BLACK	1	100	45	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
5	ADVA3347	BLACK	1	100	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
6	ADVA3335	BLACK	1	100	40	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
7	ADVA3305	BLACK	1	75	38	1/2 CIRCLE, REVERSE CUTTING	12
8	ADVA3397	BLACK	1	70	90	3/8 CIRCLE, CUTTING	12
9	ADVA164	BLACK	0	(2 X 76)	-	Non Needed	12
10	ADVA3340	BLACK	0	150	40	1/2 CIRCLE, ROUND BODIED, LOOP	12
11	ADVA3339	BLACK	0	100	45	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
12	ADVA3346	BLACK	0	100	40	1/2 CIRCLE, ROUND BODIED, HEAVY	12
13	ADVA3334	BLACK	0	100	40	1/2 CIRCLE, REVERSE CUTTING, HEAVY	12
14	ADVA3304	BLACK	0	75	38	1/2 CIRCLE, REVERSE CUTTING	12
15	ADVA3389S	BLACK	0	70	60	STRAIGHT, CUTTING	12
16	ADVA3337	BLACK	0	70	45	3/8 CIRCLE, REVERSE CUTTING	12
17	ADVA3336DK	BLACK	2-0	110	45	3/8 CIRCLE, REVERSE CUTTING	12
18	ADVA3336L	BLACK	2-0	90	45	3/8 CIRCLE, REVERSE CUTTING	12
19	ADVA3088	BLACK	2-0	75	60	STRAIGHT, CUTTING, ,	12
20	ADVA2852	BLACK	2-0	75	26	1/2 CIRCLE, ROUND BODIED	12
21	ADVA3275SL	BLACK	2-0	45	26	3/8 CIRCLE, CUTTING	12

ADVALON™ Monofilament Polyamide (Black)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
22	ADVA3333	BLACK	2-0	70	40	3/8 CIRCLE, REVERSE CUTTING	12
23	ADVA3336S	BLACK	2-0	35	45	3/8 CIRCLE, REVERSE CUTTING	12
24	ADVA3190L	BLACK	3-0	110	26	3/8 CIRCLE, REVERSE CUTTING	12
25	ADVA3328L	BLACK	3-0	90	26	3/8 CIRCLE, REVERSE CUTTING	12
26	ADVA3341	BLACK	3-0	90	25	3/8 CIRCLE, REVERSE CUTTING	12
27	ADVA3328	BLACK	3-0	70	26	3/8 CIRCLE, REVERSE CUTTING	12
28	ADVA3328G	BLACK	3-0	70	26	3/8 CIRCLE, REVERSE CUTTING	12
29	ADVA3378	BLACK	3-0	75	16	1/2 CIRCLE, TAPER POINT	12
30	ADVA3332	BLACK	3-0	70	16	3/8 CIRCLE, REVERSE CUTTING	12
31	ADVA3328S	BLACK	3-0	35	26	3/8 CIRCLE, REVERSE CUTTING	12
32	ADVA3564	BLACK	4-0	45	8	5/8 CIRCLE, TAPER POINT DA	12
33	ADVA3319	BLACK	4-0	90	16	3/8 CIRCLE, CUTTING, PRIME	12
34	ADVA3567	BLACK	4-0	70	19	3/8 CIRCLE, CUTTING, PRIME	12
35	ADVA3318	BLACK	4-0	70	16	3/8 CIRCLE, CUTTING	12
36	ADVA3357	BLACK	4-0	60	19	3/8 CIRCLE, CUTTING, PRIME	12
37	ADVA3356	BLACK	4-0	45	19	3/8 CIRCLE, CUTTING, PRIME	12
38	ADVA3326	BLACK	4-0	38	10	3/8 CIRCLE, REVERSE CUTTING	12
39	ADVA3061S1	BLACK	4-0	45	16	3/8 CIRCLE, TAPER POINT	12
40	ADVA3316	BLACK	5-0	90	12	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
41	ADVA3726S	BLACK	5-0	45	15	3/8 CIRCLE, CUTTING	12
42	ADVA3752	BLACK	5-0	70	13	3/8 CIRCLE, REVERSE CUTTING, PRIME	12

ADVALON™ Monofilament Polyamide (Black)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
43	ADVA3317	BLACK	5-0	70	12	3/8 CIRCLE, REVERSE CUTTING	12
44	ADVA3353	BLACK	5-0	60	13	3/8 CIRCLE, REVERSE CUTTING, PRIME	12
45	ADVA3315	BLACK	5-0	35	16	3/8 CIRCLE, CUTTING, PRIME	12
46	ADVA3708	BLACK	8-0	38	6	3/8 CIRCLE, ROUND BODIED, MICRO POINT	12
47	ADVA3708SP	BLACK	8-0	38	6	3/8 CIRCLE, ROUND BODIED	12
48	ADVA3709	BLACK	9-0	38	6	3/8 CIRCLE, ROUND BODIED, MICRO POINT	12
49	ADVAT3058	BLACK	10-0	30	6.19	 3/8 CIRCLE, SPATULATED MICRO POINT	12
50	ADVA1716	BLACK	8-0	30	7	 1/2 CIRCLE, SPATULATED MICRO POINT	12
51	ADVA3918	BLACK	10-0	30	5	 3/8 CIRCLE, CUTTING	12
52	ADVA3919	BLACK	10-0	75	3.8	3/8 CIRCLE, TAPER POINT	12
53	ADVA3048	BLACK	5-0	75	16	3/8 CIRCLE, CUTTING	12
54	ADVA2812	BLACK	8-0	13	5	3/8 CIRCLE, TAPER POINT	12
55	ADVA3046	BLACK	9-0	30	6	 3/8 CIRCLE, SPATULATED MICRO POINT DA	12
56	ADVA2871	BLACK	9-0	13	4	3/8 CIRCLE, TAPER POINT	12
57	ADVA2850G	BLACK	10-0	13	5	3/8 CIRCLE, TAPER POINT	12
58	ADVA3919S2	BLACK	10-0	13	3.8	3/8 CIRCLE, TAPER POINT	12
59	ADVA3188L	BLACK	0	75	36	3/8 CIRCLE, CUTTING	12
60	ADVA3058	BLACK	10-0	30	6	 3/8 CIRCLE, SPATULATED MICRO POINT	12
61	ADVA3051	BLACK	4-0	75	19	3/8 CIRCLE, CUTTING	12
62	ADVA33451	BLACK	3-0	75	30	1/2 CIRCLE, CUTTING	12
63	ADVA3780	BLACK	3-0	75	26	1/2 CIRCLE, CUTTING	12

ADVASYL™ Braided Silk (Black)



Exceptional control for confidence in gentle tissue approximation

Excellent handling combination with great strength for gentle tissue approximation

Product Description

Composition	Proteinaceous silk fibres called fibroin
Coating	Wax
Available Sizes	U.S.P. 4 to 5-0
Colour	Black
Absorption Profile	Non-Absorbable
Needle Availability	Available as Pre-cut suture and in a wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Taper Cut, Trocar, Reverse Cutting, Blunt Point, Oval Round Bodied & Diamond Point needles

Tensile Strength

Gradual loss in tensile strength over a period of time

ADVANTAGES:

- Smooth passage through tissue
- Minimal acute inflammatory reaction in tissues
- Unsurpassed strength
- Optimal elasticity and elongation properties

ADVASYL™ Braided Silk (Black)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
1	ADVA5676L	BLACK	3-0	75	30	3/8 CIRCLE, CUTTING	12
2	ADVA5326S	BLACK	3-0	75	25	1/2 CIRCLE, TAPER CUT	12
3	ADVA5275	BLACK	2-0	75	22	1/2 CIRCLE, TAPER POINT	12
4	ADVA5333A	BLACK	2-0	75	30	1/2 CIRCLE, TAPER POINT	12
5	ADVA5152	BLACK	0	75	30	3/8 CIRCLE, CUTTING	12
6	ADVA5179L	BLACK	0	75	25	1/2 CIRCLE, TAPER CUT	12
7	ADVA5142	BLACK	0	75	26	1/2 CIRCLE, TAPER POINT	12
8	ADVA5154	BLACK	3-0	75	25	3/8 CIRCLE, CUTTING	12
9	ADVA5263	BLACK	3-0	75	26	1/2 CIRCLE, TAPER POINT	12
10	ADVA5670	BLACK	2-0	76	25	1/2 CIRCLE, TAPER CUT	12
11	ADVAK959S	BLACK	2-0	75	40	3/8 CIRCLE, CUTTING	12
12	ADVAS5004	BLACK	5-0	45	10	3/8 CIRCLE SPM DA	12
13	ADVA1732	BLACK	6-0	45	6.5	3/8 CIRCLE SPM DA	12

ADVABOND™ Braided Coated Polyester (Green & White)

Ultimate security for confidence in approximation of tissues under stress



Enhancing confidence for everlasting support

Product Description

Composition	Polyethylene Terephthalate
Coating	Silicone
Available Sizes	U.S.P. 5 to 6-0
Colour	White & Green
Absorption Profile	Non-Absorbable
Tensile Strength	Permanent
Needle Availability	Round bodied (Taper point) & Taper cutting

ADVANTAGES:

- Soft and supple
- Silicon coated, designed for excellent handling and pliability
- Uniform braiding process designed to provide greater resistance for breaking
- Excellent knot security
- Low porosity PTFE felt pledgets
- Advanced intelligent dispensing system

ADVABOND™ Braided Coated Polyester (Green & White)

Sr. No.	Code	Colour	Size	Suture Lgth (cm)	Needle Lgth (mm)	Description	Units Per Box
1	ADVA4845	GREEN	5	75	55	1/2 CIRCLE, TAPER CUT, HEAVY	12
2	ADVA4842	GREEN	2	75	45	1/2 CIRCLE, TAPER CUT, HEAVY	12
3	ADVA6887S1	GREEN	2-0	45	45	3/8 CIRCLE, CUTTING	12
4	ADVA6002DA	GREEN	2-0	75	26	 1/2 CIRCLE, TAPER POINT	12
5	ADVA4856	GREEN	3-0	45	24	3/8 CIRCLE, CUTTING	12
6	ADVA6032	GREEN	3-0	75	17	1/2 CIRCLE, TAPER POINT	12
7	ADVA4857	GREEN	4-0	45	19	3/8 CIRCLE, CUTTING	12
8	ADVA6002	GREEN	2-0	75	26	1/2 CIRCLE, TAPER POINT	12
9	ADVA6355	GREEN	1	75	40	1/2 CIRCLE, CUTTING	12
10	ADVA6356	GREEN	2	75	40	1/2 CIRCLE, CUTTING	12
11	ADVA6419	GREEN	3	2 X 70	-	Non Needle	12
12	ADVA6225	GREEN	2-0	90	25	 1/2 CIRCLE, TAPER POINT	12

ADVAWAX

A sterile mixture of Natural Wax and Isopropyl Myristate



Non-Absorbable Bone Wax



Ultimate confidence for control of bleeding from bone surfaces

DESCRIPTION:

➤ ADVAWAX is a sterile mixture of bees wax and Isopropyl Myristate, a wax softening agent. It is opaque and has a waxy odor. It is available in the slab form of 2.5 gms.

Sr. No.	Code	Colour	Weight (gm)	Beeswax + IPM Combination	Paper Wrap
1	ADVAW810	WHITE	Min. 2.5	80% + 20%	62 GSM - Siliconized

Absorbable Suture Portfolio

	ADVACRYL™	ADVACRYL PLUS™	ADVACRYL RAPID™	ADVAMRYL™	ADVAPD™	ADVACAT™
Composition	90% Glycolide, 10% L-lactide Braided Coated Synthetic Absorbable Surgical Suture	90% Glycolide, 10% L-lactide Antibacterial with Triclosan Braided Coated Absorbable Surgical Suture	90% Glycolide, 10% L-lactide Braided Coated Synthetic Absorbable Surgical Suture	Poliglecaprone 25 Monofilament Synthetic Absorbable Surgical Suture	Polydioxanone Monofilament Synthetic Absorbable Surgical Suture	Collagen from Submucosal layer, Natural Absorbable Surgical Suture
Tensile Strength	70% - 14 days 50% - 21 days 20% - 28 days	75% - 14 days 40% - 55% - 21 days 15% - 28 days	42% - 5 days All of the original tensile strength is lost by 10 to 14 days post implantation	Dyed 60%, Undyed 55% - 7 days Dyed 40%, Undyed 15 - 20% - 14 days	80% - 14 days 70% - 28 days 45% - 42 days	Plain 10 to 12 days Chronic 21 to 25 days
Mass Absorption	Upto 80 days	Upto 80 days	Upto 42 days	Upto 90 days	180 - 220 days	84 days
Available Sizes	U.S.P. 2 to 6-0	USP. 2 to 4-0	USP# 1 to 5-0	U.S.P. 0 to 5-0	U.S.P. 1 to 7-0	U.S.P. 2 to 3-0
Colour	Undyed and Violet	Violet	Undyed	Undyed and Violet	Violet	Yellow and Brown

Non-Absorbable Suture Portfolio

	ADVALENE™	ADVABOND™	ADVALON™	ADVASYL™
Composition	Polypropylene Monofilament Synthetic Non-Absorbable Surgical Suture	Polyethylene Terephthalate Braided Coated Synthetic Non-Absorbable Surgical Suture	Polyamide / Nylon Monofilament Synthetic Non-Absorbable Surgical Suture	Proteinaceous Silk Fibers called Fibroin, Braided Natural Non-Absorbable Surgical Suture
Tensile Strength*	Permanent	Gradual encapsulation of suture by fibrous connective tissue	Gradual loss in tensile strength over a period of time. Gradual encapsulation of suture by fibrous connective tissue	
Mass Absorption*	Non-Absorbable	Non-Absorbable	Non-Absorbable	Non-Absorbable
Available Sizes	U.S.P. 1 to 8-0	USP. 5 to 6-0	USP# 1 to 10-0	U.S.P. 4 to 5-0
Colour	Blue	White and Green	Black	Black

* These are approximate days which may vary with size.

* Mass Absorption : Time in days required for suture to be totally absorbed by the body.



Where every suture is a seal of confidence



Corporate HQ: P-21-22, 25-26 & 34-35, GIDC Manjusar, tal - Savli, Dist. - Vadodara - 391775, Gujarat, India
T: +91 2667 271302 | E: info@amsltd.com | amsltd.com



AMSSALES/TPR/2024/00