




CATALOGUE 2021

CARDIAC CATHETERS
AND CANNULAE



Cardiovascular cannulation



Index

- 4 Arterial cannulae
- 10 Aortic catheters
- 12 Flex line venous cannulae
- 14 Flex line two stage venous
cannulae

Flex line venous cannulae

REINFORCED WITH LIGHTHOUSE TIP



PRODUCT INFORMATION


Andocor venous cannulae are intended to drain venous blood from the right atrium and caval veins.

Q Connection site: 3/8" - 36 Fr: 1/2"
Length: 39 cm
Qty/box: 20 Pcs


+ Single piece construction
Open and flexible lighthouse tip
Easy positioning
Insertion depth marker

Size	Straight tip	Right angled 90° tip
20 Fr	01V201L8	01V201L9
22 Fr	01V221L8	01V221L9
24 Fr	01V241L8	01V241L9
26 Fr	01V261L8	01V261L9
28 Fr	01V281L8	01V281L9
30 Fr	01V301L8	01V301L9
32 Fr	01V321L8	01V321L9
34 Fr	01V341L8	01V341L9
36 Fr	01V361L8	01V361L9

VENOUS CANNULA PRESSURE DROP VS FLOW



VENOUS CANNULA PRESSURE DROP VS FLOW



Flex line two stage venous cannulae

PROXIMAL AND DISTAL REINFORCED WITH LIGHTHOUSE TIP



ROUND BODY



FLAT BODY



ROUND BODY					FLAT BODY		
Size	Length	Without connector	With connector	With luer connector	Without connector	With connector	With luer connector
32/40	36	01V32L40	01V32L401	01V32L402	01V32L40FB	01V32L401FB	01V32L402FB
34/40	36	01V34L40	01V34L401	01V34L402	01V34L40FB	01V34L401FB	01V34L462FB
34/46	40	01V34L46	01V34L461	01V34L462	01V34L46FB	01V34L461FB	01V34L462FB
36/46	40	01V36L46	01V36L461	01V36L462	01V36L46FB	01V36L461FB	01V36L462FB
36/51	40	01V36L51	01V36L511	01V36L512	01V36L51FB	01V36L511FB	01V36L512FB

PRODUCT INFORMATION


Andocor two stage venous cannulae are intended for use in venous drainage via the right atrium and inferior vena cava simultaneously during cardiopulmonary bypass surgery up to six hours or less.

- 

Connection site: 1/2"
Qty/box: 20 Pcs
- 

Single piece construction
Open and flexible lighthouse tip
Easy positioning
Insertion depth marker

VENOUS CANNULA PRESSURE DROP VS FLOW




Cardioplegia

Index

- 18 Aortic root cannulae
- 20 Retrograde cardioplegia cannulae
- 21 Air aspiration needles
- 22 Ostial perfusion cannulae
- 24 Cardioplegia adapters

Aortic root cannulae



Description	Tip size	Needle size	Reference
Type A Without vent line	6 Fr	16 G	ARC07
	7 Fr	14 G	ARC04
	9 Fr	12 G	ARC01
Type B Without vent line With suture ring	6 Fr	16 G	ARC07R
	7 Fr	14 G	ARC04R
	9 Fr	12 G	ARC01R
Type C With vent line	6 Fr	16 G	ARC08
	7 Fr	14 G	ARC05
	9 Fr	12 G	ARC02
Type D With vent line	6 Fr	16 G	ARC08S
	7 Fr	14 G	ARC05S
	9 Fr	12 G	ARC02S


PRODUCT INFORMATION

Andocor aortic root cannulae are intended for use during cardiopulmonary bypass for antegrade delivery of cardioplegia solutions and venting of the heart. The cannula may also be used to aspirate air from the aorta at the conclusion of the cardiac procedure.

Q Length: 15 cm
Connection size vent line: 1/4"
Qty/box: 20 Pcs

+ Easy of use
Anti-bleeding valve helps reduce blood loss

AORTIC ROOT CANNULAE PRESSURE DROP VS FLOW



Ostial perfusion cannulae

CORONARY ARTERY OSTIAL CANNULA, BASKET TIP WITH FLANGE



Tip size	45° angled	90° angled
10 Fr	OC1045	OC1090
12 Fr	OC1245	OC1290
14 Fr	OC1445	OC1490

PRODUCT INFORMATION

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

Q Connection: Female luer lock
Length: 21 cm
Qty/box: 20 Pcs

+ Soft basket tip
Different configurations
Stainless steel shaft

CORONARY ARTERY OSTIAL CANNULA PRESSURE DROP VS FLOW

