

## Biopsy Forceps

- Standard
- Hot

Used endoscopically to obtain mucosal tissue samples from gastrointestinal tract for microscopic-histopathological examination



Art. No.		Sheath Ø (mm)	Working Length (cm)	Working Channel Ø (mm)
<b>SBF-WS-16</b>	with spike	2.4	160	2.8
<b>SBF-WS-23</b>	with spike	2.4	230	2.8
<b>SBF-16</b>	without spike	2.4	160	2.8
<b>SBF-23</b>	without spike	2.4	230	2.8
<b>SBFH-23 HOT</b>	without spike	2.4	230	2.8

## Biopsy Forceps

with PTFE Tube

Used endoscopically to obtain mucosal tissue samples from gastrointestinal tract for microscopic-histopathological examination

For easy insertion and sliding withdrawal without residues



Art. No.		Sheath Ø (mm)	Working Length (cm)	Working Channel Ø (mm)
<b>WS-2415BWH</b>	with spike	2.4	150	2.8
<b>WS-2422BWH</b>	with spike	2.4	220	2.8
<b>WS-2415BTH</b>	without spike	2.4	150	2.8
<b>WS-2422BTH</b>	without spike	2.4	220	2.8