


# LC177 PRETREATMENT PROTEK CARTRIDGE




- **General**

<b>Main Feature</b>	High Capacity Adsorption pretreatment cartridge – large size. Protects the RO membranes by removing particulates, chlorine and chloramine.
<b>Compatible with</b>	BIOPURE 60/120, CENTRA R 60/120, MEDICA PRO, PURELAB 3000, PURELAB 7000
<b>Unit of measure</b>	Individual
<b>Part Numbers of Compatible Systems</b>	CN060RDM1-115 , CN120RDM1-115 , MB060XXM1-115 , MB120XXM1-115 , MP030BBM1-115 , MP030RBM1-115 , MP030REM1-115 , MP030SBM1-115 , MP060BBM1-115 , MP060RBM1-115 , MP060REM1-115 , MP060SBM1-115 , MP120BBM1-115 , MP120RBM1-115 , MP120REM1-115 , MP120SBM1-115 , PB3060XXM1-115 , PB3120XXM1-115 , PB7060XXM1-115 , PB7120XXM1-115

## MEDPURE L1 DEIONISATION CARTRIDGE

- **General**

<b>Main Feature</b>	Removes ionic contaminants to achieve specified water quality
<b>Compatible with</b>	MEDICA PRO
<b>Unit of measure</b>	Individual
<b>Part Numbers of Compatible Systems</b>	MP030BBM1-115 , MP030RBM1-115 , MP030REM1-115 , MP030SBM1-115 , MP060BBM1-115 , MP060RBM1-115 , MP060REM1-115 , MP060SBM1-115 , MP120BBM1-115 , MP120RBM1-115 , MP120REM1-115 , MP120SBM1-115



## LC105 UV LAMP (254nm)

<b>Main Feature</b>	UV light is used as a bactericide and to break down and photo-oxidize organic contaminants to polar or ionized species. The UV bactericidal action damages DNA and RNA polymerase at low doses thereby preventing replication, while higher doses are biocidal.
<b>Compatible with</b>	BIOPURE 60/120, BIOPURE 7/15, MAXIMA MK1/MK2, MEDICA 30/60/100 Mk2 & 3, MEDICA 7/15 Mk1, MEDICA PRO, MEDICA-D 7/15, MEDICA-R 30/60/100, MEDICA-R 7/15, MEDICA-S 7/15, OPTIMA 15/30/60 & PLUS, Option 1/2/3/4/5/6, PURELAB Elect 5/10, PURELAB Option 7/15 and PS, PURELAB Option-E 25/50/80 (Elect), PURELAB Option-E 5/10, PURELAB Option-R 30/60, PURELAB Option-R 7/15, PURELAB Option-S 30/60, PURELAB Ultra Mk1, PURELAB Ultra Mk2
<b>Unit of measure</b>	Individual
<b>Part Numbers of Compatible Systems</b>	ELE025XXXX11503-1 , ELE050XXXX11503-1 , ELE080XXXX11503-1 , MB007BPM1 , MB007PK01 , MB007XXM1 , MB015BPM1 , MB015PK01 , MB015XXM1 , MB060XXM1-115 , MB120XXM1-115 , MD007BPM1 , MD007BPM2 , MD007XXM1 , MD007XXM2 , MD015BPM1 , MD015BPM2 , MD015XXM1 , MD015XXM2 , ME007BPM1 , ME007PK01 , ME007XXM1 , ME015BPM1 , ME015PK01 , ME015XXM1 , MP030BBM1-115 , MP030RBM1-115 , MP030REM1-115 , MP030SBM1-115 , MP060BBM1-115 , MP060RBM1-115 , MP060REM1-115 , MP060SBM1-115 , MP120BBM1-115 , MP120RBM1-115 , MP120REM1-115 , MP120SBM1-115 , MR007BPM1 , MR007PK01 , MR007PK02 , MR007XXM1 , MR015BPM1 , MR015PK01 , MR015PK02 , MR015XXM1 , OE005BPM2 , OE005XXM2 , OE010BPM2 , OE010XXM2 , OPT030PLXX11503-1 , OPT030XXXX11503-1 , OPT060PLXX11503-1 , OPT060XXXX11503-1 , OR007BPM1 , OR007XXM1 , OR015BPM1 , OR015XXM1 , ULXXXIOM1 ,

ULXXIOM2 , MED030XXX11503-1 , MED060XXX11503-1 ,  
MED100XXX11503-1 , MR030XXM1-115 , MR060XXM1-115 ,  
MR100XXM1-115



## LC109 UMF CARTRIDGE



- **General**

<b>Main Feature</b>	Microporous filters provide a physical barrier to the passage of particles and microorganisms in purified water systems.
<b>Compatible with</b>	MAXIMA MK1/MK2, MEDICA 30/60/100 Mk2 & 3, MEDICA PRO, MEDICA-R 30/60/100, PURELAB Chorus Analytical, PURELAB UHQ II & UHQ PS, PURELAB Ultra Mk1, PURELAB Ultra Mk2
<b>Unit of measure</b>	Individual
<b>Part Numbers of Compatible Systems</b>	MP030BBM1-115 , MP030RBM1-115 , MP030REM1-115 , MP030SBM1-115 , MP060BBM1-115 , MP060RBM1-115 , MP060REM1-115 , MP060SBM1-115 , MP120BBM1-115 , MP120RBM1-115 , MP120REM1-115 , MP120SBM1-115 , UHQ000IIXX11503-1 , UHQ000PSXX11503-1 , ULRSXANM1 , ULRSXANM2 , ULXXXANM1 , ULXXXANM2 , ULXXXIOM1 ,

ULXXIOM2 , MED030XXXX11503-1 , MED060XXXX11503-1 ,  
MED100XXXX11503-1 , MR030XXM1-115 , MR060XXM1-115 ,  
MR100XXM1-115, PC1ANRXM2

## Elga LC136M2 Composite Air Vent Filter, For Water Purification Unit

- Air vent filter prevents bacteria and other particulates from entering stored water
- Hydrophobic microporous composite filter minimizes CO2 and organic contamination
- Compatible with Elga Biopure, Centra, Maxima, Medica, and Purelab series water purification units

---

### Specifications for this item

Brand Name	Elga
Item Weight	0.010 ounces
Model Number	LC136M2
Number of Items	1
Part Number	LC136M2

UNSPSC Code

41000000

