Immunohematology




ID-System

The Widest Range of Cards and Applications

Bio-Rad's selection of card configurations eliminates the "one-size fits all" setup of other platforms. Bio-Rad not only offers more card configurations than any other supplier, we also offer the right configurations for the way you want to test.

Bio-Rad's ID-Cards allow easy results interpretation, both manually and by automated means.

Simple and flexible for maximum productivity.






Blood Grouping ABO/D

Broad range of ID-Cards complying with guidelines for ABO/D blood grouping.


- Combined forward and reverse grouping
- Confirmation ID-Cards for patients and/or donors
- Single or double determination of RhD
- Combined with DAT for newborn testing
- Different monoclonal antibodies for A, B and D (DVI- and DVI+)
- Specific protocols for newborn testing



Phenotyping Rh/K

Several ID-Card profiles covering all needs and providing different set of cell lines.

- Different combinations, including C^w
- Rh-Subgroups + K or D
- Different monoclonal antibodies



Further Investigations

Specific solutions for A subgroups and weak/partial D's investigations.


- A subgroups; anti-A₁, anti-H
- Typing of partial RhD
- Verification of weak D by IAT



Antigen Testing

Extended antigen typing of your patients and/or donors with flexible solutions: ID-Card profiles or single antigen typing cards.


- Extended antigen typing
- Antigen profiles
- Single or combined antigen testing



Direct Antiglobulin Test (DAT)

Investigation of your DAT positive samples with specific solutions to further assess the risk of hemolysis.

- With polyspecific or monospecific AHG
- Further investigations: IgG dilution and IgG subclass differentiation
- ID-Cards to assess the potential severity of a positive DAT



Antibody Screening

The comprehensive offer of red blood cell kits and ID-Cards for antibody screening.

ID-Cards

- With polyspecific or monospecific AHG
- Different options available, including enzyme technique
- For type and screen

ID-Cells

- Test cells for antibody screening
- Native and enzyme treated cells
- Kits for patients and/or donors
- Expanded kit with Lu^a and Kp^a
- Screening of maternal samples after anti-D prophylaxis





Antibody Identification

The full solution for your antibody identification process.

ID-Cards

- With polyspecific or monospecific AHG
- Combined with enzyme technique

ID-Cells

- Panel cells for antibody identification
- Additional kits for complex cases
- Native and enzyme treated cells
- Screening of maternal samples after anti-D prophylaxis


1 2 3 4 5 6 Audiquities 1 1 2 - - - 2 -



Compatibility Tests

Wide offer of compatibility testing including a complete crossmatch in one step.

- With polyspecific or monospecific AHG
- With enzyme technique
- For complete crossmatch
- Autocontrol



Quality Control for Blood Group Serology

The most comprehensive and scalable solution for your quality control of your immunohematology testing.

- For ABO forward and/or reverse grouping, RhD and Rh/K phenotyping
- Antibody screening and identification
- Crossmatch
- DAT: Differentiation between IgG and complement; 1 cell IgG coated and 1 cell C3b(c/d) coated
- Anti-D \leq 0.05 IU/ml standardised against the WHO 01/572 International Standard
- RhD typing: direct testing and IAT



Antibody Titration

Dilution medium intended to be used for the titre determination of antibodies in patient and/or donor samples with the ID-System.

ID-Cards

- LISS/Coombs (rabbit anti-IgG and monoclonal anti-C3d)
- Coombs Anti-IgG (rabbit anti-IgG)
- NaCl, Enzyme Test and Cold Agglutinins (neutral gel)

Red Blood Cells

 Native ID-Cells including reverse, screening, identification, extended panels and your own laboratory cells



Miscellaneous Solutions

Preparing, treating and storing your red blood cells.

- LISS, bromelin
- Papain
- Kit for elution
- Solution for cell freeze and thaw
- For cells stabilization suitable for ID-Cards



Visit www.bio-rad.com/immunohematology for more information.



BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

The listed applications and associated reagents are only examples. Please refer to the ID-System catalog, to your local Bio-Rad contact point or www.bio-rad.com for further information about the complete portfolio.

It can be that some of the products mentioned in this document are not available in your country.

For further information, please contact the Bio-Rad office nearest you or visit our website at www.bio-rad.com/diagnostics



Bio-Rad Laboratories, Inc.

Diagnostics Group

Website www.bio-rad.com/diagnostics Australia +61 (2) 9914 2800 Austria +43 (0) 1 877 89 01 9 Belgium +32 (0) 3 710 53 00 Brazil +55 11 3065 7550 Canada +1 514 334 4372 China +86 21 6169 8500 Czech Republic +420 241 431 660 Denmark +45 44 52 10 00 Finland +358 9 804 22 00 France +33 (0)1 47 95 60 00 Germany +49 (0) 89 31884 393 Greece +30 210 7774396 Hong Kong +85 2 2789 3300 Hungary +36 1 459 6190 India +91 124 4029300 Israel +972 03 963 6025 Italy +39 024 94 86 600 Japan +81 3 6361 7070 Korea +82 080 007 7373 Mexico +52 (55) 5488 7670 The Netherlands +31 (0)318 540 666 New Zealand +64 (9)415 2280 Norway +47 23 38 41 30 Poland +48 22 331 99 99 Portugal +351 21 47 27 700 Russia +7 495 721 1404 Singapore +65 6415 3170 South Africa +27 11 442 8508 Spain +34 91 490 6580 Sweden +46 844 98053 Switzerland +41 (0) 61 717 9555 Taiwan +886 (2) 2578-7189 Thailand (662) 651 8311 United Kingdom +44 (0)1923 471301

Bulletin 12012050 © Bio-Rad 2021/11