


# SONATA 2 Cochlear Implant

## Product Catalogue




hearLIFE



# Safety Above All Else

## SONATA 2 Cochlear Implant

With SONATA 2, your patients will have peace of mind knowing they have an implant that delivers unmatched MRI safety\*, outstanding reliability, and exceptional hearing quality.



\* Recipients with a SONATA 2 Cochlear Implant may be safely MRI scanned at 0.2, 1.0, 1.5, and 3.0 Tesla following the conditions detailed in the Medical Procedures Manual. Visit [medel.com/isi](http://medel.com/isi) for Medical Procedures Manual and other safety information for MED-EL implants.

# SONATA 2 Cochlear Implant



Actual Size

## Mi1260 SONATA 2 Cochlear Implant

Includes Mi1250 Implant Template



## Mi1250 Implant Template

Included with implant. Used to evaluate the optimal position of the CI.

36894

## with FLEX Series Electrode Array

The softest and most flexible electrode array designed for structure preservation.

FLEXSOFT	38541
FLEX28	38545
FLEX26	38544
FLEX24	38543
FLEX20	38542

## with FORM Series Electrode Array

Ideal for malformed cochleae.  
Features SEAL technology for better control of CSF leakage.

FORM24	38547
FORM19	38546

## with CLASSIC Series Electrode Array

Features wave-shaped wires and double electrode contacts.

STANDARD	38538
MEDIUM	38539
COMPRESSED	38540






# Implant Tools

These additional tools can be used for implanting SONATA 2.



## PIN Drill Guide CL

▲ ♻️  
 Central Electrode Lead Stimulator template for planning and positioning. 36878



## Processor Template

▲ ♻️  
 Used for surgical planning to position the implant at the correct distance from the ear. 01557



## Skin Flap Gauge 6

▲ ♻️  
 Used to evaluate the thickness of the skin flap covering the implant. 03543



## Surgical Claw Angled

▲ ♻️  
 Used to position and insert the electrode array into the cochlea. 00284



## Micro Forceps Right-Angled

▲ ♻️  
 Used to grip, hold, and insert the electrode. 05777



## Micro Forceps Left-Angled

▲ ♻️  
 Used to grip, hold, and insert the electrode. 05778



## Micro Forceps Set 05761

▲ ♻️  
 Includes:  
 1 Micro Forceps Right-Angled 05777  
 1 Micro Forceps Left-Angled 05778



### Mi1250 Implant Template

■ Used to evaluate the optimal position of the CI. It is included with the SONATA 2 Cochlear Implant.


36894



### Fixation Clip

■ Implantable surgical accessory that provides additional fixation of the electrode lead on the bony bridge.

09917



### FENTEXmedical Forceps

▲ ♻️ These CI electrode insertion forceps have a longitudinal groove and fits all MED-EL electrodes. This surgical tool is not a MED-EL product and may be ordered directly at your local FENTEXmedical distributor. FENTEX Item Nr.: 172400FX



### Surgical Claw Straight

▲ ♻️ Alternative tool to Surgical Claw Angled (00284). Used to position and insert the electrode array into the cochlea.

07711



### Insertion Electrodes

■ Used to establish whether the cochlear lumen is obstructed or if it is freely accessible for the preferred electrode array.

IE FLEXSOFT	08255	IE STANDARD	08254
IE FLEX28	08348	IE MEDIUM	08256
IE FLEX26	36633	IE COMPRESSED	08258
IE FLEX24	08257		



### Insertion Test Device

■ Has a stopper at 18 mm and 5 pairs of opaque markers to help determine insertion depth up to a maximum of 18 mm.

02081





# Magnet Tool Kits

SONATA 2 features a secure removable implant magnet. The Magnet Tool Kit and Magnet Replacement Kit provide all the necessary components for magnet removal and replacement.

## Magnet Tool Kit

09734



Provides the surgical tools needed for magnet removal/replacement and should be used in conjunction with the Magnet Replacement Kit.

Includes:

- 1 Magnet Insertion Tool
- 1 Magnet Removal Tool

08399  
08400



## Magnet Insertion Tool



For insertion of Non-Magnetic Spacer and Replacement Magnet.

08399



## S-Vector Magnet Replacement Kit



Provides the magnet components needed for magnet removal/insertion and should be used in conjunction with the Magnet Tool Kit.

38234



Includes:

- 1 Non-Magnetic Spacer (Purple)\*
- 1 Replacement Magnet (Silver)\*

\*Not sold separately.

## Magnet Removal Tool



For removal of implant magnet and Non-Magnetic Spacer.

08400



MED-EL recommends replacing an implant magnet always with a magnet of the original type. Nevertheless, in case of a medical reason, selecting an alternative replacement magnet may be necessary. For further information please contact your MED-EL representative.



## MED-EL Offices Worldwide

### AMERICAS

#### Argentina

medel@medel.com.ar

#### Brazil

office@br.medel.com

#### Canada

officecanada@medel.com

#### Colombia

office-colombia@medel.com

#### Mexico

office-mexico@medel.com

#### United States

implants@medelus.com

### ASIA PACIFIC

#### Australia

office@medel.com.au

#### China

office@medel.net.cn

#### Hong Kong

office@hk.medel.com

#### India

implants@medel.in

#### Indonesia

office@id.medel.com

#### Japan

office-japan@medel.com

#### Malaysia

office@my.medel.com

#### Philippines

office@ph.medel.com

#### Singapore

office@sg.medel.com

#### South Korea

office@kr.medel.com

#### Thailand

office@th.medel.com

#### Vietnam

office@vn.medel.com

### EMEA

#### Austria

office@at.medel.com

#### Belgium

office@be.medel.com

#### Denmark

nordic@medel.com

#### Finland

office@fi.medel.com

#### France

office@fr.medel.com

#### Germany

office@medel.de

#### Italy

ufficio.italia@medel.com

#### Netherlands

office@nl.medel.com

#### Norway

norge@medel.com

#### Portugal

office@pt.medel.com

#### South Africa

customerserviceZA@medel.com

#### Spain

office@es.medel.com

#### Sweden

nordic@medel.com

#### Switzerland

office-ch@medel.com

#### Turkey

medelturkey@medel.com

#### United Arab Emirates

office@ae.medel.com

#### United Kingdom

customerservice@medel.co.uk

27747 r2.0