


Which infusion set is right for you?

STEEL

Reactions to plastic cannula / History of bent cannula / Pregnant (up to 2nd trimester) / Small child

90 degree insertion

Sure-T® infusion set
Simple steel needle infusion set




PLASTIC CANNULA

Less active / Alternative site usage / Poor dexterity

90 degree insertion

Quick-set® infusion set
Our most popular infusion set



Less active / Alternative site usage

90 degree insertion


mio® infusion set
Convenient "all in one" infusion set



Active / Lean / Muscular / Pregnant (from 2nd trimester)

30-45 degree insertion

Silhouette® infusion set
Flexible insertion angle infusion set



Low/Average BMI

6mm

Wear insulin pump close to infusion site
18" or 23"

Wear insulin pump away from infusion site
32"

Average/High BMI

8mm

Wear insulin pump close to infusion site
23"

Wear insulin pump away from infusion site
32"

10mm

32"

Low/Average BMI

6mm

Wear insulin pump close to infusion site
18" or 23"

Wear insulin pump away from infusion site
32" or 43"

Average/High BMI

9mm

Wear insulin pump close to infusion site
23"

Wear insulin pump away from infusion site
32" or 43"

Low/Average BMI

6mm

Wear insulin pump close to infusion site
18" or 23"

Wear insulin pump away from infusion site
32"

Average/High BMI

9mm

32"

Low/Average BMI

13mm

Wear insulin pump close to infusion site
18" or 23"

Wear insulin pump away from infusion site
32" or 43"

Average/High BMI

17mm

Wear insulin pump close to infusion site
23"

Wear insulin pump away from infusion site
32" or 43"