



Pressure drop for all arterial HLS Cannulae 3/8"



Pressure drop for all venous HLS Cannulae 3/8"



Product order details arterial HLS cannulae

Type	Outer diameter	Insertion length	Side holes	Perforation length	Connector	Bioline Coating
PAS 1315	13 Fr (4.3 mm)	15 cm	2	1 cm	3/8" LL	BE-PAS 1315
PAS 1515	15 Fr (5.0 mm)	15 cm	2	1 cm	3/8" LL	BE-PAS 1515
PAS 1715	17 Fr (5.7 mm)	15 cm	2	1 cm	3/8" LL	BE-PAS 1715
PAS 1915	19 Fr (6.3 mm)	15 cm	2	1 cm	3/8" LL	BE-PAS 1915
PAS 2115	21 Fr (7.0 mm)	15 cm	2	1 cm	3/8" LL	BE-PAS 2115
PAS 2315	23 Fr (7.7 mm)	15 cm	2	1 cm	3/8" LL	BE-PAS 2315
PAL 1523	15 Fr (5.0 mm)	23 cm	2	1 cm	3/8" LL	BE-PAL 1523
PAL 1723	17 Fr (5.7 mm)	23 cm	2	1 cm	3/8" LL	BE-PAL 1723
PAL 1923	19 Fr (6.3 mm)	23 cm	2	1 cm	3/8" LL	BE-PAL 1923
PAL 2123	21 Fr (7.0 mm)	23 cm	2	1 cm	3/8" LL	BE-PAL 2123
PAL 2323	23 Fr (7.7 mm)	23 cm	2	1 cm	3/8" LL	BE-PAL 2323

One cannula per carton

Product order details venous HLS cannulae

Type	Outer diameter	Insertion length	Side holes	Perforation length	Connector	Bioline Coating
PVS 1938	19 Fr (6.3 mm)	38 cm	12	10 cm	3/8"	BE-PVS 1938
PVS 2138	21 Fr (7.0 mm)	38 cm	12	10 cm	3/8"	BE-PVS 2138
PVS 2338	23 Fr (7.7 mm)	38 cm	16	10 cm	3/8"	BE-PVS 2338
PVS 2538	25 Fr (8.3 mm)	38 cm	20	10 cm	3/8"	BE-PVS 2538
PVL 2155	21 Fr (7.0 mm)	55 cm	20	20 cm	3/8"	BE-PVL 2155
PVL 2355	23 Fr (7.7 mm)	55 cm	20	20 cm	3/8"	BE-PVL 2355
PVL 2555	25 Fr (8.3 mm)	55 cm	24	20 cm	3/8"	BE-PVL 2555
PVL 2955	29 Fr (9.7 mm)	55 cm	32	20 cm	3/8"	BE-PVL 2955

One cannula per carton