



PAGINA 0

**CREATA PENTRU A
PUTEA SEMNA
ELECTRONIC
DOCUMENTUL**



MEDTRONIC FULL-WALL TECHNOLOGY:
Encapsulated flat-wire braid enables thinner, more robust walls without compromising support, retention, or torque response.

BALANCED PERFORMANCE THAT'S AHEAD OF THE CURVE

The Launcher™ guide catheter offers you a blend of flexibility, support and visualisation — the capability you need to respond to your challenging cases.

Optimal Balance Between Flexibility and Support

- Flexible distal segment enables you to engage the guide for backup support
- Full-Wall Technology construction provides kink resistance and stable torque control
- Supportive secondary curve for backup support when delivering larger profile devices

Enhanced Visualisation

- Large 8 F I.D. (0.090 inches) allows for enhanced vessel visibility and device passage
- Vest-Tech Nylon outer jacket provides excellent hub to tip radiopacity
- Radiopaque marker band facilitates tip orientation for ostial engagement

Full Range of Traditional and Specialty Curves

- Extensive collection of traditional and specialty coronary curves available
- 90 cm length catheters available for your specialty cases

LAUNCHER™ 8 F GUIDE CATHETERS

LEFT CORONARY CURVES					Backup Support Left
Left Standard (JL)	Short Left (SL)	JCL	Amplatz Left	Short Amplatz Left	EBU
LA8JL30	LA8SL30	LA8JCL30	LA8AL75	LA8SAL75	LA8EBU30
LA8JL35	LA8SL35	LA8JCL35	LA8AL10	LA8SAL10	LA8EBU35
LA8JL40	LA8SL40	LA8JCL35D*	LA8AL10D*	LA8SAL15	LA8EBU35D*
LA8JL40D*	LA8SL45	LA8JCL40	LA8AL15	LA8SAL20	LA8EBU375
LA8JL45	LA8SL50	LA8JCL40D*	LA8AL20	LA8SAL25	LA8EBU40
LA8JL50	LA8SL60	LA8JCL45	LA8AL20D*	LA8SAL30	LA8EBU40D*
LA8JL60		LA8JCL45D*	LA8AL25	LA8SAL40	LA8EBU45
LA8JL70		LA8JCL50	LA8AL30		LA8EBU50
LA8JL80			LA8AL40		
			LA8AL40D*		
			LA8AL50		

RIGHT CORONARY CURVES					
Right Standard	Short Right	JCR	Amplatz Right	Shepherd's Crook Right	Backup Support Right
LA8JR30	LA8SR30	LA8JCR35	LA8AR10	LA8SCR35	LA8RBU35
LA8JR35	LA8SR35	LA8JCR40	LA8AR20	LA8SCR40	LA8RBU40
LA8JR35D	LA8SR40	LA8JCR45	LA8ALR12	LA8SCR50	LA8RBU40D
LA8JR40	LA8SR45	LA8JCR50			
LA8JR40D	LA8SR50				ECR
LA8JR45	LA8SR60				LA8ECR35
LA8JR50					LA8ECR40
LA8JR60			Short Amplatz Right		LA8ECR45
			LA8SAR10		
			LA8SAR20		
					3-Dimensional
					LA83DRC


OTHER CURVES					
Multipurpose	MAC (Multi-Aortic)	Bypass	Champ	Peripheral	NOTO
LA8HSI	LA8MAC30	LA8IMA	LA8CHAMP05	LA8CARL3D*	LA8NOTO
LA8HSII	LA8MAC35	LA8IMAD*	LA8CHAMP10	LA8CARLD*	
LA8RMSHII	LA8MAC375	LA8LCB	LA8CHAMP15	LA8CHAMP15K+	
LA8HSIII	LA8MAC40	LA8LCBD*	LA8CHAMP20	LA8CHAMP20K+	
LA8MB1	LA8MAC45	LA8RCB	LA8CHAMP25	LA8CHAMP25K+	
LA8MB1D*	LA8MAC50	LA8RCBD*	LA8CHAMP30	LA8IMAK+	
LA8MB2		LA8RCBII	LA8CHAMP35	LA8JR40K+	
LA8MB2D*			LA8CHAMP40	LA8MP1K+	EL GAMAL
LA8MP1				LA8MPHK+	LA8ELGAMAL
LA8MPST				LA8PK1W++	LA8ELGAMAL2
				LA8RDCK+	LA8ELGAMAL3
				LA8SCR40K+	
				LA8STRD*	
				LA8STRK+	
				LA8SW1	
				LA8SW1D*	

D* = 90 cm shaft length Shaft Diameter: 8 F K+ = 55 cm shaft length Inner Diameter: 0.090 in. W++ = 47 cm shaft length
 To order curves with sideholes, please add "SH" to end of model number 5 packs available in select curves

Indications for Use: The Medtronic Guiding Catheter is designed to provide a pathway through which therapeutic devices are introduced. The guiding catheter is intended to be used in the coronary or peripheral vascular system.
Important Information: Refer to *Instructions for Use* supplied with the device for indications, warnings, precautions, and adverse effects.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
 The curve recommendations in this guide are not meant to replace patient-specific clinical judgement.



Medtronic
 3576 Unocal Place
 Santa Rosa, CA 95403
 USA
 Tel: 707.525.0111

LifeLine Customer Support
 Tel: 877.526.7890
 Tel: 763.526.7890

Product Services
 Tel: 888.283.7868
 Fax: 800.838.3103

Europe
Medtronic Intl. Trading SARL
 Route du Molliau 31
 Case Postale
 CH-1131 Tolochenaz
 Switzerland
 Tel: 41.21.802.7000
 Fax: 41.21.802.7900

Canada
Medtronic of Canada Ltd.
 99 Hereford Street
 Brampton, Ontario L6Y 0R3
 Canada
 Tel: 905.460.3800
 Fax: 905.460.3998
 Toll-free: 800.268.5346

Asia Pacific
Medtronic Intl. Ltd.
 49 Changi South Avenue 2
 Nasaco Tech Centre
 Singapore 486056
 Singapore
 Tel: 65.6436.5000
 Fax: 65.6776.6335

Latin America
Medtronic
 9850 Doral Blvd., 4th Floor
 Miami, FL 33178
 USA
 Tel: 786.709.4200
 Fax: 786.709.4244

medtronic.com