



HOPKINS® Telescopes

Diameter 4 mm

For use with Hysteroscopes



26105 BA




26105 BA

HOPKINS® Forward-Oblique Telescope 30°, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red

Hysteroscopes see page 25

For use with Hysteroscopes and Resectoscopes



26105 FA



26105 FA

HOPKINS® Telescope 12°, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: black

Hysteroscopes see page 26

Resectoscopes see pages 49-52

Containers for Sterilization and Storage of Telescopes see catalog HYGIENE

Fiber Optic Light Cables

for Cold Light Fountains



Fiber Optic Light Cable with Straight Connector

- 495 NE **Fiber Optic Light Cable**, diameter 4.8 mm, length 300 cm
- 495 NCS **Fiber Optic Light Cable**, extremely heat-resistant, diameter 4.8 mm, length 250 cm
- NEW** 495 NCSC **Fiber Optic Light Cable**, extremely heat-resistant, with safety locking device, diameter 4.8 mm, length 250 cm
- 495 NB **Fiber Optic Light Cable**, diameter 4.8 mm, length 180 cm
- 495 ND **Fiber Optic Light Cable**, diameter 3.5 mm, length 300 cm
- 495 NA **Fiber Optic Light Cable**, diameter 3.5 mm, length 230 cm
- NEW** 495 NAC **Fiber Optic Light Cable**, with safety locking device, diameter 3.5 mm, length 230 cm
- 495 NL **Fiber Optic Light Cable**, diameter 3.5 mm, length 180 cm
- 495 NTA **Fiber Optic Light Cable**, diameter 2.5 mm, length 230 cm
- 495 NT **Fiber Optic Light Cable**, diameter 2.5 mm, length 180 cm
- NEW** 495 TIP **Fiber Optic Light Cable**, highly heat resistant, diameter 4.8 mm, length 300 cm, optimized for 3D TIPCAM®1 and TIPCAM®



Fiber Optic Light Cable with 90° Deflection to the Light Source

- 495 NWL **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 300 cm
- 495 NWM **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 230 cm
- 495 NW **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 3.5 mm, length 180 cm
- 495 NTX **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 2.5 mm, length 230 cm
- 495 NTW **Fiber Optic Light Cable**, with 90° deflection to the cold light fountain on the fountain side, diameter 2.5 mm, length 180 cm



495 NVC



495 NV/NVL/NVB

Fiber Optic Light Cable with 90° Deflection to the Instrument

- NEW** 495 NVC **Fiber Optic Light Cable**, with 90° deflection to the instrument, very narrow radius of curvature, diameter 4.8 mm, length 300 cm
- 495 NVB **Fiber Optic Light Cable**, with 90° deflection to the instrument, diameter 4.8 mm, length 300 cm
- 495 NVL **Fiber Optic Light Cable**, with 90° deflection to the instrument, diameter 3.5 mm, length 300 cm
- 495 NV **Fiber Optic Light Cable**, with 90° deflection to the instrument, diameter 3.5 mm, length 230 cm



Light Adaptor, angled 90°


- 495 EW **Light Adaptor**, angled 90°, diameter 4.8 mm, free rotatable, to connect with standard telescopes

Adaptor for connecting KARL STORZ fiber optic light cables with endoscopes and light sources from other manufacturers, see page TP 132

Adaptors

for Light Sources and Light Cables



HOPKINS® Telescopes with universal connector for light cables



Manufacturer	I = Instrument CF= Cold Light Fountain LC= Light Cable	KARL STORZ				
		Light Cable	Cold Light Fountain	Former Instrument	Current Instrument	Special Applications
ACMI	I	495 W	–	–	–	ACMI connector to ACMI connector 495 Y
	CF	495 PA	–	–	–	
	LC	–	495 C	–	–	
Downs	CF	495 PA	–	–	–	
Dyonics	CF	495 Q	–	–	–	
Effner	CF	495 Q	–	–	–	
Fuji	CF	495 FU	–	–	–	
KLI	LC	–	495 C	–	–	
	CF	495 PA	–	–	–	
Machida	CF	495 SM	–	–	–	
Olympus	CF	495 SON	–	–	–	Olympus Techno Fountain for KARL STORZ Cable 81495 OT, Olympus Light Cable for KARL STORZ Techno Fountain 81450/81495 OL
	LC	–	495 D	–	–	
Pentax	CF	495 T	–	–	–	
	LC	–	495 D	–	–	
HSW	CF	495 Q	–	–	–	
	LC	–	495 E	495 N	495 N	
Stryker	I	495 ST	–	–	–	
	CF	495 PA	–	–	–	
Richard Wolf	I	495 X	–	–	–	
	CF	495 Q	–	–	–	
	LC	–	495 E	495 M	–	
Olympus/ Winter & Ibe	CF	495 S	–	–	–	495 S: The use of KARL STORZ fiberscopes with integrated fiber optic light cable and KARL STORZ light cables with non-graduated distal connection (diameter 7 mm) requires an additional adaptor 495 PWS to adaptor 495 S
	LC with attached lens	–	–	495 GO	495 GO	

Home page › Online Catalog

26050SLK



Item no: 26050SLK

Resectoscope Sheath, 26 Fr.

Quantity:

–

1

+

Add to My Quote List

Scope of delivery


- 26050XA Inner Sheath, rotatable, for 26050SL (1)
- 26050SL Resectoscope Sheath, 26 Fr. (1)

Technical details	Documents & Films	Product group / field of application
-------------------	-------------------	--------------------------------------

System size	26 Fr.
profile shape	Round
Working length	185 mm
Distal end shape	Oblique beak, with ceramic insulation
For telescope direction of view	12°
For telescope size	4.0 mm
Color code	Yellow
Flow property	Continuous flow
Flow function	Inflow and outflow
Sheath rotation	Rotatable inner sheath
Connection mechanism	Rotary lock

[Home page](#) › [Online Catalog](#)

260400C



Item no: 260400C

Standard Obturator

Quantity:

—

1

+

Add to My Quote List

Technical details	Documents & Films	Product group / field of application
-------------------	-------------------	--------------------------------------

System size	26 Fr.
Distal end shape	Blunt
Color code	Yellow

Type of usage Reusable 5

STORZ
KARL STORZ — ENDOSKOPE

Customer Service

About us



Coagulating Electrode, 24 Fr.

$$-1+$$


On the offer list



Documents & Films

Product group / field of application

8



Item no: 26050NK

Coagulating Electrode, 26 Fr.

Quantity:

–

1

+

On the offer list


Technical details	Documents & Films	Product group / field of application
-------------------	-------------------	--------------------------------------

Number of poles	Unipolar
Diameter distal end	5.0 mm
Shape of distal end	Ball end
For telescope direction of view	12°
Type of usage	Reusable
Color code	Yellow
System size	26 Fr.
Packaging unit	6 pc.

Vaporization Electrodes

One-Stem Electrodes with Stabilizers, for Working Elements 27050 C/D/E

For use with 24/26 Fr. resectoscope sheaths



27050 DG

The cutting loops are delivered with a wire diameter of 0.35 mm.

Distal Tip	24/26 Fr. color code: yellow	Instrument Description
	27050 DG	Cutting Loop, diameter 0.8 mm
	27050 VG	Spike Electrode, diameter 5 mm
	27050 RG	Roller Electrode, diameter 5 mm
	27050 RK	Roller Electrode, diameter 3 mm
	27050 SG	VaporCut® Electrode, band width 1.2 mm
	27050 BG	VaporCut® Electrode, band width 0.6 mm
	27050 KG	Roller Cutting Electrode, single use recommended
	27050 WG	VaporCut® Electrode



280

280

Protection Tube, for sterilization and storage of electrodes, loops, curettes and knives

Note: Electrodes, in sterile packaging, are also available for single use. Please contact:



High Frequency Cords see page 73
High Frequency Surgical Units see chapter 15, UNITS

4-964

Configuration Planning

ENDOMAT® SELECT

The Choice is Yours



ENDOMAT® SELECT

The Choice is Yours

ENDOMAT® SELECT is a single roller pump for the irrigation or suction of fluids during surgical and diagnostic procedures.

The pump can be equipped with various software modules so that the unit can be used for multiple disciplines in various rooms according to the respective indication.

The unit offers various software packages which the user can select and combine according to individual needs.




☐ UP210

ENDOMAT® SELECT SCB, with integrated SCB module, suction and irrigation pump, incl. mains cord, power supply 100-240 VAC, 50/60 Hz

When placing an initial order, please tick the box for UP210. Furthermore, at least one software module must be included in the order. Please tick the appropriate box in the document.

Please do not tick the box if a software upgrade for an existing unit is requested.

The following software packages are available for ENDOMAT® SELECT:
Clicking on one of the fields highlighted in color takes you to the relevant page.



The following tubing sets are available for ENDOMAT® SELECT:

For pressure-controlled applications such as, for example, hysteroscopy, cystoscopy or interventions in arthroscopy or spine surgery:

031523-10* **Tubing Set**, irrigation, pressure-controlled, STERILE  sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210

Procedures that use a pressure-controlled tubing set bear the adjacent mark.



For flow-controlled applications such as, for example, laparoscopy, thoracoscopy, proctology:

- 031524-10* **Tubing Set**, irrigation, flow-controlled, STERILE sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210
- Procedures that use a flow-controlled tubing set bear the adjacent mark.



For hysteroscopic tissue ablation with the IBS® shaver system or the removal of tissue chips during resection (TUR-B/P):

- 030647-10* **Tubing Set**, suction, STERILE direct suction, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210
- Procedures that use a tubing set for direct suction bear the adjacent mark.



For use in combination with the CALCUSON lithotripsy unit in urology:

- 031647-10* **Tubing Set**, suction, STERILE via bottle, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210
- Procedures requiring a suction bottle bear the adjacent mark.



For use with the CLEARVISION® functionality in ENT/NEURO:

- 031529-10* **Tubing Set**, irrigation, CV, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210
- Procedures that use a CLEARVISION® functionality bear the adjacent mark.

SURGERY

In the daily surgical routine in the OR, the ENDOMAT® SELECT delivers a powerful yet low-noise irrigation of fluids during interventions related to the following specialties:


- Laparoscopy
- Thoracoscopy
- Proctology

Software:

- ☐ UP601 **SURGERY Software**, license, allows selection of the procedures “LAP”, “THOR” and “PROCTO”
- ☒ **Promotion Pack (031524-03*)** STERILE **Tubing Set**, irrigation, FC, sterile, for single use, package of 3

Further accessories:

- 031524-10* **Tubing Set**, irrigation, flow-controlled, STERILE sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® HAMOU® 26331120-1 and ENDOMAT® SELECT UP210



(OFFICE) HYSTEROSCOPY

ENDOMAT® SELECT provides pressure-controlled dilation of the uterine cavity for diagnostic and operative hysteroscopy in the office setting or the clinical outpatient department.


Software:

☐ UP602

HYSTEROSCOPY Software, license, allows selection of the procedure “HYS”

☒ **Promotion Pack (031523-03*)** STERILE

Tubing Set, irrigation, PC, sterile, for single use, package of 3



Further accessories:

031523-10*

Tubing Set, irrigation, pressure-regulated, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® HAMOU® 26331120-1 and ENDOMAT® SELECT UP210 STERILE

Hysteroscopic Tissue Ablation with the Intrauterine BIGATTI Shaver IBS®

In conjunction with the IBS® shaver, the ENDOMAT® SELECT can be used for the controlled suction of resected polyps, myomas or placental remnants.


Software:

☐ UP603

IBS® SHAVER, software, license, allows selection of the procedure “IBS®”

☒ **Promotion Pack (030647-03*),** STERILE

Tubing Set, suction, DS, sterile, for single use, package of 3



Mandatory accessories:

030647-10*

Tubing Set, suction, direct suction, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210 STERILE

20701070

Control Cable, connectors 2x LEMO 5-pol 0°, length 100 cm, for transmission of footswitch control signal between ENDOMAT® SELECT UP210 and UNIDRIVE® S III 20701020



UROLOGY

ENDOMAT® SELECT accommodates the wide range of applications in the field of urology, allowing seamless integration into existing systems and providing pressure-regulated irrigation for diagnostic and operative interventions.

Software:

☐ UP604

UROLOGY Software, license, allows selection of the procedures “CYST”, “RES”, “URS”, “CALCUSON”, and “PCN”

Please select the required Promotion Pack: STERILE


☒ **Promotion Pack (031523-03*) Tubing Set**, irrigation, PC, sterile, for single use, package of 3

or

☐ **Promotion Pack (030647-03*) Tubing Set**, suction, DS, sterile, for single use, package of 3

or

☐ **Promotion Pack (031647-03*) Tubing Set**, suction, BS, sterile, for single use, package of 3



Further accessories:

031523-10*

Tubing Set, irrigation, pressure-controlled, sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210 STERILE

030647-10*

Tubing Set, suction, direct suction, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210 STERILE

031647-10*

Tubing Set, suction, via bottle, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210 STERILE

20014130

One-Pedal Footswitch, digital, one-stage

20701070

Control Cable, connectors 2x LEMO 5-pol 0°, length 100 cm, for transmission of footswitch control signal between ENDOMAT® SELECT UP210 and UNIDRIVE® S III 20701020 such as between ENDOMAT® SELECT UP210 and CALCUSON 27610020



ARTHROSCOPY

In the field of arthroscopy, the ENDOMAT® SELECT regulates pressure-controlled irrigation of fluids to maintain joint distension.

Software:

☐ UP605

ARTHROSCOPY Software, license, allows selection of the procedures “KNEE”, “HIP”, “SHOULDER” and “SMALL JOINTS”,
for use with ENDOMAT® SELECT UP210

☒ **Promotion Pack (031523-03*)** STERILE

Tubing Set, irrigation, PC, sterile, for single use, package of 3

Further accessories:


- 031523-10*

Tubing Set,
irrigation, pressure-regulated, sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210

STERILE
- 20014130

One-Pedal Footswitch, digital, one-stage
- UP006

Control Cable,
connectors 2x LEMO 5-pol 0°, length 100 cm, for transmission of footswitch control signal between ENDOMAT® SELECT UP210 and UNIDRIVE® S III ARTHRO SCB 28723020-1



ENT/NEURO

ENDOMAT® SELECT includes the following software module for irrigation of the lens

Software:

☐ UP606

ENT/NEURO Software, allows selection of the procedure “CLEARVISION®”

☒ **Promotion Pack (031529-03*)** STERILE

Tubing Set, irrigation, CV, sterile, for single use, package of 3, for use with KARL STORZ ENDOMAT® SELECT UP210


Mandatory accessories:

- 031529-10*

Tubing Set,
irrigation, CV, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210

STERILE
- 20014230

One-Pedal Footswitch, digital, two-stages



SPINE SURGERY


ENDOMAT® SELECT offers the following software module for the irrigation of fluids during interventions in spine surgery:

Software:

☐ UP607 **SPINE Software**, license, allows selection of the procedures “LUMBAR” and “THORACIC”

☒ **Promotion Pack (031523-03*)** STERILE 

Tubing Set, irrigation, PC, sterile, for single use, package of 3



Further accessories:

031523-10* **Tubing Set**, STERILE 

irrigation, pressure-controlled, sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210

ADVANCED

In order to activate extended settings and additional functions, the “Advanced Package” can be selected as an additional* software package.

Software:

☐ UP610 **Advanced Package**, software, license

* The advanced package only functions when at least one other software package is available.

The following features are available with the installation of ADVANCED Software UP610:

- Adjustment of pressure (100/300/500 mmHg) during SURG applications
- Adjustment of inflow up to 3500 ml/min during SURG applications
- Adjustment of inflow (200/400/600 ml/min) during HYS, SPINE & URO applications
- Adjustment of inflow (1500/2000/2500 ml/min) during ARTHRO applications
- Duration of BOOST during ARTHRO & URO applications
- Adjustment of pump position (height) to improve pump performance
- Adjustment of volume of the button tones
- Adjustment of brightness
- SCB Integration



Bottle Cap,
for suction bottles 1.5 l and 5 l, sterilizable



Optional Accessories

for Suction and Irrigation Systems

			for use with		
			ENDOMAT® n. HAMOU® SCB 26 3311 01-1	ENDOMAT® LC SCB 20 3303 01-1	DUOMAT® 20 3210 08
	26 3100 50	Irrigation Liquid Bottle , 1 l, sterilizable	-	-	•
	20 3000 50	Suction Bottle , 5 l, sterilizable	•	-	•
	20 3000 34	Bottle Cap , for suction bottles 1.5 l and 5 l, sterilizable	•	-	•
	20 3000 32	Bottle Stand , for suction bottle 5 l	•	-	•
	20 3000 33	Bottle Stand Holder , for Bottle Stand 20 3000 32	•	-	•
	20 0142 30	One-Pedal Footswitch , digital, two-stage, for activating a higher flow for improved visibility	-	•	-
	20 0141 30	One-Pedal Footswitch , digital, one-stage, for irrigation	-	•	-
	20 3004 82	Connector Set , for ENDOMAT® LC, for use with Silicone Tubing Sets 20 3303 40, 20 3303 41 and 20 3303 43	-	•	-
	20 3303 82	Pump Tube , sterilizable, package of 25	-	•	-

SUCTION AND IRRIGATION SYSTEMS

ENDOMAT® SELECT VET

Optional Accessories

For use with program SURG

	031524-10*	Tubing Set for Irrigation , FC, sterile, for single use, package of 10, for use with HAMOU® ENDOMAT® SCB 26 3311 20-1 and ENDOMAT® SELECT	STERILE 
	UP 007	Tubing Set , irrigation, FC, reusable, sterilizable	

For use with programs ART SA and ART LA

	031523-10*	Tubing Set for Irrigation , PC, sterile, for single use, package of 10, for use with HAMOU® ENDOMAT® SCB 26 3311 20-1 and ENDOMAT® SELECT	STERILE 
	UP 008	Tubing Set , irrigation, PC, reusable, sterilizable	

For use with program SURG Direct Suction

	030647-10*	Tubing Set for Suction , DS, sterile, for single use, package of 10, for use with ENDOMAT® SELECT	STERILE 
	UP 009	Tubing Set , irrigation, DS, reusable, sterilizable	

For use with program SURG Bottle Suction








	031647-10*	Tubing Set for Suction , BS, sterile, for single use, package of 10, for use with ENDOMAT® SELECT	STERILE 
	UP 010	Tubing Set , irrigation, BS, reusable, sterilizable	



SUCTION AND IRRIGATION
SYSTEMS

UNIT SIDE

PATIENT SIDE



SUCTION AND IRRIGATION SYSTEMS

ENDOMAT® SELECT

For intraoperative cleaning of the front lens



UP220

ENDOMAT® SELECT, suction or irrigation pump

UP606

ENT Software, license, allows selection of the procedure "CLEARVISION", for use with ENDOMAT® SELECT UP210 and UP220



UF102

One-Pedal Footswitch, two-stage



031529-10

Tubing Set for Irrigation, CV, sterile, for single use, package of 10