

1.2 0.018" (0.46 mm) Guidewires

Victory™

Peripheral Guidewire

CROSSING

- Designed to break through and cross resistant lesions, offering increasing tip loads from 12g to 30g for control and a range of pushability. Shapeable distal tip. Hydrophilic distal coating. Indicated for intraluminal crossing of chronic total occlusions

| Order Number | Diameter (in) | Diameter (mm) | Length (cm) | Tip Length* (cm) | Tip Load (g) |
|-----------------|---------------|---------------|-------------|------------------|--------------|
| H74939232195120 | 0,018 | 0,46 | 195 | 2,5 | 12 |
| H74939232300120 | 0,018 | 0,46 | 300 | 2,5 | 12 |
| H74939232195180 | 0,018 | 0,46 | 195 | 2 | 18 |
| H74939232300180 | 0,018 | 0,46 | 300 | 2 | 18 |
| H74939233195250 | 0,018 | 0,46 | 195 | 2 | 25 |
| H74939233300250 | 0,018 | 0,46 | 300 | 2 | 25 |
| H74939233195300 | 0,018 | 0,46 | 195 | 2 | 30 |
| H74939233300300 | 0,018 | 0,46 | 300 | 2 | 30 |

*2 cm or 2.5 cm radiopaque tip length

ZIPWire™


Hydrophilic Guidewire

CROSSING

- Secure control when sub-selecting distal and angulated branch vessels. Excellent visibility for accurate positioning. Excellent kink-resistance. Indicated for enhanced steerability and control for maneuvering through difficult anatomy

| Order Number | Diameter (in) | Diameter (mm) | Length (cm) | Tip Taper (cm) | Tip Shape |
|--------------|---------------|---------------|-------------|----------------|-----------|
| M00146231B1 | 0,018 | 0,46 | 150 | 3 | Angled |
| M00146232B1 | 0,018 | 0,46 | 180 | 3 | Angled |
| M00146235B1 | 0,018 | 0,46 | 150 | 3 | Straight |
| M00146236B1 | 0,018 | 0,46 | 180 | 3 | Straight |

Packaged 5/box



1
Guidewires

2
Balloon
Dilatation
Catheters

3
Specialty
Balloons

4
Catheters

5
Accessories

6
Stents

7
Drug-Fluting
Technologies

8
Peripheral
Thrombectomy

9
Venous
Thrombo-
embolism

10
Peripheral
Atherectomy

11
IVUS
Intravascular
Ultrasound System

12
Embolic
Protection

13
Peripheral
Embolisations

14
SIRT

15
Non-Vascular
Interventions

16
Ablation



1.3 0.014" (0.36 mm) Guidewires

1
Guidewires

V-14™ ControlWire™

CROSSING

Guidewire with ICE™ Hydrophilic Coating

- Offers the choice of two tip loads options with an excellent balance between rail support and distal flexibility. Indicated for crossing stenosed lesions and occlusions

| Order Number | Diameter (in) (mm) | | Length (cm) | Tip Taper* | Tip Shape |
|-----------------|-----------------------|------|----------------|------------|-----------|
| H74939216718210 | 0.014 | 0.36 | 182 | Short | Straight |
| H74939216718220 | 0.014 | 0.36 | 182 | Short | Angled |
| H74939216730010 | 0.014 | 0.36 | 300 | Short | Straight |
| H74939216730020 | 0.014 | 0.36 | 300 | Short | Angled |
| H74939216918210 | 0.014 | 0.36 | 182 | Long | Straight |
| H74939216918220 | 0.014 | 0.36 | 182 | Long | Angled |
| H74939216930010 | 0.014 | 0.36 | 300 | Long | Straight |
| H74939216930020 | 0.014 | 0.36 | 300 | Long | Angled |

*2 cm radiopaque tip
Packaged 1/box

2
Balloon
Dilatation
Catheters

3
Specialty
Balloons

4
Catheters

5
Accessories

6
Stents

Victory™

CROSSING

Guidewires

- Designed to break through and cross resistant lesions, offering increasing tip loads from 12g to 30g for control and a range of pushability. Shapeable distal tip. Hydrophilic distal coating. Indicated for intraluminal crossing of chronic total occlusions

| Order Number | Diameter (in) (mm) | | Length (cm) | Tip Length* (cm) | Tip Load (g) |
|-----------------|-----------------------|------|----------------|---------------------|-----------------|
| H74939230195120 | 0.014 | 0.36 | 195 | 2.5 | 12 |
| H74939230300120 | 0.014 | 0.36 | 300 | 2.5 | 12 |
| H74939230195180 | 0.014 | 0.36 | 195 | 2.5 | 18 |
| H74939230300180 | 0.014 | 0.36 | 300 | 2.5 | 18 |
| H74939231195250 | 0.014 | 0.36 | 195 | 2 | 25 |
| H74939231300250 | 0.014 | 0.36 | 300 | 2 | 25 |
| H74939231195300 | 0.014 | 0.36 | 195 | 2 | 30 |
| H74939231300300 | 0.014 | 0.36 | 300 | 2 | 30 |

*2 cm or 2.5 cm radiopaque tip length
Packaged 1/box

7
Drug-Fluting
Technologies

8
Peripheral
Thrombectomy

9
Venous
Thrombo-
embolism

10
Peripheral
Atherectomy

Platinum Plus™

DELIVERY SUPPORT

Guidewire with Glidex™ Hydrophilic Coating

- Extra support wire for device delivery. Distal platinum coils for good visualization, 8 cm Glidex distal hydrophilic coating. Indicated for negotiation of tortuous anatomy and contralateral approaches including smaller vessels and those with short distal purchase

Platinum Plus™ 0.014" (0.36 mm) with Glidex™ Coating

| Order Number | Diameter (in) (mm) | | Length (cm) | Tip Taper (cm) | Tip Shape |
|--------------|-----------------------|------|----------------|-------------------|-----------|
| M001467400 | 0.014 | 0.36 | 180 | 9 | Straight |
| M001467120 | 0.014 | 0.36 | 300 | 7 | Straight |

*5 cm Radiopaque tip length for 180 cm
3 cm Radiopaque tip length for 300 cm
Packaged 1/box



11
IVUS
Intravascular
Ultrasound System

12
Embolic
Protection

13
Peripheral
Embolisations

14
SIRT

15
Non-Vascular
Interventions

16
Ablation



Ordering Information - 21
Peripheral Interventions

Boston
Scientific