

## MultiCap “Cup” System cap

The Cup MultiCap has one key feature. It allows the user to use standard 10 mm diameter Cup electrodes which makes this cap particularly flexible. The specially developed cup electrode mountings of the Cup MultiCap not only have the advantage that electrodes used for other monitoring systems can also be used with this cap, it also has all the benefits of the Base MultiCap and the Soft MultiCap. Particularly for new born or premature babies, this cap is a valuable aid as it eliminates the need to attach electrodes by adhesive.



### MultiCap “Cup” System cap with 21/15 electrodes holders (without electrodes)

Circumference	Comments	Item number	PU
I3 24-28 cm	brown Premature babies	MCI-3KK-14RK00-K	1 pc.
I2 28-32 cm	pink Babies	MCI-2KK-14RK00-K	1 pc.
I1 32-36 cm	lightblue Babies	MCI-1KK-14RK00-K	1 pc.
XS 36-42 cm	green Toddlers	MC-XSKK-20RK00-K	1 pc.
S 42-48 cm	yellow Children	MC-SKK-20RK00-K	1 pc.
MS 45-51 cm	red/yellow Teenager	MC-MSKK-20RK00-K	1 pc.
M 48-54 cm	red Teenager	MC-MKK-20RK00-K	1 pc.
LM 51-57 cm	blue/red Adult	MC-LMKK-20RK00-K	1 pc.
L 54-60 cm	blue Adult	MC-LKK-20RK00-K	1 pc.
XLL 57-63 cm	white/blue Adult	MC-XLLKK-20RK00-K	1 pc.
XL 60-66 cm	white Adult	MC-XLKK-20RK00-K	1 pc.

### Accessories

Item details	Comments	Item number	PU
Chest belt large L	with clip	ME-3-L	1 pc.
Chest belt small S	with clip	ME-3-S	1 pc.
Chest belt extra small	with velcro fastener	ME-3-SI	1 pc.
Chin strap large L	with velcro fastener	ME-23-L	1 pc.
Chin strap medium M	with velcro fastener	ME-23-M	1 pc.
Chin strap small S	with velcro fastener	ME-23-S	1 pc.
velcro fastener	for cable binding	ME-13	1 pc.
Electrode mounting	plastic, white	MCE-20	SET 10
Position marker	10/20 Position	ME-F-SET	SET 21



SC... (Page 26)

MCE-20



### Electrodes

You can use all stamped 10mm diameter gold or chlorinated or Ag/agCl Cup electrodes. This means you could save complete electrode sets because you can use cup electrodes in different applications. Cup electrodes are supplied as standard with a 1,5 mm DIN plug but can be supplied with a different plug on request.

### Electrode mountings

The MultiCap Cup plastic mountings are threaded for easy replacement. The cup electrodes are simply pushed into the plastic guide in the mounting. Each mounting is provided with a practical positioning tab.

### Additional electrodes

We recommend our ear electrodes as additional electrodes. These are made of sintered Ag/AgCl, tin or gold (Important – all electrodes within a single cable system must be made of the same material).



OHR-E-26