


Product-Information-Sheet



Polypectomy, EMR and ESD

Cold and Hot
Snares



- For polyp ligation prior to polypectomy
- Prevents postpolypectomy bleeding
- Preloaded loop, ready to use with handle
- Nylon-loop ligation keeps tissue damage to a minimum
- Opening diameter of 30 mm enables ligation of large polyps



| Article number | Article Name | Quantity value | Min. working channel ø value | Working length value | Loop ø | Area Compatibility |
|----------------|--------------|----------------|---|----------------------|--------|---|
| N5382130 | HX-400U-30 | 5 |  2.8 | 2300 mm | 30 mm |  |