

## Capillary blood collection system

### Capillary blood collection GK

#### Capillary blood collection GK

It consists of a prepared test vessel and a prepared plastic end-to-end collection capillary with stopper.

- Easy handling
- Plastic capillary\* with exact filling volume; complete inner surface prepared, unbreakable
- Collection vessel serves as centrifugal vessel; prepared for all common tests
- Light protected, tinted vessels for bilirubin analyses

\* with conformity certificate in accordance with the Weights and Measurement Regulations

Fill the capillary with capillary blood from a horizontal position.

After filling let the blood flow into the vessel from a vertical position (shake out remaining blood).

Remove the capillary, press on the attached stopper and mix or centrifuge.

#### 077003 EDTA 100 µl

Packing unit:

100 pcs in bag,

4,000 pcs in box

Ø 11/8 x 39 mm

