

Date (dd.mm.yyyy): 21.12.2022

**Product description:**

**NX2 ELITE, BASE CONFIGURATION**

## 1 *Product Overview*

---

**Included products are listed below:**

---

Item	Description	Qty	Rel
<b><u>NX2 ELITE, BASE CONFIGURATION</u></b>			
1.	<b>NX2 ELITE, BASE CONFIGURATION</b> 11362069 / Country of Origin: KR	1	R
2.	<b>NX2 ELITE, VA11, SOFTWARE</b> 11284501 / Country of Origin: KR	1	R
3.	<b>NX2 ELITE, VA11, INTL ENG, OPER SYS</b> 11320209 / Country of Origin: KR	1	R
4.	<b>CORDSET, 230V, GERMAN</b> 11284386 / Country of Origin: GB	1	R
5.	<b>NX2 ELITE, CLARIFY VE</b> 11284558 / Country of Origin: KR	1	R
6.	<b>NX2 ELITE, HD ZOOM</b> 11291812 / Country of Origin: KR	1	R
7.	<b>NX2 ELITE, DICOM CONNECTIVITY</b> 11291809 / Country of Origin: KR	1	R
8.	<b>NX2 ELITE, DICOM SR OB/GYN</b> 11284514 / Country of Origin: KR	1	R
9.	<b>NX2 ELITE, DICOM SR VASCULAR</b> 11284516 / Country of Origin: KR	1	R
10.	<b>NX2 ELITE, TRANSDUCER, P4-2, DL</b> 11284566 / Country of Origin: KR	1	R
11.	<b>NX2 ELITE, TRANSDUCER, EC9-4</b> 11284574 / Country of Origin: KR	1	R
12.	<b>ENDO XDCCR COVER, STERILE, 10/1</b> 11335907 / Country of Origin: US	1	R

13.	<b>NX2 ELITE, TRANSDUCER, L10-5V, DL</b> 11284582 / Country of Origin: KR	1	R
14.	<b>NX2 ELITE, TRANSDUCER, C5-2V, DL</b> 11284573 / Country of Origin: KR	1	R
15.	<b>GEL, 5L, 1/1</b> 11335940 / Country of Origin: US	2	R
16.	<b>NX2 ELITE, CONTROL PANEL COVER</b> 11291900 / Country of Origin: KR	1	R
17.	<b>NX2 ELITE, PRINTER, FRONT, BW, SONY</b> 11284587 / Country of Origin: CN	1	R
18.	<b>PRINTER PAPER BW, 84MM*12.5M, 10/1</b> 11335930 / Country of Origin: JP	1	R
19.	<b>NX2/NX2E 1.1, RUM, INST. FOR USE</b> 11361323 / Country of Origin: KR	1	R
20.	<b>NXSERIES G PREMIUM LOGISTIC SUPPORT</b> 11653945 / Country of Origin: KR	1	R
21.	<b>UPS 1.5kVA Tower 220-240V</b> 14429045 / Country of Origin:	1	R

**T O T A L quantity is for one system**

---

Optional products are listed below:

---

Item	Description	Qty	Rel
------	-------------	-----	-----

**NX2 ELITE, BASE CONFIGURATION**

**T O T A L quantity is for one system**

## 2 Technical description

Item	Description
1	<b>NX2 ELITE, BASE CONFIGURATION</b>
2	<b>NX2 ELITE, VA11, SOFTWARE</b>
3	<b>NX2 ELITE, VA11, INTL ENG, OPER SYS</b> Includes console with control panel, control panel overlay sheet, instructions for use, features and application reference, system reference and EMC Guidance Manual.
4	<b>CORDSET, 230V, GERMAN</b>
5	<b>NX2 ELITE, CLARIFY VE</b> Clarify VE technology is a real-time, adaptive, pixel-by-pixel analysis implemented through a simple, time-saving user interface that provides multiple levels of clarification to optimize tissue contrast resolution and definition of both tissue and vessels according to user preference.
6	<b>NX2 ELITE, HD ZOOM</b>
7	<b>NX2 ELITE, DICOM CONNECTIVITY</b> When installed, the ACUSON NX2 Elite ultrasound system acts as a DICOM Storage Class User and DICOM Print Class User, and can be connected to a DICOM color or grayscale printer, and/or to DICOM workstation and archive devices. In-Progress Store functionality allows individual transfer during an active exam.  Includes DICOM Region Calibration and DICOM Exchange Media (Stores exams in DICOM format to CD/DVD).  DICOM Worklist allows the user to automatically populate the "New Patient" screen with patient demographic information.  DICOM MPPS is used in conjunction with DICOM Modality Worklist to query and retrieve scheduled procedure information from a HIS and/or a RIS.  To be compatible with this functionality, the Hospital or Departmental Information System must include a DICOM Modality Performed Procedure Step server and a DICOM Modality Worklist server.
8	<b>NX2 ELITE, DICOM SR OB/GYN</b> DICOM Structured Reporting (SR) provides standardized report architecture to allow for easy transfer of measurements to offline PCs, workstations and archiving systems. DICOM Structured Reporting will automatically populate measurements to their respective fields in an external software package.
9	<b>NX2 ELITE, DICOM SR VASCULAR</b> DICOM Structured Reporting (SR) provides standardized report architecture to allow for easy

Item	Description
	transfer of measurements to offline PCs, workstations and archiving systems. DICOM Structured Reporting will automatically populate measurements to their respective fields in an external software package.
10	<b>NX2 ELITE, TRANSDUCER, P4-2, DL</b>
11	<b>NX2 ELITE, TRANSDUCER, EC9-4</b>
12	<b>ENDO XDCR COVER, STERILE, 10/1</b> The set contains: <ul style="list-style-type: none"> <li>- Sterile 3.2 tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") latex cover</li> <li>- Polyethylene cord cover 7.6 x 244 cm</li> <li>- Gel</li> <li>- Elastic bands</li> </ul>
13	<b>NX2 ELITE, TRANSDUCER, L10-5V, DL</b>
14	<b>NX2 ELITE, TRANSDUCER, C5-2V, DL</b>
15	<b>GEL, 5L, 1/1</b>
16	<b>NX2 ELITE, CONTROL PANEL COVER</b> A one-piece silicon control panel cover designed to fit the ACUSON NX2 Elite ultrasound system. The cover allows access to the DCG controls, includes separate knob covers for full rotary control, and allows easy cleaning and disinfection.
17	<b>NX2 ELITE, PRINTER, FRONT, BW, SONY</b>
18	<b>PRINTER PAPER BW, 84MM*12.5M, 10/1</b>
19	<b>NX2/NX2E 1.1, RUM, INST. FOR USE</b>
20	<b>NXSERIES G PREMIUM LOGISTIC SUPPORT</b>
21	<b>UPS 1.5kVA Tower 220-240V</b> Tower uninterruptible power supply for HPE server-based workstation/workplace supporting only 220/230/240 voltage.  UPS management software is included