



ZHONGHONG PULIN MEDICAL PRODUCTS CO., LTD
35 Songhe Ave, Luannan County, Tangshan City, Hebei Province 063500 86-10-67866277

SPECIFICATION OF NITRILE DISPOSABLE GLOVES (3.5g)
powder-free, ambidextrous, finger textured, AQL 1.5

Type: ZHPFN02
EN455 EN374 EN420

Spec	Size	X-small NM35105	Small NM35205	Medium NM35305	Larg eNM35405	X-large NM35505	Tolerance
	Total Length (mm)		240	240	240	240	240
Width of palm (mm)		75	85	95	105	115	±5
Width of cuff (mm)		85	90	95	100	105	±5
Thickness (mm)	Finger	0.08	0.08	0.08	0.08	0.08	Min
	Palm	0.06	0.06	0.06	0.06	0.06	Min
	cuff	0.05	0.05	0.05	0.05	0.05	Min
Weight per piece (gram)		2.5	3.0	3.5	4.0	4.5	±0.3
Tensile-strength (≥)min		14Mpa	14Mpa	14Mpa	14Mpa	14MPa	Min
Elongation (%)min	Before aging:	500	500	500	500	500	Min
	After aging:	400	400	400	400	400	Min

Disposable Nitrile Exam & Protective Glove

· Before use, make sure the gloves are undamaged.
· Disposable medical gloves, made of nitrile rubber, ambidextrous. Personal protective equipment of category III.

· Store in a temperature from 10°C to 40°C, in a dry and cool place in original packaging. Protect against ozone.
· After use, dispose of in accordance with applicable local regulations.



Disposable Nitrile Exam & Protective Glove

CE2777
AQL 1.5




· non-sterile
· powder-free
· ambidextrous
· latex-free

BLUE
L
100PCS

Disposable Nitrile Exam & Protective Glove

BLUE
L
100PCS




REF: ZHDFN02
GTIN: 6970405800489
LOT: ZH2208C08
2022-08
2027-07

CE
MD

BLUE
L
100PCS



Disposable Nitrile Exam & Protective Glove



MEDICAL PPE | [EU]2017/745 | [EU]2016/425 | Class I | Category III | EN 455 | EN ISO 21420, EN ISO 374



This product has been tested in accordance with EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016 and EN ISO 21420:2020, and achieved the following performance levels:

Test chemical	EN ISO 374-1:2016+A1:2018 Permeation level	EN ISO 374-4:2019 Degradation (mean value)
K Sodium Hydroxide 40%	6	-8,3%
P Hydrogen Peroxide 30%	2	34,1%
T Formaldehyde 37%	4	34,3%

EN ISO 374-1:2016+A1:2018 Type B EN ISO 374-1:2016+A1:2018 Permeation levels are based on breakthrough times as follows:

KPT	Performance level					
	1	2	3	4	5	6
Min. breakthrough times (mins)	>10	>30	>60	>120	>240	>480

EN ISO 374-4:2019 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

Disposable examination and protective gloves made of nitrile; powder-free; non-sterile.
Warning: Contact with food is allowed. Latex-free, powder-free. Compliant with MD Regulation (EU)2017/745(Class I), EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2009, PPE Regulation (EU)2016/425(CAT III), EN ISO 21420:2020, and EN ISO 374. Visit the following link to read more about and download the Declaration of Conformity: [your own link](#).

Donning & doffing gloves in a proper way is a skill that needs to be practiced by healthcare workers and others that use gloves. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. Insert the glove to the hand carefully without damaging the glove. When removing the glove, grasp the outside of the glove from the wrist area, peel the glove away from the hand, hold it in the opposite gloved hand, slide an un-gloved finger under the wrist of the remaining glove, being careful not to touch the contaminated surface of the glove, peel the remaining glove out. Avoid snapping, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby.



Resistance to bacteria and fungi - pass.
Resistance to virus - pass.
*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

CE2777

Notified Body responsible for certification and ongoing conformity:
SATRA Technology Europe Ltd,
Bracetown Business Park, Clonoe, Dublin, D15 YN2P, Ireland

Statement and Caution: This information does not reflect the actual duration of protection at the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance have been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. The result can be different if the chemical is used in a mixture. It is recommended to check whether the gloves are suitable for the intended use because the conditions (such as temperature, abrasion, and degradation) at the workplace may differ from the testing conditions. Used gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical in contact with the gloves, etc. may shorten the actual service life of gloves significantly. For corrosive chemicals, degradation can be the most important factor to consider when selecting chemical-resistant gloves. Before use, inspect the gloves for any defect or imperfections.

Zhonghong Pulin Medical Products Co.,Ltd.
West Industrial Park, Luannan County, Tangshan City,
Hebei Province, 063500 China
Phone: +86 3154167691

Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd,
2595AA, The Hague, Netherlands
peter@lotusnl.com

Importer: **Your company**

Made in China

Disposable Nitrile Exam & Protective Glove

· Disposable medical gloves, made of nitrile rubber, ambidextrous. Personal protective equipment of category III.

· Before use, make sure the gloves are undamaged.

· Store in a temperature from 10°C to 40°C, in a dry and cool place in original packaging. Protect against ozone.

· After use, dispose of in accordance with applicable local regulations.



Disposable Nitrile Exam & Protective Glove

CE2777
AQL 1.5



- non-sterile
- powder-free
- ambidextrous
- latex-free

BLUE
M
100PCS

Disposable Nitrile Exam & Protective Glove

BLUE
M
100PCS




REF: ZHDFN02
GTIN: 6970405800472
LOT: ZH2208C08
2022-08
2027-07

CE
MD

BLUE
M
100PCS



Disposable Nitrile Exam & Protective Glove



MEDICAL PPE | [EU]2017/745 | [EU]2016/425 | Class I | Category III | EN 455 | EN ISO 21420, EN ISO 374



This product has been tested in accordance with EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016 and EN ISO 21420:2020, and achieved the following performance levels:

Test chemical	EN ISO 374-1:2016+A1:2018 Permeation level	EN ISO 374-4:2019 Degradation (mean value)
K Sodium Hydroxide 40%	6	-8,3%
P Hydrogen Peroxide 30%	2	34,1%
T Formaldehyde 37%	4	34,3%

EN ISO 374-1:2016+A1:2018 Permeation levels are based on breakthrough times as follows:

Performance level	1	2	3	4	5	6
Min. breakthrough times (mins)	>10	>30	>60	>120	>240	>480

EN ISO 374-4:2019 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

Disposable examination and protective gloves made of nitrile; powder-free; non-sterile.
Warning: Contact with food is allowed. Latex-free, powder-free. Compliant with MD Regulation (EU)2017/745(Class I), EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2009, PPE Regulation (EU)2016/425(CAT III), EN ISO 21420:2020, and EN ISO 374. Visit the following link to read more about and download the Declaration of Conformity: [your own link](#).

Donning & doffing gloves in a proper way is a skill that needs to be practiced by healthcare workers and others that use gloves. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. Insert the glove to the hand carefully without damaging the glove. When removing the gloves, grasp the outside of the glove from the wrist area, peel the glove away from the hand, hold it in the opposite gloved hand, slide an un-gloved finger under the wrist of the remaining glove, being careful not to touch the contaminated surface of the glove, peel the remaining glove out. Avoid snapping, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby.



Resistance to bacteria and fungi - pass.
Resistance to virus - pass.
*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

CE2777

Notified Body responsible for certification and ongoing conformity:
SATRA Technology Europe Ltd,
Bracetown Business Park, Clonoe, Dublin, D15 YN2P, Ireland

Statement and Caution: This information does not reflect the actual duration of protection at the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance have been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. The result can be different if the chemical is used in a mixture. It is recommended to check whether the gloves are suitable for the intended use because the conditions (such as temperature, abrasion, and degradation) at the workplace may differ from the testing conditions. Used gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical in contact with the gloves, etc. may shorten the actual service life of gloves significantly. For corrosive chemicals, degradation can be the most important factor to consider when selecting chemical-resistant gloves. Before use, inspect the gloves for any defect or imperfections.

Zhonghong Pulin Medical Products Co.,Ltd.
West Industrial Park, Luannan County, Tangshan City,
Hebei Province, 063500 China
Phone: +86 3154167691

Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd,
2595AA, The Hague, Netherlands
peter@lotusnl.com

Importer: **Your company**

Made in China

Disposable Nitrile Exam & Protective Glove

· Disposable medical gloves, made of nitrile rubber, ambidextrous. Personal protective equipment of category III.

· Before use, make sure the gloves are undamaged.

· Store in a temperature from 10°C to 40°C, in a dry and cool place in original packaging. Protect against ozone.

· After use, dispose of in accordance with applicable local regulations.



Disposable Nitrile Exam & Protective Glove

CE2777
AQL 1.5


BLUE
S
100PCS

- non-sterile
- ambidextrous
- powder-free
- latex-free



Disposable Nitrile Exam & Protective Glove

BLUE
S
100PCS




6970405800465

REF: ZHDFN02
GTIN: 6970405800465
LOT: ZH2208C08
2022-08
2027-07


CE
MD

BLUE
S
100PCS



(01) 6970405800465 (17) 270701 (10) ZH2208C08

Disposable Nitrile Exam & Protective Glove



MEDICAL [EU]2017/745 Class I EN 455
PPE [EU]2016/425 Category III EN ISO 21420, EN ISO 374



This product has been tested in accordance with EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016 and EN ISO 21420:2020, and achieved the following performance levels:

Test chemical	EN ISO 374-1:2016+A1:2018 Permeation level	EN ISO 374-4:2019 Degradation (mean value)
K Sodium Hydroxide 40%	6	-8,3%
P Hydrogen Peroxide 30%	2	34,1%
T Formaldehyde 37%	4	34,3%

EN ISO 374-1:2016+A1:2018 Type B EN ISO 374-1:2016+A1:2018 Permeation levels are based on breakthrough times as follows:

Performance level	1	2	3	4	5	6
	Min. breakthrough times (mins)	>10	>30	>60	>120	>240

EN ISO 374-4:2019 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

Disposable examination and protective gloves made of nitrile; powder-free; non-sterile.
Warning: Contact with food is allowed. Latex-free, powder-free. Compliant with MD Regulation (EU)2017/745(Class I), EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2009, PPE Regulation (EU)2016/425(CAT III), EN ISO 21420:2020, and EN ISO 374. Visit the following link to read more about and download the Declaration of Conformity: [your own link](#).

Donning & doffing gloves in a proper way is a skill that needs to be practiced by healthcare workers and others that use gloves. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. Insert the glove to the hand carefully without damaging the glove. When removing the glove, grasp the outside of the glove from the wrist area, peel the glove away from the hand, hold it in the opposite gloved hand, slide an un-gloved finger under the wrist of the remaining glove, being careful not to touch the contaminated surface of the glove, peel the remaining glove out. Avoid snapping, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby.



Resistance to bacteria and fungi - pass.
Resistance to virus - pass.
*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

CE2777

Notified Body responsible for certification and ongoing conformity:
SATRA Technology Europe Ltd,
Bracetown Business Park, Clonoe, Dublin, D15 YN2P, Ireland

Statement and Caution: This information does not reflect the actual duration of protection at the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance have been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. The result can be different if the chemical is used in a mixture. It is recommended to check whether the gloves are suitable for the intended use because the conditions (such as temperature, abrasion, and degradation) at the workplace may differ from the testing conditions. Used gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical in contact with the gloves, etc. may shorten the actual service life of gloves significantly. For corrosive chemicals, degradation can be the most important factor to consider when selecting chemical-resistant gloves. Before use, inspect the gloves for any defect or imperfections.

Zhonghong Pulin Medical Products Co.,Ltd.
West Industrial Park, Luannan County, Tangshan City,
Hebei Province, 063500 China
Phone: +86 3154167691

Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd,
2595AA, The Hague, Netherlands
peter@lotusnl.com

Importer: Your company

Made in China