930D 电动吸引器

Features

- 1. Use a large-flow oil-free self-lubricating negative pressure pump with high pumping rate and no oil pollution.
- 2. The large-caliber liquid storage bottle is equipped with a bottle stopper with a sealing ring, which is convenient for users to open and close the bottle stopper and remove the dirty liquid in the bottle.
- 3. An overflow protection device is provided to prevent liquid from entering the intermediate pipe and the negative pressure pump.
- 4. The equipped air filter can prevent the negative pressure pump from being contaminated.
- 5. The manual switch and the foot switch are connected in parallel and can be selected at will.
- 6. The negative pressure value can be adjusted steplessly according to clinical needs.
- 7. The whole machine has a compact structure, stable operation and low noise.

Parameters

Ultimate negative pressure value: >0.09MPa (680mmHg Negative pressure adjustment range: 0.02MPa (150mmHg)

Vacuum rate: >40L/min

Liquid storage bottle: 2500ml (glass) x 2

(plastic PC bottle and disposable liquid absorption bag can be selected)

Power supply: ~220V, 50Hz

Input power: 250VA



Accessories

Disposable suction tip

Suction catheter

Filter