




**United**<sup>®</sup>  
Orthopedic Corporation



U2<sup>™</sup> Knee  
System

# Consistency in U2 Knee System over Time

- Share the same femoral profile of all types allow inter-changeability
- Uniform AP/ML dimensions in all tibial components enhance intraoperative flexibility






# Keep Bone Preservation in Mind...




- 13 sizes of femoral components in 2 mm increments




- Bone preserving 9 mm resection in both distal and posterior condyles



- Fully interchangeable among all sizes with the consistent 18 mm intercondylar width



- Limited intercondylar bone removal together with round shaped corner avoids the risk of intercondylar fracture



- Horizontal tibial resection surgical technique with built in 5° tibial slope





# Optimized Freedom, Better Satisfactory



- 155° flexion without implant impingement

Proven by **ENDOLAB**  
MECHANICAL ENGINEERING  
GMBH



- A matched "curve on curve" tibial post and box contact, eliminates impact and wear




- A refined Cam/Post mechanism minimizes the risk of dislocation and prosthetic loosening

# Modular Design Allow Augment and Stem for Moderate Bone Defects

Unilateral titanium augments in 5 or 10 mm thickness with screw fixation



30 mm stem made of titanium alloy to enhance stability





Distal end sealed with  
pre-assembled plug



Design without through hole  
avoids potential backside wear

# CMA<sup>TM</sup>

*Cemented Modular Augmentable*

# Baseplate

Enhanced Stability, Reduced  
Operation Time



U2 *Xuc*<sup>TM</sup>  
Ultracongruent

- Accommodates CR femoral component and the PCL sacrificed surgical technique for more bone preservation
- Up to 12.5 mm prominent anterior lip and a more conforming articulating surface provide joint stability
- Highly cross-linked UHMWPE reinforces insert longevity
- Mean flexion >150° satisfies patient's demand of activity
- Tibial Baseplate "One up / one down" size pairing

[illegible]

# Mobile Bearing Baseplate

Highly mirror-polished platform minimizes backside wear.

The pure axial rotating design based on biomechanical principles proven over decades of clinical usage.

## U2 MBC™

### Mobile Bearing Congruent

- Accommodate with CR femoral component allows more bone preservation
- Surgical flexibility: PCL either retained or sacrificed
- Central stopper enhances joint M/L stability and also allows 4.5 degrees hyperextension





# U2 MB™

## Mobile Bearing

- Accommodates with PS femoral component.
- Surgery: Both ACL and PCL sacrificed.

| MB / MBC Size Pairing | Mobile Bearing Baseplate | Insert |    |    |    |    |    |
|-----------------------|--------------------------|--------|----|----|----|----|----|
|                       |                          | #1     | #2 | #3 | #4 | #5 | #6 |
|                       | #1                       | ●      | ●  |    |    |    |    |
|                       | #2                       | ●      | ●  | ●  |    |    |    |
|                       | #3                       | ●      | ●  | ●  | ●  |    |    |
|                       | #4                       | ●      | ●  | ●  | ●  | ●  |    |
|                       | #5                       | ●      | ●  | ●  | ●  | ●  | ●  |
|                       | #6                       | ●      | ●  | ●  | ●  | ●  | ●  |

★ (MBC Insert) Size of femoral component is equal to insert.


★ (MB Insert) Size of femoral component is fully interchangeable between femoral and insert.

Defined a More Secured Stability with  
More Choices

# PSA™ Revision System

## A More Secured Stability

Enhanced locking mechanism prevents all the locking  
screws from backing out




\*No. US9044327

Tibia Trays Accommodate All U2 Knee  
Femoral Components



## Multiple Choices of Extensions

- 5 length 8 diameters of stems in straight or open slot curved from 30 mm to 200 mm
- 4 thicknesses of distal femoral augments in 4, 8, 12, 16 mm
- 2 thicknesses of posterior augments in 4, 8 mm
- 3 thickness of unilateral tibial augments in 5, 10, 15 mm
- 3 selections of offset with full range orientation in 2, 4, 6 mm






Vitamin E Stabilized Highly Cross-linked UHMWPE

# Advanced Bearing Technology


## Wear Performance



## Extraordinary Wear Performance

E-XPE insert shows 60% reduction in gravimetric wear compared to HXLPE after accelerated aging.


## Mechanical Strength



## Enhanced Mechanical Property

No heat treatment is required after cross-linking process. Therefore, E-XPE shows a 20% improvement of tensile strength than highly cross-linked polyethylene.

## Oxidative Stability



## Superior Oxidative Stability

Surface oxidative index of E-XPE shows no oxidation after in vitro accelerated aging test.


## U2 Knee Fixed Bearing System

---



**CR**

Cruciate Retaining



**XUC**

Ultracongruent



**PS**

Posterior Stabilized



## U2 Knee Mobile Bearing System

---



**MBC**

Mobile Bearing Congruent



**MB**

Mobile Bearing




## U2 Knee Revision System

---



**PSA**

Posterior Stabilized  
Augmentable



**PSA LC**

PSA Low Constrained




Reduce Table Clutter, Mitigate Infection Risk

# U2 Knee AiO™

All-in-One Block

Accommodates all 13 sizes of anterior and posterior femoral cuts in one block




\*No. US9974547



7 chamfer cutting guide sizes  
available for 13 femoral components

# U2 Knee MDT™

## Modular Disposable Trial




Single-use convenience, reduces sterilization, reprocessing costs and infection risk




When used together, the U2 Knee AiO™ and MDT™ can reduce instrument trays count from 6 to 1.5

U2 Femoral Component




| Size | #1 | #1.5 | #2 | #2.5 | #3 | #3.5 | #4 | #4.5 | #5 | #5.5 | #6 | #6.5 | #7 |
|------|----|------|----|------|----|------|----|------|----|------|----|------|----|
| AP   | 52 | 54   | 56 | 58   | 60 | 62   | 64 | 66   | 68 | 70   | 72 | 74   | 76 |
| ML   | 56 | 58   | 60 | 62   | 64 | 66   | 68 | 70   | 72 | 74   | 76 | 78   | 80 |

U2 Tibial Baseplate




| Size | #0   | #1 | #2   | #3 | #4   | #5   | #6   | #7   |
|------|------|----|------|----|------|------|------|------|
| AP   | 39.5 | 42 | 44.5 | 47 | 49.5 | 52.5 | 55.5 | 58.5 |
| ML   | 60   | 63 | 66   | 69 | 72   | 76   | 80   | 84   |

Onset Patellar Component



| Size      | XS | S  | M   | L  | XL  | XXL | EL   |
|-----------|----|----|-----|----|-----|-----|------|
| Thickness | 7  | 8  | 8.5 | 9  | 9.5 | 10  | 10.5 |
| Diameter  | 26 | 29 | 32  | 35 | 38  | 41  | 44   |

Inset Patellar Component



| Size      | S  | M  | L  | XL |
|-----------|----|----|----|----|
| Thickness | 8  | 10 | 10 | 10 |
| Diameter  | 22 | 25 | 28 | 32 |

Unit : mm

# Order Information

## U2 Femoral Component

 Special Order Items



Cemented CR

|      | Left      | Right     | Left      | Right     |
|------|-----------|-----------|-----------|-----------|
| #1   | 2103-1310 | 2103-1410 | 2103-1110 | 2103-1210 |
| #1.5 | 2103-1315 | 2103-1415 | 2103-1115 | 2103-1215 |
| #2   | 2103-1320 | 2103-1420 | 2103-1120 | 2103-1220 |
| #2.5 | 2103-1325 | 2103-1425 | 2103-1125 | 2103-1225 |
| #3   | 2103-1330 | 2103-1430 | 2103-1130 | 2103-1230 |
| #3.5 | 2103-1335 | 2103-1435 | 2103-1135 | 2103-1235 |
| #4   | 2103-1340 | 2103-1440 | 2103-1140 | 2103-1240 |
| #4.5 | 2103-1345 | 2103-1445 | 2103-1145 | 2103-1245 |
| #5   | 2103-1350 | 2103-1450 | 2103-1150 | 2103-1250 |
| #5.5 | 2103-1355 | 2103-1455 | 2103-1155 | 2103-1255 |
| #6   | 2103-1360 | 2103-1460 | 2103-1160 | 2103-1260 |
| #6.5 | 2103-1365 | 2103-1465 | 2103-1165 | 2103-1265 |
| #7   | 2103-1370 | 2103-1470 | 2103-1170 | 2103-1270 |

Porous CR



PS

|      | Left      | Right     |
|------|-----------|-----------|
| #1   | 2103-3110 | 2103-3210 |
| #1.5 | 2103-3115 | 2103-3215 |
| #2   | 2103-3120 | 2103-3220 |
| #2.5 | 2103-3125 | 2103-3225 |
| #3   | 2103-3130 | 2103-3230 |
| #3.5 | 2103-3135 | 2103-3235 |
| #4   | 2103-3140 | 2103-3240 |
| #4.5 | 2103-3145 | 2103-3245 |
| #5   | 2103-3150 | 2103-3250 |
| #5.5 | 2103-3155 | 2103-3255 |
| #6   | 2103-3160 | 2103-3260 |
| #6.5 | 2103-3165 | 2103-3265 |
| #7   | 2103-3170 | 2103-3270 |



PSA

|    | Left      | Right     |
|----|-----------|-----------|
| #1 | 2103-5110 | 2103-5210 |
| #2 | 2103-5120 | 2103-5220 |
| #3 | 2103-5130 | 2103-5230 |
| #4 | 2103-5140 | 2103-5240 |
| #5 | 2103-5150 | 2103-5250 |
| #6 | 2103-5160 | 2103-5260 |

## U2 Tibial Tray



CMA

|    |           |
|----|-----------|
| #0 | 2203-3200 |
| #1 | 2203-3210 |
| #2 | 2203-3220 |
| #3 | 2203-3230 |
| #4 | 2203-3240 |
| #5 | 2203-3250 |
| #6 | 2203-3260 |
| #7 | 2203-3270 |



Mobile Bearing

|    |           |
|----|-----------|
| #1 | 2203-7010 |
| #2 | 2203-7020 |
| #3 | 2203-7030 |
| #4 | 2203-7040 |
| #5 | 2203-7050 |
| #6 | 2203-7060 |



PSA

|    |           |
|----|-----------|
| #1 | 2203-5210 |
| #2 | 2203-5220 |
| #3 | 2203-5230 |
| #4 | 2203-5240 |
| #5 | 2203-5250 |
| #6 | 2203-5260 |

# Tibial Insert



CR Insert

|    | UHMWPE    |           |           |           |           |           |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 10 mm     | 11 mm     | 12 mm     | 13 mm     | 14 mm     | 15 mm     | 16 mm     | 17 mm     | 18 mm     |
| #0 | 2303-1201 | 2303-1206 | 2303-1202 | 2303-1207 | 2303-1203 | 2303-1208 | 2303-1204 | 2303-1209 | 2303-1200 | 2303-1205 |
| #1 | 2303-1211 | 2303-1216 | 2303-1212 | 2303-1217 | 2303-1213 | 2303-1218 | 2303-1214 | 2303-1219 | 2303-1210 | 2303-1215 |
| #2 | 2303-1221 | 2303-1226 | 2303-1222 | 2303-1227 | 2303-1223 | 2303-1228 | 2303-1224 | 2303-1229 | 2303-1220 | 2303-1225 |
| #3 | 2303-1231 | 2303-1236 | 2303-1232 | 2303-1237 | 2303-1233 | 2303-1238 | 2303-1234 | 2303-1239 | 2303-1230 | 2303-1235 |
| #4 | 2303-1241 | 2303-1246 | 2303-1242 | 2303-1247 | 2303-1243 | 2303-1248 | 2303-1244 | 2303-1249 | 2303-1240 | 2303-1245 |
| #5 | 2303-1251 | 2303-1256 | 2303-1252 | 2303-1257 | 2303-1253 | 2303-1258 | 2303-1254 | 2303-1259 | 2303-1250 | 2303-1255 |
| #6 | 2303-1261 | 2303-1266 | 2303-1262 | 2303-1267 | 2303-1263 | 2303-1268 | 2303-1264 | 2303-1269 | 2303-1260 | 2303-1265 |
| #7 | 2303-1271 | 2303-1276 | 2303-1272 | 2303-1277 | 2303-1273 | 2303-1278 | 2303-1274 | 2303-1279 | 2303-1270 | 2303-1275 |

|    | XPE       |           |           |           |           |           |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 10 mm     | 11 mm     | 12 mm     | 13 mm     | 14 mm     | 15 mm     | 16 mm     | 17 mm     | 18 mm     |
| #0 | 2303-1601 | 2303-1606 | 2303-1602 | 2303-1607 | 2303-1603 | 2303-1608 | 2303-1604 | 2303-1609 | 2303-1600 | 2303-1605 |
| #1 | 2303-1611 | 2303-1616 | 2303-1612 | 2303-1617 | 2303-1613 | 2303-1618 | 2303-1614 | 2303-1619 | 2303-1610 | 2303-1615 |
| #2 | 2303-1621 | 2303-1626 | 2303-1622 | 2303-1627 | 2303-1623 | 2303-1628 | 2303-1624 | 2303-1629 | 2303-1620 | 2303-1625 |
| #3 | 2303-1631 | 2303-1636 | 2303-1632 | 2303-1637 | 2303-1633 | 2303-1638 | 2303-1634 | 2303-1639 | 2303-1630 | 2303-1635 |
| #4 | 2303-1641 | 2303-1646 | 2303-1642 | 2303-1647 | 2303-1643 | 2303-1648 | 2303-1644 | 2303-1649 | 2303-1640 | 2303-1645 |
| #5 | 2303-1651 | 2303-1656 | 2303-1652 | 2303-1657 | 2303-1653 | 2303-1658 | 2303-1654 | 2303-1659 | 2303-1650 | 2303-1655 |
| #6 | 2303-1661 | 2303-1666 | 2303-1662 | 2303-1667 | 2303-1663 | 2303-1668 | 2303-1664 | 2303-1669 | 2303-1660 | 2303-1665 |
| #7 | 2303-1671 | 2303-1676 | 2303-1672 | 2303-1677 | 2303-1673 | 2303-1678 | 2303-1674 | 2303-1679 | 2303-1670 | 2303-1675 |



CR Insert

| E-XPE |           |           |           |           |           |           |           |           |           |           |
|-------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|       | 9 mm      | 10 mm     | 11 mm     | 12 mm     | 13 mm     | 14 mm     | 15 mm     | 16 mm     | 17 mm     | 18 mm     |
| #0    | 2303-1801 | 2303-1806 | 2303-1802 | 2303-1807 | 2303-1803 | 2303-1808 | 2303-1804 | 2303-1809 | 2303-1800 | 2303-1805 |
| #1    | 2303-1811 | 2303-1816 | 2303-1812 | 2303-1817 | 2303-1813 | 2303-1818 | 2303-1814 | 2303-1819 | 2303-1810 | 2303-1815 |
| #2    | 2303-1821 | 2303-1826 | 2303-1822 | 2303-1827 | 2303-1823 | 2303-1828 | 2303-1824 | 2303-1829 | 2303-1820 | 2303-1825 |
| #3    | 2303-1831 | 2303-1836 | 2303-1832 | 2303-1837 | 2303-1833 | 2303-1838 | 2303-1834 | 2303-1839 | 2303-1830 | 2303-1835 |
| #4    | 2303-1841 | 2303-1846 | 2303-1842 | 2303-1847 | 2303-1843 | 2303-1848 | 2303-1844 | 2303-1849 | 2303-1840 | 2303-1845 |
| #5    | 2303-1851 | 2303-1856 | 2303-1852 | 2303-1857 | 2303-1853 | 2303-1858 | 2303-1854 | 2303-1859 | 2303-1850 | 2303-1855 |
| #6    | 2303-1861 | 2303-1866 | 2303-1862 | 2303-1867 | 2303-1863 | 2303-1868 | 2303-1864 | 2303-1869 | 2303-1860 | 2303-1865 |
| #7    | 2303-1871 | 2303-1876 | 2303-1872 | 2303-1877 | 2303-1873 | 2303-1878 | 2303-1874 | 2303-1879 | 2303-1870 | 2303-1875 |

# Tibial Insert



PS Insert

| UHMWPE |           |           |           |           |           |           |           |           |           |           |
|--------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|        | 9 mm      | 10 mm     | 11 mm     | 12 mm     | 13 mm     | 14 mm     | 15 mm     | 16 mm     | 17 mm     | 18 mm     |
| #0     | 2303-3001 | 2303-3006 | 2303-3002 | 2303-3007 | 2303-3003 | 2303-3008 | 2303-3004 | 2303-3009 | 2303-3000 | N/A       |
| #1     | 2303-3011 | 2303-3016 | 2303-3012 | 2303-3017 | 2303-3013 | 2303-3018 | 2303-3014 | 2303-3019 | 2303-3010 | 2303-3015 |
| #2     | 2303-3021 | 2303-3026 | 2303-3022 | 2303-3027 | 2303-3023 | 2303-3028 | 2303-3024 | 2303-3029 | 2303-3020 | 2303-3025 |
| #3     | 2303-3031 | 2303-3036 | 2303-3032 | 2303-3037 | 2303-3033 | 2303-3038 | 2303-3034 | 2303-3039 | 2303-3030 | 2303-3035 |
| #4     | 2303-3041 | 2303-3046 | 2303-3042 | 2303-3047 | 2303-3043 | 2303-3048 | 2303-3044 | 2303-3049 | 2303-3040 | 2303-3045 |
| #5     | 2303-3051 | 2303-3056 | 2303-3052 | 2303-3057 | 2303-3053 | 2303-3058 | 2303-3054 | 2303-3059 | 2303-3050 | 2303-3055 |
| #6     | 2303-3061 | 2303-3066 | 2303-3062 | 2303-3067 | 2303-3063 | 2303-3068 | 2303-3064 | 2303-3069 | 2303-3060 | 2303-3065 |
| #7     | 2303-3071 | 2303-3076 | 2303-3072 | 2303-3077 | 2303-3073 | 2303-3078 | 2303-3074 | 2303-3079 | 2303-3070 | 2303-3075 |

| XPE |           |           |           |           |           |           |           |           |           |           |
|-----|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|     | 9 mm      | 10 mm     | 11 mm     | 12 mm     | 13 mm     | 14 mm     | 15 mm     | 16 mm     | 17 mm     | 18 mm     |
| #0  | 2303-3601 | 2303-3606 | 2303-3602 | 2303-3607 | 2303-3603 | 2303-3608 | 2303-3604 | 2303-3609 | 2303-3600 | N/A       |
| #1  | 2303-3611 | 2303-3616 | 2303-3612 | 2303-3617 | 2303-3613 | 2303-3618 | 2303-3614 | 2303-3619 | 2303-3610 | 2303-3615 |
| #2  | 2303-3621 | 2303-3626 | 2303-3622 | 2303-3627 | 2303-3623 | 2303-3628 | 2303-3624 | 2303-3629 | 2303-3620 | 2303-3625 |
| #3  | 2303-3631 | 2303-3636 | 2303-3632 | 2303-3637 | 2303-3633 | 2303-3638 | 2303-3634 | 2303-3639 | 2303-3630 | 2303-3635 |
| #4  | 2303-3641 | 2303-3646 | 2303-3642 | 2303-3647 | 2303-3643 | 2303-3648 | 2303-3644 | 2303-3649 | 2303-3640 | 2303-3645 |
| #5  | 2303-3651 | 2303-3656 | 2303-3652 | 2303-3657 | 2303-3653 | 2303-3658 | 2303-3654 | 2303-3659 | 2303-3650 | 2303-3655 |
| #6  | 2303-3661 | 2303-3666 | 2303-3662 | 2303-3667 | 2303-3663 | 2303-3668 | 2303-3664 | 2303-3669 | 2303-3660 | 2303-3665 |
| #7  | 2303-3671 | 2303-3676 | 2303-3672 | 2303-3677 | 2303-3673 | 2303-3678 | 2303-3674 | 2303-3679 | 2303-3670 | 2303-3675 |



PS Insert

| E-XPE |           |           |           |           |           |           |           |           |           |           |
|-------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|       | 9 mm      | 10 mm     | 11 mm     | 12 mm     | 13 mm     | 14 mm     | 15 mm     | 16 mm     | 17 mm     | 18 mm     |
| #0    | 2303-3801 | 2303-3806 | 2303-3802 | 2303-3807 | 2303-3803 | 2303-3808 | 2303-3804 | 2303-3809 | 2303-3800 | N/A       |
| #1    | 2303-3811 | 2303-3816 | 2303-3812 | 2303-3817 | 2303-3813 | 2303-3818 | 2303-3814 | 2303-3819 | 2303-3810 | 2303-3815 |
| #2    | 2303-3821 | 2303-3826 | 2303-3822 | 2303-3827 | 2303-3823 | 2303-3828 | 2303-3824 | 2303-3829 | 2303-3820 | 2303-3825 |
| #3    | 2303-3831 | 2303-3836 | 2303-3832 | 2303-3837 | 2303-3833 | 2303-3838 | 2303-3834 | 2303-3839 | 2303-3830 | 2303-3835 |
| #4    | 2303-3841 | 2303-3846 | 2303-3842 | 2303-3847 | 2303-3843 | 2303-3848 | 2303-3844 | 2303-3849 | 2303-3840 | 2303-3845 |
| #5    | 2303-3851 | 2303-3856 | 2303-3852 | 2303-3857 | 2303-3853 | 2303-3858 | 2303-3854 | 2303-3859 | 2303-3850 | 2303-3855 |
| #6    | 2303-3861 | 2303-3866 | 2303-3862 | 2303-3867 | 2303-3863 | 2303-3868 | 2303-3864 | 2303-3869 | 2303-3860 | 2303-3865 |
| #7    | 2303-3871 | 2303-3876 | 2303-3872 | 2303-3877 | 2303-3873 | 2303-3878 | 2303-3874 | 2303-3879 | 2303-3870 | 2303-3875 |

# Tibial Insert



XUC Insert

|    | XPE       |           |           |           |           |           |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 10 mm     | 11 mm     | 12 mm     | 13 mm     | 14 mm     | 15 mm     | 16 mm     | 17 mm     | 18 mm     |
| #0 | 2303-1401 | 2303-1406 | 2303-1402 | 2303-1407 | 2303-1403 | 2303-1408 | 2303-1404 | 2303-1409 | 2303-1400 | 2303-1405 |
| #1 | 2303-1411 | 2303-1416 | 2303-1412 | 2303-1417 | 2303-1413 | 2303-1418 | 2303-1414 | 2303-1419 | 2303-1410 | 2303-1415 |
| #2 | 2303-1421 | 2303-1426 | 2303-1422 | 2303-1427 | 2303-1423 | 2303-1428 | 2303-1424 | 2303-1429 | 2303-1420 | 2303-1425 |
| #3 | 2303-1431 | 2303-1436 | 2303-1432 | 2303-1437 | 2303-1433 | 2303-1438 | 2303-1434 | 2303-1439 | 2303-1430 | 2303-1435 |
| #4 | 2303-1441 | 2303-1446 | 2303-1442 | 2303-1447 | 2303-1443 | 2303-1448 | 2303-1444 | 2303-1449 | 2303-1440 | 2303-1445 |
| #5 | 2303-1451 | 2303-1456 | 2303-1452 | 2303-1457 | 2303-1453 | 2303-1458 | 2303-1454 | 2303-1459 | 2303-1450 | 2303-1455 |
| #6 | 2303-1461 | 2303-1466 | 2303-1462 | 2303-1467 | 2303-1463 | 2303-1468 | 2303-1464 | 2303-1469 | 2303-1460 | 2303-1465 |
| #7 | 2303-1471 | 2303-1476 | 2303-1472 | 2303-1477 | 2303-1473 | 2303-1478 | 2303-1474 | 2303-1479 | 2303-1470 | 2303-1475 |



XUC Insert

|    | E-XPE     |           |           |           |           |           |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 10 mm     | 11 mm     | 12 mm     | 13 mm     | 14 mm     | 15 mm     | 16 mm     | 17 mm     | 18 mm     |
| #0 | 2303-1701 | 2303-1706 | 2303-1702 | 2303-1707 | 2303-1703 | 2303-1708 | 2303-1704 | 2303-1709 | 2303-1700 | 2303-1705 |
| #1 | 2303-1711 | 2303-1716 | 2303-1712 | 2303-1717 | 2303-1713 | 2303-1718 | 2303-1714 | 2303-1719 | 2303-1710 | 2303-1715 |
| #2 | 2303-1721 | 2303-1726 | 2303-1722 | 2303-1727 | 2303-1723 | 2303-1728 | 2303-1724 | 2303-1729 | 2303-1720 | 2303-1725 |
| #3 | 2303-1731 | 2303-1736 | 2303-1732 | 2303-1737 | 2303-1733 | 2303-1738 | 2303-1734 | 2303-1739 | 2303-1730 | 2303-1735 |
| #4 | 2303-1741 | 2303-1746 | 2303-1742 | 2303-1747 | 2303-1743 | 2303-1748 | 2303-1744 | 2303-1749 | 2303-1740 | 2303-1745 |
| #5 | 2303-1751 | 2303-1756 | 2303-1752 | 2303-1757 | 2303-1753 | 2303-1758 | 2303-1754 | 2303-1759 | 2303-1750 | 2303-1755 |
| #6 | 2303-1761 | 2303-1766 | 2303-1762 | 2303-1767 | 2303-1763 | 2303-1768 | 2303-1764 | 2303-1769 | 2303-1760 | 2303-1765 |
| #7 | 2303-1771 | 2303-1776 | 2303-1772 | 2303-1777 | 2303-1773 | 2303-1778 | 2303-1774 | 2303-1779 | 2303-1770 | 2303-1775 |



MBC Insert

|    | UHMWPE    |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 11 mm     | 13 mm     | 15 mm     | 18 mm     |
| #1 | 2303-7211 | 2303-7212 | 2303-7213 | 2303-7214 | 2303-7215 |
| #2 | 2303-7221 | 2303-7222 | 2303-7223 | 2303-7224 | 2303-7225 |
| #3 | 2303-7231 | 2303-7232 | 2303-7233 | 2303-7234 | 2303-7235 |
| #4 | 2303-7241 | 2303-7242 | 2303-7243 | 2303-7244 | 2303-7245 |
| #5 | 2303-7251 | 2303-7252 | 2303-7253 | 2303-7254 | 2303-7255 |
| #6 | 2303-7261 | 2303-7262 | 2303-7263 | 2303-7264 | 2303-7265 |



MB Insert

|    | UHMWPE    |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 11 mm     | 13 mm     | 15 mm     | 18 mm     |
| #1 | 2303-7011 | 2303-7012 | 2303-7013 | 2303-7014 | 2303-7015 |
| #2 | 2303-7021 | 2303-7022 | 2303-7023 | 2303-7024 | 2303-7025 |
| #3 | 2303-7031 | 2303-7032 | 2303-7033 | 2303-7034 | 2303-7035 |
| #4 | 2303-7041 | 2303-7042 | 2303-7043 | 2303-7044 | 2303-7045 |
| #5 | 2303-7051 | 2303-7052 | 2303-7053 | 2303-7054 | 2303-7055 |
| #6 | 2303-7061 | 2303-7062 | 2303-7063 | 2303-7064 | 2303-7065 |

# Tibial Insert



PSA Insert

|    | UHMWPE    |           |           |           |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 11 mm     | 13 mm     | 15 mm     | 18 mm     | 21 mm     | 25 mm     | 30 mm     |
| #1 | 2303-5011 | 2303-5012 | 2303-5013 | 2303-5014 | 2303-5015 | 2303-5016 | 2303-5017 | 2303-5018 |
| #2 | 2303-5021 | 2303-5022 | 2303-5023 | 2303-5024 | 2303-5025 | 2303-5026 | 2303-5027 | 2303-5028 |
| #3 | 2303-5031 | 2303-5032 | 2303-5033 | 2303-5034 | 2303-5035 | 2303-5036 | 2303-5037 | 2303-5038 |
| #4 | 2303-5041 | 2303-5042 | 2303-5043 | 2303-5044 | 2303-5045 | 2303-5046 | 2303-5047 | 2303-5048 |
| #5 | 2303-5051 | 2303-5052 | 2303-5053 | 2303-5054 | 2303-5055 | 2303-5056 | 2303-5057 | 2303-5058 |
| #6 | 2303-5061 | 2303-5062 | 2303-5063 | 2303-5064 | 2303-5065 | 2303-5066 | 2303-5067 | 2303-5068 |

|    | XPE       |           |           |           |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 11 mm     | 13 mm     | 15 mm     | 18 mm     | 21 mm     | 25 mm     | 30 mm     |
| #1 | 2303-5611 | 2303-5612 | 2303-5613 | 2303-5614 | 2303-5615 | 2303-5616 | 2303-5617 | 2303-5618 |
| #2 | 2303-5621 | 2303-5622 | 2303-5623 | 2303-5624 | 2303-5625 | 2303-5626 | 2303-5627 | 2303-5628 |
| #3 | 2303-5631 | 2303-5632 | 2303-5633 | 2303-5634 | 2303-5635 | 2303-5636 | 2303-5637 | 2303-5638 |
| #4 | 2303-5641 | 2303-5642 | 2303-5643 | 2303-5644 | 2303-5645 | 2303-5646 | 2303-5647 | 2303-5648 |
| #5 | 2303-5651 | 2303-5652 | 2303-5653 | 2303-5654 | 2303-5655 | 2303-5656 | 2303-5657 | 2303-5658 |
| #6 | 2303-5661 | 2303-5662 | 2303-5663 | 2303-5664 | 2303-5665 | 2303-5666 | 2303-5667 | 2303-5668 |



PSA Insert

|    | E-XPE     |           |           |           |           |           |           |           |
|----|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm      | 11 mm     | 13 mm     | 15 mm     | 18 mm     | 21 mm     | 25 mm     | 30 mm     |
| #1 | 2303-5811 | 2303-5812 | 2303-5813 | 2303-5814 | 2303-5815 | 2303-5816 | 2303-5817 | 2303-5818 |
| #2 | 2303-5821 | 2303-5822 | 2303-5823 | 2303-5824 | 2303-5825 | 2303-5826 | 2303-5827 | 2303-5828 |
| #3 | 2303-5831 | 2303-5832 | 2303-5833 | 2303-5834 | 2303-5835 | 2303-5836 | 2303-5837 | 2303-5838 |
| #4 | 2303-5841 | 2303-5842 | 2303-5843 | 2303-5844 | 2303-5845 | 2303-5846 | 2303-5847 | 2303-5848 |
| #5 | 2303-5851 | 2303-5852 | 2303-5853 | 2303-5854 | 2303-5855 | 2303-5856 | 2303-5857 | 2303-5858 |
| #6 | 2303-5861 | 2303-5862 | 2303-5863 | 2303-5864 | 2303-5865 | 2303-5866 | 2303-5867 | 2303-5868 |



PSA Insert Low Constrained

|    | XPE, LC type |           |           |           |           |           |           |           |
|----|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|    | 9 mm         | 11 mm     | 13 mm     | 15 mm     | 18 mm     | 21 mm     | 25 mm     | 30 mm     |
| #1 | 2303-5211    | 2303-5212 | 2303-5213 | 2303-5214 | 2303-5215 | 2303-5216 | 2303-5217 | 2303-5218 |
| #2 | 2303-5221    | 2303-5222 | 2303-5223 | 2303-5224 | 2303-5225 | 2303-5226 | 2303-5227 | 2303-5228 |
| #3 | 2303-5231    | 2303-5232 | 2303-5233 | 2303-5234 | 2303-5235 | 2303-5236 | 2303-5237 | 2303-5238 |
| #4 | 2303-5241    | 2303-5242 | 2303-5243 | 2303-5244 | 2303-5245 | 2303-5246 | 2303-5247 | 2303-5248 |
| #5 | 2303-5251    | 2303-5252 | 2303-5253 | 2303-5254 | 2303-5255 | 2303-5256 | 2303-5257 | 2303-5258 |
| #6 | 2303-5261    | 2303-5262 | 2303-5263 | 2303-5264 | 2303-5265 | 2303-5266 | 2303-5267 | 2303-5268 |

# Patella Component

 Special Order Items



Inset Patella

| Size | UHMWPE    | XPE       | E-XPE     |
|------|-----------|-----------|-----------|
| S    | 2401-1010 | 2403-3010 | 2403-5010 |
| M    | 2401-1020 | 2403-3020 | 2403-5020 |
| L    | 2401-1030 | 2403-3030 | 2403-5030 |
| XL   | 2401-1040 | 2403-3040 | 2403-5040 |




Onset Patella

| Size | UHMWPE    | XPE       | E-XPE     |
|------|-----------|-----------|-----------|
| XS   | 2403-1010 | 2403-3210 | 2403-5210 |
| S    | 2403-1020 | 2403-3220 | 2403-5220 |
| M    | 2403-1030 | 2403-3230 | 2403-5230 |
| L    | 2403-1040 | 2403-3240 | 2403-5240 |
| XL   | 2403-1050 | 2403-3250 | 2403-5250 |
| XXL  | 2403-1060 | 2403-3260 | 2403-5260 |
| EL   | 2403-1070 | 2403-3270 | 2403-5270 |

# Extensions

## Femoral Augment Set



|    | Distal     |            |            |            | Distal    |           | Posterior |           |
|----|------------|------------|------------|------------|-----------|-----------|-----------|-----------|
|    | 4 mm LM/RL | 4 mm LL/RM | 8 mm LM/RL | 8 mm LL/RM | 12 mm     | 16 mm     | 4 mm      | 8 mm      |
| #1 | 2603-5111  | 2603-5211  | 2603-5112  | 2603-5212  | 2603-5313 | 2603-5314 | 2603-5011 | 2603-5012 |
| #2 | 2603-5121  | 2603-5221  | 2603-5122  | 2603-5222  | 2603-5323 | 2603-5324 | 2603-5021 | 2603-5022 |
| #3 | 2603-5131  | 2603-5231  | 2603-5132  | 2603-5232  | 2603-5333 | 2603-5334 | 2603-5031 | 2603-5032 |
| #4 | 2603-5141  | 2603-5241  | 2603-5142  | 2603-5242  | 2603-5343 | 2603-5344 | 2603-5041 | 2603-5042 |
| #5 | 2603-5151  | 2603-5251  | 2603-5152  | 2603-5252  | 2603-5353 | 2603-5354 | 2603-5051 | 2603-5052 |
| #6 | 2603-5161  | 2603-5261  | 2603-5162  | 2603-5262  | 2603-5363 | 2603-5364 | 2603-5061 | 2603-5062 |

# Extensions

 Special Order Items

## Tibial Augment



|    | 5 mm      | 10 mm     | 15 mm LM/RL | 15 mm LL/RM |
|----|-----------|-----------|-------------|-------------|
| #0 | 2803-5201 | 2803-5202 | N/A         | N/A         |
| #1 | 2803-5211 | 2803-5212 | 2803-5113   | 2803-5213   |
| #2 | 2803-5221 | 2803-5222 | 2803-5123   | 2803-5223   |
| #3 | 2803-5231 | 2803-5232 | 2803-5133   | 2803-5233   |
| #4 | 2803-5241 | 2803-5242 | 2803-5143   | 2803-5243   |
| #5 | 2803-5251 | 2803-5252 | 2803-5153   | 2803-5253   |
| #6 | 2803-5261 | 2803-5262 | 2803-5163   | 2803-5263   |
| #7 | 2803-5271 | 2803-5272 | N/A         | N/A         |

## Offset Adaptor



|      |           |
|------|-----------|
| 2 mm | 2903-3010 |
| 4 mm | 2903-3020 |
| 6 mm | 2903-3030 |

## PSA Femoral Screw



|            |           |
|------------|-----------|
| M5 x 14 mm | 2903-1014 |
|------------|-----------|

## Extension Stem



|     | Straight Stem |           |           |           |           | Curved Stem |           |
|-----|---------------|-----------|-----------|-----------|-----------|-------------|-----------|
|     | 30 mm         | 75 mm     | 100 mm    | 150 mm    | 200 mm    | 150 mm      | 200 mm    |
| Ø10 | N/A           | 2703-5011 | 2703-5021 | 2703-5051 | 2703-5061 | 2703-5031   | 2703-5041 |
| Ø12 | N/A           | 2703-5012 | 2703-5022 | 2703-5052 | 2703-5062 | 2703-5032   | 2703-5042 |
| Ø14 | 2703-5003     | 2703-5013 | 2703-5023 | 2703-5053 | 2703-5063 | 2703-5033   | 2703-5043 |
| Ø16 | N/A           | 2703-5014 | 2703-5024 | 2703-5054 | 2703-5064 | 2703-5034   | 2703-5044 |
| Ø18 | N/A           | 2703-5016 | 2703-5026 | 2703-5055 | 2703-5066 | 2703-5036   | 2703-5046 |
| Ø20 | N/A           | 2703-5016 | 2703-5026 | 2703-5056 | 2703-5066 | 2703-5036   | 2703-5046 |

Note:

5 mm / 10 mm augments and 30 mm straight stem are compatible with CMA Baseplate #0, #7 augments are not available with PSA baseplate



*Each Step* **We Care**

#### Contact Us

■ **Taiwan**

**United Orthopedic Corporation**

12F, No.80, Sec.1, Chenggong Rd., Yonghe Dist.,  
New Taipei City 23452, Taiwan

Tel: +886 2 2929 4567 Fax: +886 2 2922 4567

■ **China**

**United Medical Instrument Co., Ltd.**

Room 402, Building 7, No.697, Ling Shi Road,  
ShangHai 200072, China

Tel: +86 21 3661 8055 Fax: +86 21 6628 3031

■ **Japan**

**United Biomech Japan Inc.**

Ginyo Bldg 5th floor, 2-9-40 Kitasaiwai, Nishi-Ku  
Yokohama City, Kanagawa, 220-0004, Japan

Tel: +81 45 620 0741 Fax: +81 45 620 0742

■ **USA**

**UOC USA Inc.**

20, Fairbanks, Suite 173, Irvine 92618, CA

Tel: +1 949 328 3366 Fax: +1 949 328 3367

■ **EU Representative**

**mdi Europa GmbH**

Langenhagener Strasse 71, Langenhagen 30855, Germany

Tel: +49 511 3908 9530 Fax: +49 511 3908 9539

■ **Europe**

**UOC (Suisse) SA**

Avenue Général Guisan 60A, 1009 Pully, Switzerland

Tel: +41 21 634 70 70 Fax: +41 21 634 70 76

**UOC (France) SAS**

7 Allée des Peupliers, 54180 Houdemont, France

Tel: +33 3 83 23 39 72 Fax: +33 3 83 23 39 10

**UOC (UK) Limited**

Baileys Barn, Pimlico Farm, Austrey Lane,  
No Mans Heath, Warwickshire, B79 0PF, UK

Tel: +44 1827 214773



[www.uoc.com.tw](http://www.uoc.com.tw)