


Vascular Intervention // **Peripheral**
PTA Balloon Catheter / 0.035" / OTW

Passeo[®]-35 Xeo

Low profile PTA balloon with excellent deliverability in a wide range of sizes.



Improved crossability




Excellent deliverability



Low profile, wide range of sizes

Passeo[®]-35 Xeo

Vascular
Intervention
Peripheral



The Passeo-35 Xeo is indicated to dilate stenosis in the iliac*, femoral, popliteal and infrapopliteal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. Passeo-35 Xeo is also recommended for post-dilation of balloon expandable and self-expanding stents in the peripheral vasculature.**

Technical Data - Balloon Catheter

Catheter type	OTW	Balloon markers	2 swaged markers
Recommended guide wire	0.035"	Sizes	ø 3.0-12.0 mm; L: 20-250 mm
Tip	Low entry profile, colored	Shaft	5.1-5.4F, dual-lumen, hydrophobic coating
Balloon material	SCP (Semi-Crystalline Polymer), controlled compliance	Usable length	90, 130 and 170 cm
Balloon coating	Hydrophobic patchwork coating	Guide wire lumen	Hydrophobic coating

	Balloon ø (mm)	NP ^a	RBP ^a	Catheter length 90 cm									
				Balloon length (mm)									
				20	40	60	80	100	120	150	170	200	250
5F	3.0	9	21	428777	428786	428795	428804	-	-	-	-	-	-
	4.0	9	20	428778	428787	428796	428805	428814	428823	-	-	-	-
	5.0	9	19	428779	428788	428797	428806	428815	428824	428833	-	428843	428848
	6.0	9	16	428780	428789	428798	428807	428816	428825	428834	-	428844	428849
	7.0	9	15	428781	428790	428799	428808	428817	-	-	-	-	-
6F	8.0	9	14	428782	428791	428800	428809	428818	-	-	-	-	-
	9.0	9	12	428783	428792	428801	428810	428819	-	-	-	-	-
	10.0	8	11	428784	428793	428802	428811	428820	-	-	-	-	-
7F	12.0	8	10	-	428794	428803	428812	428821	-	-	-	-	-

	Balloon ø (mm)	NP ^a	RBP ^a	Catheter length 130 cm									
				Balloon length (mm)									
				20	40	60	80	100	120	150	170	200	250
5F	3.0	9	21	428851	428860	428869	428878	428887	428896	428905	428910	428915	428920
	4.0	9	20	428852	428861	428870	428879	428888	428897	428906	428911	428916	428921
	5.0	9	19	428853	428862	428871	428880	428889	428898	428907	428912	428917	428922
	6.0	9	16	428854	428863	428872	428881	428890	428899	428908	428913	428918	428923
	7.0	9	15	428855	428864	428873	428882	428891	428900	428909	428914	428919	428924
6F	8.0	9	14	428856	428865	428874	428883	428892	428901	-	-	-	-
	9.0	9	12	428857	428866	428875	428884	428893	428902 ^b	-	-	-	-
	10.0	8	11	428858	428867	428876	428885	428894	428903 ^b	-	-	-	-
7F	12.0	8	10	-	428868	428877	428886	428895	428904 ^b	-	-	-	-

	Balloon ø (mm)	NP ^a	RBP ^a	Catheter length 170 cm									
				Balloon length (mm)									
				20	40	60	80	100	120	150	170	200	250
5F	4.0	9	20	-	428935 ^b	428944 ^b	428953 ^b	-	428971 ^b	428980 ^b	-	428990 ^b	428995 ^b
	5.0	9	19	428927 ^b	428936 ^b	428945 ^b	428954 ^b	-	428972 ^b	428981 ^b	-	428991 ^b	428996 ^b
	6.0	9	16	428928 ^b	428937 ^b	428946 ^b	428955 ^b	-	428973 ^b	428982 ^b	-	428992 ^b	428997 ^b
	7.0	9	15	428929 ^b	428938 ^b	428947 ^b	428956 ^b	-	-	-	-	-	-
	8.0	9	14	428930 ^b	428939 ^b	428948 ^b	428957 ^b	-	-	-	-	-	-
6F	9.0	9	12	-	428940 ^b	428949 ^b	428958 ^b	-	-	-	-	-	-
	10.0	8	11	-	428941 ^b	428950 ^b	428959 ^b	-	-	-	-	-	-

^aNP = Nominal Pressure and RBP = Rated Burst Pressure in atm (1 atm = 1.013 bar)
^b8 weeks pre-order only

*Note for Australia: Passeo-35 Xeo not approved by Therapeutic Goods Administration for use in the common iliac arteries. **Indication as per IFU.

Passeo and Xeo are trademarks or registered trademarks of the BIOTRONIK Group of Companies.

BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com

© 2022 BIOTRONIK AG – All rights reserved.
Specifications are subject to modification, revision and improvement.

BIOTRONIK
excellence for life