

Lot 21



Sodium Hyaluronate 1.8% w/v (Can also be supplied in 1.0%)

Salient Features :

- Supplied Sterile in a 1ml Pre Filled Syringe with a 27G cannula
- This is a sterile, non-pyrogenic, viscoelastic, preparation, of sodium hyaluronate. It contains 18mg/ml of sodium hyaluronate dissolved in physiological sodium chloride phosphate buffer. (Other concentrations like 1.0% also available)
- Sodium Hyaluronate is derived from bacterial fermentation. Its viscosity is responsible for maintenance of space in the anterior chamber and allows the manipulation of tissue during the surgery without escape through the incision.
- The 1.8% solution is transparent . It does not interfere with epithelization and normal wound healing.
- Uses : It is indicated for use as a surgical aid in cataract extraction (intra and extra capsular extraction), IOL implant, corneal transplant, glaucoma filtration and retinal detachment surgery.
- Contraindication : At present there are no known contraindications of the use of this sodium hyaluronate.
- Cost economical.

Another Quality Product from :

Action

Action Medical Mktg Pvt Ltd.

(An ISO 9001:2008 & ISO 13485:2003 Certified Company)

Gat No 528, Koregaon-Bhima, Near SBI, Pune-Nagar Road

Taluka Shirur, Dist Pune Pin 412216, INDIA

Tel : 0091-2137-253617 / 18. Fax : 0091-2137-253620.

E-mail : action@actionmedihealth.com

Website : www.actionmedihealth.com

Lot 22



ActionVisc

2%w/v Hydroxypropyl Methylcellulose OPHTHALMIC VISCOSURGICAL DEVICE

Description

The Ophthalmic Viscosurgical Device (OVD) is a sterile, isotonic, apyrogenic viscoelastic solution of highly purified, non inflammatory Hydroxy Propyl Methyl Cellulose (HPMC) solution with a high molecular weight greater than 80,000 daltons for intraocular use during the surgery of the anterior segment of the eye. For a 2% w/v HPMC solution each ml of solution provides 20mg of HPMC dissolved in balanced salt solution. The viscosity of this solution is 4000 to 5500 cps

Indications

The OVD is indicated for use as an ophthalmic surgical aid in anterior segment surgical procedures, such as cataract extraction and intraocular lens implantation. The OVD maintains a deep chamber during anterior segment surgery and thereby allowing for more efficient manipulation with significantly less trauma to the corneal endothelium and other ocular tissues

Contraindications

At present, there are no known contraindications, when used by a trained Ophthalmic Surgeon, as indicated.

Warnings

Do not use the package if the integrity of the sterile package is compromised.
The OVD is for intraocular use only.
The OVD is for single use only.
The device is not to be placed in the direct sunlight or temperature exceeding 40°C.
Do not resterilize the device

Storage

Store at a temperature between 15° Celsius and 35° Celsius and a relative humidity less than 65%. Protect from freezing.
Do not expose containers to direct sunlight

Supply is a sterile solution of 2%w/v Hydroxypropyl Methylcellulose is available in 2 ml ,3 ml PFS in pouch pack, 5 ml Vial with pouch in pack of ten.

FOR YOUR SUPPLIES CONTACT :

Action Medical Mktg Pvt Ltd.
Gat No 528, Korvegaon, Bhima, Tal Shirur, Dist Pune, INDIA.
Tel : 0091-2137-253617 / 18

www.actionmedilifealth.com

Lot 23

ActionBLU :

(Trypan Blue Ophthalmic Solution 0.06% w/v)



Description:

ActionBlu 0.06% is a dye which safely stains the anterior lens capsule & is mainly used in cataract surgery. Special feature of Optiblu is that it stains the anterior capsule without affecting the corneal endothelium.

Indication:

ActionBlu is indicated for facilitating surgical procedure for cataract extraction (IOL Surgery and phacosurgery) and as a staining agent.

Contraindication:

ActionBlu is contraindicated in patients with hypersensitivity to the dye or any of its components.

Warning & Precautions:

- Do not use the product after the Expiry date
- Do not use on children
- For single use only
- Not for injection

Storage:

Store at a temperature between 15° Celsius and 35° Celsius and relative humidity less than 65% protect from freezing.

Do not expose containers to direct sunlight

Supply:

ActionBlu is a sterile solution is available in 1ml with pouch in pack of Ten and ActionBlu PFS 1ml in sterile glass syringe.

Another Quality Product from:



Action Medical Pvt Ltd.

(An ISO 9001:2008 & ISO 13485:2003 Certified Company)

Gat No 528, Koregaon-Bhima, Near SBI, Pune-Nagar Road

Taluka Shirur, Dist Pune Pin 412216, INDIA

Tel: 0091-2137-253617 / 18. Fax: 0091-2137-253620.

E-mail : action@actionmedihealth.com

Website: www.actionmedihealth.com