


SV600/650/800/850 Ventilator

Service Manual

Intellectual Property Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this product and this manual. This manual may refer to information protected by copyrights or patents and does not convey any license under the patent rights of Mindray, nor the rights of others. Mindray does not assume any liability arising out of any infringements of patents or other rights of third parties.

Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden. Release, amendment, reproduction, distribution, rent, adaption and translation of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

mindray、 and **MINDRAY** are the registered trademarks or trademarks owned by Mindray in China and other countries. All other trademarks that appear in this manual are used only for editorial purposes without the intention of improperly using them. They are the property of their respective owners.

Contents of this manual are subject to changes without prior notice.


Revision History

This manual has a revision number. This revision number changes whenever the manual is updated due to software or technical specification change. Contents of this manual are subject to change without prior notice. Revision 1.0 is the initial release of the document.

- Revision number: 5.0
- Release time: 2021-09

© Copyright 2017-2021 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

11.3 Display Assembly



11.4 Parts List

Item No.	Order Number	Part Descriptions	QTY	Remarks
1	043-004516-00	EV20 Knob	1	/
2	115-049775-00	SV600 Front cover FRU (light column)	1	refer to the machine's configuration
	115-049776-00	SV650 Front cover FRU (light column)		
	115-049777-00	SV800 Front cover FRU (light column)		
	115-049778-00	SV850 Front cover FRU (light column)		
3	051-002524-00	VT70 alarm light PCBA	1	/
4	049-001311-00	Silence button (VT70-SY)	1	/
5	051-002526-00	VT70 button and encoder PCBA	1	/
6	049-001310-00	Power Button (VT70-SY)	1	/
7	051-002525-00	VT70 switch touch screen control PCBA	1	/
8	021-000313-00	18.5 inch FIVE WIRE RESISTIVE TOUCH PANEL (ELO)	1	refer to the machine's configuration
	021-000310-00	15.6 inch FIVE WIRE RESISTIVE TOUCH PANEL (ELO)		
9	115-049780-00	18.5 inch screen FRU (cable)	1	refer to the machine's configuration
	115-049779-00	15.6 inch screen FRU (cable)		
10	020-000014-01	Speaker 2W 4ohm 500Hz	1	/
11	115-065359-00	VT70 main control PCBA (software/battery)	1	Need to provide machine's SN and function configuration
12	043-012117-00	Back Bezel of 15.6 Inch Monitor (VT70-SY)	1	refer to the machine's configuration
	043-012120-00	Back Bezel of 18.5 Inch Monitor (VT70-SY)		
13	115-038569-00	screen arm(VT70)	1	/