



Additional staining reagents and solutions

Fuchsin Acid Van Gieson reagent

For use in special staining kits.
Synonym: Picrofuchsin reagent.

100 mL	FAG-OT-100
--------	------------

Fuchsin Basic reagent

For use in bacteriology.

100 mL	FBR-OT-100
500 mL	FBR-OT-500
1000 mL	FBR-OT-1L

Fuchsin Basic, solution

For use with H.B.F.P. kit.

100 mL	FBO-OT-100
--------	------------

Gelatin, gel in transport jar

For preparing developing solution with Warthin Starry staining kit.

12x8 mL	GEL-128
---------	---------



Glutaraldehyde, 4% buffered solution

Solution used with Sudan Black B kit.

100 mL	G4-OT-100
--------	-----------



Gold chloride solutions

Aqueous solutions of gold chloride (aurichloric acid, HAuCl_4), component of many special staining kits.

Gold chloride, 0.2% solution	100 mL	ZK02-OT-100
Gold chloride, 0.6% solution	100 mL	ZK06-OT-100
Gold chloride, 1% solution	50 mL	ZK1-OT-50
Gold chloride, 2% solution	50 mL	ZK2-OT-50

Hematoxylin, Verhoeff A

For use with Verhoeff and Movat kits.

100 mL	HEMV-OT-100
--------	-------------



Hematoxylin, Weigert A

Hematoxylin solution in ethanol, component of Hematoxylin W kit.
Part of several special staining kits.

100 mL	HEMA-OT-100
500 mL	HEMA-OT-500
1000 mL	HEMA-OT-1L



Hydroquinone, solution

For preparation of developing solution in Warthin Starry kit.

50 mL	HQ-OT-50
-------	----------



Hydroquinone-sodium sulphite

For use with Grimelius kit

15 x 3,5 g	HNS-15-35
------------	-----------

Iodine solution, Verhoeff C

For use with Verhoeff and Movat kits.

100 mL	JODV-OT-100
--------	-------------

Iron chloride, solution

For use with Fouchet-Van Gieson kit.

100 mL	FK-OT-100
--------	-----------

Kovac indole reagent

Reagent for identification of indole-positive and indole-negative microorganisms.

100 mL	KOV-OT-100
--------	------------



Ladewig's reagent

For use with Ladewig Trichrome kit, contains Aniline Blue, Acid Fuchsin and Orange G dyes.

100 mL	LDWR-OT-100
--------	-------------

Light Green S.F. reagent

Counterstain used with trichrome special staining kits.

100 mL	LGS-OT-100
--------	------------

Lithium carbonate, Luxol

Lithium carbonate solution used with Luxol Fast Blue kit.

100 mL	LKL-OT-100
--------	------------

Lithium carbonate, solution

Lithium carbonate solution used with Von-Kossa staining kit.

100 mL	LK-OT-100
--------	-----------

Lugol's solution, M.I.F.

For use with M.I.F. kit.

100 mL	LUGM-OT-100
--------	-------------

Luxol Fast Blue, solution

For use with Luxol Fast Blue kit.

100 mL	LFB-OT-100
--------	------------



M.F. reagent

For use with M.I.F. kit.

100 mL	MF-OT-100
--------	-----------

Martius Yellow, solution

For use with MSB kit.

100 mL	MAY-OT-100
--------	------------

Mayer's mucicarmine

For use with Mucicarmine kit.

100 mL	MMU-OT-100
--------	------------



Metanil Yellow, solution

For use with Mucicarmine kit.

100 mL	MY-OT-100
--------	-----------