Prestige Plus Complex Xtra Pack

Prestige Plus Helical Xtra Pack

Model No.	Implant size
PRES0202CXPPL	2mm x 2cm
PRES0204CXPPL	2mm x 4cm
PRES0206CXPPL	2mm x 6cm
PRES0307CXPPL	3mm x 7cm
PRESO310CXPPL	3mm x 10cm

Plus Complex Pack

Model No.	Implant size
PRESO204HXPPL	2mm × 4cm
PRESO206HXPPL	2mm x 6cm
PRESO208HXPPL	2mm x 8cm
PRESO310HXPPL	3mm × 10cm
PRESO315HXPPL	3mm x 15cm

Prestige

Model No.	Implant size
PRESO407CPKPL	4mm x 7cm
PRESO410CPKPL	4mm x 10cm
PRESO413CPKPL	4mm x 13cm
PRESO510CPKPL	5mm x 10cm
PRESO517CPKPL	5mm x 17cm
PRESO620CPKPL	6mm x 20cm
PRESO720CPKPL	7mm x 20cm
PRESO830CPKPL	8mm x 30cm
PRESO930CPKPL	9mm x 30cm
PRES1030CPKPL	10mm x 30cm
PRES1040CPKPL	10mm x 40cm

Complex

Model No.	Implant size
PRES1240COM18	12mm x 40cm
PRES1447COM18	14mm x 47cm
PRES1655COM18	16mm x 55cm
PRES1860COM18	18mm x 60cm
PRES2065COM18	20mm x 65cm
PRES2465COM18	24mm x 65cm

Prestige Plus Helical Pack

Model No.	Implant size
PRESO410HPKPL	4mm x 10cm
PRESO420HPKPL	4mm x 20cm
PRESO515HPKPL	5mm x 15cm
PRES0520HPKPL	5mm x 20cm
PRESO620HPKPL	6mm x 20cm
PRESO720HPKPL	7mm x 20cm
PRESO830HPKPL	8mm x 30cm
PRESO930HPKPL	9mm x 30cm
PRES1040HPKPL	10mm x 40cm

ULTRA Detachment Controller

Model No.	Box content
ULTRA	1 detachment controller

The Prestige Coil System is intended for use in the peripheral vasculature to endovascularly obstruct or occlude blood flow in vascular abnormalities of the peripheral vassels. Class III CE0297 in compliance with Medical Device Directive (MDD 93/42/ EEC amended by 2007/47/EC). Manufactured by BALT USA LLC, 29 Parker, Invine, CA 92618. Carefully read the instructions for use before use. First CE marking: 2021. For a list of indications, warnings, precautions and contraindications, refer to IFU-060.


The content of this document, in particular data, information, trademarks and logos is BALT S.A.S and diffiliates' sole property. Consequently, all representation and/ or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT S.A.S and affiliates' copyrights and other intellectual proprietary rights. 2020? BALT S.A.S and affiliates all rights reserved. This document with associated prictures is non-contractual and is solely dedictated to healthcare professionals and BALT S.A.S and affiliates' distributors. The products commercialized by BALT S.A.S and affiliates shall exclusively be used in accordance with the package inserts included in the boxes. DC05768 (05/2021)





prestige plus

Intended for use in the peripheral vasculature to endovascularly obstruct or occlude blood flow in vascular abnormalities of the peripheral vessels



entire size range deliverable through

.0165"-.021/.022" catheters

Design

Premium coil softness through variations of filar and primary wind for deliverability and conformability

Body coil/hypotube pusher design provides increased ease of delivery, even in tortuous anatomy

Smooth, precise embolization thanks to implant softness and enhanced pusher technology

Optimal detachment

Instantaneous detachment in under 1 second* allows for rapid treatment & efficient case times

Reliable detachment* system provides confidence during coiling cases

Complete offer

With a full range of Complex and Helical coils to accommodate a variety of clinical application

Unique large volume coils*
(Complex 1 8 and Complex Pack)
for dense packing






Implant design

- ✓ Platinum, non fibered coil → fully visible
- ✓ Complete set of configurations: Helical & Complex
- ✓ Deliverable through microcatheters with ID = 0.0165" 0.027"
- ✓ Second radiopaque mark: possibility to use it with double or single mark microcatheters.
- ✓ New 185 cm total lenght.




helical



The content of this document, in particular data, information, trademarks and logos is BALT SAS and affiliates' sole property. ©2018 BALT and affiliates, all rights reserved. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT SAS and affiliates' copyrights and other intellectual proprietary rights. This document with associated pictures is non-contractual and is solely dedicated to BALTs distributors (BALT supplier's distributors). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. The distribution, the sale and the use of BALTs products is under distributor's liability.






The content of this document, in particular data, information, trademarks and logos is BALT SAS and affiliates' sole property. @2018 BALT and affiliates, all rights reserved. Consequently, all representation and/or reproduction, whether in part or in fulf, is forbidden and would be considered a violation of BALT SAS and affiliates' copyrights and other intellectual proprietary rights. This document with associated pictures is non-contractual and is solely dedicated to BALT's distributors (BALT's supplier's distributors). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. The distribution, the sale and the use of BALT's products is under distributor's liability.

EXTERNAL USE -CONFIDENTIAL

Components



Prestige coil system

To use with (not included in the package):

- ✓ ULTRA Detachment Controller
- ✓ Compatible guidewire
- ✓ Delivery microcatheter I.D. between .0165" and .027"
- ✓ Compatible sheat / guiding catheter

The content of this document, in particular data, information, trademarks and logos is BALT SAS and affiliates' sole property. ©2018 BALT and affiliates, all rights reserved. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT SAS and affiliates' copyrights and other intellectual proprietary rights. This document with associated pictures is non-contractual and is solely dedicated to BALT's distributors (BALT's supplier's distributions). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. The distribution, the sale and the use of BALT's products is under distributor's liability.

