

Absorbent underpad

TD-03-I.2.c-1.1

edition I

date: 2021-05-26 amendment of: N/A

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com

CLASSIFICATION

Class I medical device compliant with Medical Devices Regulation (EU) 2017/745 (MDR)

INTENDED PURPOSE

 disposable device is used to absorb urine and other body fluids. It is used on the medical bed to directly protect the patient from contamination and infection

FEATURES

- quilted absorbent insert
- composed of layers: polypropylene non-woven, cellulose wadding, cellulose pulp, cellulose wadding, PE foil
- white absorbent insert
- blue foil layer
- single-use
- non-sterile
- shelf life: 5 years



| CATALOGUE NUMBER | DIMENSIONS | COLOUR | TYPE OF PACKAGING | SALES PACKAGING | Transport Packaging |
|---------------------|------------|------------|-------------------|--------------------|------------------------|
| BT-PC4060 | 40 x 60 cm | white/blue | foil | 25 pcs | 4 x 25 pcs |
| BT-PC6090 | 60 x 60 cm | white/blue | foil | 25 pcs | 4 x 25 pcs |
| BT-PC6090 | 60 x 90 cm | white/blue | foil | 25 pcs | 4 x 25 pcs |