



174200

174430

174700

174800

174200

COTTLE Metal Mallet,
length 18 cm

174420

Metal Mallet, small model,
diameter 20 mm, length 19.5 cm

174426

Metal Mallet, medium model,
diameter 26 mm, length 19.5 cm

174430


Metal Mallet, large model,
diameter 30 mm, length 19.5 cm

174700

Metal Mallet, with plastic
replacement head, small,
lightweight model, length 18 cm

174800

Metal Mallet, with plastic
replacement head, length
22.5 cm



175840 –
175804

176040 –
176001



175840

SPRATT Mastoid Curette,
size 0000, length 16.5 cm



175830

Same, size 000



175820

Same, size 00



175810

Same, size 0



175801

SPRATT Mastoid Curette,
size 1, length 17 cm



175802

Same, size 2



175803

Same, size 3



175804

Same, size 4



176040

LEMPERT Mastoid Curette,
size 0000, length 21 cm



176030

Same, size 000



176020

Same, size 00



176010

Same, size 0





176001

Same, size 1

Components/Spare Parts see chapter 16



Malleus Nippers, working length 8 cm, dismantling



222800

HOUSE-DIETER Malleus Nipper, dismantling, upbiting, working length 8 cm



222900

Same, downbiting



222801

FISCH-DIETER Malleus Nipper, dismantling, extra delicate, upbiting, width 0.8 mm, working length 8 cm



223000

Instrument Guard, Teflon, autoclavable, for micro forceps and micro scissors with working length 8 cm, length 11.5 cm



Scissors, working lengths 8 cm and 9 cm



with cleaning connector



222401

WULLSTEIN Scissors, delicate, straight, with cleaning connector, working length 9 cm



222402

Same, curved upwards



222403

Same, curved to left



222404



Same, curved to right



222405

WULLSTEIN Scissors, extra delicate, straight, with cleaning connector, working length 8 cm

Scissors and Micro Scissors, working length 8 cm, dismantling



dismantling



222500

BELLUCCI Scissors, dismantling, cutting length 8 mm, working length 8 cm



222501

FISCH-BELLUCCI Scissors, dismantling, extra fine blade, cutting length 7 mm, working length 8 cm



222600

BELLUCCI Scissors, dismantling, delicate standard model, cutting length 8 mm, working length 8 cm



222601

FISCH-BELLUCCI Micro Ear Scissors, dismantling, beaked scissors blades, cutting length 7 mm, working length 8 cm



222602

HOUSE-BELLUCCI Scissors, dismantling, extra delicate, working length 8 cm




222604 R

BELLUCCI Scissors, dismantling, delicate, curved to right, rounded, working length 8 cm



222605 L

Same, curved to left



221100

HARTMANN Ear Forceps, extra delicate, serrated, 1 x 4.5 mm, working length 8 cm



221150

Same, working length 12.5 cm



221105

FISCH-HARTMANN Ear Forceps, delicate, serrated, 1.2 x 8 mm, working length 8 cm

221107

Same, 1.2 x 12 mm



221110

FISCH-HARTMANN Ear Forceps, extra delicate, smooth, 1 x 4.5 mm, working length 8 cm



221111

FISCH Ear Forceps, extra delicate, smooth, 0.4 x 3.5 mm, working length 8 cm



221111 FN

FISCH Ear Forceps, very fine, smooth, with antireflection coating, 0.4 x 3.5 mm, working length 8 cm



221112

Ear Forceps, straight, with spring, smooth, 4 mm, working length 8 cm



221200

FISCH Ear Forceps, extra delicate, pointed, serrated, 1 x 4.5 mm, working length 8 cm



221201

FISCH Ear Forceps, extra delicate, serrated, 0.4 x 3.5 mm, working length 8 cm



221304

Ear Forceps, extra delicate, serrated, curved to right, working length 8 cm




221305

Same, curved to left



221307

Same, curved upwards




221800

HEERMANN Ear Forceps, straight, with cleaning connector, working length 10 cm



222000

Same, extra delicate, working length 8 cm



Item no: 451501B

RHINOFORCE II Nasal Forceps, size 1

Quantity:

–

1


+

On the offer list

Technical details	Documents & Films	Product group / field of application
-------------------	-------------------	--------------------------------------

Author	GRÜNWALD-HENKE
Width of distal end	3.5 mm
Cutting property	Through-cutting, smooth cut
Instrument shape	BLAKESLEY shape
Working length	13 cm
Brand	RHINOFORCE®
Size	1
Jaw curvature	Upturned
Jaw angle	45°
Cleaning connector	Yes

STAMMBERGER Antrum Punch, backwards cutting



459010

STAMMBERGER RHINOFORCE® II Antrum Punch,
upside backward cutting, with cleaning connector,
working length 10 cm




459011

Same, right side backward cutting



459012

Same, left side backward cutting



459016

STAMMBERGER Antrum Punch, backward
cutting, sheath 360° rotatable, with fixing screw,
dismantling, working length 10 cm, for use with
Cleaning Adaptor 459015 LL



459015 LL

Cleaning Adaptor

[Home page](#) > [Online Catalog](#)



Item no: 452802

HEUWIESER Nasal Forceps, 45°

Quantity:

–

1

+

Add to My Quote List 



Product details

- Particularly for frontal sinuses and frontal recess

Technical details		Documents & Films	Product group / field of application
Author	HEUWIESER		
Cutting property	Through-cutting		
Model	Extra delicate		
Working length	12.5 cm		
Brand	RHINOFORCE®		
Jaw opening	Backward opening		
Cutting width	1.5 mm		
Cleaning connector	Yes		
Bending angle sheath	45°		
Sheath curvature distal	Upturned		

Standard Instruments for Rhinology and Rhinoplasty

Cotton Applicators, Gauze Applicator, Balloon Catheters,
Nasal Dressing Forceps, Suction Tube




- 417013 **Cotton Applicator**, standard model, triangular, serrated, diameter 1.3 mm, length 15 cm
- 417017 **Same**, diameter 1.7 mm
- 417217 **ZIEGLER Cotton Applicator**, with thread, diameter 1.7 mm, length 17.5 cm
- 417600 **FROHN Cotton Applicator**, with thread, length 16 cm
- 420000 **SIEBENMANN Gauze Applicator**, length 18 cm
- 423700 **REUTER Epistaxis Balloon Catheter**, sterile, for single use
- 426516 JANSSEN Nasal Dressing Forceps**, bayonet-shaped, length 16.5 cm
- 426620 **GRÜNWALD Nasal Dressing Forceps**, bayonet-shaped, length 20 cm
- 426900 **TRÖLTSCH Nasal Dressing Forceps**, angular, working length 10 cm
- 203835 **Suction Tube**, conical, LUER, outer diameter 3.5 mm, working length 11 cm

Suction and Irrigation Tubes

Size 10 mm

Operating instruments, **length 36 cm**,
for use **with trocars size 11 mm**




26120 S



26120 S

Ascites Suction Tube, size 10 mm,
length 36 cm



28175 YR



28175 YR

OLSEN Suction or Irrigation Tube,
with protective basket, with trumpet
valve, size 10 mm, length 36 cm


Please note:

Units for use with Handles for Suction and Irrigation see chapter 20, UNITS

Trocars

Size 6 mm

For use with telescopes and instruments size 5 mm



Multifunctional Valve



Automatic Valve

Size:		6 mm			
Working length:		8.5 cm	10.5 cm	8.5 cm	10.5 cm
Color code:		black-white	black	black-white	black
		Multifunctional Valve		Automatic Valve	
	Trocar, with conical tip including: Cannula Trocar only Valve	30160 WY	30160 MC	30160 FY	30160 AC
		30160 H1	30160 H2	30160 H1	30160 H2
		30160 Y	30160 C	30160 Y	30160 C
		30160 M1	30160 M1	30160 A1	30160 A1
	Trocar, with pyramidal tip including: Cannula Trocar only Valve	30160 WZ	30160 MP	30160 FZ	30160 AP
		30160 H1	30160 H2	30160 H1	30160 H2
		30160 Z	30160 P	30160 Z	30160 P
		30160 M1	30160 M1	30160 A1	30160 A1
	Trocar, with blunt tip including: Cannula Trocar only Valve	30160 WX	30160 MA	30160 FX	30160 AA
		30160 H1	30160 H2	30160 H1	30160 H2
		30160 X	30160 A	30160 X	30160 A
		30160 M1	30160 M1	30160 A1	30160 A1
	Trocar including: Cannula Trocar only Valve	—	30160 MB	—	30160 AB
			30160 H2		30160 H2
			30160 B		30160 B
			30160 M1		30160 A1


Telescopes for use with trocars size 6 mm see page 33

HEINKEL-SEMM Dilation Set see page 78, Accessories for Trocars see pages 76-77

Trocars for Adipose Patients

Size 6, 11 and 13.5 mm

XL



Size:		6 mm	11 mm	13.5 mm
Working length:		15 cm	15 cm	15 cm
Color code:		black-red	green-red	blue-red
	Trocar, with conical tip including: Cannula Trocar only Multifunctional Valve	31160 MC 31160 H1 31160 C 30160 M1	31103 MC 31103 H1 31103 C 30103 M1	—
	Trocar, with pyramidal tip including: Cannula Trocar only Multifunctional Valve	31160 MP 31160 H1 31160 P 30160 M1	31103 MP 31103 H1 31103 P 30103 M1	31108 MP 31108 H1 31108 P 30108 M1
	Trocar, with blunt tip including: Cannula Trocar only Multifunctional Valve	31160 MA 31160 H1 31160 A 30160 M1	31103 MA 31103 H1 31103 A 30103 M1	—

Telescopes for use with trocars size 6 mm see page 33



Telescopes for use with trocars size 11 mm see pages 31, 35, 38-40

HEINKEL-SEMM Dilation Set see page 78, Accessories for Trocars see pages 76-77

Trocars

Size 11 mm

For use with telescopes and instruments size 11 mm



Multifunctional Valve



Automatic Valve


Size:		11 mm			
Working length:		8.5 cm	10.5 cm	8.5 cm	10.5 cm
Color code:		green-white	green	green-white	green
		Multifunctional Valve		Automatic Valve	
	Trocar, with conical tip including: Cannula Trocar only Valve	30103 WY	30103 MC	30103 FY	30103 AC
		30103 H1	30103 H2	30103 H1	30103 H2
		30103 Y	30103 C	30103 Y	30103 C
		30103 M1	30103 M1	30103 A1	30103 A1
	Trocar, with pyramidal tip including: Cannula Trocar only Valve	30103 WZ	30103 MP	30103 FZ	30103 AP
		30103 H1	30103 H2	30103 H1	30103 H2
		30103 Z	30103 P	30103 Z	30103 P
		30103 M1	30103 M1	30103 A1	30103 A1
	Trocar, with blunt tip including: Cannula Trocar only Valve	30103 WX	30103 MA	30103 FX	30103 AA
		30103 H1	30103 H2	30103 H1	30103 H2
		30103 X	30103 A	30103 X	30103 A
		30103 M1	30103 M1	30103 A1	30103 A1
	MOTSON Trocar, with conical, blunt tip, including: Cannula Trocar only Valve	—	30103 MMO	—	30103 AMO
			30103 H2		30103 H2
			30103 MO		30103 MO
			30103 M1		30103 A1
	Trocar including: Cannula Trocar only Valve	—	30103 WSB	—	30103 FSB
			30103 H2		30103 H2
			30103 SB		30103 SB
			30103 M1		30103 A1

Telescopes for use with trocars size 11 mm see pages 31, 35, 38-40
HEINKEL-SEMM Dilation Set see page 78, Accessories for Trocars see pages 76-77

Fan Retractors

Size 5 mm

Operating instrument, **length 36 cm**,
for use **with trocars size 6 mm**




30623 F

Fan Retractor, dismantling, retractable,
size 5 mm, length 36 cm, opening angle in
extended position approx. 16°

Size 10 mm

Operating instrument, **length 36 cm**,
for use **with trocars size 11 mm**



30623 FP

Fan Retractor, dismantling, retractable,
size 10 mm, length 36 cm, opening angle in
extended position approx. 60°

Size 12 mm

Operating instrument, **length 36 cm**,
for use **with trocars size 13 mm**



26173 FP

Fan Retractor, dismantling, retractable, with
proximal rotating wheel for expansion, size 12 mm,
length 36 cm, opening angle in extended position
approx. 60°

Components/Spare Parts see chapter 18

Dissecting and Grasping Forceps

CLICK'line® – rotating, dismantling, insulated,
with connector pin for unipolar coagulation

Size 5 mm

Operating instruments, **lengths 30 and 36 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**, for
use with telescopes with inbuilt working
channel and trocars size 6 mm



Length	Handle					
	33151	33152	33153	33156	33121	33149
30 cm						
36 cm						
43 cm						


Double action jaws

Working Insert	Complete Instrument					
33210 MD	33251 MD	33252 MD	33253 MD	33256 MD	33221 MD	33249 MD
33310 MD	33351 MD	33352 MD	33353 MD	33356 MD	33321 MD	33349 MD
33410 MD	33451 MD	33452 MD	33453 MD	33456 MD	33421 MD	33449 MD



CLICK'line® KELLY Dissecting and Grasping Forceps

33210 ML	33251 ML	33252 ML	33253 ML	33256 ML	33221 ML	33249 ML
33310 ML	33351 ML	33352 ML	33353 ML	33356 ML	33321 ML	33349 ML
33410 ML	33451 ML	33452 ML	33453 ML	33456 ML	33421 ML	33449 ML



CLICK'line® KELLY Dissecting and Grasping Forceps, long

33310 DK	33351 DK	33352 DK	33353 DK	33356 DK	33321 DK	33349 DK
----------	----------	----------	----------	----------	----------	----------




CLICK'line® DeBAKEY Grasping Forceps,
curved and slender jaws, atraumatic

33310 UY	33351 UY	33352 UY	33353 UY	33356 UY	33321 UY	33349 UY
----------	----------	----------	----------	----------	----------	----------



CLICK'line® UYAMA Dissecting and Grasping Forceps,
long jaws, with reservoir, fine, for dissection

33310 NO	33351 NO	33352 NO	33353 NO	33356 NO	33321 NO	33349 NO
----------	----------	----------	----------	----------	----------	----------



CLICK'line® NOSHIRO Dissecting and Grasping Forceps,
short jaws, fine, for dissection

Scissors

CLICK'linē – rotating, dismantling,
with and without connector pin for unipolar coagulation

Size 5 mm

Operating instruments, **length 30, 36 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm

Length	Handle					
	33151	33121	33149	33161	33131	33147
30 cm						
36 cm						
43 cm						

Double action jaws

Working Insert	Complete Instrument					
34210 MS	34251 MS	34221 MS	34249 MS	34261 MS	34231 MS	34247 MS
34310 MS	34351 MS	34321 MS	34349 MS	34361 MS	34331 MS	34347 MS
34410 MS	34451 MS	34421 MS	34449 MS	34461 MS	34431 MS	34447 MS



CLICK'linē METZENBAUM Scissors, curved

34210 MA	34251 MA	34221 MA	34249 MA	34261 MA	34231 MA	34247 MA
34310 MA	34351 MA	34321 MA	34349 MA	34361 MA	34331 MA	34347 MA
34410 MA	34451 MA	34421 MA	34449 MA	34461 MA	34431 MA	34447 MA



CLICK'linē Scissors, serrated blades, spoon-shaped, curved

34210 MW	34251 MW	34221 MW	34249 MW	34261 MW	34231 MW	34247 MW
34310 MW	34351 MW	34321 MW	34349 MW	34361 MW	34331 MW	34347 MW
34410 MW	34451 MW	34421 MW	34449 MW	34461 MW	34431 MW	34447 MW



CLICK'linē Scissors, serrated, curved, conical

34210 MD	34251 MD	34221 MD	34249 MD	34261 MD	34231 MD	34247 MD
34310 MD	34351 MD	34321 MD	34349 MD	34361 MD	34331 MD	34347 MD



CLICK'linē Scissors, straight



Item no: 26596CL

Knot Tier

Quantity:

–

1

+

On the offer list

Technical details	Documents & Films	Product group / field of application
Author	CICE	
Model	CLERMONT-FERRAND	
Size	5 mm	
Working length	36 cm	
Application	For extracorporeal knotting	

Dissecting and Grasping Forceps

CLICK^{line} – rotating, dismantling, insulated,
with connector pin for unipolar coagulation

Size 5 mm

Operating instruments, lengths 30 and 36 cm,
for use with trocars size 6 mm

Operating instruments, length 43 cm,
for use with telescopes with inbuilt
working channel and trocars size 6 mm



Length	Handle						
	33151	33152	33153	33156	33121	33125	33149
30 cm							
36 cm							
43 cm							

Double action jaws

Working Insert	Complete Instrument						
33210 R	33251 R	33252 R	33253 R	33256 R	33221 R	33225 R	33249 R
33310 R	33351 R	33352 R	33353 R	33356 R	33321 R	33325 R	33349 R
33410 R	33451 R	33452 R	33453 R	33456 R	33421 R	33425 R	33449 R



CLICK^{line} Dissecting and Grasping Forceps, right-angled

33210 UL	33251 UL	33252 UL	33253 UL	33256 UL	33221 UL	33225 UL	33249 UL
33310 UL	33351 UL	33352 UL	33353 UL	33356 UL	33321 UL	33325 UL	33349 UL
33410 UL	33451 UL	33452 UL	33453 UL	33456 UL	33421 UL	33425 UL	33449 UL



CLICK^{line} REDDICK-OLSEN Dissecting and Grasping Forceps, robust

33210 HM	33251 HM	33252 HM	33253 HM	33256 HM	33221 HM	33225 HM	33249 HM
33310 HM	33351 HM	33352 HM	33353 HM	33356 HM	33321 HM	33325 HM	33349 HM
33410 HM	33451 HM	33452 HM	33453 HM	33456 HM	33421 HM	33425 HM	33449 HM



CLICK^{line} Dissecting and Grasping Forceps, “alligator jaws”

Accessories

Reduction Sleeves/Extractors and Reducers

Reduction Sleeves/Extractors



30140 DB

Instrument Diameter	Outer Diameter of Trocar Cannula – Color Code					
	6 mm black	11 mm green	12 mm white	13 mm black	15 mm	22 mm
3 mm	30140 KA	30140 DA	–	–	–	–
5 mm	–	30140 DB	30140 EB	30140 HB	30140 FB	30140 GB
8 mm	–	30140 DD	–	–	–	–
10 mm	–	–	30140 EE	30140 HE	30140 FE	30140 GE

Reducers, for use with trocars with automatic valves and TERNAMIAN EndoTIP cannulas with automatic valves



30141 AA



30141 DB

30141 AA Reducer, 6/3 mm

30141 DB Same, 11/5 mm

30141 HB Same, 13/5 mm and 13.5/5 mm

30141 HE Same, 13/10 mm

Double Reducer



30142 HB

30142 HB Double Reducer, 13/10 mm, 13.5/10 mm, 13/5 mm and 13.5/5 mm

Reducers, for use with trocars with silicone leaflet valves and TERNAMIAN EndoTIP cannulas with silicone leaflet valves



30160 RE

30160 RE Reducer, 6/3 mm

30103 RE Same, 11/5 mm

Coagulation and Dissection Electrodes

without suction channel, insulated sheath,
with connector pin for unipolar coagulation

Size 5 mm

Operating instruments, **lengths 30 and 36 cm**,
for use with **trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm



Length	Instrument
30 cm	
36 cm	
43 cm	

Distal Tip	Instrument	
	26772 UF	Coagulation and Dissection Electrode, L-shaped
	26775 UF	
	26778 UF	
	26775 UE	Coagulation and Dissection Electrode, spatula-shaped, blunt
	26778 UE	
	26775 S	Coagulation and Dissection Electrode, knife-shaped
	26775 R	MANGESHIKAR Coagulation Electrode, with ball end

Special Features:

- The hook electrode is suitable for resection, insulating structures as well as dissection and coagulation.
 - The distal tip is semicircular: The thick outer curvature enables safe dissection.
- The thin inner curvature enables fast work with minimal coagulation.
 - The striated handle allows easier gripping.
 - The hook length enables ergonomic work.

Length	Instrument
36 cm	
43 cm	

Distal Tip	Instrument	
	26775 C	CADIERE Coagulation and Dissection Electrode, L-shaped, with cm-marking, distal tip tapered
	26775 CL	

Units for use with Coagulation and Dissection Electrodes see chapter 20, UNITS

ROBI® Bipolar Scissors

rotating, dismantling, with connector pin for bipolar coagulation,
Modell CLERMONT-FERRAND

Size 5 mm



Operating instruments, **length 36 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm

Outer Sheath	Handle
	38151
Length 36 cm	
Length 43 cm	

Single action jaws

Working Insert	Complete Instrument
38610 MT	38651 MT
 23	ROBI® Scissors , CLERMONT-FERRAND model, straight jaws, for cutting and bipolar coagulation of vessels and tissue layers
38610 MZ	38651 MZ
 23	ROBI® Scissors , CLERMONT-FERRAND model, straight jaws, serrated, for cutting and bipolar coagulation

Double action jaws

38610 MW	38651 MW
38710 MW	38751 MW
 20	ROBI® METZENBAUM Scissors , CLERMONT-FERRAND model, curved jaws, slender scissor blades, for cutting and bipolar coagulation

Please note:
For the instrument only the **individual component parts** are numbered. The catalog number for the **complete instrument** is not on the instrument. Please take this number from the numbers indicated in the **red** background of the table above. The color **green** indicates the working inserts.
Bipolar High Frequency Cords see page 176

ROBI® Bipolar Grasping Forceps ^{NEW}

rotating, dismantling, with connector pin for bipolar coagulation,
CLERMONT-FERRAND model

STORZ
KARL STORZ — ENDOSKOPE

Size 5 mm



Operating instruments, **length 36 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm

Outer Sheath	Handle
	38151
Length 36 cm	
Length 43 cm	

Double action jaws

Working insert	Complete Instrument
38610 ON	38651 ON
38710 ON	38751 ON
 18	ROBI® Grasping Forceps , CLERMONT-FERRAND model, fenestrated, with especially fine atraumatic serration
38610 WA	38651 WA
 19	ROBI® WATTIEZ Grasping Forceps , with multiple teeth, for atraumatic and accurate grasping
38610 MD	38651 MD
38710 MD	38751 MD
 19	ROBI® KELLY Dissecting and Grasping Forceps , CLERMONT-FERRAND model, especially suitable for dissection
38610 ML	38651 ML
 24	ROBI® KELLY Dissecting and Grasping Forceps , CLERMONT-FERRAND model, especially suitable for dissection
38610 MA	38651 MA
 24	ROBI® MALZONI Dissecting and Grasping Forceps , especially suitable for dissection

Please note:

For the instrument only the **individual component parts** are numbered. The catalog number for the **complete instrument** is not on the instrument. Please take this number from the numbers indicated in the **red** background of the table above. The color **green** indicates the working inserts.

Bipolar High Frequency Cords see accessories page 158

2-15