



Stat Fax® 4700

Analizator imunologic/Cititor de microstripuri

The Stat Fax® 4700 este un cititor de microstripuri compact și autonom. Designul său simplificat oferă o interfață cu ecran tactil, o optică excelentă, software de încărcare curbă și imprimantă încorporată pentru a satisface cerințele din laboratoare moderne.



AWARENESS TECHNOLOGY, Inc. PO Box 1679, Palm City, FL 34991 USA (772) 283 6540 fax (772) 283 8020
web: <http://www.awaretech.com> email: info@awaretech.com

Cost Effective by Design®

Stat Fax[®] 4700

Microstrip Reader Specificatii

Fotometric

Interval de măsurare liniară:	0,0 până la 3,0 Unități de absorbție (A).
Precizia fotometrică:	+/- (1% din lectură + 0,005A).
Stabilitate:	Drift de cel mult 0,005A în 8 ore / bicromatic.
Sursă de lumină:	Lampa de wolfram cu funcție de salvare a lămpilor.
Lungimi de undă standard:	Standard: 405, 450, 492, și 630nm. 6-VIS: 405, 450, 492, 545, 600, și 630nm. 6-UV: 340, 405, 450, 492, 545, și 630nm. (filtre alternative disponibile de la 340 la 700 nm).
Tip filtru:	Interferență IAD hardcoat, 10nm jum bandpass.
Recipient:	Benzi simple, duble sau despărțitoare până la 12 godeuri. Capacitate de încărcare cu 3 benzi cu 2 opțiuni: <input checked="" type="checkbox"/> 3x12 careta pentru benzi cu o singură distrugere sau 12 benzi. <input checked="" type="checkbox"/> 3x8 purtător pentru benzi care nu se desprind (rigide) sau 2x8.

Electronic

Afișa:	Ecran tactil interactiv de 3.5 "LCD, afișaj grafic color.
Printer:	Matricea termică a punctelor, cu capacitate grafică.
Cerinte de putere:	115V sau 230V AC, 1.5A, 50-60Hz (intrare universală).
Interfata:	Mouse USB.

Software-ul

Viteza:	Citește, calculează și tipărește rezultatele, câte 3 secunde pe tub.
Modul de calculare:	Calibrarea unui singur punct în funcție de standard sau factor, multipoint calibrare cu potrivire curbă punct la punct, regresii liniare cu jurnal, lin, și selecții log-logit, regresii neliniare cu 4 parametri și selecții spline cubice.
Test Meniu:	Peste 50 de canale deschise pentru stocarea testelor. Stochează toți parametrii, inclusiv lungimile de undă, calcule, coduri de unitate, intervale, interpretări, valori calibrator, nume de test și curba anterioară.

Altele

Carcasă:	Acoperă și bază din plastic ABS ignifugă vopsită.
Dimensiuni:	Aproximativ. 9x13.5x5in. (24x34x13cm) cântărește 10 kg (4,5 kg).
Certificarea:	NRTL listat, marca CE.

Dri-Dye[®] benzi de verificare

Pentru verificarea performanței Stat Fax[®], ChroMate[®],

și alte cititoare de microstrip și microplata



Awareness Technology products are designed and manufactured in the U.S.A.

Awareness Technology reserves the right to change specifications without notice due to design improvements.

CERTIFICATE OF REGISTRATION

This is to certify that the management system of:

Awareness Technology, Inc.

Main Site: 1935 SW Martin Highway

Palm City, Florida 34990 USA

Additional site: 2325 SW Martin Highway, Palm City, Florida 34990 USA

has been assessed by Intertek as conforming to the requirements of:

ISO 13485:2016

The quality management system is applicable to:

The design, development, manufacture, distribution, installation and service of IVDD General Laboratory Instruments.

Main site: Management, QA, Design, Sales, Service

Additional site: Manufacturing, Quality Control, Distribution, Shipping, Installation and Service.

Certificate Number:

9362

Revision Level: 09

Initial Certification Date:

March 28, 2012

Date of Certification Decision:

December 19, 2023

Issuing Date:

December 19, 2023

Valid Until:

December 18, 2026



SCC Accredited
CB-MS



OCSM
Accredited CCN

intertek

The SCC Accredited Symbol is an official symbol of the Standards Council of Canada, used under license.

Calin Moldovean
President

Intertek Testing Services NA Inc. dba Intertek,
900 Chelmsford Street,
Lowell, MA, USA





Awareness Technology, Inc.

Declaration of Conformity

Product identification

Product name : Stat Fax® Microstrip Reader
 Model/Type : 4700 Series
 EDMS Code : 21 01 10 01
 Class : Other IVD Devices – Self-declared

Manufacturer

Name : Awareness Technology, Inc.
 Address : P.O. Box 1679
 Palm City, Florida 34991
 Country : USA

Authorized Representative in Europe

Representative:
 Name : Emergo Europe
 Address : Molenstraat 15
 2513 BH The Hague
 Country : The Netherlands
 tel: +31 70 345 8570
 fax: +31 70 346 7299

Means of Conformity

Awareness Technology, Inc. declares that the product listed is in conformity with the Annex III, essential requirements and provisions of Council Directive:

98/79/EC

and is in conformance with the following standards:

IEC 61326-1 / IEC 61326-2-6
 EN 61010-1 / EN 61010-2-101

Signature

Place and Date: Awareness Technology, Inc.

December 12, 2012

Signature :

Name :

Steve Andrus
 Quality Manager

Doug Shields
 Compliance Engineer



MEDICAL
DEVICES
ISO 13485:2016
NSAI Certified

DECLARATION OF CONFORMITY

Product Family TOTAL AND FREE PROSTATE SPECIFIC ANTIGEN (PSA and FPSA)

Specific Product Details						
Product Description	Item # ELISA	Item # CLIA	EDMS Code	GMDN ELISA Code	GMDN CLIA Code	Risk Class
Total PSA Test System	2125-300A 2125-300B	2175-300A 2175-300B	12.03.01.32.00	54664	54665	High/ List B
Total PSA Extra Sensitive Test System	8725-300A 8725-300B	8775-300A 8775-300B	12.03.01.32.00	54664	54665	High/ List B
Free PSA Test System	2325-300A 2325-300B	2375-300A 2375-300B	12.03.01.33.00	54668	54669	High/ List B
Cancer VAST Test System	8425-300B 8425-300D 8425-300E	8475-300B 8475-300D 8475-300E	12.03.01.32.00	54664	54665	High/ List B
Multi Ligand Control	ML-300B	ML-300B	12.03.01.32.00	38207	38207	High/ List B

Manufacturer

Name Monobind Inc.
Address 100 North Pointe, Lake Forest, CA 92630
Country United States

Representative

Name CEpartner4U BV,
Address Esdoornlaan 13, 3951DB Maarn
Country The Netherlands
Telephone +31 (0)6 – 516.536.26

Notified Body

Name NSAI
Body ID Number 0050
CE Cert # 304.1006
Registration # NL-CA002-2011-23306

Means of Conformity

Monobind Inc. declares that the product listed is in conformity with the Annex IV, IVD Type List B essential requirements and provisions of Council Directive: 98/79/EC

And is in conformance with the following standards:

EN 13612:2002	EN 15223-1:2016	EN ISO 14971:2019
EN ISO 18113:2011	EN 13641:2002	EN ISO 23640:2015
Under the principles of	EN ISO 13485:2016	

Signature

Place and effective date

Monobind Inc.

January 30, 2021 revision 04

Signature

Ashatola

Name

Tony Shatola

Title

QA Director




Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HCV Ab CODES: CVAB.CE (192 tests) CVAB.CE.96 (96 tests) CVAB.CE.480 (480 tests) CVAB.CE.960 (960 tests) CVAB.CE.DB (192 tests)
CLASSIFICATION	ANNEX II – LIST A
CONFORMITY ASSESSMENT ROUTE	ANNEX IV

**WE HEREBY DECLARE THAT THE ABOVE MENTIONED
PRODUCT MEETS THE PROVISIONS OF THE COUNCIL
DIRECTIVE 98/79/EC
FOR IN VITRO DIAGNOSTIC DEVICES.**

NOTIFIED BODY	AEMPS – n° 0318
(EC) CERTIFICATE(S)	<ul style="list-style-type: none">• FULL QUALITY ASSURANCE SYSTEM N° 2003 12 0388 CT (in accordance with Annex IV – except Section IV) of the Directive 98/79/EC), RELEASED BY EC NOTIFIED BODY N° 0318• DESIGN CERTIFICATE N° 2003 12 0392 ED RELEASED BY EC NOTIFIED BODY N° 0318• UNE EN ISO 13485 N° 2013 11 0039 EN, RELEASED BY EC NOTIFIED BODY N° 0318

PLACE & DATE OF FIRST ISSUE	MILANO – JANUARY 2004
PLACE & DATE OF CURRENT EMISSION	SESTO SAN GIOVANNI (MI) – MARCH 2018
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	 DIA.PRO DIAGNOSTIC BIOPROBES s.r.l.

Rev: 0318



Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HBs Ab CODE: SAB.CE (96 tests)
CLASSIFICATION	ANNEX II – LIST A
CONFORMITY ASSESSMENT ROUTE	ANNEX IV

**WE HEREBY DECLARE THAT THE ABOVE MENTIONED
PRODUCT MEETS THE PROVISIONS OF THE COUNCIL
DIRECTIVE 98/79/EC
FOR IN VITRO DIAGNOSTIC DEVICES.**

NOTIFIED BODY	AEMPS – n° 0318
(EC) CERTIFICATE(S)	<ul style="list-style-type: none">• FULL QUALITY ASSURANCE SYSTEM N° 2003 12 0388 CT (in accordance with Annex IV – except Section IV) of the Directive 98/79/EC), RELEASED BY EC NOTIFIED BODY N° 0318• DESIGN CERTIFICATE N° 2003 12 0390 ED RELEASED BY EC NOTIFIED BODY N° 0318• UNE EN ISO 13485 N° 2013 11 0039 EN, RELEASED BY EC NOTIFIED BODY N° 0318

PLACE & DATE OF FIRST ISSUE	MILANO – JANUARY 2004
PLACE & DATE OF CURRENT EMISSION	SESTO SAN GIOVANNI (MI) – MAY 2018
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	 DIA.PRO DIAGNOSTIC BIOPROBES S.R.L.

Rev: 05/2018



Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HBs Ag one Version ULTRA CODES: SAG1ULTRA.CE (192 tests) SAG1ULTRA.CE.96 (96 tests) SAG1ULTRA.CE.480 (480 tests) SAG1ULTRA.CE.960 (960 tests) SAG1ULTRA.CE.DB (192 tests)
CLASSIFICATION	ANNEX II – LIST A
CONFORMITY ASSESSMENT ROUTE	ANNEX IV

**WE HEREBY DECLARE THAT THE ABOVE MENTIONED
PRODUCT MEETS THE PROVISIONS OF THE COUNCIL
DIRECTIVE 98/79/EC
FOR IN VITRO DIAGNOSTIC DEVICES.**

NOTIFIED BODY	AEMPS – n° 0318
(EC) CERTIFICATE(S)	<ul style="list-style-type: none">• FULL QUALITY ASSURANCE SYSTEM N° 2003 12 0388 CT (in accordance with Annex IV – except Section IV) of the Directive 98/79/EC), RELEASED BY EC NOTIFIED BODY N° 0318• DESIGN CERTIFICATE N° 2008 12 0588 ED RELEASED BY EC NOTIFIED BODY N° 0318• UNE EN ISO 13485 N° 2013 11 0039 EN, RELEASED BY EC NOTIFIED BODY N° 0318

PLACE & DATE OF FIRST ISSUE	MILANO – DECEMBER 2008
PLACE & DATE OF CURRENT EMISSION	SESTO SAN GIOVANNI (MI) – MARCH 2018
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 03/2018

DECLARATION OF CONFORMITY

1) Manufacturer (Name, department): **Monobind Inc.**

Address: **100 North Pointe, LAKE FOREST, CA 92630. UNITED STATES**

and

2) European authorized representative: **CEpartner4U BV,**

Address: **ESDOORNLAAN 13, 3951DB MAARN, THE NETHERLANDS;**

(on product labels printed as:

CEpartner4U , ESDOORNLAAN 13, 3951DB MAARN, THE NETHERLANDS. www.cepartner4u.com)

3) Product(s) (name, type or model/batch number, etc.):

Immunoassay products;

AccuBind® ELISA,

AccuLite® CLIA,

QSure® Control,

Instruments

see appendix

4) The product(s) described above is in conformity with:

<u>Document No.</u>	<u>Title</u>
98/79/EC	<i>In vitro</i> Diagnostic Medical Devices Directive

5) Additional information (Conformity procedure, Notified Body, CE certificate, Registration nr., etc.):

Conformity assessment procedure for CE marking: *In vitro* Diagnostic Medical Device Directive, Annex III

Registration nr. : **NL- CA002-22758 and NL- CA002-22762**

Lake Forest, USA; 2021-09-20

(Place & date of issue (yyyy-mm-dd))



Tony Shatola; QA Director, Monobind Inc.
(name, function and signature of manufacturer)

Appendix

Date: 2021-09-20

List of devices.

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
Allergy & Anemia							
Ferritin Test System	2825-300A 2825-300B	2875-300A 2875-300B			12.07.01.02.00	Low	2005-11-11
Folate Test System	7525-300A 7525-300B	7575-300A 7575-300B			12.07.01.03.00	Low	2010-06-29
Immunoglobulin E (IgE) Test System	2525-300A 2525-300B	2575-300A 2575-300B			12.02.01.02.00	Low	2005-11-11
Transferrin Soluble Receptor (sTfR) Test System	8625-300A 8625-300B	8675-300A 8675-300B			12.07.01.06.00	Low	2010-06-29
Vitamin B-12 (Vit B12) Test System	7625-300A 7625-300B	7675-300A 7675-300B			12.07.02.04.00	Low	2011-09-26
Folate, Vitamin B-12 (Anemia Panel VAST) Test System	7825-300A 7825-300B	7875-300A 7875-300B			12.07.01.00.00	Low	2013-09-16
Autoimmune							
Anti-Cyclic Citrullinated Peptide IgG (Anti-CCP IgG) Test System	12725-300A 12725-300B	12775-300A 12775-300B			12.11.01.90.00	Low	2019-04-03
Anti-Thyroglobulin (Anti-Tg) Test System	1025-300A 1025-300B	1075-300A 1075-300B			12.10.03.04.00	Low	2005-11-11
Anti-Thyropoxidase (Anti-TPO) Test System	1125-300A 1125-300B	1175-300A 1175-300B			12.10.03.01.00	Low	2005-11-11
Bone Metabolism & Growth							
Calcitonin Test System	9325-300A 9325-300B	9375-300A 9375-300B			12.06.03.02.00	Low	2019-04-03
Growth Hormone (hGH) Test System	1725-300A 1725-300B	1775-300A 1775-300B			12.06.04.02.00	Low	2005-11-11
Parathyroid Hormone (PTH) Test System	9025-300A 9025-300B	9075-300A 9075-300B			12.06.03.13.00	Low	2011-09-26
Parathyroid Hormone (PTH) 3rd & 2nd Gen (VAST) Test System	10025-300A 10025-300B	10075-300A 10075-300B			12.06.03.13.00	Low	2019-04-03
25(OH) Vitamin D Total Direct (Vit D-Direct) Test System	7725-300A 7725-300B	7775-300A 7775-300B			12.06.03.10.00	Low	2017-07-05
Cancer Markers							
Alpha-Fetoprotein (AFP) Test System	1925-300A 1925-300B	1975-300A 1975-300B			12.03.90.01.00	Low	2005-11-11
CA-125 Test System	3025-300A 3025-300B	3075-300A 3075-300B			12.03.01.06.00	Low	2005-11-11
CA 15-3 Test System	5625-300A 5625-300B	5675-300A 5675-300B			12.03.01.02.00	Low	2010-06-29
CA 19-9 Test System	3925-300A 3925-300B	3975-300A 3975-300B			12.03.01.03.00	Low	2005-11-11
Carcinoembryonic Antigen (CEA) Test System	1825-300A 1825-300B	1875-300A 1875-300B			12.03.01.31.00	Low	2005-11-11
Next Generation Carcinoembryonic Antigen	4625-300A	4675-300A			12.03.01.31.00	Low	2010-06-29

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
(CEA-Next Gen) Test System	4625-300B	4675-300B					
Free β-Subunit Human Chorionic Gonadotropin (Free Beta hCG) Test System	2025-300A 2025-300B	2075-300A 2075-300B			12.03.01.90.00	Low	2005-11-11
Cardiac Markers							
CK-MB Test System	2925-300A 2925-300B	2975-300A 2975-300B			12.13.01.02.00	Low	2005-11-11
Digoxin (DIG) Test System	925-300A 925-300B	975-300A 975-300B			12.08.01.01.00	Low	2005-11-11
High Sensitivity CRP (hs-CRP) Test System	3125-300A 3125-300B	3175-300A 3175-300B			12.13.01.90.00	Low	2005-11-11
Myoglobin Test System	3225-300A 3225-300B	3275-300A 3275-300B			12.13.01.05.00	Low	2005-11-11
Troponin I (cTnI) Test System	3825-300A 3825-300B	3875-300A 3875-300B			12.13.01.07.00	Low	2005-11-11
Diabetes							
C-Peptide Test System	2725-300A 2725-300B	2775-300A 2775-300B			12.06.01.01.00	Low	2005-11-11
Insulin Test System	2425-300A 2425-300B	2475-300A 2475-300B			12.06.01.03.00	Low	2005-11-11
Rapid Insulin Test System	5825-300A 5825-300B				12.06.01.03.00	Low	2010-06-29
Insulin - C-Peptide (Diabetes Panel VAST)	7325-300A 7325-300B	7375-300A 7375-300B			12.06.01.03.00	Low	2005-11-11
Endocrine							
ACTH Test System	10625-300	10675-300			12.06.04.01.00	Low	2019-04-03
Aldosterone Test System	10125-300	10175-300			12.06.02.01.00	Low	2019-04-03
Leptin Test System	10925-300	10975-300			12.06.90.17.00	Low	2019-04-03
Fertility & Prenatal							
Anti-Müllerian Hormone (AMH) Test System	9725-300A 9725-300B	9775-300A 9775-300B			12.05.02.16.00	Low	2019-04-03
Folicle Stimulating Hormone (FSH) Test System	425-300A 425-300B	475-300A 475-300B			12.05.01.04.00	Low	2005-11-11
B-Human Chorionic Gonadotropin (hCG) Test System	825-300A 825-300B	875-300A 875-300B			12.05.02.05.00	Low	2005-11-11
B-Human Chorionic Gonadotropin Extended Range (hCG-XR) Test System	8825-300A 8825-300B	8875-300A 8875-300B			12.05.02.05.00	Low	2013-09-16
Rapid B-Human Chorionic Gonadotropin (Rapid hCG) Test System	3325-300A 3325-300B				12.05.02.05.00	Low	2005-11-11
Inhibin A Test System	9525-300A 9525-300B	9575-300A 9575-300B			12.05.01.90.00	Low	2019-04-03
Inhibin B Test System	9625-300A 9625-300B	9675-300A 9675-300B			12.05.01.90.00	Low	2019-04-03
Luteinizing Hormone (LH) Test System	625-300A 625-300B	675-300A 675-300B			12.05.01.05.00	Low	2005-11-11
Pregnancy Associated Plasma Protein – A Mass Units (PAPP-A Mass Units) Test System	12625-300A 12625-300B	12675-300A 12675-300B			12.05.02.10.00	Low	2017-07-05
Prolactin Hormone (PRL) Test System	725-300A 725-300B	775-300A 775-300B			12.05.01.08.00	Low	2005-11-11

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
Prolactin Hormone Sequential (PRLs) Test System	4425-300A 4425-300B	4475-300A 4475-300B			12.05.01.08.00	Low	2005-11-11
Human Chorionic Gonadotropin (hCG) , Human Prolactin (hPRL), Human Luteinizing Hormone (hLH), Follicle Stimulating Hormone (FSH) (Fertility Panel VAST) Test System	8325-300B 8325-300D 8325-300E	8375-300B 8375-300D 8375-300E			12.05.01.90.00	Low	2006-08-24
Alpha-Fetoprotein (AFP), Human Chorionic Gonadotropin (hCG), Unconjugated Estiol (u-E3) Triple Screen (Triple Screen Panel VAST) Test System	8525-300A 8525-300B	8575-300A 8575-300B			12.05.01.90.00	Low	2010-06-29
Infectious Diseases							
Anti-H. Pylori IgG (H. Pylori Ab IgG) Test System	1425-300A 1425-300B	1475-300A 1475-300B			15.01.04.03.00	Low	2005-11-11
Anti-H. Pylori IgM (H. Pylori Ab IgM) Test System	1525-300A 1525-300B	1575-300A 1575-300B			15.01.04.03.00	Low	2005-11-11
Anti-H. Pylori IgA (H. Pylori Ab IgA) Test System	1625-300A 1625-300B	1675-300A 1675-300B			15.01.04.03.00	Low	2005-11-11
Anti-SARS-CoV-2 (COVID-19) IgG Test System	11925-300A 11925-300B	11975-300A 11975-300B			15.04.80.90.00	Low	2020-08-25
Anti-SARS-CoV-2 (COVID-19) IgM Test System	11725-300A 11725-300B	11775-300A 11775-300B			15.04.80.90.00	Low	2020-08-25
Anti-SARS-CoV-2 (COVID-19) IgA Test System	11825-300A 11825-300B	11875-300A 11875-300B			15.04.80.90.00	Low	2020-08-25
Anti-SARS-CoV-2 (COVID-19) S1-RBD IgG Test System	12025-300A 12025-300B	12075-300A 12075-300B			15.04.80.90.00	Low	2021-09-20
D-Dimer Test System	9225-300A 9225-300B	9275-300A 9275-300B			13.02.05.03.00	Low	2020-08-25
Procalcitonin (PCT) Test System	1425-300A 1425-300B	1475-300A 1475-300B			12.06.90.16.00	Low	2017-07-05
Neonatal							
Neonatal 17OHP (N-17OHP) Test System	5525-300A 5525-300B				12.05.01.07.00	Low	2008-02-01
Neonatal (N-T4) Thyroxine Test System	2625-300A 2625-300B				12.04.01.12.00	Low	2005-11-11
Neonatal TBG (N-TBG) Test System	8925-300A 8925-300B				12.04.01.09.00	Low	2013-09-16
Neonatal TSH (N-TSH) Test System	3425-300A 3425-300B 3425-300D 3425-300E				12.04.01.90.00	Low	2005-11-11
Steroid							
Androstenedione (ANST) Test System	12425-300A 12425-300B	12475-300A 12475-300B			12.05.01.01.00	Low	2021-09-20
Cortisol Test System	3625-300A 3625-300B	3675-300A 3675-300B			12.06.02.04.00	Low	2005-11-11
Dehydroepiandrosterone (DHEA) Test System	7425-300A 7425-300B	7475-300A 7475-300B			12.05.01.02.00	Low	2011-09-26
Dehydroepiandrosterone Sulfate (DHEA-S) Test System	5125-300A 5125-300B	5175-300A 5175-300B			12.05.01.02.00	Low	2010-06-29
Estrone (E1) Test System	10325-300A 10325-300B	10375-300A 10375-300B			12.05.02.04.00	Low	2019-04-03

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
Estradiol (E2) Test System	4925-300A 4925-300B	4975-300A 4975-300B			12.05.01.03.00	Low	2010-06-29
Unconjugated Estiol (u-E3) Test System	5025-300A 5025-300B	5075-300A 5075-300B			12.05.02.02.00	Low	2010-06-29
Progesterone Test System	4825-300A 4825-300B	4875-300A 4875-300B			12.05.01.06.00	Low	2010-06-29
17-OH Progesterone (17-OHP) Test System	5225-300A 5225-300B	5275-300A 5275-300B			12.05.01.07.00	Low	2010-06-29
17-OH Progesterone SI (17-OHP-SI) Test System	9925-300A 9925-300B	9975-300A 9975-300B			12.05.01.07.00	Low	2010-10-18
Sex Hormone Binding Globulin (SHBG) Test System	9125-300A 9125-300B	9175-300A 9175-300B			12.05.01.09.00	Low	2013-09-16
Testosterone Test System	3725-300A 3725-300B	3775-300A 3775-300B			12.05.01.10.00	Low	2007-11-01
Free Testosterone Test System	5325-300A 5325-300B	5375-300A 5375-300B			12.05.01.10.00	Low	2010-06-29
Thyroid							
Total Triiodothyronine (tT3) Test System	125-300A 125-300B 125-300D 125-300E	175-300A 175-300B 175-300D 175-300E			12.04.01.05.00	Low	2005-11-11
Free Triiodothyronine (fT3) Test System	1325-300A 1325-300B 1325-300A 1325-300B	1375-300A 1375-300B 1375-300D 1375-300E			12.04.01.01.00	Low	2005-11-11
Total Triiodothyronine (tT3 SBS) Test System	8125-300A 8125-300B	8175-300A 8175-300B			12.04.01.01.00	Low	2010-06-29
Rapid Total Triiodothyronine (Rapid -tT3) Test System	11225-300A 11225-300B				12.04.01.01.00	Low	2017-07-05
T3-Uptake (T3U) Test System	525-300A 525-300B	575-300A 575-300B			12.04.01.06.00	Low	2005-11-11
Thyroxine (tT4) Test System	225-300A 225-300B 225-300D 225-300E	275-300A 275-300B 275-300D 275-300E			12.04.01.07.00	Low	2005-11-11
Free Thyroxine (fT4) Test System	1225-300A 1225-300B 1225-300D 1225-300E	1275-300A 1275-300B 1275-300D 1275-300E			12.04.01.02.00	Low	2005-11-11
Total Thyroxine (tT4 SBS) Test System	8225-300A 8225-300B	8275-300A 8275-300B			12.04.01.01.00	Low	2010-06-29
Rapid Total Thyroxine (Rapid -tT4) Test System	11125-300A 11125-300B				12.04.01.01.00	Low	2017-07-05
Thyrotropin (TSH) Test System	325-300A 325-300B 325-300D 325-300E	375-300A 375-300B 375-300D 375-300E			12.04.01.11.00	Low	2005-11-11
Rapid TSH Test System	6025-300A 6025-300B	6075-300A 6075-300B			12.04.01.11.00	Low	2010-06-29
Thyroxine-Binding Globulin (TBG) Test System	3525-300A 3525-300B	3575-300A 3575-300B			12.04.01.09.00	Low	2005-11-11
Thyroglobulin (Tg) Test System	2225-300A	2275-300A			12.04.01.08.00	Low	2005-11-11

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
	2225-300B	2275-300B					
Total Thyroxine (tT4), Total Triiodothyronine (tT3) & Thyroid Stimulating Hormone (TSH) (Thyroid Panel VAST) Test System	8025-300B 8025-300D 8025-300E	8075-300B 8075-300D 8075-300E			12.04.01.01.00	Low	2005-11-11
Free Thyroxine (fT4), Free Triiodothyronine (fT3) & Thyroid Stimulating Hormone (TSH) (Free Thyroid Panel VAST) Test System	7025-300B 7025-300D 7025-300E	7075-300B 7075-300D 7075-300E			12.04.01.01.00	Low	2010-06-29

Miscellaneous Controls							
Anti-H. Pylori Control (IgA, IgG, IgM) – Positive & Negative			HPC-300		12.50.01.16.00	Low	2013-09-16
Anti-Tg & Anti-TPO Control – Positive & Negative			AIT-101		12.50.01.16.00	Low	2010-06-29
Maternal Control – (AFP, uE3, hCG, Free beta hCG) Tri Level			MC-300		12.50.01.16.00	Low	2010-06-29
TBG Control – Tri-Level			TBG-300		12.50.01.16.00	Low	2013-09-16
Tg Control – Tri-Level			TG-300		12.50.01.16.00	Low	2010-06-29
Tumor Marker Control – (CA 125, CA 15-3, CA 19-9) Tri-Level			TMC-300		12.50.01.16.00	Low	2013-09-16

Miscellaneous Instruments							
Autoplex® ELISA & CLIA Analyzer				IN006	21.02.10.01	Low	2010-06-29
Autoplex® G2 ELISA & CLIA Analyzer				IN006-2	21.02.10.01	Low	2013-09-16
Autoplex® G3 ELISA & CLIA Analyzer				IN006-3	21.02.10.01	Low	2017-07-05
NeoEldex® ELISA Analyzer				IN009	21.02.10.01	Low	2011-09-26
Impulse® 3 CLIA Analyzer				IN007	21.02.10.01	Low	2010-06-29
NeoLumax® CLIA Analyzer				IN010	21.02.10.01	Low	2011-09-26
LuMatic® CLIA Analyzer				IN008	21.02.10.01	Low	2011-09-26
PrisMatic® ELISA Analyzer				IN013	21.02.10.01	Low	2013-09-16
PlateWash - Immunoassay Washer				IN002	21.02.10.01	Low	2010-06-29
TITIN® ELISA & CLIA Analyzer				IN015-EC	21.02.10.01	Low	2017-07-05
TITIN® ELISA Analyzer				IN015-E	21.02.10.01	Low	2017-07-05
TITIN-s® ELISA & CLIA Analyzer				IN016-EC	21.02.10.01	Low	2017-07-05
TITIN-s® ELISA Analyzer				IN016-E	21.02.10.01	Low	2017-07-05