

 Copy	Copy Image	Copies images.
 Target E.I.	Target E.I.	Sets Target E.I for each patient size.
 None	Grid Suppression (None)	Not applying the grid correction.
 Auto	Grid Suppression (Auto)	Detects and corrects grid lines automatically.
 VLine	Grid Suppression (Vertical)	Corrects vertical grid lines.
 HLine	Grid Suppression (Horizontal)	Corrects horizontal grid lines.
 None	None	Saves the left/right patient information indicated on the image.
 Left	Left	Saves the left/right patient information indicated on the image.
 Right	Right	Saves the left/right patient information indicated on the image.
 Both	Both	Saves the left/right patient information indicated on the image.
 Union	Union	Saves the left/right patient information indicated on the image.
 Orienta...	Patient Orientation	Sets patient orientation of the image.
 Up	Up	Moves the selected original image to upward.
 Down	Down	Moves the selected original image to downward.
 PS Image	Image Pixel Spacing	Measures the actual size
 Full View	Open Full View	Displays the current image as a full view.
 Suspend	Suspend Study	Suspends and closes the current study.
 Reset	Reset	Resets the image to the step setting value.
 Change Step	Change Step	Changes the step of image that is currently selected.
 Guide	Guide (Grid)	Displays grid guides on the screen.