

Date: 05/Jan/2023

STATEMENT

We, Atlas Medical having a registered office at Ludwig-Erhard-Ring 3, 15827 Blankenfelde-Mahlow, Berlin, Germany assign SRL Sanmedico having a registered office at A. Corobceanu Street 7A, apt.9, Chisinau MD-2012, Moldova, as authorized representative in correspondence with the conditions of directive 98/79/EEC.

We declare that the company mentioned above is authorized to register, notify, renew or modify the registration of medical devices on the territory of the Republic of Moldova.

On Behalf of Manufacturer:

General Manager

Haya Amawi

Signature: _____

Date: _____

Atlas Medical GmbH
Ludwig - Erhard Ring 3
15827 Blankenfelde - Mahlow
Tel. (0049) 33708 - 355030

Atlas Medical: Ludwig-Erhard-Ring 3, 15827 Blankenfelde-Mahlow, Berlin, Germany,
Tel: +4933708355030

Regulatory Office: William James House, Cowley Rd, Cambridge, CB4 0WX, United Kingdom
Tel: +44 (0) 1223 858 910

Middle East Site: P.O Box 204, King Abdullah II Industrial Estate, Amman, 11512, Jordan
Tel: +962 6 4026468


CE Declaration of Conformity

We,
Atlas Medical GmbH
 Head office: Ludwig-Erhard-Ring 3
 15827 Blankenfelde-Mahlow Germany
 Tel: +49(0)33708355030
 Email: info@atlas-medical.com

Middle East Site: Sahab Industrial Zone Area, King Abdullah II Industrial City
 Amman 11512, Jordan
 Tel.: +962 6 4026468
 Fax: +962 6 4022588
 Email: info@atlas-medical.com

Declare our responsibility that the following product:

Product Code	Product Name	Class	GMDN code
8.00.18.0.0005	RPR Carbon Antigen Reagent, 5 ml/vial	General-IVD	32450
8.00.18.2.1000	RPR Carbon Antigen 1000ml/bottle	General-IVD	32450
8.00.18.0.0050	RPR Carbon Antigen Kit, 50 Tests	General-IVD	32450
8.00.18.1.0050	RPR Carbon Antigen Kit, 50 Tests, White Glass Slide.	General-IVD	32450
8.00.18.2.0500	RPR Carbon Antigen Kit, 500 Tests (2ml latex, 2x0.5 ml control) Without card.	General-IVD	32450
8.00.18.3.0500	RPR Carbon Antigen Kit, 500 Tests (10ml latex, 2x0.5 ml control) Without card, stirring sticks.	General-IVD	32450
8.00.18.0.0100	RPR Carbon Antigen Kit, 100 Tests (2ml latex, 2x0.5 ml control)	General-IVD	32450
8.00.18.2.0100	RPR Carbon Antigen Kit, 100 Tests (2ml latex, 2x0.5 ml control +White Glass slide stirring sticks)	General-IVD	32450
8.00.18.0.0025	RPR Carbon Antigen Kit, 25 Tests (0.5ml latex, 2x0.5 ml control)	General-IVD	32450
8.00.18.0.0150	RPR Carbon Antigen Kit, 150 Tests	General-IVD	32450
8.00.18.0.0200	RPR Carbon Antigen Kit, 200 Tests	General-IVD	32450
8.00.18.0.0250	RPR Carbon Antigen Kit, 250 Tests	General-IVD	32450

Atlas Medical	First issue date	Date of review	Management approval	MRXDO10F.10 08.02.2011
	September.2021	06.09.2021		

Amin Al-Habashat
 RA Manager

Declaration Ref No: DC21-0193

8.00.18.0.0500	RPR Carbon Antigen Kit,500 Tests	General-IVD	32450
8.00.18.0.1000	RPR Carbon Antigen Kit, 1000 Tests	General-IVD	32450
8.00.18.4.0500	RPR Carbon Antigen Kit,500 Tests (3x3.4ml reagent,2x1 controls)	General-IVD	32450
8.00.18.5.0500	RPR Carbon Antigen Kit, 500 Tests, (3x3.4ml reagent,2x1 controls)	General-IVD	32450
8.00.18.8.0500	RPR Carbon Antigen 500 Test (10ml reagent) without Control's.	General-IVD	32450
8.00.18.9.0050	RPR Carbon Antigen Kit, (5x10ml Reagent,2x2ml Control) , white glass Slide, Stirring Stick.	General-IVD	32450
8.33.04.0.0001	RPR Positive control	General-IVD	32450
8.33.04.1.0001	RPR Positive control ,Bulk	General-IVD	32450
8.33.04.0.0100	RPR Positive control(100ml/vial)	General-IVD	32450
8.33.04.0.0500	RPR Positive control(500ml/bottle)	General-IVD	32450
8.33.08.0.0001	RPR Negative control	General-IVD	32450

Is produced under Atlas quality system (ISO13485: 2016) supported by GMED certificate:

Certificate N^o.: 36655 rev 1

Expiry Date: October 8th.2023

and complies with the essential requirements of
In Vitro Diagnostic Medical Devices Directive 98/79/EC Annex I
And

EN ISO 18113-1, -2 :2011, EN ISO 15223:2016

**EN ISO 14971:2019, EN ISO 23640 :2015 , ISO 2859 :2017,
EN 13612:2002, EN 13641:2002 , EN 13975:2003, ISO 13485:2016**

And

Intended for In-Vitro Professional use only.

This Declaration includes the batches produced beyond this day according to the product Lot Log.


Manufacturer

Atlas Medical GmbH

Ludwig-Erhard-Ring 3

15827 Blankenefelde-Mahlow Germany.



Atlas Medical	First issue date	Date of review	Management approval	MRXDO10F.10 08.02.2011
	September.2021	06.09.2021		

Anwar Al-Hadad
RA Manager

GMED certifie que le système de management de la qualité développé par
GMED certifies that the quality management system developed by

ATLAS MEDICAL GmbH
Ludwig-Erhard-Ring 3
15827 Blankenfelde-Mahlow GERMANY

pour les activités
for the activities

Conception et développement, fabrication et vente de dispositifs médicaux de diagnostic in vitro .

Design and Development, Manufacturing and Sales of in vitro diagnostic medical devices.

réalisées sur le(s) site(s) de
performed on the location(s) of

Voir addendum

See addendum

est conforme aux exigences des normes internationales
complies with the requirements of the international standards

ISO 13485: 2016

Début de validité / Effective date October 9th, 2020 (included)

Valable jusqu'au / Expiry date : October 8th, 2023 (included)

Etabli le / Issued on : October 8th, 2020




**CERTIFICATION
DE SYSTEMES
DE MANAGEMENT**
Accréditation n°4-0608
Liste des sites accrédités
et portée disponible sur
www.cofrac.fr

GMED N° 36655-1

Ce certificat est délivré selon les règles de certification GMED / This certificate is issued according to the rules of GMED certification

Renouvelle le certificat 36655-0



On behalf of the President
Béatrice LYS
Technical Director

Ce certificat couvre les activités et les sites suivants :

This certificate covers the following activities and sites:

French version :

Conception et développement, fabrication et vente de dispositifs médicaux de diagnostic *in vitro* à usage professionnel et/ ou d'autodiagnostic, dans les domaines du groupage sanguin, de la microbiologie, de la biochimie, de la toxicologie, de l'oncologie, de la cardiologie, de l'histologie, de l'endocrinologie et des maladies infectieuses, dans les techniques d'Agglutination/ ELISA/ Tests rapides/ Colorimétrie/ Disques antibiotiques.

English version:

Design and Development, Manufacturing and Sales of in vitro diagnostic medical devices for professional use and/or for self-testing, in the field of Immunohematology, Microbiology, Biochemistry, Toxicology, Oncology, Cardiology, Histology, Endocrinology Biosensors and Infectious diseases, in techniques of Agglutination/ ELISA/ Rapid tests/ Colorimetry/Antibiotic disks.

**ATLAS MEDICAL GmbH
Ludwig-Erhard-Ring 3
15827 Blankenfelde-Mahlow
GERMANY**

French version:

Siège social, responsable de la mise sur le marché

English version:

Headquarter, legal manufacturer

**Sahab Industrial Zone Area
King Abdullah II Industrial City
Amman 11512
JORDAN**

French version:

Conception, fabrication et contrôle final

English version:

Design, manufacture and final control

**William James House
Cowley Road,
Cambridge, CB OWX
United Kingdom**

French version:

Contact réglementaire

English version:

Regulatory Administration

3 sites / 3 sites

DocuSigned by:

Beatrice Lys
EF33BDA9BAA04A3...


**On behalf of the President
Béatrice LYS
Technical Director**

RPR SYPHILIS CARD TEST

IVD For In-Vitro diagnostic and professional use only

2°C 8°C
Store at 2 to 8 °C

INTRODUCTION

Syphilis is a disease caused by infection with the spirochete *Treponema pallidum*. The infection is systemic and the disease is characterized by periods of latency. These features, together with the fact that *T. pallidum* cannot be isolated in culture, mean that serologic techniques play a major role in the diagnosis and follow-up of treatment for syphilis.

Syphilis is categorized by an early primary infection in which patients may have non-specific symptoms, and potentially, genital lesions. Patients tested by serology during the primary phase may be negative for antibodies, especially if testing is performed during the first 1 to 2 weeks after symptom onset. As the disease progresses into the secondary phase, antibodies to *T. pallidum* reach peak titers, and may persist indefinitely regardless of the disease state or prior therapy. Therefore, detection of antibodies to nontreponemal antigens, such as cardiolipin (a lipoidal antigen released by host cells damaged by *T. pallidum*) may help to differentiate between active and past syphilis infection. Nontreponemal antibodies are detected by the rapid plasma reagin (RPR) assay, which is typically positive during current infection and negative following treatment or during late/latent forms of syphilis.

PRINCIPLE

RPR utilises carbon particles coated with cardiolipin antigen to detect reagin antibodies present in serum or plasma of syphilitic persons.

Specimens that contain reagin cause aggregation of the carbon particles which appear as dark clumps against a white background. The aggregation can be read macroscopically. Non-reactive samples typically appear as a smooth non-aggregated pattern which may form buttons in the centre of the test area.

MATERIALS

MATERIALS PROVIDED

- **RPR carbon antigen reagent:** Contains less than 0.1% sodium azide.
- **Positive Control :** Contains less than 0.1% sodium azide.
- **Negative control:** Contains less than 0.1% sodium azide

- **RPR test cards (Optional).**
- **Plastic sticks.**
- **Package insert.**

NOTE: This package insert is also used for individually packed reagent.

MATERIALS NEEDED BUT NOT PROVIDED

- Rotator (100rpm).
- Timer.
- Pipettes.

SAMPLES

Fresh serum or plasma. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Always use a fresh pipette tip for every test.
- Handle all negative and positive in the manner as patient specimens .
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- The used test should be discarded according to local regulations.
- Components of different human origin have been tested and found to be negative for the presence of antibodies anti- HIV 1+2 and anti-HCV, as well as for HBsAg. However, the controls should be handled cautiously as potentially infectious.

STORAGE AND STABILITY

All components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C.

PROCEDURES

QUALITATIVE PROCEDURE

- **Mix well the RPR reagent before use.**
- 1. Bring the reagents and samples to room temperature.
- 2. Dispense **50 µL of each sample** into a separate circle on the card. Use a separate tip for each sample.
- 3. Dispense **1 drop of each of positive and negative controls** into two additional circles.
- 4. Gently shake the dispensing vial and slightly press to remove air bubbles from the needle and the drop obtained is correct.

5. Dispense **1 drop (17.5 µl) of RPR antigen** to each circle next to the sample to be tested.
6. Place the card on a mechanical rotator and rotate at 100 r.p.m. for 8 minutes.
7. Observe macroscopically for agglutination within a minute after removing the card from the rotator.

SEMI-QUANTITATIVE PROCEDURE


- **Mix well the RPR reagent before use.**
- 1. Make doubling dilutions from Undiluted to 1:16 normal saline.
- 2. Place 50 µl of each dilution in to a separate circle on the test card.
- 3. Spread each dilution evenly over the test circle.
- 4. Continue as from Qualitative procedure .
The titer of the sample is expressed as the final dilution which shows aggregation of the carbon particles.

PERFORMANCE CHARACTERISTICS


1. **Sensitivity:** 100%.
2. **Specificity:** 100%.

INTERPRETATION OF TEST RESULTS

1. **Strong Reactive:** Large clumps of carbon particles with a clear background.




2. **Reactive:** Large clumps of carbon particles somewhat more disperse than Strong Reactive pattern.



3. **Weak Reactive:** Small clumps of carbon particles with light grey background.



4. **Trace Reactive:** Slight clumping of carbon particles typically seen as a button of aggregates in the centre of the test circle or dispersed around the edge of the test circle.



5. **Non-Reactive:** Typically a smooth grey pattern or a button of non-aggregated carbon particles in the centre of the test circle.



















REFERENCES

- Falcone V.H., Stout G.W. and Moore M.B. Jr., PHR 79: 491-495, 1964.

 **ATLAS MEDICAL**
Ludwig-Erhard Ring 3
15827 Blankenfelde-Mahlow
Germany
Tel: +49 - 33708 – 3550 30
Email: Info@atlas-medical.com
Website: www.atlas-medical.com

PPI2074A01
Rev B (15.03.2021)

	Catalogue Number		Temperature limit
	<i>In Vitro</i> diagnostic medical device		Caution
	Contains sufficient for <n> tests and Relative size		Consult instructions for use (IFU)
	Batch code		Manufacturer
	Fragile, handle with care		Use-by date
	Manufacturer fax number		Do not use if package is damaged
	Manufacturer telephone number		Date of Manufacture
	Keep away from sunlight		Keep dry