

CV Series Accessory List

1

Cataract

2

Vitrectomy

3

Tubing

4

Diathermy

5

Laser

6

Other Optional Accessories

NIDEK

- ... Applicable
- ... Standard Configuration
- A... Standard Accessory Set A (CV-30000 A, AP)*
- P... Standard Accessory Set P (CV-30000 AP)*

* To complete the standard configuration of CV-30000,
paid-for accessory sets are needed.






(Type A --> Standard Accessory Set A
Type AP --> Standard Accessory Set A & P)

Cataract

Phaco Tip / Sleeve 20G									
Product name	Product code	Reusable		Remarks	30000	9000	9000R	24000	7000
US tip Straight 0°	18254 0290			6pcs/box	○	○	○	○	○
	18252 0290	○		AC (≤135°C) 10 times of autoclaving	○	○	○	○	○
US tip Straight 15°	18254 0210			6pcs/box	○	○	○	○	○
	18252 0210	○		AC (≤135°C) 10 times of autoclaving	○	○	○	○	○
US tip Straight 30°	18254 0230			6pcs/box	○	○	○	○	○
	18252 0230	○		AC (≤135°C) 10 times of autoclaving	○	○	○	○	○
US tip Straight 30° Flared	18254 0235			6pcs/box	○	○	○	○	○
	18252 0235	○		AC (≤135°C) 10 times of autoclaving	○	●	●	○	○
US tip Curved 30° Flared	18254 0237			6pcs/box	○	○	○	○	○
	18252 0237	○		AC (≤135°C) 10 times of autoclaving	○	○	○	○	○
US tip Straight 45°	18254 0270			6pcs/box	○	○	○	○	○
	18252 0270	○		AC (≤135°C) 10 times of autoclaving	○	○	○	○	○
I/A tip Straight (φ0.3)	18258 0110	○		AC (≤135°C) Usable life: 1 year or 200 times of autoclaving	○	●	●	○	○
I/A tip Curved (φ0.3)	18258 0120	○		AC (≤135°C) Usable life: 1 year or 200 times of autoclaving	○	○	○	○	○
Silicone Sleeve Set type S	18241 3020			2.75-2.8mm 10pcs/box (Silicone Sleeve x2, Test Chambar, Wrench)	○	○	○	○	○
Silicone Sleeve Set type S *	18241 3020A			2.75-2.8mm 10pcs/box Made & sterilized in Thailand.	○	○	○	○	○
Silicone Sleeve type S	18241 M012L	○		2.75-2.8mm AC (≤135°C) 5 times of autoclaving	○	●	●	○	○













20G (Incisin Size 2.75 - 2.8mm)

Cataract

		Handpiece					30000	9000	9000R	24000	7000
Product name	Product code	Reusable		Remarks							
Handpiece	US Handpiece	18241 1600	○		AC ($\leq 135^{\circ}\text{C}$) Usable life: 1 year or 500 times of autoclaving	●	●	●	○	○	
	I/A Handpiece	18164 0000	○		Standard configuration (x2) AC ($\leq 135^{\circ}\text{C}$) Usable life: 1 year	●	●	●	○	○	
	Irrigation Handpiece (Long)	18026	○		AC ($\leq 135^{\circ}\text{C}$) Usable life: 1 year	○	○	○	○	○	
	Irrigation Handpiece (Short)	18027	○		Standard configuration (x2) AC ($\leq 135^{\circ}\text{C}$) Usable life: 1 year	○	○	○	○	○	
	Irrigation Handpiece (Very short)	18028	○		AC ($\leq 135^{\circ}\text{C}$) Usable life: 1 year	○	○	○	○	○	
	Aspiration Handpiece (Very short)	18074	○		AC ($\leq 135^{\circ}\text{C}$) Usable life: 1 year	○	○	○	○	○	



		Phaco Accessories					30000	9000	9000R	24000	7000
Product name	Product code	Reusable		Remarks							
Phaco Accessories	Wrench for Tip	18241 M042	○			●	●	●	○	○	
	Test Chamber	18084 M012	○			○	●	●	○	○	
	Sterilization Case (For US-handpiece)	18186 M018	○			●	●	●	○	○	
	Sterilization Case (For I/A and others)	18186 M019	○			●	○	○	○	○	

Tubing

Tubing										
Product name	Product code	Reusable		Remarks	30000	9000	9000R	24000	7000	
Phaco Surgery Cassette Pack	18601 7510			10pcs/box	●A	○				
Phaco Surgery Cassette w/ Infusion Tube *	18601 7540A			Made & sterilized in Thailand.	○	○				
Vit Surgery Cassette Pack	18601 7530			10pcs/box with Male to Male Connector	●P					
Phaco Surgery Cassette w/o Tube Pack	18602 7610			6pcs/box without I/A Tube	○	●				
I/A Tube	18602 7620	○		AC ($\leq 135^{\circ}\text{C}$) 5 times of autoclaving	○	●				
Infusion Tube	18241 M026				●A	○	●	○	○	
Vitreotomy Irrigation Tube	18271 2101P			10pcs/box	●P			○		
Air Irrigation Tubing	18271 2102P			10pcs/box	●P			○		
Injector	18271 2103P			10pcs/box	●P					
Extention Tube (180cm)	18271 1440P			25pcs/box ($\phi 2.9 \times \phi 4.9 \times 1800\text{mm}$)	○					
Irrigation Sleeve for Vitreous Cutter	18271 7140			Standard configuration For 20G	●A	●	●	○		
I/A Tube Kit	18603 7511	○		AC ($\leq 135^{\circ}\text{C}$) 10 times of autoclaving			●			
Drainage Bag	18603 M120P			10pcs/box			●			

* Please check the clearance of the regulatory procedure in your country or region.

Diathermy

Diathermy										
Product name	Product code	Reusable		Remarks	30000	9000	9000R	24000	7000	
Diathermy Cable	18241 M981	○		2.0m	●A	●	●	○	○	
Diathermy Cable	18241 0991			10pcs/box	○	○	○	○	○	
Diathermy Pencil (Straight 18G)	18241 0976			10pcs/box	○	○	○	○	○	
Diathermy Pencil (Straight 18G)	18241 M866	○		Nonsterile 20 times of autoclaving	●A	●	●	○	○	
Diathermy Pencil (Curved 18G)	18241 M867	○		Nonsterile 20 times of autoclaving	○	○	○	○	○	
Diathermy Pencil (Straight 20G)	18241 M868	○		Nonsterile 20 times of autoclaving	○	○	○	○	○	
Diathermy Pencil (23G/ Tapered)	18241 0978			5pcs/box	●P	○	○	○	○	
Diathermy Pencil (Straight 23G/ Tapered)	18271 M969	○		Nonsterile 20 times of autoclaving	○	○	○	○	○	
Diathermy Pencil (25G/ Tapered)	18241 0979			5pcs/box	○	○	○	○	○	
Diathermy Forceps (Coaptation)	18241 M861	○		Nonsterile 20 times of autoclaving	○	○	○	○	○	
Diathermy Forceps (Straight, Insulated)	18241 M862	○		Nonsterile 20 times of autoclaving	○	○	○	○	○	
Diathermy Forceps (Angled)	18241 M864	○		Nonsterile 20 times of autoclaving	○	○	○	○	○	
Diathermy Forceps (Straight)	18241 M865	○		Nonsterile 20 times of autoclaving	●A	○	○	○	○	

Diathermy



DECLARATION OF CONFORMITY

Manufacturer's name NIDEK Co. Ltd.
 Manufacturer's address 34-14 Maehama, Hiroishi-cho, Gamagori, Aichi 443-0038, Japan
NIDEK s.a.
 European Representative Europarc, 13 rue Auguste Perret, 94042 Créteil, France
 Identification of device Ophthalmic Surgical System
 Model No. CV-9000
 Classification(Annex IX, MDD) II b
 Category (for RoHS) 8

We herewith declare under our sole responsibility that the above mentioned products meet the provisions of the following EC Council Directives and Standards. All supporting documentation is retained under the premises of the manufacturer. The supporting technical documentation is also retained at NIDEK s.a., Europarc, 13 rue Auguste Perret, 94042 Créteil, France.

General applicable directives	Notified Body	Date CE Mark was affixed
COUNCIL DIRECTIVE 93/42/EEC of 14 June 1993 concerning medical devices.	TÜV SÜD Product Service GmbH, Ridlerstr. 65, 80339 München, Germany Certificate : G1 18 04 23653 195 (Annex II , Section 3 of MDD)	February 6, 2012
COUNCIL DIRECTIVE 2011/65/EU of 8 June 2011 concerning restriction of the use of certain hazardous substances.	N/A	June 30, 2014

Place: Aichi, Japan

Effective date : July 24, 2018

Signed by

H. Torii

Date of signature : June 1, 2018

Hiroyuki Torii
Director,
Top Management,
NIDEK Co., Ltd.