## ordering information


Silk Vista

Reference	Nominal stent		Unconstrained stent		Delivery
	Ø (mm)	Length (mm)	Ø (mm)	Length (mm)	catheter I.D.
SILK_V_3,50X15	3,50	15	3,75	11,50	021"
SILK_V_3,50X20		20		15,50	
SILK_V_3,50X25		25		19,50	
SILK_V_3,50X30		30		22,50	
SILK_V_3,50X35		35		26,50	
SILK_V_3,50X40		40		30	
SILK_V_3,75X15	3,75	15	4,00	12	
SILK_V_3,75X20		20		15	
SILK_V_3,75X25		25		19	
SILK_V_3,75X30		30		22	
SILK_V_4,00X15	4,00	15	4,25	11,50	
SILK_V_4,00X20		20		15	
SILK_V_4,00X25		25		18,50	
SILK_V_4,00X30		30		22	
SILK_V_4,25X15	4,25	15	- 4,50	11,50	
SILK_V_4,25X20		20		15	
SILK_V_4,25X25		25		18	
SILK_V_4,25X30		30		22,50	
SILK_V_4,50X15	4,50	15	4,75	11,50	
SILK_V_4,50X20		20		15,50	
SILK_V_4,50X25		25		19	
SILK_V_4,50X30		30		22,50	
SILK_V_4,75X15	4,75	15	5,00	11,50	
SILK_V_4,75X20		20		15	
SILK_V_4,75X25		25		18	
SILK_V_4,75X30		30		21	

The self-expandable SIk vista stents are designed for the treatment of intractantal aneutysms. Class III CE0297 in compliance with Medical Device Directive (MDD 93/42/EEC amended by 2007/47/EC). Manufactured by BAIT Extrusion S.A.S. Carefully read the instructions for use before use. First CE marking: 2020. The content of this document, in particular data, information, trademarks and logos is BAIT SAS and affiliates (see propery, © 2020 BAIT SAS and affiliates) and affiliate instruction and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BAIT SAS and offiliates (copyrights and other intellectual proprietary rights. This document with associated private is non-contractual and is solely dedicated to healthcare professionals and BAITs distributors (BAITs augments). The products commercialized by BAIT shall exclusively be used in accordance with the instructions for use included in the baxes. DC056GB (06/2020)

www.balt-corp.com







Advanced flow diversion, enhanced navigability\*\*

Self-expandable flow diverter designed for the treatmen of intracranial aneurysms.



.021"

delivery catheter for  $\varnothing$  up to 5,00mm unconstrained

## Leading design, towards more control

Bringing a high wall apposition force\* for a similar chronic outward force To improve wall apposition while ensuring an atraumatic & smooth opening of the stent

> Homogenous stent surface Improved compared to previous generation of ED\*

Escalating flow diversion\*
Thanks to an increased mesh density
of the 48 braided wires

## Upgraded navigability widening opportunitie

Compatibility

All references are delivered through a .021" catheter for Ø up to 5,00mm unconstrained

Advanced drawn filled tube technolog

Entirely visible device

To enhance the radiopacity of each & every wire of the stent









high wall apposition force\*

with the courtesy of Dr Pereira, Canada