

# SOFIA® EX

Intracranial Support Catheter

EMEA

## EXcellent Trackability. EXtra Stability.

### EXcellent Trackability for distal navigation and positioning.<sup>1</sup>

- Soft distal tip designed for easy navigation to the desired target location




### EXtra Stability for reliable device delivery.<sup>2</sup>

- Exhibits minimal tip movement when deploying and re-sheathing treatment devices<sup>2</sup>
- Kink resistant construction provides consistent support for device delivery in tortuous anatomy


Distal Soft Length: 9 cm

Proximal Length: 96 cm



Distal Soft Length: 9 cm

Proximal Length: 106 cm



## SOFIA® EX Intracranial Support Catheter

Catalog Number	Distal/Proximal ID		Distal OD			Proximal OD			Usable Length (cm)	Distal Soft Length (cm)	Tip Shape
	(in)	(mm)	(F)	(in)	(mm)	(F)	(in)	(mm)			
ISC5105ST	0.058	1.47	5.2	0.068	1.73	5.4	0.071	1.80	105	9	Straight
ISC5115ST									115		

1. TR19-145A; n=3
2. TR19-073A; n=3

### INDICATIONS FOR USE

The SOFIA® EX Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA® EX Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA® EX Catheter is not intended for use in coronary arteries.

### For Professional Use Only

MICROVENTION and SOFIA® are registered trademarks of MicroVention, Inc.  
© 2019 MicroVention, Inc. MM872(i) EMEA 06/19

### MicroVention Worldwide Innovation Center

PH +1.714.247.8000

35 Enterprise

Aliso Viejo, CA 92656, USA

MicroVention UK Limited

MicroVention Europe, S.A.R.L.

MicroVention Deutschland GmbH

Web

PH +44 (0) 191 258 6777

PH +33 (1) 39 21 77 46

PH +49 211 210 798-0

microvention.com