

MEDIVATORS™ ISA™

Endoscope Reprocessor
Washer-Disinfecter



EVOLUTIONARY PROCESS, REVOLUTIONARY SYSTEM

The MEDIVATORS ISA Endoscope Reprocessor is a rapid, efficient and flexible reprocessing system designed for endoscopy units demanding the most advanced technology and the highest quality standards.



FAST & SIMPLE OPERATION FOR IMPROVED EFFICIENCY

- Full cycle in 20 minutes with channel monitoring
- Simple, user-friendly handling and operation
- Insertion of endoscopes and operator information through RFID system - helps to avoid errors and saves time



MEDIVATORS ISA Endoscope Reprocessor The Fast And Flexible Reprocessing System

COMPACT & ERGONOMIC DESIGN FOR MAXIMUM USER COMFORT

- Large and easy to load basin accommodates all types of flexible endoscopes
- Intuitive interface assures easy and fast access to all the cycle information
- Compact size contributes to easy installation, operation and maintenance



SAFETY FOR OPERATORS, PATIENTS & ENDOSCOPES

- Hands-free operation and RFID-based technology for improved infection control
- Full traceability through printed and electronic cycle documentation
- Compliant to EN ISO 15883 parts 1, 4 and 5. Tested with multiple test soils and in accordance with German cleaning standards



Foot pedal allows hands-free door lid opening and closing

THE RIGHT FORMULA FOR THE BEST PERFORMANCE & RESULT


- Dedicated room temperature chemicals compatible with endoscope material
- Provides low temperature reprocessing at 20° C to 40° C
- Improved water and chemical consumption
- Designed specifically for use with MEDIVATORS ISA Endoscope Reprocessor



PART OF OUR COMPLETE CIRCLE OF GI REPROCESSING PRODUCTS

Efficient and effective high-level disinfection.

Incorporating STERIS products into every phase of your endoscope reprocessing will help improve department efficiency, reduce the risk of cross-contamination and provide reliable outcomes.



PRODUCT SPECIFICATIONS

Dimensions

700 mm (l) x 1400 mm (h) x 650 mm (d) - monitor included

Weight

75 kg

Electrical Requirements

100 V-240 V, 50-60Hz

Operating Room Temperature

5° - 40° C

Humidity

20 - 80%

Reprocessing Time

20 minutes (full cycle), 12 minutes only (HLD cycle),
25 minutes (double cleaning cycle)

Process Capacity

1 endoscope per cycle

High Level Disinfectant/Sterilant

ISASPOR™ SINGLE SHOT High-Level Disinfectant and Sterilant,
two parts, Peracetic acid and ISAZONE™ Molecule

Detergent

ISACLEAN™ AER Detergent and Decontaminant, Multienzymatic
PROTEAZONE™ PLUS Detergent and Cleaner, Enzymatic

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	QUANTITY
MEDISA3	MEDIVATORS ISA Endoscope Reprocessor Washer-Disinfectant using ISASPOR SINGLE SHOT High-Level Disinfectant and Sterilant 5L and ISACLEAN AER Detergent 1L	1 x AER
MEDISA4	MEDIVATORS ISA Endoscope Reprocessor Washer-Disinfectant using ISASPOR SINGLE SHOT High-Level Disinfectant and Sterilant 10L and ISACLEAN AER Detergent 1L	1 x AER
MEDISA6	MEDIVATORS ISA Endoscope Reprocessor Washer-Disinfectant using ISASPOR SINGLE SHOT High-Level Disinfectant and Sterilant 5L and ISACLEAN AER Detergent 5L	1 x AER
MEDISA7	MEDIVATORS ISA Endoscope Reprocessor Washer-Disinfectant using ISASPOR SINGLE SHOT High-Level Disinfectant and Sterilant 10L and ISACLEAN AER Detergent 5L	1 x AER
MEDISA8	MEDIVATORS ISA Endoscope Reprocessor Washer-Disinfectant using ISASPOR SINGLE SHOT High-Level Disinfectant and Sterilant 10L and PROTEAZONE PLUS Detergent 5L	1 x AER
MEDISA9	MEDIVATORS ISA Endoscope Reprocessor Washer-Disinfectant using ISASPOR SINGLE SHOT High-Level Disinfectant and Sterilant 5L and PROTEAZONE PLUS Detergent 5L	1 x AER
500187	ISASPOR SINGLE SHOT High-Level Disinfectant and Sterilant	10L Part A + 10L Part B
500188	ISASPOR SINGLE SHOT High-Level Disinfectant and Sterilant	5L Part A + 5L Part B
4005635	ISACLEAN AER Detergent and Decontaminant (EMEA 1)	12 x 1L
4005636	ISACLEAN AER Detergent and Decontaminant (EMEA 2)	12 x 1L
4005637	ISACLEAN AER Detergent and Decontaminant (EMEA 1)	2 x 5L
4006082	ISACLEAN AER Detergent and Decontaminant (EMEA 2)	2 x 5L
500199	PROTEAZONE PLUS Detergent and Cleaner	2 x 5L
11061031	MEDIVATORS ISA Water Filters Kit (0,1µm + 0,45µm)	1
ISA4400	ISA Air Filter	1
4008833	MEDIVATORS ISA Water Filter (0,1µm)	1
4008832	MEDIVATORS ISA Water Filter (0,45µm)	1
ISAROTS2	ISA Printer Roll	1

Contact your STERIS representative for more information or for the detergent and disinfectant combination available in your country.