

Surgical Needles

MEDITÜ

KFT. Needle Technology for Surgery

IMPORTANT!

Dear Customers!

In as much as you would like to print the catalogue
and place an order by it, please make certain that you adjust the document to 100 % when printing,
you can thus avoid minor discrepancies of the potential automatic settings,
which can lead to improper sizes and inappropriate order.

CATALOGUE



MEDITÚ Surgical needles
KFT. Needle Technology for Surgery

MEDITÚ KFT. Needle Technology for Surgery

For 30 years the Meditú Kft. has been **manufacturing surgical needles** and with its decades of experience it is sure to guarantee quality to its customers. Our production facility in Makó in south-eastern Hungary offers the following range of products:

- **Surgical needles**
- **Atraumatic needles**
- **Needles for Drains**
- **Special, unique medical needles**

The success of our company is based on **paying attention to fulfilling our customers' needs** and **ensuring the highest quality**. Our needles are manufactured from **high-quality stainless steel**. We place great emphasis on following all environmental compliances during the manufacturing process as well as with our products themselves.

Meditú Ltd. exports needles to over **40 countries** throughout the world.

Our quality is ensured by fulfilling the requirements of the standards of ISO 9001:2000 and we circulate our products with CE conformity marking.

Diameters and Lengths of Surgical Needles / mm /

Type		Diameter	Length
A	000	1,80	100
B	00	1,60	92
G	0	1,60	85
AR	1	1,40	80
BR	2	1,40	75
GR	3	1,30	70
BT	4	1,30	65
GT	5	1,20	60
	6	1,20	55
	7	1,10	50
	8	1,00	45
	9	1,00	40
	10	0,90	36
	11	0,80	32
	12	0,80	28
	13	0,70	25
	14	0,70	22
	15	0,60	20
	16	0,60	18
GA	000	2,00	90
GAR	00	1,80	80
GAT	0	1,80	70
	1	1,60	60
	2	1,60	55
	3	1,50	50
	4	1,40	45
	5	1,30	40
	6	1,20	35
	7	1,10	30
	8	1,00	25
	9	0,90	22
	10	0,80	19
	11	0,80	16
E	000	1,30	70
ET	00	1,20	60
	0	1,20	55
	1	1,00	47
	2	1,00	40
	3	0,80	32
	4	0,80	28
	5	0,70	22
	6	0,60	18
	7	0,50	15
	8	0,50	12
OB	0	0,60	36
OG	1	0,60	30
OA	2	0,50	25
	3	0,50	21
	4	0,50	18
	5	0,50	15
	6	0,50	12
H	00	1,40	75
HR	0	1,20	60
	1	1,10	52
	2	1,00	45
	3	0,90	38
	4	0,80	30
	5	0,70	22

Type		Diameter	Length
PB	00	0,80	75
PC	0	0,80	65
PD	1	0,70	55
PH	2	0,70	45
PA	3	0,60	36
	4	0,60	30
	5	0,60	25
	6	0,60	22
	7	0,50	20
	8	0,50	18
	9	0,50	15
	10	0,50	12
AV	1	2,00	100
	2	2,00	90
	3	1,80	80
	4	1,80	70
	5	1,60	60
	6	1,60	55
	7	1,50	50
	8	1,40	45
	9	1,30	40
	10	1,20	35
	11	1,10	30
	12	1,00	25
	13	0,90	22
BV	1	2,00	100
	2	2,00	90
	3	1,80	80
	4	1,80	70
	5	1,60	60
	6	1,60	55
	7	1,50	50
	8	1,40	45
	9	1,30	40
	10	1,20	35
	11	1,10	30
	12	1,00	25
	13	0,90	22
BE	1	1,40	80
GE	2	1,40	75
BER	3	1,30	70
GER	4	1,30	65
	5	1,20	60
	6	1,20	55
	7	1,10	52
	8	1,10	50
	9	1,00	47
	10	1,00	45
	11	1,00	40
	12	0,90	38
	13	0,90	36
	14	0,80	34
	15	0,80	32
	16	0,80	30
	17	0,80	28
	18	0,70	25
	19	0,70	23
	20	0,70	22
	21	0,60	20
	22	0,60	18

Type		Diameter	Length
Keith	0	1,00	80
"Keith"R"	1	0,90	75
	2	0,90	70
	3	0,90	65
	4	0,80	60
	5	0,80	55
	6	0,80	50
	7	0,70	45
König	1	1,20	55
	2	1,10	50
	3	1,00	45
	4	1,00	40
	5	0,90	36
	6	0,80	34
	7	0,80	32
	8	0,80	30
	9	0,80	26
	10	0,70	23
	11	0,60	20
	12	0,60	18
Neumann	1	0,70	23
	2	0,60	18
Kalt	00	0,35	30
	0	0,35	25
	1	0,35	22
	2	0,35	18
	3	0,35	16
	4	0,35	12
Mayo	1	1,50	50
Mayo T	2	1,40	45
	3	1,40	42
	4	1,30	40
PM/S	1	2,00	145
PM/A	2	2,00	130
	3	2,00	115
	4	2,00	100
	5	1,80	85
	6	1,80	75
Payr U	00	0,35	30
	0	0,35	25
	1	0,35	22
	2	0,35	18
	3	0,35	16
	4	0,35	12
Payr ger.	1	0,35	30
	2	0,35	25
	3	0,35	22
	4	0,35	20
	5	0,35	16
	6	0,35	12
Ferguson	1	1,40	74
	2	1,30	63
	3	1,10	50
	4	1,00	40
	5	0,80	34
	6	0,70	23
Hernia	1	1,40	45
	2	1,30	40
	3	1,20	35
	4	1,10	30

The illustrations represent the natural size of the needles.

Eye



regular eye

odd numbers



spring eye

even numbers

Curvature



straight



half curved



$\frac{3}{8}$ circle



$\frac{1}{2}$ circle



$\frac{5}{8}$ circle



hook-shaped



s-shaped

Point



reverse cutting



round



trocar



cutting

Ref.-No.

ordering example

<i>B 502</i>	<i>11</i>	<i>-</i>
type of needle	size	non-sterile













Packaging


1 doz. / PE-pouch (non-sterile)
individual packaging on request









Special items


We produce also on request. Please contact us.

G / GR / GT ½ -circle












Ref.-No.:	Eye	Profile
G 503		
G 504		
GR 513		
GR 514		
GT 515		
GT 516		


G / GR / GT

E / ET ½-circle













Ref.-No.:	Eye	Profile
E 529		
E 530		
ET 527		
ET 528		


E / ET


GA / GAR / GAT ½-circle

Ref.-No.:	Eye	Profile
GA 507		
GA 508		
GAR 497		
GAR 498		
GAT 499		
GAT 500		

GA / GAR / GAT

B / BR / BT ¾-circle

Ref.-No.:	Eye	Profile
B 501		
B 502		
BR 509		
BR 510		
BT 511		
BT 512		


B / BR / BT


Hernia

½-circle


Ref.-No.:	Eye	Profile
Hernia 571		
Hernia 572		

Hernia




Keith / Keith R

straight









Ref.-No.:	Eye	Profile
Keith 517		
Keith 518		
KeithR 417		
KeithR 418		

Keith / Keith R




H / HR

5⁄8-circle

Ref.-No.:	Eye	Profile
H 521		
H 522		
HR 421		
HR 422		

H / HR



A / AR

half curved

Ref.-No.:

A 505


A 506

AR 495


AR 496

Eye

Profile



A / AR



König

hook-shaped

Ref.-No.:

König 525


König 526

Eye









Profile




König







BE / BER ¾-circle

Ref.-No.:	Eye	Profile
BE 559		
BE 560		
BER 561		
BER 562		


BE / BER











Ferguson ½-circle


Ref.-No.:	Eye	Profile
Ferguson 545		
Ferguson 546		

Ferguson




GE / GER ½-circle

Ref.-No.:	Eye	Profile
GE 547		
GE 548		
GER 565		
GER 566		



GE / GER

Mayo / Mayo T ½-circle

Ref.-No.:	Eye	Profile
Mayo 541		
Mayo 542		
MayoT 441		
MayoT 442		


Mayo / Mayo T


PD

$\frac{3}{8}$ -circle



Ref.-No.:	Eye	Profile
PD 533		●
PD 534		●

PD




PB

$\frac{1}{2}$ -circle



Ref.-No.:	Eye	Profile
PB 537		●
PB 538		●

PB




PC

straight



Ref.-No.:	Eye	Profile
PC 535		●
PC 536		●

PC




PA

half curved



Ref.-No.:	Eye	Profile
PA 531		●
PA 532		●

PA




PH

5/8-circle





Ref.-No.:	Eye	Profile
PH 543		●
PH 544		●

PH




OA

half curved





Ref.-No.:	Eye	Profile
OA 549		
OA 550		

OA




OB

$\frac{3}{8}$ -circle





Ref.-No.:	Eye	Profile
OB 551		
OB 552		

OB




OG

$\frac{1}{2}$ -circle

Ref.-No.:	Eye	Profile
OG 553		
OG 554		

OG




Kalt

$\frac{3}{8}$ -circle



Ref.-No.:	Eye	Profile
Kalt 557		

Kalt




Neumann

straight

Ref.-No.:	Eye	Profile
Neumann 573		
Neumann 574		

Neumann




Payr C

straight

Ref.-No.: Eye Profile
PayrC 555



Payr C




Payr G

½-circle

Ref.-No.: Eye Profile
PayrG 655



Payr G



PM/A

half curved

Ref.-No.:


Eye

Profile

PM/A 623



PM/A



PM/S

s-shaped

Ref.-No.:


Eye

Profile

PM/S 627



PM/S




BV

$\frac{3}{8}$ -circle

Ref.-No.:
BV 580



BV




AV

half curve

Ref.-No.:
AV 590



AV



Sales and transportation conditions

Transportation occurs with the following conditions. We also accept special requests for transportation based on prior agreement, which only refers to that particular business or business partner. The validation of the contract expires upon completion of the agreement.

Prices and offers

Our prices and offers do not require engagement on our part without a written delivery agreement. For international customers, the price of delivery is always determined by the actual exchange rate at the time of the order.

Orders: occur when you give a model number, the type of the needle eye's formation and the amount.

Special manufacturing

We accept orders in small amounts as well. When submitting the offer please present a master pattern or at least one dimensional drawing or sketch. In the case of the latter two, the placer of the order should check the manufacture piece/s. As may be easily understood we are not going to start investigation if special designs are protected in some way and therefore we decline any responsibility.

Catalogue types

We reserve the right to change the construction type whenever necessary.

Delivery time

We do our best to deliver precisely, and according to the stipulated time agreed upon. Partial deliveries are also accepted by prior agreement. We reserve the right to withhold delivery or retire from the contract without granting indemnity in cases of improper business behavior on behalf of our partners.

Dispatch

We deliver at our expense regardless of the value of the delivered products using the most favorable packaging and transport methods.

Claims for faculty goods can only be accepted in writing at the latest, two weeks after arrival of delivery.

Payment: must be effected at the term and time agreed upon based on the amount shown on the invoice. In case of delay of payment, interest rates will be applied based on actual regulations. The buyer is not permitted to withhold payments or to reckon up any service.

Property: The delivered product- even if they have been manufactured- remain our property and can not be pledged or sold until fully paid for.

Exportation of our products requires our special agreement.