

## Electrodes

### Two-Stem Electrodes with Stabilizers, for Working Elements 27040DB/EB

For use with 24/26 Fr. resectoscope sheaths

The cutting loops are delivered with a wire diameter of 0.35 mm. Loops with 30 or 40 as the last digit of the order number indicate a wire diameter of 0.30 mm or 0.40 mm.



27040GP1

#### Cutting Loops

Distal Tip	24/26 Fr. color code: yellow	Instrument Description
	27040GP1	<b>Cutting Loop</b> , bipolar
	27040GP140	<b>Cutting Loop</b> , bipolar
	27040GD1	<b>Cutting Loop</b> , bipolar, small
	27040BL1	<b>Cutting Loop</b> , bipolar, pointed

#### Cutting Loops

Distal Tip	24/26 Fr. color code: yellow/orange	Instrument Description
	27040GP130	<b>Cutting Loop</b> , bipolar, diameter 0.30 mm
	27040JB1	<b>Cutting Loop</b> , bipolar, longitudinal
	27040JB130	<b>Cutting Loop</b> , bipolar, longitudinal, diameter 0.30 mm
	27040JBE130	<b>Cutting Loop</b> , bipolar, rectangular, longitudinal, diameter 0.30 mm

#### VapoEnucleation Electrode

Distal Tip	24/26 Fr. color code: yellow	Instrument Description
	27040VE	<b>VapoEnucleation Electrode</b> , bipolar, hemispherical