

PRODUCT INFORMATION

10X TBE Electrophoresis Buffer

#_ 1 liter

Lot: _ Expiry Date: _

Store at room temperature.

www.thermoscientific.com/onebio

1X Buffer Composition

89 mM Tris, 89 mM boric acid, 2 mM EDTA.

Applications

- ☐ Nucleic acid agarose and polyacrylamide gel electrophoresis;
- ☐ Agarose and polyacrylamide gel preparation.

Note

- ☐ Buffer concentrate should be diluted to a working concentration of 1X before use.
- ☐ For each electrophoresis fresh 1X buffer should be used.

CERTIFICATE OF ANALYSIS

Tested for the absence of deoxyribonucleases, ribonucleases and proteases.

Quality authorized by:



Jurgita Zilinskiene

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Rev.4.