ordering information

Reference	Vessel Ø (mm)	Silk+ length at	Unconstrained stent		Delivery
		nominal ø (mm)	Ø (mm)	Length (mm)	catheter
SILK 2,0×15	1,50 to 2,25	15	2,5	10	VASCO+21
SILK 2,0×20	1,30 to 2,23	20	2,3	11	
SILK 2,5×15		15		10	VASCO+21
SILK 2,5x20	2,00 to 2,75	20	3,0	13	
SILK 2,5x25		25		16	
SILK 3,0×15		15		10	VASCO+21
SILK 3,0×20	0.50. 0.05	20	3,5	13	
SILK 3,0×25	2,50 to 3,25	25		16	
SILK 3,0x30		30		19	
SILK 3,5x15		15	4,0	11	VASCO+21
SILK 3,5×20		20		12	
SILK 3,5×25	3,00 to 3,75	25		14	
SILK 3,5×30		30		15	
SILK 3,5×35		35		1 <i>7</i>	
SILK 4,0x15	- 3,50 to 4,25	15	4,5	10	- VASCO+21
SILK 4,0x20		20		12	
SILK 4,0×25		25		14	
SILK 4,0×30		30		16	
SILK 4,0x35		35		18	
SILK 4,0x40		40		20	
SILK 4,5x15		15		8	VASCO+21
SILK 4,5×20	4,00 to 4,75	20	5,0	10	
SILK 4,5×25		25		12	
SILK 4,5x30		30		15	
SILK 4,5×35		35		17	
SILK 4,5x40		40		19	
SILK 5,0x25		25		11	
SILK 5,0×30	4,50 to 5,25	30	5,5	13	VASCO+25
SILK 5,0×40		40		17	
SILK 5,5×25		25		11	
SILK 5,5x30	5,00 to 5,75	30	6,0	13	VASCO+25
SILK 5,5×40		40		1 <i>7</i>	

				nominal ø (mm)			
		prox	dist	prox	dist	Ø (mm)	Length (mr
pered	SILKP4,0D3,0x30	3,50 to 4,25	2,50 to 3,25	12	10	4,5/3,5	15
k+	SILKP 4,5D3,0x25	4,00 to 4,75	2,50 to 3,25	12	10	5,0/3,5	20
	SILKP 4,5D3,5x30	4,00 to 4,75	3,00 to 3,75	11	15	5,0/4,0	16

Vessel Ø (mm)

Reference

Silk+ length at Unconstrained stent

FI. Francis turprian.


SIIK- is a self-expanding intracranial stent indicated for use in the treatment of intracranial aneurysms in association with embolization coils. The content of this document, in particular data, information, trademarks and logos are BALLS.A.S and offiliates' sole property. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALLS.A.S and affiliates' copyrights and other intellectual proprietory rights @2017 BALLS.A.S and affiliates all rights reserved. This document with associated pictures are non-contractual and are solely dedicated to healthcare professionals and BALLS.A.S and affiliates' distributors. The products commercialized by BALLS.A.S and affiliates shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes. SIKEr is closs III.C.E marked [INEX/C-Mode (ED459)] according to the Medical Device Directive 93/42/EEC since 2012. DC025GB [06/2017]

www.balt.fr



Delivery catheter

VASCO+21





Combining reliability & effectiveness¹

82%12 months

Adequate behaviour even in

Flexibility thanks to the nickel titanium materials

Conformability given by the sliding-cell technology

Smooth wall apposition

Navigability with VASCO+21 low profile delivery microcatheter

Aneurysm exclusion & vessel

Resheathability







BALT SILK BROCHURE 2017 indd 3-4 09/06/2017 15:15