



**CILITA® Bimanual Irrigation/Aspiration handpieces** come in multiple options to provide best outcomes in AC irrigation and cortex removal.

### Color coding

helps differentiate between the instruments:

- Gold handle indicates aspiration;
- Blue handle indicates irrigation.

The **ribbed surface** of the handle ensures firm grip of the instrument.



### Regular tip

The regular smooth tip is a default option for most handpieces.



### Polishing tip

The polishing sandblasted tip allows to polish the capsule without additional instruments.

### Features of CILITA® bimanual handpieces:

- Bullet-shaped tip for easier introduction of the instrument into the AC
- Aspiration port  $D=0.35/0.3$  mm increases safety during aspiration and provides for optimum stability of the AC
- Port near the top of the tip and its specially beveled edges prevent closure of the port by the capsule wall, thus avoiding possible rupture of the capsule
- $90^\circ$  angle between the 2 fluid jets results in optimal irrigation
- Available in 2 options: with 20 mm and 18 mm long tubes (cannulas)
- Compatible with your existing Phaco or I/A System by standard Luer-connectors

# Product range

## Irrigation, 20 mm cannula

21 Ga	P-081	2 x 0.35 ports, Smooth handle
	P-081R	2 x 0.35 ports, Ribbed handle
22 Ga	P-081-22	2 x 0.35 ports, Smooth handle
	P-081R-22	2 x 0.35 ports, Ribbed handle
23 Ga	P-081-23	2 x 0.35 ports, Smooth handle
	P-081R-23	2 x 0.35 ports, Ribbed handle

## Irrigation, 18 mm cannula

21 Ga	P-083-21	2 x 0.4 ports, Ribbed handle
22 Ga	P-083-22	2 x 0.35 ports, Ribbed handle
	P-0831-22	2 x 0.35 ports, polishing tip, Ribbed handle
23 Ga	P-083-23	2 x 0.35 ports, Ribbed handle
	P-0831-22	2 x 0.35 ports, polishing tip, Ribbed handle

## Aspiration, 20 mm cannula

21 Ga	P-082	0.35 port, Smooth handle
	P-082R	0.35 port, Ribbed handle
	P-0821-21	0.35 port, polishing tip, Smooth handle
22 Ga	P-0821	0.35 port, polishing tip, Smooth handle
	P-0821R	0.35 port, polishing tip, Ribbed handle
23 Ga	P-0821-23	0.3 port, polishing tip, Smooth handle
	P-0821R-23	0.3 port, polishing tip, Ribbed handle

## Aspiration, 18 mm cannula

21 Ga	P-084-21	0.35 port, Ribbed handle
22 Ga	P-0841-22	0.3 port, polishing tip, Ribbed handle
23 Ga	P-0841-23	0.3 port, polishing tip, Ribbed handle



20 Moskovskoye shosse,  
390044 Ryazan, Russia  
+7 (4912) 779-678  
sales@cilita.com

EXCELLENT INSTRUMENTS  
FOR PROFESSIONAL HANDS

[www.cilita.com](http://www.cilita.com)