

14 Accessories

The accessory material that contacts the patients has undertaken the bio-compatibility test and is verified to be in compliance with ISO 10993-1.

WARNING

- **Use accessories specified in this chapter. Using other accessories may cause damage to the equipment or not meet the claimed specifications.**
- **Single-use accessories are not designed to be reused. Reuse may cause a risk of contamination and affect the measurement accuracy.**

CAUTION

- **The accessories may not meet the performance specifications if stored or used outside the specified temperature and humidity ranges. If accessory performance is degraded due to aging or environmental conditions, contact your service personnel.**
- **Check the accessories and their packages for any sign of damage. Do not use them if any damage is detected.**
- **Use the accessories before the expiry date if their expiry date is indicated.**
- **The disposable accessories shall be disposed of according to hospital's regulations.**

14.1 ECG Accessories

| Model | PN | Description | Applicable Patient |
|-------------------|----------------------|--|--------------------|
| EC6408** | 040-001642-00 | 12-pin integrative patient cable, reusable, AHA, Φ 4 banana plugs, defibrillation-proof | Adult |
| EC6409** | 040-001643-00 | 12-pin integrative patient cable, reusable, AHA, clip (MR), defibrillation-proof | Adult, Pediatric |
| EC6410** | 040-001644-00 | 12-pin integrative patient cable, reusable, IEC, Φ 4 banana plugs, defibrillation-proof | Adult |
| EC6411** | 040-001645-00 | 12-pin integrative patient cable, reusable, IEC, clip (MR), defibrillation-proof | Adult, Pediatric |
| ECG-FD09X4 | 040-007423-00 | 12-pin integrative patient cable, reusable, AHA, Φ 4 banana plugs, defibrillation-proof | Adult |
| ECG-FD10X4 | 040-007643-00 | 12-pin integrative patient cable, reusable, IEC, Φ 4 banana plugs, defibrillation-proof | Adult |
| ECG-FJ01 | 040-007640-00 | Multifunction limb electrode, reusable, IEC, copper alloy (MR) | Adult |
| ECG-FJ01 | 040-007641-00 | Multifunction limb electrode, reusable, IEC, Ag/AgCl (MR) | Adult |
| ECG-FJ02 | 040-007639-00 | Multifunction limb electrode, reusable, AHA, copper alloy (MR) | Adult |
| ECG-FJ02 | 040-007642-00 | Multifunction limb electrode, reusable, Ag/AgCl (MR) | Adult |
| ECG-FQX41 | 040-007428-00 | Multifunction chest electrode, reusable | Adult |
| ECG-FQX42 | 040-007429-00 | Multifunction chest electrode, reusable | Adult |
| SF06 | 040-002711-00 | Disposable electrode, 5 pcs/package | Adult |

Cablu ECG cu 10 fire pentru electrocardiograful BeneHeart R3, Mindray:

Model: ECG-FD10X4

PN: 040-007643-00





Name: ECG cable
 产品名称: 心电图线
 Model specification/ 原机型号: ECG-0310X4
 2511017
 2025-12
 2028-12
 4.8 Barona Plug EC
 Quantity/数量: 1
 其他附件详见说明书

mindray
 MPN:040-007643-00

生产许可证号: 粤监备证字第20150003号
 产品名称: 心电图线
 生产许可证号: 粤监备证字第20150004号

温州科晖医疗电子有限公司
 Wenzhou Kehui Medical Electronics Co., Ltd.
 地址: 温州市瓯海区
 温州科晖医疗电子有限公司 研发中心
 温州市瓯海区金康路1号 2号楼3层92号
 TEL: +86-755-27772721 Email: sales@kh13.com

WALKANG Ltd (www.CE-marking.eu)
 Category: Iah, IIV, Business Complete, Saraghemore
 A.S. Dany, 8748 RSE, Northern Ireland, UK.

