


Product-Information-Sheet

Polypectomy, EMR and ESD

Cold and Hot
Snares





- For polyp ligation prior to polypectomy
- Prevents postpolypectomy bleeding
- Preloaded loop, ready to use with handle
- Nylon-loop ligation keeps tissue damage to a minimum
- Opening diameter of 30 mm enables ligation of large polyps





Article number	Article Name	Quantity value	Min. working channel \varnothing value	Working length value	Loop \varnothing	Area Compatibility
N5382130	HX-400U-30	5	 2.8	2300 mm	30 mm	

Endoscopic Ancillaries

ENDOCUFF
VISION™



Article number	Article Name	Quantity value	Description	Endoscope Compatibility	Area Compatibility
027797	MH-438	1	Air/water valve, lubricant-free	Gastro-, colono-, and duodenoscopes: OES 40, EVIS 140, 145, 160, 165, 180, 185, 170, 190, 240, 260, 290, 1100, 1200, 1500, TGF-UC180J	



Article number	Article Name	Quantity value	Description	Endoscope Compatibility	Area Compatibility
027798	MH-443	1	Suction valve, lubricant-free	Gastro-, colono-, duodenoscopes: OES 40, EVIS 140, 145, 160, 165, 180, 185, 190, 240, 260, 290, 1100, 1200, 1500, TGF-UC180J	



Article number	Article Name	Quantity value	For Use With	Description	Area Compatibility
029421	MH-481	10	SSU-2; KV-2, -4, -5	Suction-line filter	

