


5F Bipolar Pacing Catheter™ w/o Balloon

Intracardiac Catheter for
Temporary Cardiac Pacing



Ordering Information

Model	Curve	Spacing	Shaft Diameter	Length	Order Number
5F Bipolar Pacing Catheter w/o Balloon BPX2502-10	J-Curve	10 mm	5 F	110 cm	445697

Cables	Length	Order Number
PK-67-S	0.8 m	128085
PK-67-L	2.6 m	123672

Adapters	Order Number
PA-1-C	349723
PK-1-B	123751

Product Highlights

Quality and Reliability

Allows safe, secure and accurate bipolar pacing to support critical procedures.

5F Bipolar Pacing Catheter™ w/o Balloon









Technical Data

Catheter Specification	
Product packaging contains	Temporary catheter
Tip electrode length	3 mm
Ring electrode length	3 mm
Number of electrodes	2
Electrode spacing	10 mm
Electrode material	Stainless steel
Shaft diameter	5F (1.67 mm)
Curve type	J-Curve
Curve shape	Fixed
Effective catheter length	110 cm
Connection plug	2 mm lead connector with touch protection
Recommended introducer size	5F or 6F

Transport and Storage Conditions	
Storage	Please store in a cool dark place to avoid exposure to fluorescent or sunlight
Light	Keep away from sunlight

Patient Cable	Cable Length	Suitable for Connection to:	TC-Adapt (2 x 2mm Pins)	Order Number	
Cable					
PK-67-S	sterile	0.8m	Reocor D, Reocor S, EDP 20/B, EDP 30/B, EDP 30/BP	-	128085
PK-67-L	sterile	2.6m	Reocor D, Reocor S, EDP 20/B, EDP 30/B, EDP 30/BP	-	123672
Adapter					
PA-1-C	sterile		PK-67-S / -L	1	349723
PA-1-B	sterile		PK-67-S / -L	2	123751

Temporary Catheter with Shrouded 2mm Pins

Patient Side (Electrode)	Adapter	BIOTRONIK Cable	Adapter	Device Side
	 PA-1-B	 PK-67*		Reocor S
	 PA-1-C	 PK-67*		Reocor D

Please attend the instructions for use regarding the connection and numbers of adapters and catheters!

* Using a dual chamber cable for a single chamber device, do always connect at the ventricular cable pin.

Using a single chamber cable for a dual chamber device, do always connect to the ventricular output of the device.

5F Bipolar Pacing Catheter is a trademark of Bioptimal International Pte. Ltd.

© BIOTRONIK SE&Co. KG
All rights reserved. Specifications are subject to modification, revision and improvement.

BIOTRONIK SE&Co. KG
Woermannkehe 1
12359 Berlin / Germany
www.biotronik.com

Manufacturer
Bioptimal International Pte. Ltd.
36 Jalan Tukang
Singapore 619266 / Singapore

 **BIOTRONIK**
excellence for life