

ORDIN DE PLATA NR.: 76 TIP.DOC. 1 :  
DATA EMITERII:30 mai 2022 :  
===== :  
PLATITI: 42000-00 LEI: Patruzeci si Doua Mii lei 00 b :  
ani :  
===== :  
PLATITOR: (R) S.C. "OXIVI CONTUL DE PLATI/CODUL IBAN :  
T-MED" S.R.L. MD44ML000000002251729503 :  
CODUL FISCAL :1007600044280 / :  
===== :  
PRESTATORUL PLATITOR CODUL BANCII: :  
BC"Moldindconbank"S.A. fil."Invest" Chisinau :MOLDMD2X329: :  
===== :  
BENEFICIAR (R) IMSP Spitalu CONTUL DE PLATI/CODUL IBAN :  
l Clinic Republican "Timofei MD32ML000000002251502448 :  
Mosnea CODUL FISCAL :1003600150783 / :  
===== :  
PRESTATORUL BENEFICIAR CODUL BANCII: :  
BC"Moldindconbank"S.A. :MOLDMD2X :  
===== :  
DESTINATIA PLATII: Pentru garantia pentru: TIPUL TRANSFERULUI :  
oferta la procedura de achizitie public: NORMAL/URGENT :N: :  
a nr. ocds-b3wdpl-MD-1652184471176 din 3: :  
1.05.2022 :  
: :  
: :  
: :  
===== :  
CODUL TRANZACTIEI:001: :  
DATA PRIMIRII:30/05/2022 : SEMNATURILE :  
DATA EXECUTARII: : EMITENTULUI :  
:----- :  
CONDUCTOR:Web Kojevnikov Dmitrii :  
MIIGfAYJKoZIhvcNAQcCoIIIGbTCCBmkCAQExCzAJBgUrDgMCGGUAMAsGCSqGSIb3: :  
DQEHAACCBIUwggSBMIIDaaADAgECAhNHAACEjCA/4xcrKCbfAAAAAISMMA0GCSqG: :  
SIb3DQEBCwUAMCIXIDAEbGNVBAMTF0NFULQxLUNBLU1vbGRpbmRjb25iYW5rMB4X: :  
DTIwMDMxNjA4NDUwM1oXDTIzMDMxNjA4NTUwM1owgbgxSzAJBgNVBAYTAk1EMRow: :  
YDVQIExeFSZXB1YmxyY2EgTW9sZG92YTERMA8GA1UEBxMIQ2hpc2luYXUxZzAV :  
\_\_\_\_\_  
(semnatura electronica)  
CONTABIL-SEF:Web Kojevnikov Dmitrii :  
MIIGfAYJKoZIhvcNAQcCoIIIGbTCCBmkCAQExCzAJBgUrDgMCGGUAMAsGCSqGSIb3: :  
DQEHAACCBIUwggSBMIIDaaADAgECAhNHAACEjCA/4xcrKCbfAAAAAISMMA0GCSqG: :  
SIb3DQEBCwUAMCIXIDAEbGNVBAMTF0NFULQxLUNBLU1vbGRpbmRjb25iYW5rMB4X: :  
DTIwMDMxNjA4NDUwM1oXDTIzMDMxNjA4NTUwM1owgbgxSzAJBgNVBAYTAk1EMRow: :  
YDVQIExeFSZXB1YmxyY2EgTW9sZG92YTERMA8GA1UEBxMIQ2hpc2luYXUxZzAV :  
\_\_\_\_\_  
L.S. (semnatura electronica)  
CONDUCTOR: \_\_\_\_\_  
(semnatura manuala)  
CONTABIL-SEF: \_\_\_\_\_  
(semnatura manuala)  
SEMNATURA PRESTATORUL L.S. :  
:----- :  
MOTIVUL REFUZULUI : L.S.



REPUBLICA



MOLDOVA

# CERTIFICAT DE ÎNREGISTRARE

**Societatea Comercială "OXIVIT-MED" S.R.L.**  
**ESTE ÎNREGISTRATĂ LA CAMERA ÎNREGISTRĂRII DE STAT**

*Numărul de identificare de stat - codul fiscal*  
**1007600044280**

*Data înregistrării*

**30.07.2007**

*Data eliberării*

**30.07.2007**

**Bordeianu Tatiana, registrator de stat**

*Funcția, numele, prenumele persoanei  
care a eliberat certificatul*

*semnătura*

**MD 0067985**





**I.P. "AGENȚIA SERVICII PUBLICE"**

Departamentul înregistrare și licențiere a unităților de drept

**EXTRAS**  
**din Registrul de stat al persoanelor juridice**

nr. 8871 din 05.05.2021

Denumirea completă: **Societatea Comercială «OXIVIT-MED» S.R.L.**

Denumirea prescurtată: **S.C. «OXIVIT-MED» S.R.L.**

Forma juridică de organizare: **Societate cu Răspundere Limitată.**

Numărul de identificare de stat și codul fiscal: **1007600044280.**

Data înregistrării de stat: **30.07.2007.**

Sediul: **MD-2032, bd. Decebal, 82, ap.(of.) 90, mun. Chișinău, Republica Moldova.**

Modul de constituire: **nou creată.**

Obiectul principal de activitate:

- 1 Importul, fabricarea, comercializarea, asistența tehnică și (sau) reparația dispozitivelor medicale și (sau) a opticii;**
- 2 Comerțul cu ridicata al parfumurilor și produselor cosmetice;**
- 3 Comerțul cu amănuntul al produselor cosmetice și de parfumerie, articolelor de toaletă;**
- 4 Intermedieri pentru vânzarea unui asortiment larg de mărfuri;**
- 5 Alte tipuri de comerț cu amănuntul în magazine nespecializate;**
- 6 Alte tipuri de comerț cu ridicata;**
- 7 Închirierea altor mașini și echipamente.**

Capitalul social: **5400 lei.**

**Administrator: KOJEVNIKOV DMITRII, IDNP 0972305012362,**

**Asociați:**

**1. KOJEVNIKOV DMITRII , IDNP 0972305012362**

**cota 5400.00 lei, ce constituie 100 %.**

Prezentul extras este eliberat în temeiul art. 34 al Legii nr. 220-XVI din 19 octombrie 2007 privind înregistrarea de stat a persoanelor juridice și a întreprinzătorilor individuali și confirmă datele din Registrul de stat la data de: 05.05.2021.

Specialist coordonator  
tel. 022-207-840

Lazari Aliona



EEI 0358094

# OXIVIT MED

c/f: 1007600044280; adresa: str. Decebal 82-90, or. Chișinău, Republica Moldova

telefon: + 373 22 808002; fax: + 373 22 808003

web: [www.oxivit-med.com](http://www.oxivit-med.com); e-mail: [info@oxivit-med.com](mailto:info@oxivit-med.com)

## **Lista fondatorilor companiei SRL „Oxivit-Med”**

| Nr. | Numele, Prenumele  | Codul Personal |
|-----|--------------------|----------------|
| 1   | Kojevnikov Dmitrii | 0972305012362  |



Benannt durch/Designated by  
Zentralstelle der Länder  
für Gesundheitsschutz  
bei Arzneimitteln und  
Medizinprodukten  
www.zlg.de  
ZLG-BS-244.10.08



Product Service

# EC Certificate

Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)  
(Devices in Class IIa, IIb or III)

**No. G1 077608 0079 Rev. 00**

**Manufacturer:**

**Covidien llc**

15 Hampshire Street  
Mansfield, MA 02048  
USA

**Product Category(ies):** Medical Instruments, Surgical Products  
and Hemostatic Materials:

- Surgical Suture Products, Pledgets and Retention Tapes
- Endoscopy Instruments and  
Accessories including Lubricant
- Surgical Staple, Clip Products and Accessories
- Manual Surgical Instruments
- Implantable Wound Dressing Materials
- Ultrasonic Surgical Devices and Accessories
- Suction / Irrigation Devices and Accessories
- Arthroscopy Implants, Instruments and Accessories
- Bone Wax
- Temporary Cardiac Pacing Lead
- Powered Stapling Systems

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

**Report No.:**

713164286

**Valid from:**

2019-09-13

**Valid until:**

2024-05-26

**Date,**

2019-09-13

Stefan Preiß  
Head of Certification/Notified Body

Page 1 of 2

TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany

TÜV®

ZERTIFIKAT • CERTIFICATE • 認證證書 • CERTIFICADO • CERTIFICAT



Product Service

## EC Certificate

### Full Quality Assurance System

**Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)  
(Devices in Class IIa, IIb or III)**

**No. G1 077608 0079 Rev. 00**

**Facility(ies):**

**Covidien Inc**  
15 Hampshire Street,  
Mansfield, MA 02048, USA

4

# CERTIFICATE

Number: 2090418

The management system of the organization(s) and locations mentioned on the addendum belonging to:

## Medtronic EMEA Medtronic B.V.

Earl Bakkenstraat 10  
6422 PJ Heerlen  
The Netherlands

including the implementation meets the requirements of the standard:

## ISO 9001:2015 EN ISO 13485:2016


### Scope:

Sales, order management, warehousing and distribution of medical devices.  
Including regulatory affairs, post market surveillance, technical service, customer education and spine  
loaner operations

Certificate expiry date: 1 July 2024  
Certificate effective date: 1 July 2021  
Certified since: 1 July 2006

This certificate is valid for the organization(s) and/or locations mentioned on the addendum.

DEKRA Certification B.V.



B.T.M. Holtus  
Managing Director



J.A. van Vugt  
Certification Manager

© Integral publication of this certificate and adjoining reports is allowed



# ADDENDUM

To certificate: 2090418

The management system of the organization(s) and/or location(s) of:

## Medtronic EMEA Medtronic B.V.

Earl Bakkenstraat 10  
6422 PJ Heerlen

Certified organization(s) and/or locations:

Different scope

Medtronic Trading NL B.V.  
Larixplein 4  
5616 VB Eindhoven  
The Netherlands

Sales, order management and distribution of medical devices.  
Including customer education

Medtronic Italia S.p.A.  
Via Varesina 162  
20156 Milano  
Italy

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Danmark A/S.  
Arne Jacobsens Alle 17  
2300 Copenhagen  
Denmark

Sales, order management and distribution of medical devices.  
Including customer education

Medtronic Finland Oy  
Lentajantie 3  
01530 Vantaa  
Finland

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic AB  
P.O. Box 1034  
164 21 Kista  
Sweden

Sales, order management and distribution of medical devices.  
Including customer education

Medtronic Norge AS  
Martin Linges vei 25  
1364 Fornebu  
Norway

Sales, order management and distribution of medical devices.  
Including customer education.

# ADDENDUM

To certificate: 2090418

The management system of the organization(s) and/or location(s) of:

## Medtronic EMEA Medtronic B.V.

Earl Bakkenstraat 10  
6422 PJ Heerlen

Medtronic Africa (Pty) Ltd.  
Waterfall Distribution Campus CNR  
K101 and Bridal Veil Road Waterfall  
Midrand  
1685 Gauteng  
South Africa

Sales, order management, warehousing and distribution of medical devices. Including customer education and spine loaner operations.

Medtronic Medikal Teknoloji Ticaret Ltd  
Sti  
Saray Mah. Esnaf Sk. Akkom Ofis Park  
Laodik Plaza Sitesi B Blok Apt: 2/8  
34764 Umraniye - Istanbul  
Turkey

Sales, order management and distribution of medical devices.  
Including customer education

Medtronic Ibérica S.A.  
Calle de Maria de Portugal, 11  
28050 Madrid  
Spain

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Ibérica S.A.  
WTC Almeda Park Placa de la Pau, s/n.  
Edificio 7, 3 piso Cornellà de Llobregat  
08940 Barcelona  
Spain

Sales, order management and distribution of medical devices.

Medtronic Portugal LDA-  
Rua Tomas da Fonseca Torre E, 11  
 piso  
1600 Lisboa  
Portugal

Sales, Order Management and distribution of medical devices  
including customer education.

Warehousing and distribution of medical devices, including spine  
loaner operations

# ADDENDUM

To certificate: 2090418

The management system of the organization(s) and/or location(s) of:

## Medtronic EMEA Medtronic B.V.

Earl Bakkenstraat 10  
6422 PJ Heerlen

Medtronic Portugal, LDA-  
Avenida Gomes Pereira 61B  
Benfica  
1600 Lisboa  
Portugal

Sales, Order Management and distribution of medical devices.  
Including customer education.

Warehousing and distribution of medical devices, including spine  
loaner operations.

Medtronic GmbH  
Earl-Bakken-Platz 1  
40670 Meerbusch  
Germany

Scope for EN ISO 13485:2016: Sales, order management and  
distribution of medical devices. Including customer education.  
ISO 9001:2015 excluded

Medtronic GmbH  
Mollsfeld 12  
40670 Meerbusch  
Germany

Scope for EN ISO 13485:2016: Sales, order management and  
distribution of medical devices. Including customer education.  
ISO 9001:2015 excluded

Medtronic Österreich GmbH  
Millennium Tower, 20th floor Handelskai  
94-96  
1200 Wien  
Austria

Sales, order management, warehousing and distribution of  
medical devices. Including customer education

Medtronic (Schweiz) AG  
Talstrasse 9  
3053 Munchenbuchsee  
Switzerland

Sales, order management, warehousing and distribution of  
medical devices. Including customer education

Medtronic France SAS  
9, boulevard Romain Rolland  
75014 Paris  
France

Sales, order management and distribution of medical devices.  
Including customer education

# ADDENDUM

To certificate: 2090418

The management system of the organization(s) and/or location(s) of:

## Medtronic EMEA Medtronic B.V.

Earl Bakkenstraat 10  
6422 PJ Heerlen

Medtronic Hellas S.A.  
Avenue Kifisias 24 Building B  
151 25 Marousi Pref. Attica  
Greece

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Hellas S.A. Diabetes Shop  
Mesogeion Avenue 2-4  
115 27 Athens  
Greece

Sales, order management and distribution of diabetes medical  
devices. Including customer education.

Medtronic Romania SRL  
Ploiesti 42-44, Building B, B2 Wing, 2nd  
floor, district 1 Baneasa Business &  
Technology Park  
013696 Bucharest  
Romania

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Hungária Kft.  
Bocskai ut 134-146 Cepuleit 3. emelet  
1113 Budapest  
Hungary

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Serbia Ltd.  
Bulevar Zorana Djindjica, 64a  
11070 Belgrade  
Serbia

Sales, order management and distribution of medical devices.  
Including customer education.

# ADDENDUM

To certificate: 2090418

The management system of the organization(s) and/or location(s) of:

## Medtronic EMEA Medtronic B.V.

Earl Bakkenstraat 10  
6422 PJ Heerlen

Medtronic Poland Sp.z o.o Medtronic  
Customer Care Center of Experience  
Warsaw  
Polna 11  
00-633 Warszawa  
Poland

Order management of medical devices.

Medtronic Trading Ltd.  
10 Hamada Street  
4673344 Herzliya  
Israel

Import, sales, order management and distribution of medical  
devices.  
Including customer education

Medtronic Czechia s.r.o.  
Prosek Point, Budova B, Prosecka  
852/66  
852 66 Praha  
Czech Republic

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Bulgaria EOOD  
48 Sitnyakovo blvd., R-N OBORISHT  
DISTR., floor 7  
1505 Sofia  
Bulgaria

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Limited  
Building 9, Croxley ParkHatters Ln  
WD18 8WW Watford  
United Kingdom

Sales, order management and distribution of medical devices.  
Including customer education.

# ADDENDUM

To certificate: 2090418

The management system of the organization(s) and/or location(s) of:

## Medtronic EMEA Medtronic B.V.

Earl Bakkenstraat 10  
6422 PJ Heerlen

Medtronic Ireland Limited  
Block 3090-3094 Lake Drive, Citywest  
Business Campus  
D24 NW2F Dublin  
Ireland

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic B.V.  
Medtronic Service & Repair EMEA  
Jan Campertstraat 21-A  
6416 SG Heerlen

Order management, warehousing and technical service of  
medical devices including field service EMEA.

Medtronic Slovakia s.r.o.  
CBC III, Karadzicova 12  
821 08 Bratislava  
Slovak Republic

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Belgium  
Burgemeester E. Demunterlaan 5  
1090 Brussel  
Belgium

Sales, Order Management and distribution of medical devices.  
Including customer education

Medtronic Croatia  
Folnegoviceva 1c  
10000 Zagreb  
Croatia

Sales, order management and distribution of medical devices.  
Including customer education.

Medtronic Slovenia  
Ameriska Ulica 8  
1000 Ljubljana  
Slovenia

Sales, order management and distribution of medical devices

Addendum expiry date: 1 July 2024  
Addendum effective date: 1 July 2021

# STAPLING SOLUTIONS FOR CHALLENGING APPLICATIONS


Stapling  
Product  
Catalogue  
2016



**Medtronic**  
Further. Together

# Stapling solutions for challenging tissue and applications

Surgical stapling products from Medtronic enable surgeons to handle the broadest range of tissues and applications with outstanding clinical performance.



# Table of Contents

## Endoscopic Staplers

|                                                                                                       |            |
|-------------------------------------------------------------------------------------------------------|------------|
| <b>Endo GIA™ reloads with Tri-Staple™ technology and Endo GIA™ Ultra universal staplers . . . . .</b> | <b>.05</b> |
| Endo GIA™ Ultra universal handle . . . . .                                                            | .05        |
| Endo GIA™ loading units . . . . .                                                                     | .05        |
| Endo GIA™ curved tip reload with Tri-Staple™ technology . . . . .                                     | .06        |
| Endo GIA™ radial reload with Tri-Staple™ technology . . . . .                                         | .06        |
| <b>iDrive™ Ultra powered stapling system . . . . .</b>                                                | <b>.07</b> |
| iDrive™ Ultra powered stapling system . . . . .                                                       | .07        |
| <b>Endo GIA™ universal stapling system . . . . .</b>                                                  | <b>.08</b> |
| Endo GIA™ universal stapling system. . . . .                                                          | .08        |
| Endo GIA™ universal Reticulator™ loading units . . . . .                                              | .08        |
| Endo GIA™ universal straight loading units . . . . .                                                  | .09        |
| <b>MultiFire Endo GIA™ 30 staplers and reloads . . . . .</b>                                          | <b>.10</b> |
| MultiFire Endo GIA™ 30 staplers . . . . .                                                             | .10        |
| MultiFire Endo GIA™ 30 reloads. . . . .                                                               | .10        |
| <b>MultiFire Endo TA™ 30 staplers and reloads . . . . .</b>                                           | <b>.11</b> |
| Endo TA™ 30 MultiFire stapling system . . . . .                                                       | .11        |
| Endo TA™ 30 MultiFire loading units. . . . .                                                          | .11        |

## Staplers for Open Surgery

|                                                              |            |
|--------------------------------------------------------------|------------|
| <b>Circular Staplers . . . . .</b>                           | <b>.12</b> |
| EEA™ haemorrhoid and prolapse stapler set . . . . .          | .12        |
| EEA™ staplers with DST Series™ technology . . . . .          | .12        |
| EEA™ XL staplers with DST Series™ technology . . . . .       | .13        |
| OrVil™ device . . . . .                                      | .13        |
| EEA™ introducer device with DST Series™ technology . . . . . | .14        |
| Purstring™ device . . . . .                                  | .14        |
| EEA™ sizers . . . . .                                        | .14        |
| <b>Staplers for Open Surgery . . . . .</b>                   | <b>.15</b> |
| GIA™ staplers with DST Series™ technology . . . . .          | .15        |
| GIA™ loading units . . . . .                                 | .15        |
| Poly GIA™ 75 stapler – 0.060 mm . . . . .                    | .15        |
| TA™ stapler with DST Series™ technology. . . . .             | .16        |
| TA™ loading units . . . . .                                  | .16        |
| Premium MultiFire TA™ stapler . . . . .                      | .17        |
| Premium MultiFire TA™ loading units . . . . .                | .17        |
| Reticulator™ staplers . . . . .                              | .17        |
| Reticulator 55 Poly™ staplers . . . . .                      | .17        |
| Premium Poly CS™ 57 staplers . . . . .                       | .18        |

## Table of Contents

### Staplers for Open Surgery

|                                                |           |
|------------------------------------------------|-----------|
| <b>Reusable Stapling Instruments</b> . . . . . | <b>19</b> |
| GIA 50 Premium™ stapler . . . . .              | 19        |
| GIA 50 Premium™ loading units. . . . .         | 19        |
| GIA 90 Premium™ loading units. . . . .         | 19        |
| ILA™ stapler. . . . .                          | 19        |
| ILA™ loading units. . . . .                    | 20        |
| TA Premium™ stapler . . . . .                  | 20        |
| TA Premium™ loading units . . . . .            | 20        |
| Premium Polysorb™ 55 loading units . . . . .   | 21        |
| TA™ 90 B loading units. . . . .                | 21        |
| <b>Skin Staplers</b> . . . . .                 | <b>22</b> |
| MultiFire Premium™ skin staplers . . . . .     | 22        |
| MultiFire Premium™ loading units . . . . .     | 22        |
| Appose™ ULC 35 staplers . . . . .              | 22        |
| Premium™ skin staple remover . . . . .         | 22        |
| DFS™ staplers . . . . .                        | 23        |



# Endo GIA™ reloads with Tri-Staple™ technology and Endo GIA™ Ultra universal staplers




## Endo GIA™ Ultra universal handle

45° of articulation, approved for 25 firings. Compatible with Tri-Staple™ technology and Endo GIA™ Ultra universal single use loading units (SULUs).

| Order Code | Description                                                      | Box Qty |
|------------|------------------------------------------------------------------|---------|
| EGIAUSHORT | Endo GIA™ Ultra universal short stapler<br>6 cm shaft length     | 3       |
| EGIAUSTND  | Endo GIA™ Ultra universal standard stapler<br>16 cm shaft length | 3       |
| EGIAUXL    | Endo GIA™ Ultra universal XL stapler<br>26 cm shaft length       | 3       |

## Endo GIA™ loading units

The reload design with Tri-Staple™ technology, intended to be used over a wider range of tissue thicknesses. They are compatible with the Endo GIA™ Ultra and the standard Endo GIA™ stapling system.



| Order Code | Color Code                                                                          | Description                                                                         | Box Qty |
|------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------|
| EGIA30AV   |  | Endo GIA™ 30 mm articulating vascular reload                                        | 6       |
| EGIA30AVM  |  | Endo GIA™ 30 mm articulating vascular/medium reloads<br>With Tri-Staple™ technology | 6       |
| EGIA30AMT  |  | Endo GIA™ 30 mm articulating medium/thick reloads<br>With Tri-Staple™ technology    | 6       |
| EGIA45AV   |  | Endo GIA™ 45 mm articulating vascular reloads*                                      | 6       |
| EGIA45AVM  |  | Endo GIA™ 45 mm articulating vascular/medium reloads<br>With Tri-Staple™ technology | 6       |
| EGIA45AMT  |  | Endo GIA™ 45 mm articulating medium/thick reloads<br>With Tri-Staple™ technology    | 6       |
| EGIA45AXT  |  | Endo GIA™ 45 mm articulating x-tra thick reloads<br>With Tri-Staple™ technology     | 6       |
| EGIA60AVM  |  | Endo GIA™ 60 mm articulating vascular medium reloads<br>With Tri-Staple™ technology | 6       |
| EGIA60AMT  |  | Endo GIA™ 60 mm articulating medium/thick reloads<br>With Tri-Staple™ technology    | 6       |

\*Code EGIA30AV and EGIA45AV do not come with Tri-Staple™ technology, have a slimmer fixed anvil

# Endo GIA™ reloads with Tri-Staple™ technology

## Endo GIA™ curved tip reload with Tri-Staple™ technology



| Order Code  | Color Code                                                                          | Description                                                                                 | Box Qty |
|-------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------|
| EGIA30CTAV  |    | Endo GIA™ curved tip 30mm reload articulating vascular reload*                              | 6       |
| EGIA30CTAVM |    | Endo GIA™ curved tip 30 mm articulating vascular/medium reloads With Tri-Staple™ technology | 6       |
| EGIA45CTAV  |    | Endo GIA™ curved tip 45 mm articulating vascular reload*                                    | 6       |
| EGIA45CTAVM |    | Endo GIA™ curved tip 45 mm articulating vascular/medium reloads With Tri-Staple™ technology | 6       |
| EGIA60CTAVM |    | Endo GIA™ curved tip 60 mm articulating vascular/medium reloads With Tri-Staple™ technology | 6       |
| EGIA60CTAMT |  | Endo GIA™ curved tip 60 mm articulating medium/thick reloads With Tri-Staple™ technology    | 6       |

\* Code EGIA30CTAV and EGIA45CTAV vascular reloads do not come with Tri-Staple™ technology, have a slimmer fixed anvil

## Endo GIA™ radial reload with Tri-Staple™ technology



| Order Code | Color Code                                                                          | Description                                                              | Box Qty |
|------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------|---------|
| EGIARADVM  |  | Endo GIA™ radial reload – medium tissue With Tri-Staple™ technology      | 6       |
| EGIARADMT  |  | Endo GIA™ radial reload – thick tissue With Tri-Staple™ technology       | 6       |
| EGIARADXT  |  | Endo GIA™ radial reload – extra-thick tissue With Tri-Staple™ technology | 6       |

# iDrive™ Ultra powered stapling system

## iDrive™ Ultra powered stapling system

The iDrive™ Ultra powered stapling system provides surgeons with improved precision,<sup>1</sup> enabling them to position and keep the stapler exactly where they intend, completely one-handed and at a push of a button.



| Order Code  | Description                              | Box Qty |
|-------------|------------------------------------------|---------|
| IDRVULTRA1  | iDrive™ Ultra powered handle             | 1       |
| EGIAADAPT   | Endo GIA™ adapter                        | 1       |
| EGIAADAPTXL | Endo GIA™ extra long adapter             | 1       |
| IDRVBIG     | iDrive™ battery insertion guide          | 1       |
| INTB100     | iDrive™ battery pack                     | 1       |
| INTBIC1NC   | iDrive™ battery charger and power supply | 1       |
| IDRVTRAY    | iDrive™ sterilization tray               | 1       |
| IDRVRET     | iDrive™ Ultra manual adapter tool        | 1       |

Average 61% reduction in reload tip travel during firing when compared to Ethicon Endo-Surgery Echelon Flex™ in indicated media, n=10 surgeons, 172 total trials, p<0.0005. Seils D, Tantawy T, Peterson D. Final report on results, University of Connecticut Biodynamics study. 2012; p.7, table 3. Internally funded study.

# Endo GIA™ universal stapling system



## Endo GIA™ universal stapling system

It accommodates straight and articulating loading units, approved for 25 firings.

| Order Code | Description                                                                                     | Box Qty |
|------------|-------------------------------------------------------------------------------------------------|---------|
| 030403     | GIA™ universal 12 mm<br>Single use instrument for open surgery, 6 cm shaft length               | 3       |
| 030449     | Endo GIA™ universal 12 mm<br>Single use instrument, 16 cm shaft length                          | 3       |
| EGIAUNIVXL | Endo GIA™ universal XL 12 mm<br>Single use instrument for bariatric surgery, 26 cm shaft length | 3       |

## Endo GIA™ universal Reticulator™ loading units

Single use loading units (SULUs) with titanium staples for use with the Endo GIA™ universal (030449), GIA™ universal (030403), and Endo GIA™ universal XL (EGIAUNIVXL) articulation up to 45°. Also compatible with the Endo GIA™ Ultra stapling system.



**NOT AVAILABLE any more  
USE BLACK LS page 5**

**NOT AVAILABLE any more  
USE BLACK LS page 5**

| Order Code | Color Code                                                                          | Description                                                           | Box Qty |
|------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------|---------|
| 030450     |  | Endo GIA™ universal Reticulator™ 30 – 2.0 mm                          | 6       |
| 030451     |  | Endo GIA™ universal Reticulator™ 30 – 2.5 mm                          | 6       |
| 030452     |  | Endo GIA™ universal Reticulator™ 30 – 3.5 mm                          | 6       |
| 030453     |  | Endo GIA™ universal Reticulator™ 45 – 4.8 mm<br>Requires 15 mm trocar | 6       |
| 030454     |  | Endo GIA™ universal Reticulator™ 45 – 2.5 mm                          | 6       |
| 030455     |  | Endo GIA™ universal Reticulator™ 45 – 3.5 mm                          | 6       |
| 030456     |  | Endo GIA™ universal Reticulator™ 45 – 4.8 mm<br>Requires 15 mm trocar | 6       |
| 030457     |  | Endo GIA™ universal Reticulator™ 60 – 2.5 mm                          | 6       |
| 030458     |  | Endo GIA™ universal Reticulator™ 60 – 3.5 mm                          | 6       |
| 030459     |  | Endo GIA™ universal Reticulator™ 60 – 4.8 mm<br>Requires 15 mm trocar | 6       |

# Endo GIA™ universal stapling system

## Endo GIA™ universal straight loading units

Single use loading units (SULUs) with titanium staples for use with Endo GIA™ universal (030449), GIA™ Universal (030403) and Endo GIA™ universal XL (EGIAUNIVXL) single use instruments. Also compatible with the Endo GIA™ Ultra stapling system.



**NOT AVAILABLE any more  
USE BLACK LS page 5**

**NOT AVAILABLE any more  
USE BLACK LS page 5**

| Order Code | Color Code                                                                          | Description                                                       | Box Qty |
|------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------|---------|
| 030416     |    | Endo GIA™ universal straight 30 – 2.0 mm                          | 6       |
| 030418     |    | Endo GIA™ universal straight 30 – 2.5 mm                          | 6       |
| 030419     |    | Endo GIA™ universal straight 30 – 3.5 mm                          | 6       |
| 030426     |    | Endo GIA™ universal straight 45 – 2.0 mm                          | 6       |
| 030425     |    | Endo GIA™ universal straight 45 – 2.5 mm                          | 6       |
| 030422     |    | Endo GIA™ universal straight 45 – 3.5 mm                          | 6       |
| 030423     |  | Endo GIA™ universal straight 45 – 4.8 mm<br>Requires 15 mm trocar | 6       |
| 030412     |  | Endo GIA™ universal straight 60 – 2.5 mm                          | 6       |
| 030414     |  | Endo GIA™ universal straight 60 – 3.5 mm                          | 6       |
| 030415     |  | Endo GIA™ universal straight 60 – 4.8 mm<br>Requires 15 mm trocar | 6       |

## Multifire Endo GIA™ 30 staplers and reloads

### Multifire Endo GIA™ 30 staplers

Single use staplers with titanium staples, approved for 8 firings.



| Order Code | Color Code                                                                        | Description                                      | Box Qty |
|------------|-----------------------------------------------------------------------------------|--------------------------------------------------|---------|
| 030811     |  | MultiFire Endo GIA™ 30 – 2.5 mm<br>12 mm stapler | 3       |
| 030813     |  | MultiFire Endo GIA™ 30 – 3.5 mm<br>12 mm stapler | 3       |



### Multifire Endo GIA™ 30 reloads

Single use loading units (SULU's) with titanium staples for use with Multifire Endo GIA™ 30 – 12mm stapler.

| Order Code | Color code                                                                          | Description                     | Box Qty |
|------------|-------------------------------------------------------------------------------------|---------------------------------|---------|
| 030805L    |  | MultiFire Endo GIA™ 30 – 2.5 mm | 6       |
| 030807L    |  | MultiFire Endo GIA™ 30 – 3.5 mm | 6       |

## MultiFire Endo TA™ 30 staplers and reloads



### Endo TA™ 30 multifire stapling system

Single use staplers with titanium staples, approved for 8 firings.

| Order Code | Color Code               | Description                                                | Box Qty |
|------------|--------------------------|------------------------------------------------------------|---------|
| 010901     | <input type="checkbox"/> | MultiFire Endo TA™ 30 – 2.5 mm<br>12 mm single use stapler | 3       |




### Endo TA™ 30 multiFire loading units

Single use loading units (SULU's) with titanium staples for use with MultiFire Endo TA™ 30 staplers.

| Order Code | Color Code               | Description                    | Box Qty |
|------------|--------------------------|--------------------------------|---------|
| 010911L    | <input type="checkbox"/> | MultiFire Endo TA™ 30 – 2.5 mm | 6       |

## Circular Staplers



### EEA™ haemorrhoid and prolapse stapler with DST Series™ technology

Shaft length 12 cm, 33 mm diameter, 32 DST Series™ titanium staples, 3.5 and 4.8 mm size. Shell volume 20 cc. Comes with a transparent anoscope, a detachable anvil and a conical dilator.

| Order Code | Color Code                                                                        | Description                                                                                 | Box Qty |
|------------|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------|
| HEM3335    |  | EEA™ haemorrhoid and prolapse stapler set 33 mm with DST Series™ technology. 3.5 mm staples | 3       |
| HEM3348    |  | EEA™ haemorrhoid and prolapse stapler set 33 mm with DST Series™ technology. 4.8 mm staples | 3       |

### EEA™ staplers with DST Series™ technology

Circular staplers, available in 22 and 35 cm shaft length and in 3.5 and 4.8 mm staple heights.

| Order Code | Color Code                                                                          | Description                                                                              | Box Qty |
|------------|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|---------|
| EEA21      |  | EEA™ 21mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples  | 3       |
| EEA2135    |  | EEA™ 21 mm stapler with DST Series™ technology<br>Single use stapler with 3.5 mm staples | 3       |
| EEA25      |  | EEA™ 25 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |
| EEA2535    |  | EEA™ 25 mm stapler with DST Series™ technology<br>Single use stapler with 3.5 mm staples | 3       |
| EEA28      |  | EEA™ 28 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |
| EEA2835    |  | EEA™ 28 mm stapler with DST Series™ technology<br>Single use stapler with 3.5 mm staples | 3       |
| EEA31      |  | EEA™ 31 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |
| EEA33      |  | EEA™ 33 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |



Directional  
stapling  
technology




# Circular Staplers

## EEA™ XL staplers with DST Series™ technology



| Order Code | Color Code                                                                          | Description                                                                                 | Box Qty |
|------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------|
| EEAXL21    |    | EEA™ XL 21 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |
| EEAXL2135  |    | EEA™ XL 21 mm stapler with DST Series™ technology<br>Single use stapler with 3.5 mm staples | 3       |
| EEAXL25    |    | EEA™ XL 25 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |
| EEAXL2535  |    | EEA™ XL 25 mm stapler with DST Series™ technology<br>Single use stapler with 3.5 mm staples | 3       |
| EEAXL28    |    | EEA™ XL 28 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |
| EEAXL2835  |    | EEA™ XL 28 mm stapler with DST Series™ technology<br>Single use stapler with 3.5 mm staples | 3       |
| EEAXL31    |  | EEA™ XL 31 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |
| EEAXL33    |  | EEA™ XL 33 mm stapler with DST Series™ technology<br>Single use stapler with 4.8 mm staples | 3       |



## OrVil™ device

Pre-tilted anvil trans-oral delivery device, mounted on a 90 cm long PVC nasogastric tube  
Attention: Only compatible with XL staplers 21 and 25.

| Order Code | Color Code                                                                          | Description                                                                                | Box Qty |
|------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|---------|
| EEAORVIL21 |  | EEA™ OrVil™ 21 mm device with DST Series™ technology<br>Compatible with XL stapler 21 only | 3       |
| EEAORVIL25 |  | EEA™ OrVil™ 21 mm device with DST Series™ technology<br>Compatible with XL stapler 25 only | 3       |

## Circular Staplers



### EEA™ introducer device with DST Series™ technology

Facilitate the smooth insertion of the EEA™ circular stapler during bariatric and colorectal procedures.

| Order Code | Color Code                                                                        | Description                                      | Box Qty |
|------------|-----------------------------------------------------------------------------------|--------------------------------------------------|---------|
| EEATAID21D |  | EEA™ 21 mm introducer                            | 6       |
| Endoscopic |  | EEA™ 25 mm introducer                            | 6       |
| EEATAID28D |                                                                                   | EEA28 disposable trans-anal/abdominal introducer | 6       |
| EEATAID31D |                                                                                   | EEA31 disposable trans-anal/abdominal introducer | 6       |
| EEATAID33D |                                                                                   | EEA33 disposable trans-anal/abdominal introducer | 6       |

### Purstring™ device

Single use instruments with stainless steel staples.



| Order Code | Description                                                                                                            | Box Qty |
|------------|------------------------------------------------------------------------------------------------------------------------|---------|
| 020242     | <b>Purstring™ 65</b><br>Single use instrument with Monosof™ 2 – 0 mm non absorbable monofilament nylon suture          | 3       |
| 020730     | <b>Purstring™ 45</b><br>Single use instrument with Surgidac™ 2 – 0 mm non absorbable braided polyester surgical suture | 3       |

### EEA™ sizers




| Order Code | Description                             | Box Qty |
|------------|-----------------------------------------|---------|
| 020250     | EEA™ reusable sizer set (25, 28, 31 mm) | 1       |

# Staplers for Open Surgery



Directional  
stapling  
technology



## GIA™ staplers with DST Series™ technology

Single use reloadable stapler with titanium staples, approved for 8 firings.

| Order Code | Color code | Description                 | Box Qty |
|------------|------------|-----------------------------|---------|
| GIA6025S   |            | GIA™ 60 mm – 2.5 mm         | 3       |
| GIA6038S   |            | GIA™ 60 mm – 3.8 mm         | 3       |
| GIA6048S   |            | GIA™ 60 mm – 4.8 mm         | 3       |
| GIA8038S   |            | GIA™ 80 mm – 3.8 mm         | 3       |
| GIA8048S   |            | GIA™ 80 mm – 4.8 mm         | 3       |
| GIA10038S  |            | GIA™ 100 mm – 3.8 mm        | 3       |
| GIA10048S  |            | GIA™ 100 mm – 4.8 mm        | 3       |
| SGIA6038S  |            | SGIA™ 60 – 3.8 mm Knifeless | 3       |

## GIA™ loading units

Single use loading units (SULUs) with titanium staples for use with GIA™ DST Series™ reloadable staplers.



| Order Code | Color code | Description       | Box Qty |
|------------|------------|-------------------|---------|
| GIA6025L   |            | GIA™ 60 – 2.5 mm  | 6       |
| GIA6038L   |            | GIA™ 60 – 3.8 mm  | 6       |
| GIA6048L   |            | GIA™ 60 – 4.8 mm  | 6       |
| GIA8038L   |            | GIA™ 80 – 3.8 mm  | 6       |
| GIA8048L   |            | GIA™ 80 – 4.8 mm  | 6       |
| GIA10038L  |            | GIA™ 100 – 3.8 mm | 6       |
| GIA10048L  |            | GIA™ 100 – 4.8 mm | 6       |

## Poly GIA™ 75 stapler – 0.060 mm

Single use stapler with Lactomer™ absorbable staples.



| Order Code | Color code | Description                                   | Box Qty |
|------------|------------|-----------------------------------------------|---------|
| 030775     |            | Poly GIA™ 75 – 0.060 mm<br>Closure gap 1.5 mm | 6       |

## Staplers for Open Surgery



Directional  
stapling  
technology



### TA™ stapler with DST Series™ technology

Single use reloadable stapler with titanium staples, approved for 8 firings.

| Order Code | Color Code                                                                          | Description               | Box Qty |
|------------|-------------------------------------------------------------------------------------|---------------------------|---------|
| TA30V3S    |    | TA™ 30 mm – V3 (Vascular) | 3       |
| TA3035S    |    | TA™ 30 – 3.5 mm           | 3       |
| TA3048S    |    | TA™ 30 – 4.8 mm           | 3       |
| TA4535S    |    | TA™ 45 – 3.5 mm           | 3       |
| TA4548S    |    | TA™ 45 – 4.8 mm           | 3       |
| TA6035S    |    | TA™ 60 – 3.5 mm           | 3       |
| TA6048S    |    | TA™ 60 – 4.8 mm           | 3       |
| TA9035S    |    | TA™ 90 – 3.5 mm           | 3       |
| TA9048S    |  | TA™ 90 – 4.8 mm           | 3       |

### TA™ loading units

Single use loading units (SULUs) with titanium staples for use with TA™ staplers with DST Series™ technology reloadable staplers.

| Order Code | Color Code                                                                          | Description               | Box Qty |
|------------|-------------------------------------------------------------------------------------|---------------------------|---------|
| TA30V3L    |  | TA™ 30 mm – V3 (Vascular) | 6       |
| TA3035L    |  | TA™ 30 – 3.5 mm           | 6       |
| TA3048L    |  | TA™ 30 – 4.8 mm           | 6       |
| TA4535L    |  | TA™ 45 – 3.5 mm           | 6       |
| TA4548L    |  | TA™ 45 – 4.8 mm           | 6       |
| TA6035L    |  | TA™ 60 – 3.5 mm           | 6       |
| TA6048L    |  | TA™ 60 – 4.8 mm           | 6       |
| TA9035L    |  | TA™ 90 – 3.5 mm           | 6       |
| TA9048L    |  | TA™ 90 – 4.8 mm           | 6       |



# Staplers for Open Surgery



## Premium MultiFire TA™ stapler

Single use reloadable stapler with titanium staples, approved for 8 firings.

| Order Code | Color Code                                                                        | Description                                                                  | Box Qty |
|------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------|---------|
| 010315     |  | Premium MultiFire TA™ 30 mm – V3 (Vascular)<br>Single use reloadable stapler | 3       |



## Premium MultiFire TA™ loading units

Single use loading units (SULUs) with titanium staples for use with Premium MultiFire TA™ single use reloadable stapler.

| Order Code | Color Code                                                                        | Description                                                            | Box Qty |
|------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------|---------|
| 010316L    |  | Premium MultiFire TA™ 30 mm – V3<br>(Vascular) single use loading unit | 6       |

## Rotulator™ staplers

Single use staplers with titanium staples.



| Order Code | Color Code                                                                          | Description                   | Box Qty |
|------------|-------------------------------------------------------------------------------------|-------------------------------|---------|
| 017619     |  | Rotulator™ 30 – V3 (Vascular) | 3       |
| 017615     |  | TA™ Rotulator™ 30 – 3.5 mm    | 3       |
| 017617     |  | TA™ Rotulator™ 30 – 4.8 mm    | 3       |
| 017612     |  | TA™ Rotulator™ 55 – 3.5 mm    | 3       |
| 017614     |  | TA™ Rotulator™ 55 – 4.8 mm    | 3       |

## Rotulator 55 Poly™ staplers

Single use staplers with absorbable Lactomer™ staples.



| Order Code | Color Code                                                                          | Description                   | Box Qty |
|------------|-------------------------------------------------------------------------------------|-------------------------------|---------|
| 013602     |  | Rotulator 55 Poly™ – 0.170 mm | 3       |
| 013601     |  | Rotulator 55 Poly™ – 0.200 mm | 3       |

## Staplers for Open Surgery



### Premium Poly™ CS 57 staplers

Single use staplers with absorbable Lactomer™ staples (2 instruments per kit).

| Order Code | Color Code                                                                        | Description                    | Box Qty |
|------------|-----------------------------------------------------------------------------------|--------------------------------|---------|
| 015140     |  | Premium Poly™ CS 57 – 0.140 mm | 3       |
| 015170     |  | Premium Poly™ CS 57 – 0.170 mm | 3       |

## EC CERTIFICATE

According to Annex II of the Directive 93/42/EEC on Medical Devices

### Full Quality Assurance System

Certificate Number: 2195-MED-2012001

**Manufacturer:**

**MERIL ENDO SURGERY PVT. LTD.**

Head Office Third Floor, E1-E3, Meril Park, Survey No. 135/2/B & 174/2,  
& Factory I: Muktanand Marg, Chala, Vapi – 396 191, Gujarat, INDIA  
Factory II: Plot. No. 688 / 10 & 11, Siddhivinayak Industrial Estate,  
Somnath Road, Daman – 396210, INDIA

**Product(s):**

**Sterile Single Use Endo Mechanical Products**

**Model(s):**

1. MIRUS™ Endoscopic Linear Cutter and Reloads
2. MIRUS™ Linear Stapler and Reloads
3. MIRUS™ Linear Cutter and Reloads
4. MIRUS™ Circular Stapler
5. MIRUS™ Hemorrhoids Stapler
6. MIRUS™ Disposable Skin Stapler
7. MONIK™ Disposable Endoscopic Trocar

**Reference Report No:** MM0755-P002-R01, MM0755-P002-R02

Szutest, Notified Body 2195, declares that the aforementioned manufacturer has implemented a quality assurance system according to Annex II (excluding section 4), Section 3 of the directive 93/42/EEC on medical devices. This quality assurance system covers those aspects of manufacturing concerned with securing and maintaining safe conditions of the respective product(s) and conforms to the provisions of this Directive. The approved quality system is subject to surveillance pursuant to Annex II, Section 5 of Directive 93/42/EEC and unannounced audits.

Szutest must be informed of any significant changes in the design and/or construction of the product(s). For class I devices with sterile conditions the quality management system evaluation is restricted to the aspects of manufacture concerned with securing and maintaining sterile conditions. For class I devices with measuring function the quality management system evaluation is restricted to the aspects of manufacture concerned with the conformity of the devices with metrological requirements.

*This EC certificate is valid till 2024-05-26.*

Issue Date:

2020-04-29



Rukiye BALKAN  
Deputy General Manager



# Certificate

No. Q5 105557 0001 Rev. 00

**Holder of Certificate:** **Meril Endo Surgery Pvt. Ltd.**  
Third floor, E1- E3, Meril Park  
Survey No 135/2/B & 174/2, Muktanand Marg, Chala  
396191 Vapi, Gujarat  
INDIA

**Facility(ies):** Meril Endo Surgery Pvt. Ltd.  
Third floor, E1- E3, Meril Park , Survey No 135/2/B & 174/2,  
Muktanand Marg, Chala, 396191 Vapi, Gujarat, INDIA

Meril Endo Surgery Pvt. Ltd  
Type A-2, Shed No. 11, Survey No. 725/P, Phase – I, GIDC,  
396195 Vapi, Gujarat, INDIA

Meril Endo Surgery Pvt. Ltd  
Plot No 688/10 & 11, Siddivinayak Industrial Estate, Somnath  
Road, 396210 Daman, INDIA

## Certification Mark:



**Scope of Certificate:** Design, Development, Production, Testing, Storage, Sales and Distribution of Sterile / Non-sterile Endo Surgical Products e.g., Bulk and Finished Surgical Sutures, Polytetrafluoroethylene Pledgets, Sutures with Polytetrafluoroethylene Pledgets, Contraceptive Devices, Surgical Meshes, Mesh Fixation Devices, Bone Wax, Umbilical Cotton Tape, Surgical Kits, Surgical Staplers and Staples, Tourniquet Devices, Disposable Endoscopic Trocars, Ligating Clips, Surgical Haemostats, Skin Adhesives, Surgical Adhesives, Sealants & Surgical Needles”

**Applied Standard(s):** EN ISO 13485:2016  
Medical devices - Quality management systems -  
Requirements for regulatory purposes  
(ISO 13485:2016)  
DIN EN ISO 13485:2016


The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

**Report No.:** IND2019074  
**Valid from:** 2020-03-30  
**Valid until:** 2023-03-29

**Date,** 2020-03-30

Christoph Dicks  
Head of Certification/Notified Body

# MIRUS™ Circular Stapler



II Rows for better and secure hemostasis




Stainless steel anvil shaft with longer purse tying area



Also available in III Rows for superior hemostasis

Staple compression gauge to ensure accurate staple formation that suits variable tissue thicknesses



Ergonomic anvil design

Safety lock

Wing nut to assist control and dexterity while opening and closing anvil

| Rows of Staples | Stapler Code | Close Leg Length | Head Diameter | Staple Quantity |
|-----------------|--------------|------------------|---------------|-----------------|
| Two Rows        | MCS-19       | 1.0 - 2.5 mm     | 19 mm         | 16              |
|                 | MCS-21       | 1.0 - 2.5 mm     | 21 mm         | 16              |
|                 | MCS-24       | 1.0 - 2.5 mm     | 24 mm         | 18              |
|                 | MCS-25       | 1.0 - 2.5 mm     | 25 mm         | 20              |
|                 | MCS-26       | 1.0 - 2.5 mm     | 26 mm         | 20              |
|                 | MCS-29       | 1.0 - 2.5 mm     | 29 mm         | 24              |
|                 | MCS-32       | 1.0 - 2.5 mm     | 32 mm         | 26              |
| Three Rows      | MCS-21R3     | 1.0 - 2.5 mm     | 21 mm         | 21              |
|                 | MCS-24R3     | 1.0 - 2.5 mm     | 24 mm         | 24              |
|                 | MCS-25R3     | 1.0 - 2.5 mm     | 25 mm         | 27              |
|                 | MCS-26R3     | 1.0 - 2.5 mm     | 26 mm         | 33              |
|                 | MCS-29R3     | 1.0 - 2.5 mm     | 29 mm         | 36              |
|                 | MCS-31R3     | 1.0 - 2.5 mm     | 31 mm         | 36              |
|                 | MCS-32R3     | 1.0 - 2.5 mm     | 32 mm         | 42              |

# In anorectal surgeries, why compromise on hemostasis & security?

## OPTIMIZED STAPLE FORMATION TECHNIQUE

- Depending upon the thickness of prolapsed tissue mass, surgeon can control the level of compression which needs to be applied with Mirus Hemorrhoids Stapler, as it delivers **Adjustable Height Staples**
- The tissue compression indicator adjusts accordingly so that, when fired, Mirus Hemorrhoids Stapler delivers right compression, proper staple line with optimum closed height and superior hemostasis

## Benefits of MIPH

- Less pain
  - Faster regaining of normal bowel function
  - Shorter hospital stay
  - Less post operative care after discharge
- Shorter wound healing time
  - Shorter operating time
  - Enhanced patient satisfaction
  - Faster return to normal work

## Product Specification Chart

| Meril Code | External Diameter (mm) | No. of Rows | Staple Quantity | Closed Staple Height (mm) | Housing capacity (CC) |
|------------|------------------------|-------------|-----------------|---------------------------|-----------------------|
| Two Rows   |                        |             |                 |                           |                       |
| MPPH32     | 32                     | 2           | 32              | 0.75 -1.6                 | 14                    |
| MPPH34     | 34                     | 2           | 32              | 0.75 - 1.6                | 14                    |
| Three Rows |                        |             |                 |                           |                       |
| MPPH34R3   | 34                     | 3           | 48              | 0.75 - 1.6                | 14                    |
| MPPH36R3   | 36                     | 3           | 48              | 0.75 - 1.6                | 14                    |



MERIL ENDO-SURGERY PVT. LTD.  
Corporate Headquarters, Vapi  
E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg,  
Chala, Vapi- 396 191, Gujarat, India.  
Tel.: +91 260 3052 100, Toll Free No.: 1800-4194-433

India Sales and Marketing, Mumbai  
601-612, Midas, Sahar Plaza Complex, M.V. Road, J.B Nagar,  
Andheri (East), Mumbai 400 059, India. Tel.: +91 22 307 32 400  
Email: endosurgery@merillife.com, Web: www.merillife.com



# Revolutionary 3 Rows Approach for MIPH & STARR


## MIRUS™ Hemorrhoids Stapler




## Mirus Hemorrhoids Stapler with Three Circular Staggered Staple Lines (3 Rows)

15,000+  
successful procedures registered and counting....





- Transparent Purse String Suture Anoscope with the Markings



- Circular Anal Dilator (CAD)




- Obturator



- Suture Threader




- Filaprop for Purse String Suturing (Polypropylene)




# MIRUS™ Linear Cutter




Cuts and separates tissue between two double staggered rows



New knife with every reload



Every reload has safety lock to prevent firing over empty reload




| Stapler Code | Reload Code               | Number of Firings | Number of Staples | Staple Line Length | Cut Length |
|--------------|---------------------------|-------------------|-------------------|--------------------|------------|
| MLC60        | MLCR60-3.8   MLCR60-4.8   | 8                 | 64                | 64 mm              | 60 mm      |
| MLC80        | MLCR80-3.8   MLCR80-4.8   | 8                 | 84                | 84 mm              | 80 mm      |
| MLC100       | MLCR100-3.8   MLCR100-4.8 | 8                 | 104               | 104 mm             | 100 mm     |


# MIRUS™ Linear Stapler(Auto)



One-hand operation to close, fire and release



Long slim shaft for excellent access & visibility



| Stapler Code | Reload Code             | Number of Firings | Number of Staples | Rows of Staples |
|--------------|-------------------------|-------------------|-------------------|-----------------|
| MALS30       | MLSR30-3.5   MLSR30-4.8 | 8                 | 11                | 2               |
| MALS45       | MLSR45-3.5   MLSR45-4.8 | 8                 | 15                | 2               |
| MALS60       | MLSR60-3.5   MLSR60-4.8 | 8                 | 21                | 2               |
| MALS90       | MLSR90-3.5   MLSR90-4.8 | 8                 | 33                | 2               |

# MIRUS™ Skin Stapler

## Outstanding performance across wider skin wounds

Skin Staplers are typically used at the end of a surgical procedure for routine wound closure. The staples are then removed with staple extractor after the skin has had an opportunity to heal.

The MIRUS™ Skin Stapler's handle design is characterized by good texture for an user friendly experience.


The angled head provides a clear view to ensure precise staple placement and allows the pins to staple the tissue easily.

- Staple material : 316L Stainless Steel
- Sterilization : Ethylene Oxide
- Packaging : 8 pcs/Box
- Code : MSS35 & MSSP35
- Size : 35
- Staples : Regular

Angled head contains 35 staples



## Staple Dimension:



MSSP35

Angled head contains 35 staples  
Rectangle staple formation



**Larger diameter wire, greater width, and extra staple height leads to a secure wound closure**

|                       |               |
|-----------------------|---------------|
| Stapler               | Skin Stapler  |
| Procedure Application | Wound closure |