



DECLARATION OF CONFORMITY

MANUFACTURER NAME : ZİMED MEDİKAL SANAYİ VE TİCARET LİMİTED ŞİRKETİ

ADRESS : Aydınlar Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-
Türkiye

DEVICE NAME : UPPER-LOWER EXTREMITY BONE PLATES AND SCREWS

NO	PRODUCT NO	PRODUCT NAME	RAW MATERIAL	REF NO	GMDN	STERILIZATION
1.	8680941231903	LOCKING ANTEROLATERAL DISTAL TIBIA PLATE 3.5 / 4.0 mm , 5 HOLES LEFT	ISO 5832-1	864.005	46647	NON-STERILE
2.	8680941231910	LOCKING ANTEROLATERAL DISTAL TIBIA PLATE 3.5 / 4.0 mm , 7 HOLES LEFT	ISO 5832-1	864.007	46647	NON-STERILE
3.	8680941231927	LOCKING ANTEROLATERAL DISTAL TIBIA PLATE 3.5 / 4.0 mm , 9 HOLES LEFT	ISO 5832-1	864.009	46647	NON-STERILE
4.	8680941231934	LOCKING ANTEROLATERAL DISTAL TIBIA PLATE 3.5 / 4.0 mm , 11 HOLES LEFT	ISO 5832-1	864.011	46647	NON-STERILE
5.	8680941231941	LOCKING ANTEROLATERAL DISTAL TIBIA PLATE 3.5 / 4.0 mm , 5 HOLES RIGHT	ISO 5832-1	864.105	46647	NON-STERILE
6.	8680941231958	LOCKING ANTEROLATERAL DISTAL TIBIA PLATE 3.5 / 4.0 mm , 7 HOLES RIGHT	ISO 5832-1	864.107	46647	NON-STERILE
7.	8680941231965	LOCKING ANTEROLATERAL DISTAL TIBIA PLATE 3.5 / 4.0 mm , 9 HOLES RIGHT	ISO 5832-1	864.109	46647	NON-STERILE
8.	8680941231972	LOCKING ANTEROLATERAL DISTAL TIBIA PLATE 3.5 / 4.0 mm , 11 HOLES RIGHT	ISO 5832-1	864.111	46647	NON-STERILE
9.	8680941231989	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 4 HOLES	ISO 5832-1	157.204	46647	NON-STERILE
10.	8680941231996	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 5 HOLES	ISO 5832-1	157.205	46647	NON-STERILE
11.	8680941232009	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 6 HOLES	ISO 5832-1	157.206	46647	NON-STERILE
12.	8680941232016	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 7 HOLES	ISO 5832-1	157.207	46647	NON-STERILE
13.	8680941232023	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 8 HOLES	ISO 5832-1	157.208	46647	NON-STERILE
14.	8680941232030	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 9 HOLES	ISO 5832-1	157.209	46647	NON-STERILE
15.	8680941232047	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 10 HOLES	ISO 5832-1	157.210	46647	NON-STERILE
16.	8680941232054	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 12 HOLES	ISO 5832-1	157.212	46647	NON-STERILE
17.	8680941232061	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 14 HOLES	ISO 5832-1	157.214	46647	NON-STERILE
18.	8680941232078	NARROW LOCKING COMPRESSION PLATE 4.5 / 5.0 mm , 16 HOLES	ISO 5832-1	157.216	46647	NON-STERILE
19.	8680941231996	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 6 HOLES	ISO 5832-1	158.206	46647	NON-STERILE
20.	8699131670903	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 7 HOLES	ISO 5832-1	158.207	46647	NON-STERILE

21.	8699131670941	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 8 HOLES	ISO 5832-1	158.208	46647	NON-STERILE
22.	8699131670958	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 9 HOLES	ISO 5832-1	158.209	46647	NON-STERILE
23.	8699131600641	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 10 HOLES	ISO 5832-1	158.210	46647	NON-STERILE
24.	8699131600658	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 12 HOLES	ISO 5832-1	158.212	46647	NON-STERILE
25.	8699131600665	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 14 HOLES	ISO 5832-1	158.214	46647	NON-STERILE
26.	8699131600672	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 16 HOLES	ISO 5832-1	158.216	46647	NON-STERILE
27.	8699131600696	BROAD LOCKING COMPRESSION PLATE 4.5 / 5.0 mm 18 HOLES	ISO 5832-1	158.218	46647	NON-STERILE
28.	8699131600702	CORTICAL SCREW 4.5 mm Self Tapping 14 mm	ISO 5832-1	106.214	56642	NON-STERILE
29.	8699131600719	CORTICAL SCREW 4.5 mm Self Tapping 16 mm	ISO 5832-1	106.216	56642	NON-STERILE
30.	8699131600689	CORTICAL SCREW 4.5 mm Self Tapping 18 mm	ISO 5832-1	106.218	56642	NON-STERILE
31.	8699131600597	CORTICAL SCREW 4.5 mm Self Tapping 20 mm	ISO 5832-1	106.220	56642	NON-STERILE
32.	8699131600566	CORTICAL SCREW 4.5 mm Self Tapping 22 mm	ISO 5832-1	106.222	56642	NON-STERILE
33.	8680941219543	CORTICAL SCREW 4.5 mm Self Tapping 24 mm	ISO 5832-1	106.224	56642	NON-STERILE
34.	8680941219550	CORTICAL SCREW 4.5 mm Self Tapping 26 mm	ISO 5832-1	106.226	56642	NON-STERILE
35.	8699131629253	CORTICAL SCREW 4.5 mm Self Tapping 28 mm	ISO 5832-1	106.228	56642	NON-STERILE
36.	8680941227548	CORTICAL SCREW 4.5 mm Self Tapping 30 mm	ISO 5832-1	106.230	56642	NON-STERILE
37.	8680941227555	CORTICAL SCREW 4.5 mm Self Tapping 32 mm	ISO 5832-1	106.232	56642	NON-STERILE
38.	8680941227562	CORTICAL SCREW 4.5 mm Self Tapping 34 mm	ISO 5832-1	106.234	56642	NON-STERILE
39.	8680941293413	CORTICAL SCREW 4.5 mm Self Tapping 36 mm	ISO 5832-1	106.236	56642	NON-STERILE
40.	8680941293406	CORTICAL SCREW 4.5 mm Self Tapping 38 mm	ISO 5832-1	106.238	56642	NON-STERILE
41.	8680941293390	CORTICAL SCREW 4.5 mm Self Tapping 40 mm	ISO 5832-1	106.240	56642	NON-STERILE
42.	8680941293383	CORTICAL SCREW 4.5 mm Self Tapping 42 mm	ISO 5832-1	106.242	56642	NON-STERILE
43.	8680941293352	CORTICAL SCREW 4.5 mm Self Tapping 44 mm	ISO 5832-1	106.244	56642	NON-STERILE
44.	8680941293369	CORTICAL SCREW 4.5 mm Self Tapping 46 mm	ISO 5832-1	106.246	56642	NON-STERILE
45.	8680941293376	CORTICAL SCREW 4.5 mm Self Tapping 48 mm	ISO 5832-1	106.248	56642	NON-STERILE
46.	8699131674093	CORTICAL SCREW 4.5 mm Self Tapping 50 mm	ISO 5832-1	106.250	56642	NON-STERILE
47.	8699131674109	CORTICAL SCREW 4.5 mm Self Tapping 52 mm	ISO 5832-1	106.252	56642	NON-STERILE
48.	8699131674116	CORTICAL SCREW 4.5 mm Self Tapping 54 mm	ISO 5832-1	106.254	56642	NON-STERILE
49.	8699131699591	CORTICAL SCREW 4.5 mm Self Tapping 56 mm	ISO 5832-1	106.256	56642	NON-STERILE
50.	8699131674130	CORTICAL SCREW 4.5 mm Self Tapping 58 mm	ISO 5832-1	106.258	56642	NON-STERILE
51.	8699131674161	CORTICAL SCREW 4.5 mm Self Tapping 60 mm	ISO 5832-1	106.260	56642	NON-STERILE
52.	8699131674178	CORTICAL SCREW 4.5 mm Self Tapping 62 mm	ISO 5832-1	106.262	56642	NON-STERILE
53.	8699131674185	CORTICAL SCREW 4.5 mm Self Tapping 64 mm	ISO 5832-1	106.264	56642	NON-STERILE
54.	8699131674192	CORTICAL SCREW 4.5 mm Self Tapping 66 mm	ISO 5832-1	106.266	56642	NON-STERILE
55.	8699131674208	CORTICAL SCREW 4.5 mm Self Tapping 68 mm	ISO 5832-1	106.268	56642	NON-STERILE
56.	8699131661802	CORTICAL SCREW 4.5 mm Self Tapping 70 mm	ISO 5832-1	106.270	56642	NON-STERILE
57.	8699131661833	CORTICAL SCREW 4.5 mm Self Tapping 72 mm	ISO 5832-1	106.272	56642	NON-STERILE
58.	8699131661888	CORTICAL SCREW 4.5 mm Self Tapping 74 mm	ISO 5832-1	106.274	56642	NON-STERILE
59.	8699131661956	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 25 mm	ISO 5832-1	110.725	56642	NON-STERILE
60.	8699131661994	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 30 mm	ISO 5832-1	110.730	56642	NON-STERILE
61.	8699131662021	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 35 mm	ISO 5832-1	110.735	56642	NON-STERILE
62.	8699131661819	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 40 mm	ISO 5832-1	110.740	56642	NON-STERILE
63.	8699131661840	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 45 mm	ISO 5832-1	110.745	56642	NON-STERILE
64.	8699131661918	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 50 mm	ISO 5832-1	110.750	56642	NON-STERILE

65.	8699131661963	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 55 mm	ISO 5832-1	110.755	56642	NON-STERILE
66.	8699131662007	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 60 mm	ISO 5832-1	110.760	56642	NON-STERILE
67.	8699131662038	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 65 mm	ISO 5832-1	110.765	56642	NON-STERILE
68.	8699131661826	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 70 mm	ISO 5832-1	110.770	56642	NON-STERILE
69.	8699131661864	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 75 mm	ISO 5832-1	110.775	56642	NON-STERILE
70.	8699131661932	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 80 mm	ISO 5832-1	110.780	56642	NON-STERILE
71.	8699131661987	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 85 mm	ISO 5832-1	110.785	56642	NON-STERILE
72.	8699131662014	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 90 mm	ISO 5832-1	110.790	56642	NON-STERILE
73.	8699131662045	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 95 mm	ISO 5832-1	110.795	56642	NON-STERILE
74.	8699131601334	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 100 mm	ISO 5832-1	110.800	56642	NON-STERILE
75.	8699131607039	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 105 mm	ISO 5832-1	110.805	56642	NON-STERILE
76.	8699131625743	5.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD – SELF TAPPING , 110 mm	ISO 5832-1	110.810	56642	NON-STERILE
77.	8699131627587	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 10 mm	ISO 5832-1	102.310	56642	NON-STERILE
78.	8699131601372	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 12 mm	ISO 5832-1	102.312	56642	NON-STERILE
79.	8699131627709	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 14 mm	ISO 5832-1	102.314	56642	NON-STERILE
80.	8699131656501	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 16 mm	ISO 5832-1	102.316	56642	NON-STERILE
81.	8699131601419	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 18 mm	ISO 5832-1	102.318	56642	NON-STERILE
82.	8699131601433	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 20 mm	ISO 5832-1	102.320	56642	NON-STERILE
83.	8699131601457	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 22 mm	ISO 5832-1	102.322	56642	NON-STERILE
84.	8699131601471	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 24 mm	ISO 5832-1	102.324	56642	NON-STERILE
85.	8699131601495	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 26 mm	ISO 5832-1	102.326	56642	NON-STERILE
86.	8699131603116	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 28 mm	ISO 5832-1	102.328	56642	NON-STERILE
87.	8699131603123	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 30 mm	ISO 5832-1	102.330	56642	NON-STERILE
88.	8699131603130	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 32 mm	ISO 5832-1	102.332	56642	NON-STERILE
89.	8699131603147	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 34 mm	ISO 5832-1	102.334	56642	NON-STERILE
90.	8699131603154	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 36 mm	ISO 5832-1	102.336	56642	NON-STERILE
91.	8699131603161	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 38 mm	ISO 5832-1	102.338	56642	NON-STERILE
92.	8699131661543	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 40 mm	ISO 5832-1	102.340	56642	NON-STERILE
93.	8699131661581	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 42 mm	ISO 5832-1	102.342	56642	NON-STERILE
94.	8699131661642	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 44 mm	ISO 5832-1	102.344	56642	NON-STERILE
95.	8699131661703	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 46 mm	ISO 5832-1	102.346	56642	NON-STERILE
96.	8699131661727	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 48 mm	ISO 5832-1	102.348	56642	NON-STERILE
97.	8699131661758	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 50 mm	ISO 5832-1	102.350	56642	NON-STERILE
98.	8699131661550	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 55 mm	ISO 5832-1	102.355	56642	NON-STERILE
99.	8699131661611	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 60 mm	ISO 5832-1	102.360	56642	NON-STERILE
100.	8699131661673	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 65 mm	ISO 5832-1	102.365	56642	NON-STERILE
101.	8699131661710	3.5 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 70 mm	ISO 5832-1	102.370	56642	NON-STERILE
102.	8699131661789	ONE THIRD TUBULER PLATE 3.5 mm 2 HOLES	ISO 5832-1	138.102	46647	NON-STERILE

103.	8699131600603	ONE THIRD TUBULER PLATE 3.5 mm 3 HOLES	ISO 5832-1	138.103	46647	NON-STERILE
104.	8699131600832	ONE THIRD TUBULER PLATE 3.5 mm 4 HOLES	ISO 5832-1	138.104	46647	NON-STERILE
105.	8699131604199	ONE THIRD TUBULER PLATE 3.5 mm 5 HOLES	ISO 5832-1	138.105	46647	NON-STERILE
106.	8699131604205	ONE THIRD TUBULER PLATE 3.5 mm 6 HOLES	ISO 5832-1	138.106	46647	NON-STERILE
107.	8699131604212	ONE THIRD TUBULER PLATE 3.5 mm 7 HOLES	ISO 5832-1	138.107	46647	NON-STERILE
108.	8699131604229	ONE THIRD TUBULER PLATE 3.5 mm 8 HOLES	ISO 5832-1	18.108	46647	NON-STERILE
109.	8699131604236	ONE THIRD TUBULER PLATE 3.5 mm 9 HOLES	ISO 5832-1	138.109	46647	NON-STERILE
110.	8699131604243	ONE THIRD TUBULER PLATE 3.5 mm 10 HOLES	ISO 5832-1	138.110	46647	NON-STERILE
111.	8699131604250	ONE THIRD TUBULER PLATE 3.5 mm 12 HOLES	ISO 5832-1	138.112	46647	NON-STERILE
112.	8699131604267	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 5 HOLES	ISO 5832-1	144.105	46647	NON-STERILE
113.	8699131604274	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 6 HOLES	ISO 5832-1	144.106	46647	NON-STERILE
114.	8699131604281	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 7 HOLES	ISO 5832-1	144.107	46647	NON-STERILE
115.	8699131604298	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 8 HOLES	ISO 5832-1	144.108	46647	NON-STERILE
116.	8699131604304	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 9 HOLES	ISO 5832-1	144.109	46647	NON-STERILE
117.	8699131604311	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 10 HOLES	ISO 5832-1	144.110	46647	NON-STERILE
118.	8699131604328	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 12 HOLES	ISO 5832-1	144.112	46647	NON-STERILE
119.	8699131604335	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 14 HOLES	ISO 5832-1	144.114	46647	NON-STERILE
120.	8699131604342	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 16 HOLES	ISO 5832-1	144.116	46647	NON-STERILE
121.	8699131604359	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 18 HOLES	ISO 5832-1	144.118	46647	NON-STERILE
122.	8699131604366	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 20 HOLES	ISO 5832-1	144.120	46647	NON-STERILE
123.	8699131604373	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 22 HOLES	ISO 5832-1	144.122	46647	NON-STERILE
124.	8699131604380	STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 mm 24 HOLES	ISO 5832-1	144.124	46647	NON-STERILE
125.	8699131604397	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 10 mm	ISO 5832-1	102.210	56642	NON-STERILE
126.	8699131604403	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 12 mm	ISO 5832-1	102.212	56642	NON-STERILE
127.	8699131604410	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 14 mm	ISO 5832-1	102.214	56642	NON-STERILE
128.	8699131604427	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 16 mm	ISO 5832-1	102.216	56642	NON-STERILE
129.	8699131604434	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 18 mm	ISO 5832-1	102.218	56642	NON-STERILE
130.	8699131604441	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 20 mm	ISO 5832-1	102.220	56642	NON-STERILE
131.	8699131604458	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 22 mm	ISO 5832-1	102.222	56642	NON-STERILE
132.	8699131604465	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 24 mm	ISO 5832-1	102.224	56642	NON-STERILE
133.	8699131604472	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 26 mm	ISO 5832-1	102.226	56642	NON-STERILE
134.	8699131604489	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 28 mm	ISO 5832-1	102.228	56642	NON-STERILE
135.	8699131604496	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 30 mm	ISO 5832-1	102.230	56642	NON-STERILE
136.	8699131604502	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 32 mm	ISO 5832-1	102.232	56642	NON-STERILE
137.	8699131604519	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 34 mm	ISO 5832-1	102.234	56642	NON-STERILE
138.	8699131604526	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 36 mm	ISO 5832-1	102.236	56642	NON-STERILE
139.	8699131604533	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 38 mm	ISO 5832-1	102.238	56642	NON-STERILE
140.	8699131604540	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 40 mm	ISO 5832-1	102.240	56642	NON-STERILE
141.	8699131604557	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 42 mm	ISO 5832-1	102.242	56642	NON-STERILE

142.	8699131604564	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 44 mm	ISO 5832-1	102.244	56642	NON-STERILE
143.	8699131604571	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 46 mm	ISO 5832-1	102.246	56642	NON-STERILE
144.	8699131604588	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 48 mm	ISO 5832-1	102.248	56642	NON-STERILE
145.	8699131604595	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 50 mm	ISO 5832-1	102.250	56642	NON-STERILE
146.	8699131604601	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 55 mm	ISO 5832-1	102.255	56642	NON-STERILE
147.	8699131604618	3.5 mm CORTICAL SCREW – SELF TAPPING (HEXDRIVE) 60 mm	ISO 5832-1	102.260	56642	NON-STERILE
148.	8699131604625	CORTICAL SCREW 2.7 X 10 MM	ISO 5832-1	2212-2710	56642	NON-STERILE
149.	8699131604632	CORTICAL SCREW 2.7 X 12 MM	ISO 5832-1	2212-2712	56642	NON-STERILE
150.	8699131604649	CORTICAL SCREW 2.7 X 14 MM	ISO 5832-1	2212-2714	56642	NON-STERILE
151.	8699131604656	CORTICAL SCREW 2.7 X 16 MM	ISO 5832-1	2212-2716	56642	NON-STERILE
152.	8699131604663	CORTICAL SCREW 2.7 X 18 MM	ISO 5832-1	2212-2718	56642	NON-STERILE
153.	8699131655948	CORTICAL SCREW 2.7 X 20 MM	ISO 5832-1	2212-2720	56642	NON-STERILE
154.	8699131655986	CORTICAL SCREW 2.7 X 22 MM	ISO 5832-1	2212-2722	56642	NON-STERILE
155.	8699131656037	CORTICAL SCREW 2.7 X 24 MM	ISO 5832-1	2212-2724	56642	NON-STERILE
156.	8699131663103	CORTICAL SCREW 2.7 X 26 MM	ISO 5832-1	2212-2726	56642	NON-STERILE
157.	8699131648162	CORTICAL SCREW 2.7 X 28 MM	ISO 5832-1	2212-2728	56642	NON-STERILE
158.	8699131648285	CORTICAL SCREW 2.7 X 30 MM	ISO 5832-1	2212-2730	56642	NON-STERILE
159.	8699131655955	CORTICAL SCREW 2.7 X 32 MM	ISO 5832-1	2212-2732	56642	NON-STERILE
160.	8699131655993	CORTICAL SCREW 2.7 X 34 MM	ISO 5832-1	2212-2734	56642	NON-STERILE
161.	8699131656020	CORTICAL SCREW 2.7 X 36 MM	ISO 5832-1	2212-2736	56642	NON-STERILE
162.	8699131663110	CORTICAL SCREW 2.7 X 38 MM	ISO 5832-1	2212-2738	56642	NON-STERILE
163.	8699131662878	CORTICAL SCREW 2.7 X 40 MM	ISO 5832-1	2212-2740	56642	NON-STERILE
164.	8699131662830	CORTICAL SCREW 2.7 X 42 MM	ISO 5832-1	2212-2742	56642	NON-STERILE
165.	8680941239060	CORTICAL SCREW 2.7 X 44 MM	ISO 5832-1	2212-2744	56642	NON-STERILE
166.	8699131604861	CORTICAL SCREW 2.7 X 46 MM	ISO 5832-1	2212-2746	56642	NON-STERILE
167.	8699131604878	CORTICAL SCREW 2.7 X 48 MM	ISO 5832-1	2212-2748	56642	NON-STERILE
168.	8699131604885	CORTICAL SCREW 2.7 X 50 MM	ISO 5832-1	2212-2750	56642	NON-STERILE
169.	8699131604892	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 12 mm	ISO 5832-1	105.512	56642	NON-STERILE
170.	8699131604908	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 14 mm	ISO 5832-1	105.514	56642	NON-STERILE
171.	8699131604915	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 16 mm	ISO 5832-1	105.516	56642	NON-STERILE
172.	8699131604922	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 18 mm	ISO 5832-1	105.518	56642	NON-STERILE
173.	8680941239077	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 20 mm	ISO 5832-1	105.520	56642	NON-STERILE
174.	8699131604946	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 22 mm	ISO 5832-1	105.522	56642	NON-STERILE
175.	8699131604960	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 24 mm	ISO 5832-1	105.524	56642	NON-STERILE
176.	8699131604984	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 26 mm	ISO 5832-1	105.526	56642	NON-STERILE
177.	8699131605004	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 28 mm	ISO 5832-1	105.528	56642	NON-STERILE
178.	8699131605028	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 30 mm	ISO 5832-1	105.530	56642	NON-STERILE

179.	8699131605042	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 32 mm	ISO 5832-1	105.532	56642	NON-STERILE
180.	8699131605059	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 34 mm	ISO 5832-1	105.534	56642	NON-STERILE
181.	8699131605066	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 36 mm	ISO 5832-1	105.536	56642	NON-STERILE
182.	8699131605073	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 38 mm	ISO 5832-1	105.538	56642	NON-STERILE
183.	8699131605080	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 40 mm	ISO 5832-1	105.540	56642	NON-STERILE
184.	8699131605097	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 42 mm	ISO 5832-1	105.542	56642	NON-STERILE
185.	8699131605103	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 44 mm	ISO 5832-1	105.544	56642	NON-STERILE
186.	8699131605110	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 46 mm	ISO 5832-1	105.546	56642	NON-STERILE
187.	8699131605127	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 48 mm	ISO 5832-1	105.548	56642	NON-STERILE
188.	8699131648384	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 50 mm	ISO 5832-1	105.550	56642	NON-STERILE
189.	8699131648391	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 55 mm	ISO 5832-1	105.555	56642	NON-STERILE
190.	8699131633441	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 60 mm	ISO 5832-1	105.560	56642	NON-STERILE
191.	8699131633458	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 65 mm	ISO 5832-1	105.565	56642	NON-STERILE
192.	8699131665671	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 70 mm	ISO 5832-1	105.570	56642	NON-STERILE
193.	8699131633465	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 75 mm	ISO 5832-1	105.575	56642	NON-STERILE
194.	8699131623404	4.00 mm LOCKING HEAD CANCELLOUS SCREW FULL THREAD SELF TAPPING (HEXDRIVE) 80 mm	ISO 5832-1	105.580	56642	NON-STERILE
195.	8699131623411	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 40	ISO 5832-1	2952-0040	56642	NON-STERILE
196.	8699131637920	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 45	ISO 5832-1	2952-0045	56642	NON-STERILE
197.	8699131633472	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 50	ISO 5832-1	2952-0050	56642	NON-STERILE
198.	8699131665756	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 55	ISO 5832-1	2952-0055	56642	NON-STERILE
199.	8699131633489	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 60	ISO 5832-1	2952-0060	56642	NON-STERILE
200.	8699131665800	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 65	ISO 5832-1	2952-0065	56642	NON-STERILE
201.	8699131605288	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 70	ISO 5832-1	2952-0070	56642	NON-STERILE
202.	8699131623428	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 75	ISO 5832-1	2952-0075	56642	NON-STERILE
203.	8699131623442	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 80	ISO 5832-1	2952-0080	56642	NON-STERILE
204.	8699131637944	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 85	ISO 5832-1	2952-0085	56642	NON-STERILE
205.	8699131633724	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 90	ISO 5832-1	2952-0090	56642	NON-STERILE
206.	8699131633700	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 95	ISO 5832-1	2952-0095	56642	NON-STERILE
207.	8699131633687	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 100	ISO 5832-1	2952-0100	56642	NON-STERILE
208.	8699131633670	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 105	ISO 5832-1	2952-0105	56642	NON-STERILE
209.	8699131633649	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 110	ISO 5832-1	2952-0110	56642	NON-STERILE
210.	8699131633625	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 115	ISO 5832-1	2952-0115	56642	NON-STERILE

211.	8699131616673	CANCELLOUS SCREW 32 MM THREAD. 6.5 X 120	ISO 5832-1	2952-0120	56642	NON-STERILE
212.	8680941228804	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 4 HOLES	ISO 5832-1	140.204	46647	NON-STERILE
213.	8699131617731	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 5 HOLES	ISO 5832-1	140.205	46647	NON-STERILE
214.	8680941228811	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 6 HOLES	ISO 5832-1	140.206	46647	NON-STERILE
215.	8699131617748	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 7 HOLES	ISO 5832-1	140.207	46647	NON-STERILE
216.	8680941292591	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 8 HOLES	ISO 5832-1	140.208	46647	NON-STERILE
217.	8699131617861	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 9 HOLES	ISO 5832-1	140.209	46647	NON-STERILE
218.	8699131605646	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 10 HOLES	ISO 5832-1	140.210	46647	NON-STERILE
219.	8680941228828	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 11 HOLES	ISO 5832-1	140.211	46647	NON-STERILE
220.	8699131605653	SMALL LOCKING COMPRESSION PLATE 3.5 mm , 12 HOLES	ISO 5832-1	140.212	46647	NON-STERILE
221.	8680941228835	LOCKING OLECRANON PLATE 3.5 / 4.00 mm , 3 HOLES LEFT	ISO 5832-1	865.003	46647	NON-STERILE
222.	8699131605660	LOCKING OLECRANON PLATE 3.5 / 4.00 mm , 5 HOLES LEFT	ISO 5832-1	865.005	46647	NON-STERILE
223.	8680941292607	LOCKING OLECRANON PLATE 3.5 / 4.00 mm , 7 HOLES LEFT	ISO 5832-1	865.007	46647	NON-STERILE
224.	8699131605677	LOCKING OLECRANON PLATE 3.5 / 4.00 mm , 9 HOLES LEFT	ISO 5832-1	865.009	46647	NON-STERILE
225.	8699131605691	LOCKING OLECRANON PLATE 3.5 / 4.00 mm , 3 HOLES RIGHT	ISO 5832-1	865.103	46647	NON-STERILE
226.	8699131605714	LOCKING OLECRANON PLATE 3.5 / 4.00 mm , 5 HOLES RIGHT	ISO 5832-1	865.105	46647	NON-STERILE
227.	8699131605738	LOCKING OLECRANON PLATE 3.5 / 4.00 mm , 7 HOLES RIGHT	ISO 5832-1	865.107	46647	NON-STERILE
228.	8699131605752	LOCKING OLECRANON PLATE 3.5 / 4.00 mm , 9 HOLES RIGHT	ISO 5832-1	865.109	46647	NON-STERILE
229.	8680941292614	PROXIMAL HUMERAL MULTI ANGLE LOCKING PLATE 3.5 / 4.00 mm -LONG , 5 HOLES	ISO 5832-1	861.205	46647	NON-STERILE
230.	8680941239503	PROXIMAL HUMERAL MULTI ANGLE LOCKING PLATE 3.5 / 4.00 mm -LONG , 6 HOLES	ISO 5832-1	861.206	46647	NON-STERILE
231.	8680941239510	PROXIMAL HUMERAL MULTI ANGLE LOCKING PLATE 3.5 / 4.00 mm -LONG , 8 HOLES	ISO 5832-1	861.208	46647	NON-STERILE
232.	8680941292621	PROXIMAL HUMERAL MULTI ANGLE LOCKING PLATE 3.5 / 4.00 mm -LONG , 10 HOLES	ISO 5832-1	861.210	46647	NON-STERILE
233.	8699131605707	PROXIMAL HUMERAL MULTI ANGLE LOCKING PLATE 3.5 / 4.00 mm -LONG , 12 HOLES	ISO 5832-1	861.212	46647	NON-STERILE
234.	8699131605721	LOCKING ANTERIOR CLAVICLE PLATE S SHAPED 3.5 mm , 6 HOLES LEFT	ISO 5832-1	877.006	46647	NON-STERILE
235.	8699131605745	LOCKING ANTERIOR CLAVICLE PLATE S SHAPED 3.5 mm , 7 HOLES LEFT	ISO 5832-1	877.007	46647	NON-STERILE
236.	8699131605769	LOCKING ANTERIOR CLAVICLE PLATE S SHAPED 3.5 mm , 8 HOLES LEFT	ISO 5832-1	877.008	46647	NON-STERILE
237.	8680941292638	LOCKING ANTERIOR CLAVICLE PLATE S SHAPED 3.5 mm , 6 HOLES RIGHT	ISO 5832-1	877.106	46647	NON-STERILE
238.	8680941239527	LOCKING ANTERIOR CLAVICLE PLATE S SHAPED 3.5 mm , 7 HOLES RIGHT	ISO 5832-1	877.107	46647	NON-STERILE
239.	8680941239534	LOCKING ANTERIOR CLAVICLE PLATE S SHAPED 3.5 mm , 8 HOLES RIGHT	ISO 5832-1	877.108	46647	NON-STERILE
240.	8680941292645	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 5 HOLES LEFT	ISO 5832-1	853.215	46647	NON-STERILE
241.	8699131605899	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 7 HOLES LEFT	ISO 5832-1	853.218	46647	NON-STERILE
242.	8699131605912	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 9 HOLES LEFT	ISO 5832-1	853.219	46647	NON-STERILE
243.	8699131605929	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 11 HOLES LEFT	ISO 5832-1	853.221	46647	NON-STERILE
244.	8699131637296	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 13 HOLES LEFT	ISO 5832-1	853.223	46647	NON-STERILE
245.	8699131637302	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 5 HOLES RIGHT	ISO 5832-1	853.205	46647	NON-STERILE
246.	8699131637319	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 7 HOLES RIGHT	ISO 5832-1	853.207	46647	NON-STERILE
247.	8699131637326	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 9 HOLES RIGHT	ISO 5832-1	853.209	46647	NON-STERILE
248.	8699131664896	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 11 HOLES RIGHT	ISO 5832-1	853.211	46647	NON-STERILE

249.	8699131664926	LOCKING DISTAL FEMUR PLATE 4.5 / 5.0 mm , 13 HOLES RIGHT	ISO 5832-1	853.213	46647	NON-STERILE
250.	8699131664940	LOCKING CALCANEAL PLATE 3.5 mm , LEFT	ISO 5832-1	137.265	46647	NON-STERILE
251.	8699131664971	LOCKING CALCANEAL PLATE 3.5 mm , RIGHT	ISO 5832-1	137.260	46647	NON-STERILE
252.	8699131664988	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 4 HOLES LEFT	ISO 5832-1	160.304	46647	NON-STERILE
253.	8699131665008	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 5 HOLES LEFT	ISO 5832-1	160.305	46647	NON-STERILE
254.	8699131665039	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 6 HOLES LEFT	ISO 5832-1	160.306	46647	NON-STERILE
255.	8699131665046	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 8 HOLES LEFT	ISO 5832-1	160.308	46647	NON-STERILE
256.	8699131665060	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 10 HOLES LEFT	ISO 5832-1	160.310	46647	NON-STERILE
257.	8699131665084	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 12 HOLES LEFT	ISO 5832-1	160.312	46647	NON-STERILE
258.	8699131665091	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 4 HOLES RIGHT	ISO 5832-1	160.314	46647	NON-STERILE
259.	8699131665114	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 5 HOLES RIGHT	ISO 5832-1	160.315	46647	NON-STERILE
260.	8699131605943	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 6 HOLES RIGHT	ISO 5832-1	160.316	46647	NON-STERILE
261.	8699131605967	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 8 HOLES RIGHT	ISO 5832-1	160.318	46647	NON-STERILE
262.	8699131605981	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 10 HOLES RIGHT	ISO 5832-1	160.320	46647	NON-STERILE
263.	8699131606001	LOCKING PROXIMAL MEDIAL TIBIA PLATE 4.5 / 5.00 mm , 12 HOLES RIGHT	ISO 5832-1	160.322	46647	NON-STERILE
264.	8699131606025	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 4 HOLES LEFT	ISO 5832-1	854.004	46647	NON-STERILE
265.	8699131605950	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 6 HOLES LEFT	ISO 5832-1	854.006	46647	NON-STERILE
266.	8699131605974	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 8 HOLES LEFT	ISO 5832-1	854.008	46647	NON-STERILE
267.	8699131605998	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 10 HOLES LEFT	ISO 5832-1	854.010	46647	NON-STERILE
268.	8699131606018	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 12 HOLES LEFT	ISO 5832-1	854.012	46647	NON-STERILE
269.	8699131606032	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 14 HOLES LEFT	ISO 5832-1	854.014	46647	NON-STERILE
270.	8699131606049	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 4 HOLES RIGHT	ISO 5832-1	854.104	46647	NON-STERILE
271.	8699131606063	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 6 HOLES RIGHT	ISO 5832-1	854.106	46647	NON-STERILE
272.	8699131606087	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 8 HOLES RIGHT	ISO 5832-1	854.108	46647	NON-STERILE
273.	8699131606100	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 10 HOLES RIGHT	ISO 5832-1	854.110	46647	NON-STERILE
274.	8699131606124	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 12 HOLES RIGHT	ISO 5832-1	854.112	46647	NON-STERILE
275.	8699131606056	LOCKING DISTAL MEDIAL TIBIAL PLATE 3.5 / 4.00 mm , 14 HOLES RIGHT	ISO 5832-1	854.114	46647	NON-STERILE
276.	8699131606070	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 5 HOLES LEFT	ISO 5832-1	878.005	46647	NON-STERILE
277.	8699131606094	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 7 HOLES LEFT	ISO 5832-1	878.007	46647	NON-STERILE
278.	8699131606117	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 9 HOLES LEFT	ISO 5832-1	878.009	46647	NON-STERILE
279.	8699131606131	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 11 HOLES LEFT	ISO 5832-1	878.011	46647	NON-STERILE
280.	8699131606247	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 13 HOLES LEFT	ISO 5832-1	878.013	46647	NON-STERILE
281.	8699131606261	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 5 HOLES RIGHT	ISO 5832-1	878.105	46647	NON-STERILE
282.	8699131606285	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 7 HOLES RIGHT	ISO 5832-1	878.107	46647	NON-STERILE
283.	8699131606308	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 9 HOLES RIGHT	ISO 5832-1	878.109	46647	NON-STERILE

284.	8699131648209	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 11 HOLES RIGHT	ISO 5832-1	878.111	46647	NON-STERILE
285.	8699131606254	LOCKING PERIARTICULAR PROXIMAL FEMUR LATERAL PLATE 4.5 / 5.00 mm , 13 HOLES RIGHT	ISO 5832-1	878.113	46647	NON-STERILE
286.	8699131606278	CURVED RECONSTRUCTION PLATE 3 HOLE	ISO 5832-1	1932-0003	46647	NON-STERILE
287.	8699131606292	CURVED RECONSTRUCTION PLATE 4 HOLE	ISO 5832-1	1932-0004	46647	NON-STERILE
288.	8699131606315	CURVED RECONSTRUCTION PLATE 5 HOLE	ISO 5832-1	1932-0005	46647	NON-STERILE
289.	8699131648216	CURVED RECONSTRUCTION PLATE 6 HOLE	ISO 5832-1	1932-0006	46647	NON-STERILE
290.	8699131606346	CURVED RECONSTRUCTION PLATE 7 HOLE	ISO 5832-1	1932-0007	46647	NON-STERILE
291.	8699131606360	CURVED RECONSTRUCTION PLATE 8 HOLE	ISO 5832-1	1932-0008	46647	NON-STERILE
292.	8699131606384	CURVED RECONSTRUCTION PLATE 9 HOLE	ISO 5832-1	1932-0009	46647	NON-STERILE
293.	8699131606407	CURVED RECONSTRUCTION PLATE 10 HOLE	ISO 5832-1	1932-0010	46647	NON-STERILE
294.	8699131606421	CURVED RECONSTRUCTION PLATE 11 HOLE	ISO 5832-1	1932-0011	46647	NON-STERILE
295.	8699131606353	CURVED RECONSTRUCTION PLATE 12 HOLE	ISO 5832-1	1932-0012	46647	NON-STERILE
296.	8699131606377	CURVED RECONSTRUCTION PLATE 13 HOLE	ISO 5832-1	1932-0013	46647	NON-STERILE
297.	8699131606391	CURVED RECONSTRUCTION PLATE 14 HOLE	ISO 5832-1	1932-0014	46647	NON-STERILE
298.	8699131606414	CURVED RECONSTRUCTION PLATE 15 HOLE	ISO 5832-1	1932-0015	46647	NON-STERILE
299.	8699131606438	CURVED RECONSTRUCTION PLATE 16 HOLE	ISO 5832-1	1932-0016	46647	NON-STERILE
300.	8699131606483	CURVED RECONSTRUCTION PLATE 18 HOLE	ISO 5832-1	1932-0018	46647	NON-STERILE
301.	8699131606506	CURVED RECONSTRUCTION PLATE 20 HOLE	ISO 5832-1	1932-0020	46647	NON-STERILE
302.	8699131606520	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 14 mm	ISO 5832-1	106.314	56642	NON-STERILE
303.	8699131606544	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 16 mm	ISO 5832-1	106.316	56642	NON-STERILE
304.	8699131606568	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 18 mm	ISO 5832-1	106.318	56642	NON-STERILE
305.	8699131606490	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 20 mm	ISO 5832-1	106.320	56642	NON-STERILE
306.	8699131606513	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 22 mm	ISO 5832-1	106.322	56642	NON-STERILE
307.	8699131606537	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 24 mm	ISO 5832-1	106.324	56642	NON-STERILE
308.	8699131606551	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 26 mm	ISO 5832-1	106.326	56642	NON-STERILE
309.	8699131606575	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 28 mm	ISO 5832-1	106.328	56642	NON-STERILE
310.	8699131636725	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 30 mm	ISO 5832-1	106.330	56642	NON-STERILE
311.	8699131636732	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 32 mm	ISO 5832-1	106.332	56642	NON-STERILE
312.	8699131636756	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 34 mm	ISO 5832-1	106.334	56642	NON-STERILE
313.	8699131663264	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 36 mm	ISO 5832-1	106.336	56642	NON-STERILE
314.	8699131663271	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 38 mm	ISO 5832-1	106.338	56642	NON-STERILE
315.	8699131663288	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 40 mm	ISO 5832-1	106.340	56642	NON-STERILE
316.	8699131606643	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 42 mm	ISO 5832-1	106.342	56642	NON-STERILE
317.	8699131606667	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 44 mm	ISO 5832-1	106.344	56642	NON-STERILE
318.	8699131606650	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 46 mm	ISO 5832-1	106.346	56642	NON-STERILE
319.	8699131606674	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 48 mm	ISO 5832-1	106.348	56642	NON-STERILE
320.	8699131626870	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 50 mm	ISO 5832-1	106.350	56642	NON-STERILE
321.	8699131606681	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 55 mm	ISO 5832-1	106.355	56642	NON-STERILE
322.	8699131606704	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 60 mm	ISO 5832-1	106.360	56642	NON-STERILE
323.	8699131626887	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 65 mm	ISO 5832-1	106.365	56642	NON-STERILE

324.	8699131606698	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 70 mm	ISO 5832-1	106.370	56642	NON-STERILE
325.	8699131606711	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 75 mm	ISO 5832-1	106.375	56642	NON-STERILE
326.	8699131606728	5.00 mm LOCKING HEAD SCREW – SELF TAPPING (HEXDRIVE) 80 mm	ISO 5832-1	106.380	56642	NON-STERILE
327.	8699131606742	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 50 mm	ISO 5832-1	108.550	56642	NON-STERILE
328.	8699131606766	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 55 mm	ISO 5832-1	108.555	56642	NON-STERILE
329.	8699131606780	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 60 mm	ISO 5832-1	108.560	56642	NON-STERILE
330.	8699131606735	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 65 mm	ISO 5832-1	108.565	56642	NON-STERILE
331.	8699131606759	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 70 mm	ISO 5832-1	108.570	56642	NON-STERILE
332.	8699131606773	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 75 mm	ISO 5832-1	108.575	56642	NON-STERILE
333.	8699131606797	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 80 mm	ISO 5832-1	108.580	56642	NON-STERILE
334.	8699131673454	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 85 mm	ISO 5832-1	108.585	56642	NON-STERILE
335.	8699131643501	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 90 mm	ISO 5832-1	108.590	56642	NON-STERILE
336.	8699131643525	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 95 mm	ISO 5832-1	108.595	56642	NON-STERILE
337.	8699131643549	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 100 mm	ISO 5832-1	108.600	56642	NON-STERILE
338.	8699131643563	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 105 mm	ISO 5832-1	108.605	56642	NON-STERILE
339.	8699131673447	5.00 mm LOCKING HEAD CANCELLOUS SCREW 16 mm THREAD – SELF TAPPING HEX DRIVE – 110 mm	ISO 5832-1	108.310	56642	NON-STERILE
340.	8699131643594	LOCKING PROXIMAL TIBIA RAFT PLATE , 3 HOLES LEFT	ISO 5832-1	863.003	46647	NON-STERILE
341.	8699131643617	LOCKING PROXIMAL TIBIA RAFT PLATE , 5 HOLES LEFT	ISO 5832-1	863.005	46647	NON-STERILE
342.	8699131643631	LOCKING PROXIMAL TIBIA RAFT PLATE , 7 HOLES LEFT	ISO 5832-1	863.007	46647	NON-STERILE
343.	8699131643655	LOCKING PROXIMAL TIBIA RAFT PLATE , 9 HOLES LEFT	ISO 5832-1	863.009	46647	NON-STERILE
344.	8699131663622	LOCKING PROXIMAL TIBIA RAFT PLATE , 11 HOLES LEFT	ISO 5832-1	863.011	46647	NON-STERILE
345.	8699131607138	LOCKING PROXIMAL TIBIA RAFT PLATE , 3 HOLES RIGHT	ISO 5832-1	863.103	46647	NON-STERILE
346.	8699131633502	LOCKING PROXIMAL TIBIA RAFT PLATE , 5 HOLES RIGHT	ISO 5832-1	863.105	46647	NON-STERILE
347.	8699131633519	LOCKING PROXIMAL TIBIA RAFT PLATE , 7 HOLES RIGHT	ISO 5832-1	863.107	46647	NON-STERILE
348.	8699131633526	LOCKING PROXIMAL TIBIA RAFT PLATE , 9 HOLES RIGHT	ISO 5832-1	863.109	46647	NON-STERILE
349.	8699131607145	LOCKING PROXIMAL TIBIA RAFT PLATE , 11 HOLES RIGHT	ISO 5832-1	863.111	46647	NON-STERILE

PRODUCT CLASS

: CLASS IIb (According to Annex 9 Rule 8)

GMDN CODE

: 56642, 46647

CONF. ASSESSMENT ROUTE

: ANNEX-II.3

We here by declare that the above mentioned products meet the provisions of the Council Directive 93/42/EEC by annex directive 2007/47/EC for medical devices. All supporting documentation is retained under the premises of the manufacturer.

STANDARDS APPLIED

: EN ISO 13485:2016, TS EN ISO 9001:2015, EN ISO 14971:2012, EN ISO 62366:2015, EN ISO 11737-1 : 2006, EN ISO 14630:2010, EN ISO 14602:2010, EN ISO 10993-1:2009, ISO 5832-3:2016, ASTM F136, ISO 5832-2:2018, ASTM F67, EN ISO 15223-1:2016, EN 1041:2008, MEDDEV 2.7.1 REV4

NOTIFIED BODY

: **KIWA BELGELENDİRME HİZMETLERİ A.Ş., İTOSB 9.**
Cad. No:15 Tepeören Tuzla İstanbul Türkiye , NB1984.

EC CERTIFICATE NUMBER

: **1984-MDD-11-129**

START OF CE MARKING DATE

: **11 KASIM 2011**

PLACE, DATE OF ISSUE

: **17.04.2020, GAZİANTEP**

M.SERHAT DİRGAR
General Manager

MEDICAL DEVICE TIC. LTD. ŞTİ.
Aydınlı Mh. 03070 M. 03070 4 Şehitkamil/Gaziantep
V.D. ŞEHİTKAMİL 998 070 3416
TEL. +90 342 238 43 43 +90 342 238 44 11