

# VESSEL SEALING DEVICE COMPARISON

## Open Devices

MEDTRONIC

## Open Devices

COMPETITIVE



NAME	LigaSure Impact™ Curved, Large Jaw, Sealer/Divider Nano-Coated	LigaSure Atlas™ Tissue Fusion, Sealer/Divider	LigaSure™ Small Jaw, Sealer/Divider	LigaSure™ Exact Dissector, Sealer/Divider	LigaSure™ Maryland Jaw Sealer/Divider Nano-Coated	LigaSure™ Blunt Tip Sealer/Divider Nano-Coated	Sonicision™ Cordless Ultrasonic Dissector
PART NO.	LF4418	LS1020	LF1212A	LF2019	LF1923	LF1823	SCD13, SCD26
SEAL LENGTH	36 mm	22 mm	17.3 mm	20.6 mm	20.3 mm	18.3 mm	14.5 mm
CUT LENGTH	34 mm	20 mm	15.2 mm	19.8 mm	18.5 mm	16.8 mm	14.5 mm

	Caiman™ 12	Enseal™ G2 Super Jaw	Enseal™ X1 Large Jaw	Harmonic Focus™, Focus™ +, Focus™ Long Shears
SEAL LENGTH	50 mm	40 mm	38 mm	16 mm
CUT LENGTH	47 mm	37 mm	30 mm	16 mm

JAW PROFILE



**Medtronic**  
Further, Together

# VESSEL SEALING DEVICE COMPARISON

## Laparoscopic Devices

MEDTRONIC



NAME	LigaSure™ Retractable L-Hook Sealer/Divider	LigaSure™ Maryland Jaw Sealer/Divider Nano-Coated	LigaSure™ Maryland Jaw Thoracic Sealer/Divider Nano-Coated	LigaSure™ Blunt Tip Sealer/Divider Nano-Coated	LigaSure Atlas™ Tissue Fusion, Sealer/Divider	Sonicision™ Cordless Ultrasonic Dissector
PART NO.	LF5637, LF5644	LF1937, LF1944	LF1930T	LF1837, LF1844	LS1037	SCD396, SCD48
SEAL LENGTH	20 mm	20.3 mm	20.3 mm	18.3 mm	22 mm	14.5 mm
CUT LENGTH	18.5 mm	18.5 mm	18.5 mm	16.8 mm	20 mm	14.5 mm



Product images not actual size

**Resources:**

Medtronic product measurements taken from: [medtronic.com/covidien](http://medtronic.com/covidien).

Caiman™ measurements taken from: [https://www.aesculapusa.com/assets/base/doc/DOC1098\\_Rev\\_D\\_Caiman\\_Vessel\\_Sealers\\_Brochure.pdf](https://www.aesculapusa.com/assets/base/doc/DOC1098_Rev_D_Caiman_Vessel_Sealers_Brochure.pdf)

Harmonic™ Focus measurements taken from: [http://www.accessdata.fda.gov/cdrh\\_docs/pdf10/K100597.pdf](http://www.accessdata.fda.gov/cdrh_docs/pdf10/K100597.pdf)

Enseal™ G2 Super Jaw measurements taken from:

[http://gb.ethicon.com/sites/default/files/EES\\_ENSEAL\\_G2SuperJawTissueSealer\\_ComparativeSummary\\_BR544.pdf](http://gb.ethicon.com/sites/default/files/EES_ENSEAL_G2SuperJawTissueSealer_ComparativeSummary_BR544.pdf)

Medtronic device measurements based on internal report #R0050236\_A.

Enseal™ X1 Large Jaw measurements taken from [http://www.ethicon.com/sites/default/files/managed-documents/061968-161020\\_enseal\\_x1\\_productsellsheet\\_120916\\_v17.pdf](http://www.ethicon.com/sites/default/files/managed-documents/061968-161020_enseal_x1_productsellsheet_120916_v17.pdf), as well as internal data for the seal length measurements.

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. ™\* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 05/2021 - US-SE-2100124 - [WF# 2135923]

5920 Longbow Drive  
Boulder, CO 80301 800.722.8772 [medtronic.com/covidien](http://medtronic.com/covidien)

## Laparoscopic Devices

COMPETITIVE



NAME	Caiman™ 5	Enseal™ G2 Lap	Gyrus™ PKS 5 mm	Thunderbeat™	Harmonic™ Ace+
PART NO.					
SEAL LENGTH	26.5 mm	15.69 mm	13.7 mm	19.24 mm	13.52 mm
CUT LENGTH	23.5 mm	14.65 mm	11.8 mm	19.24 mm	13.52 mm

