

Multiband Ligator

Barrel with 6 bands

Used for endoscopic ligation of Esophageal varices at or above gastroesophageal junction

Fits on wide range of endoscopes
Second last Band light colour

Standard Band Ligator

requires 9.5 to 11.2 mm
endoscope diameter


Paediatric Band Ligator

requires 8.0 to 9.5 mm
endoscope diameter



Sclerotherapy Injection Needles

Used for endoscopic injection into esophageal varices above esophageal junction



Sharp-Bevel-Needle with JET-Pressure effect
metal protection

Especially developed needle to provide
secure firm grip with inner tube and prevent
possible leakage from the joint of inner tube
and needle.

All tubes are made of PTFE (Teflon)

| | Art. No. | Product name |
|------------|----------|--|
| Standard | SGK-6 | Set incl. Applicator barrel, thread puller and flushing needle |
| | SGB-6 | Barrel and thread puller only |
| Paediatric | SGK-6-S | Set incl. Applicator barrel, thread puller and flushing needle |
| | SGB-6-S | Barrel and thread puller only |

| Art. No. | Product name | Needle gauge | Needle length | Working length | Working Channel Ø |
|------------|----------------------|--------------|---------------|----------------|-------------------|
| SGN-2321 | Sclerotherapy Needle | 21G | 7 mm | 170 cm | 2.8 mm |
| SGN-2321-C | Sclerotherapy Needle | 21G | 7 mm | 230 cm | 2.8 mm |
| SGN-2323 | Sclerotherapy Needle | 23G | 7 mm | 170 cm | 2.8 mm |
| SGN-2323-C | Sclerotherapy Needle | 23G | 7 mm | 230 cm | 2.8 mm |
| SGN-2324 | Sclerotherapy Needle | 24G | 4 mm | 170 cm | 2.8 mm |
| SGN-2324-C | Sclerotherapy Needle | 24G | 4 mm | 230 cm | 2.8 mm |

C = Colono

Loop Basket

Used for endoscopic retrieval of foreign bodies from gastrointestinal tract and to retrieve polyps through colonoscope

Shape: oval with flush port
Sheath: PTFE



| Art. No. | Loop Size (mm) | Working Length (cm) | Working Channel Ø (mm) |
|---------------------|----------------|---------------------|------------------------|
| SFBH-LB | 60 | 170 | 2.8 |
| SFBH-LB-C | 60 | 230 | 2.8 |
| SFBH-LB-40 | 40 | 170 | 2.8 |
| SFBH-LB-C-40 | 40 | 230 | 2.8 |
| SFBH-LB-90 | 90 | 170 | 2.8 |
| SFBH-LB-C-90 | 90 | 230 | 2.8 |



Flat Basket available

Class IIa device – non sterile – for single use only

Polypectomy Snares

Used for endoscopic excision of polyps using monopolar electrosurgical unit and retrieval of foreign bodies from gastrointestinal tract

Sheath: PTFE
Loop Diameter: 6/13/15/20/25/30/35/40/45/50/55mm
Loop Diameter for Duckbill Snare is 25mm only



| Art. No. | Loop Shape | | Working Length (cm) | Working Channel Ø (mm) |
|--------------------|--------------|-----------|------------------------|---------------------------|
| SPSO | Oval | | 220 | 2.8 |
| SPSH | Hexagonal | | 220 | 2.8 |
| SPSC | Crescent | | 220 | 2.8 |
| SPSD-25 | Duckbill | | 220 | 2.8 |
| SPSM | Monofilament | | 220 | 2.8 |
| | | | | |
| SPSO-R | Oval | Rotatable | 220 | 2.8 |
| SPSH-R | Hexagonal | Rotatable | 220 | 2.8 |
| SPSC-R | Crescent | Rotatable | 220 | 2.8 |
| SPSM-R | Monofilament | Rotatable | 220 | 2.8 |
| SPSO-30-DBE | Oval | Rotatable | 280 | 2.8 |

Class IIb device – sterile – for single use only